Product Introduction


Thank you for choosing YUNMAI 2 Smart Scale



- A ITO conductive coating x 4
- B Concealed Display Screen
- C Reset Kev
- D USB Type-C charge port
- E Speaker

1 Activation

Place YUNMAI 2 on a flat floor and press the surface until the screen shows "YUNMAI" to activate it.



2 App Installation

Search "YUNMAI" in the iOS App Store or Android Application Market or scan the following QR Code to download, install it and it will get connected to the scale once turned on.



Scan QR Code Download <YUNMAI App>

3 Get Connected

Please lay the scale at a place with strong WiFi signal or near the router when it first get connected. WiFi signal will influence connecting time and stability.



Place in the scope of WiFi in first connection

3-1 Connect to the scale via Bluetooth

Open YUNMAI App and Bluetooth and it will get connected automatically, and your real-time data will sync to the App.



iOS: Slide from the bottom of the screen to open the Control Center and click on the [Bluetooth] Logo



Android: Click on [Allow] to turn on Bluetooth when the system reminds you when open YUNMAI App

3-2 Connect via WiFi

Under WiFi mode, as long as the scale and the phone are placed in the scope of the WiFi, then the data will be synchronized automatically to the App even without opening it. Open WiFi option on the phone and YUNMAI App and get connected as the following instructions:







2 Click on [Device]



(3) Click on [Device Connection]



4 Finish the process as indicated

4 Start Using

Set each foot on each side of the central point naturally as shown below:





Feet go out of the verge may affect the accuracy



Stand in wrong direction will shield the screen



Feet not in a line may affect the accuracy

Note: To ensure the accuracy please stand on the scale barefoot

YUNMAI2

Display and Status



WiFi Connection



WiFi Signal Error

Low Battery Power





Alarm Clock



Full Battery Power

Weight gained



Weight lost

Upgrading

Body Fat



lb (Unit)







Body Statistics Table

Male	Age	Body Fat	Muscle Ratio	BMI
	18–39	0-9.9% Low 10-20.9% Normal 21-25.9% Slightly obese 26-60% Obese	<51% Low >51% Normal	25-29.9 Overweight
	40-59	0-10.9% Low 11-21.9% Normal 22-26.9% Slightly obese 27-60% Obese		
	60-120	0-12.9% Low 13-23.9% Normal 24-28.9% Slightly obese 29-60% Obese		

	Age	Body Fat	Muscle Ratio	BMI
Female	18–39	0-9.9% Low 10-33.9% Normal 34-38.9% Slightly obese 39-60% Obese	<46% Low >46% Normal	
	40-59	0-20.9% Low 21-34.9% Normal 35-39.9% Slightly obese 40-60% Obese		< 18.5 Underweight 18.5–24.9 Normal 25–29.9 Overweight 30–34.9 Obese 3 5 Grossly Overweight
	60-120	0-21.9% Low 22-35.9% Normal 36-40.9% Slightly obese 41-60% Obese		

Note: Juveniles have rapid physical developments and their body statistics, such as body fat, BMI and muscle rate may fluctuate largely, so it is difficult to get an accurate result except the weight.

Mobile Device Support

iOS 8.0 and above

iPhone 4S/5/5S/6/6 Plus/SE/6S/6S Plus/7/7 Plus iPod touch /iPad Android 4.3 and above with Bluetooth 4.0

Important Note

Measurement results of YUNMAI 2 is only for reference, but professional medical advice. Measurement Range: 3-150 Kg.

Safe Power Charging

Do not use low-quality plugs and cables. Refrain from long-time charging in case of short circuit.

Do not charge in damp environment in case of electric

Maintenance

Use a wet cloth to clean the scale. Do not use any detergent and don't use water to wash it directly.

Warranty

YUNMAI 2 Smart Scale guarantees a year of technical support and hardware warranty.

Product Specification

Product Name: YUNMAI 2 Smart Scale Measurement Unit: Kg /Lb Connection: Bluetooth 4.0, WiFi Size: $300 \times 300 \times 23 \text{ mm}$ Executive Standard: Q / HC 1 - 2015 Net Weight: 1.9 Kg Weighing Range: 3 - 180 Kg Accuracy: ±50g Working Temperature: -5 °C - 40 °C Exterior Material: Eco-friendly ABS + PC Battery Type : Li-Polymer Battery

FCC ID: 2ADEB-M1601

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- -- Reorient or relocate the receiving antenna.
- -- Increase the separation between the equipment and receiver.
- -- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- -- Consult the dealer or an experienced radio/TV technician for help. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
- (1)This device may not cause harmful interference. (2) This device must accept any interference received, including interference that may cause undesired operation.
- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

For more information please visit www.iyunmai.com Hotline: +86-755-8671-6230 Facebook: www.facebook.com/iyunma twitter: twitter.com/YunmaiGlobal Manufacturer: Shenzhen YUNMAI Technology Co., Ltd Address: Room 10D, 10F, 4C Building, Software Industry Base, Nanshan District, Shenzhen, China

⊕YUNMAI

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.