YUNMAI S

Smart Scale Instructions



Do not disassmeble, squeeze, puncture the product or put it into water, fire, or exposes it to temperature higher than 60°C (140°F)

Important Notice:

YUNMAI Smart Scale should not be used to diagnose or treat any medical situation.

Notice before Use

- * Please ingore the first weighing reading at your first time use or after move the scale.
- * After scale power on or reset, please wait for the LED display flash "0.0" for several seconds, and then start your weight-in.
- * After each move of the scale, please wait for 5 seconds for the sensor to self-recalibarate with the floor, and then weight for accurate reading.

Statement:

Following groups may not get accurate results:

- Juveniles
- Professional athletes
- Pregnant women
- Users with Cardiac Pacemaker device

Weighing Notice


- * Do not weight when charging battery.
- * The scale surface may become slippery in humid places, please watch your steps and don't slip during weighing.
- * "weighing before meal or before bath" is the most recommend time for weighing, because eating and bathing may affect the body moisture and temperature, which might affect the sensor accuracy.
- * Recommend to weigh yourself at the same time and same condition in each day for more accurate data record.

Maintenance

- * Use a damp cloth to wipe the surface of the scale
- * Do not use any detergent to clean the scale, because detergent might damage the electrodes which may affect the accuracy of detecting body fat data.

Products Instruction and Installation

Thank you for choosing YUNMAI Smart Scale.



- (1) Note: Charging port. Plug the charge cable into the charging port to recharge the battery.
- (2) Note: Power Switch. Press the button down to power on.

1. Install the YUNMAI APP

Download method

Search "YUNMAI" in the Android Application Market or the iOS App Store or scan the following QR Code to download, install it and it will get connected to the scale once turned on.

App down loaded via other channels may be unusable. Before pairing with your equipment, make sure your equipment meets the following requirements:

- Compatible with iPhone/iPod touch with iOS 9.0 and above.
- · Android 4.3 and above with Bluetooth 4.0





Connect to the scale via Bluetooth

Open YUNMAI App and Bluetooth and it will get connected automatically, and your real-time data will sync to the app.



iOS: Slide from the bottom of the screen to open the Control Center and click on the [Bluetooth] Logo.




Android: Click on [Allow] to turn on Bluetooth when the system reminds you when open YUNMAI App.

2. Start Using


Correct Standing Posture

Make sure your feet placement and standing posture is correct as shown below:



To ensure the accuracy, please stand on the scale barefoot with the right posture.

Display



(1) Remarks : LED on: Charging LED flash: Low Power

Body Mass Index

BMI (for men and women)					
Age	Underweight	Normal	Overweight	Obese	Grossly Overweight
18–120	< 18.5	18.5–24.9	25–29.9	30-34.9	≥ 35

Body mass index (BMI) = weight (kg) ÷ height ^ 2 (m)

For example: BMI = $70\text{kg} \div (1.75\text{m} \times 1.75\text{m}) = 22.86$, belongs to the normal range.

Note: BMI Index is a standard that is widely used to measure a person's health condition.

3. YUNMAI APP

Download YUNMAI APP to use together with YUNMAI S scale, you can record, analyze and set goals for your personal health.

YUNMAI saves all your information to cloud and provides brief health reports and improvement recommendations.



The picture is for reference only. For more information, please download the latest version of YUNMAI App in App Store.

FCC Statement

FCC ID: 2ADEB-M1805

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference

in a residential installation. This equipment generates uses and can radiate

radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However,

there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Product Specifications

Product Name: YUNMAI S Smart Scale

Model: M1805

Measurement Unit: Kg/Lb Connection: Bluetooth 4.0 Size: 300 × 300 × 20 mm

Executive Standard: Q / YM 004 - 2017

Net Weight: 1.74 Kg

Weighing Range: 5 Kg-180 Kg

Accuracy: ± 50g

Working Temperature: 23 °F-104 °F (-5 °C - 40 °C)

Battery Capacity: 1000mAh

Input Voltage: DC 5V

Battery type: Lithium Battery Charging Input: DC 5V 0.5A

Input Current: 0.5A

Material: Tempered Glass (Panel)+ ABS&PS (Bottom Cover)

For more information please visit www.iyunmai.com

Hotline: +86-755-8671-6230

Facebook: www.facebook.com/iyunmai Twitter: twitter.com/iyunmaiGlobal

Manufacturer: ShenZhen YUNMAI Technology Co., Ltd

Address: Room 2201, Block B, Sunshine TDI Building, Nanshan District,

Shenzhen 518000, China