

#whdentalwerk  
f @ in >  
iodent.wh.com

Now at  
your dealer  
or **wh.com**



## Quick Installation Guide – ioDent® system



1

Connect the ioDent® gateway mini to power supply



2

Plug in the ioDent® Wi-Fi dongle into the device



3

Connect your device to power supply



4

Follow the link sent to your e-mail address or follow-up ioDent.com



Enter e-mail address and password

5

Run the initial set up on the ioDent.com platform



6

You have successfully connected your device

W&H Dentalwerk Bürmoos GmbH  
Ignaz-Glaser-Straße 53, Postfach 1  
5111 Bürmoos, Austria

t +43 6274 6236-0  
f +43 6274 6236-55  
office@wh.com **wh.com**



51026 AEN Rev. 000 / 26.03.2019  
Subject to alterations