# INVACARE® TDX®SI GROUP 2 POWER WHEELCHAIRS





The Invacare® TDX® SI-2 and Invacare® TDX® SI-HD Power Wheelchairs are the first TDX Power Wheelchairs ever introduced for the Group 2 consumer. Available in both 300 and 450 lb. weight capacities, the TDX SI Wheelchair delivers legendary TDX power and performance to meet the driving needs of almost any consumer. Center-wheel drive and a narrow footprint provide exceptional indoor maneuverability, while its brightly painted frame and integrated rear suspension combine for a simple yet elegant design. The end result is an eye-catching power wheelchair with a go-anywhere attitude.

## Invacare® TDX® SI-2 Wheelchair Features

Available with captain's seat, captain's seat with solid seat pan or rehab seating

Captain's seat available in 16", 18", 20" and 22" widths Invacare® MK5 75 amp controller coupled with Invacare® SPJ™+ joystick deliver optimal power to the motors Invacare® SureStep® Suspension with Stability Assist climbs 2" with ease while maintaining six wheels on the ground Integrated rear suspension provides a non-jarring, smooth ride when tackling uneven terrain or thresholds

Removable battery access door with slide-out battery tray makes servicing a snap

Advanced 2-pole motors provide high torque and a 4.9 mph top speed, while 14" drive wheels provide ideal traction Center-wheel drive and narrow base combine for tight turning radius and easy access through doors and hallways Available in Electric Blue and Electric Red Transport tie-down option\*

# HCPCS Code: K0822, K0823







- 1. Invacare® SureStep® Suspension with Stability Assist
- 2. Integrated rear suspension
- 3. Superior serviceability



\*TRBKTS has been crash-tested in accordance with ANSI/RESNA WC Vol 1. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving whicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

The TDX SI-2 Wheelchair for Group 2 is available in both Standard Configurations and Custom Configurations.

Standard Configurations: Pre-configured with the most popular options. Just select the seat size and frame color.

Custom Configurations: For those who need ASBA rehab seating, extra wide seating or other more customizable features.

Model #	Code	Seating	Configuration Type
TDXSI-2-S	K0822	Captain's w/ Solid Pan	Standard
TDXSI-2	K0822	Rehab Seat	Custom
TDXSIV-2-S	K0823	Captain's Seat	Standard
TDXSIV-2	K0823	Captain's Seat	Custom

## Invacare® TDX® SI-HD Wheelchair Features

Available with captain's seat, captain's seat with solid seat pan or rehab seating

Captain's seat available in 18", 20", 22" and 24" widths Invacare® MK6i™ 90 amp controller coupled with Invacare® SPJ™+ joystick delivers optimal power to the motors Invacare® SureStep® Suspension with Stability Assist climbs 2" with ease while maintaining six wheels on the ground Integrated rear suspension provides a non-jarring, smooth ride when tackling uneven terrain or thresholds

Removable battery access door with slide-out battery tray makes servicing a snap

Advanced 4-pole motors provide high torque and a 5 mph top speed, while 14" drive wheels provide ideal traction Center-wheel drive and narrow base combine for tight turning radius and easy access through doors and hallways Available in Electric Blue and Electric Red Transport tie-down option\*

# HCPCS Code: K0824, K0825







- 1. Invacare® SureStep® Suspension with Stability Assist
- 2. Integrated rear suspension
- 3. Superior serviceability



\*TRBKTS has been crash-tested in accordance with ANSI/RESNA WC Vol 1. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

The TDX SI-HD Wheelchair for Group 2 is available in both Standard Configurations and Custom Configurations.

Standard Configurations: Pre-configured with the most popular options. Just select the seat size and frame color.

Custom Configurations: For those who need ASBA rehab seating, extra wide seating or other more customizable features.

Model #	Code	Seating	Configuration Type
TDXSI-HD-S	K0824	Captain's w/ Solid Pan	Standard
TDXSI-HD	K0824	Rehab Seat	Custom
TDXSIV-HD-S	K0825	Captain's Seat	Standard
TDXSIV-HD	K0825	Captain's Seat	Custom















口	
Op.	

Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
TDXSI-2, TDXSI-2-S, TDXSIV-2, TDXSIV-2-S	38.5" L* 24" W 38" H**	16"-22" W 16"-22" D	16.5"-20.5"	14"	3"	24"	9°
TDXSI-HD, TDXSI-HD-S TDXSIV-HD, TDXSIV-HD-S	35.2" L 26.2" W 37" H	18"-24" W 16"-22" D	16.5"-20.5"	14"	3"	20.5"	6°

	$\overline{\mathbf{C}}$
$\angle$	lb.\













3 Years Seat Frame 1 Year Electronics/

Motor/Gearbox

				006			
	Total Weight	Weight Capacity	Motor	Max Speed	Battery	Battery Range	Limited Warranty
TDXSI-2, TDXSI-2-S, TDXSIV-2, TDXSIV-2-S	254 lb.	300 lb.	2-Pole	4.9 mph	22NF	Up To 17 Miles	5 Years Base Frame 3 Years Seat Frame 1 Year Electronics/ Motor/Gearbox
TDXSI-HD, TDXSI-HD-S TDXSIV-HD,	254 lb.	450 lb.	4-Pole HD	4.8 mph	22NF	Up To 12.57 Miles	5 Years Base Frame

\*With front riggings

TDXSIV-HD-S

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# **Invacare Corporation**

www.invacare.com

# USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

## Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

©2012 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols <sup>™</sup> and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 08-060 Rev. 09/12 120359



 $<sup>{\</sup>tt **Captain's \ seat \ without \ headrest, \ for \ other \ seating \ option \ specifications \ refer \ to \ the \ Invacare \ Visual \ Configurator}$ Do not operate on roads, streets or highways