3M Health Care Solutions for Vascular Access Catalogue





Health care professionals around the globe look to 3M for innovative solutions that help to reduce the risk of infection, provide better patient outcomes and reduce health care costs.

Our trusted brands continue to set the industry standard and offer the finest in skin health and infection prevention products, supported by educational services, professional and technical support, and ongoing research and development.

What sets us apart from our competition is our ability to respond to customers' needs in a variety of ways.

- Experienced sales team, including specially-trained representatives who provide in-service training, continuing education and conversion programs
- Clinical specialists provide support using a variety of methods:
- Professional continuing education programs
- Protocol and formulary development to promote clinical outcomes, standardization and utilization within your departments
- Solutions-based website for 24 hours-per-day, 7 days-per-week access to information on a wide range of skin and wound clinical conditions. Visit www.3M.com/skinhealth
- Toll-free helpline that offers experts to assist customers with technical and product information
- Customer service representatives who are available to handle daily calls for order entry, and special shipment requests to ensure customer satisfaction

Table of Contents

Hand Hygiene

3M™ Quick Foam Hand Antiseptic	4
3M™ Hand Antiseptic Gel	
3M™ Avagard™ D Instant Hand Antiseptic	6
3M [™] Avagard [™] CHG Surgical Scrub	7

Clip

3M™ Surgical Clippers

Prep

Solu I.V.™ Maxi Swabsticks	10
Solu I.V.™ Swabsticks	11
Solu I.V.™ Swabs	12

Drape

```
3M<sup>™</sup> Surgical Drapes.....
```

Protect

3M[™] Cavilon[™] No Sting Barrier Film......14

Secure

3M [™] Tegaderm [™] CHG I.V. Securement Dressings	. 15
3M™ Tegaderm™ I.V. Transparent Film Dressings	
3M [™] Transparent Film Dressings	. 17
3M [™] Tegaderm [™] HP Transparent Film HP Dressings	.18
3M™ Transparent Film Roll	.19



.13

3M[™] Quick Foam Hand Antiseptic NPN 80008794

- 4 Non-aerosol foaming hand antiseptic containing 70% v/v ethyl alcohol with moisturizers which contain emollients to moisturize the skin.
 - Foam hand antiseptic kills harmful bacteria without soap or water
 - Dye and Fragrance free with no sticky residue after use.
 - Spreads easily and evenly while drying quickly
 - This product is easy to use for facilities as it is drip and splatter resistant and has versatile point-of-care and dispenser options

Suggested applications:

Apply to clean dry hands. Thoroughly wet hands and fingernails. Rub and let dry.

ORDERING INFORMATION:

Catalogue Number	Description	Items/Case
9320C	3M [™] Quick Foam Hand Antiseptic, 50 mL	25
9321	3M [™] Quick Foam Hand Antiseptic, 500 mL	12
9322	3M [™] Quick Foam Hand Antiseptic, 1 L refill	6
9330	Belt Clip for 3M $^{\rm \tiny M}$ Quick Foam Hand Antiseptic, 50 mL	500
9323	3M [™] Wall Dispenser Standard for 1 L	1
9328	$3M^{\mbox{\tiny TM}}$ Wall Dispenser UX10C Lockable for 1 L	1
9329	3M [™] Wall Dispenser UX10T Touchless for 1 L	1
9332	3M [™] Wire Wall Bracket for 500 mL	1

3M[™] Hand Antiseptic Gel NPN 80007158

Instant hand antiseptic with emollients containing 70% v/v ethyl alcohol which contain emollients to moisturize the skin

- Gel hand antiseptic kills harmful bacteria without soap or water
- Latex and fragrance free
- Spreads easily and evenly while drying quickly
- This product is easy to use for facilities as it is drip and splatter resistant and has versatile point-of-care and dispenser options

Suggested applications:

Apply to clean dry hands. Thoroughly wet hands and fingernails. Rub and let dry.

5

Catalogue Number	Description	Items/Case
BELT CLIP	Belt Clip for 3M [™] Hand Antiseptic Gel, 115 mL	EACH
9325	3M [™] Hand Antiseptic Gel, 115 mL	48
9326	3M [™] Hand Antiseptic Gel, 500 mL	12
9327	3M [™] Hand Antiseptic Gel, 1 L refill	6
9323	3M [™] Wall Dispenser Standard for 1 L	1
9328	$3M^{\scriptscriptstyle\rm TM}$ Wall Dispenser UX10C Lockable for 1 L	1
9329	3M [™] Wall Dispenser UX10T Touchless for 1 L	1
9332	3M [™] Wire Wall Bracket for 500 mL	1



7

3M[™] Avagard[™] D Instant Hand Antiseptic with Moisturizer NPN 02242852

Rapid bactericidal activity while helping to maintain the skin's integrity.

- Provides rapid broad spectrum kill while helping to maintain skin integrity and enhance skin hydration
- CHG Compatible

6

 Compatible with latex and non-latex gloves. Latex glove tested: Safeskin. Non-latex gloves tested: Triflex PVC, Elastylon

Suggested applications:

- Avagard D instant hand antiseptic can be used to disinfect hands when hands are not visibly soiled
- Apply to clean dry hands. Wet hands thoroughly with product and allow to dry

3M[™] Avagard[™] CHG Surgical Scrub (Chlorhexidine Gluconate 1% Solution and Ethyl Alcohol 61% w/w) Surgical and Healthcare Personnel Hand Antiseptic with Moisturizer

DIN 02246888

- Contains two active ingredients; alcohol for fast kill and chlorhexidine gluconate for persistent and cummulative activity
- Provides rapid antimicrobial kill of a broad spectrum of microorganisms, including antibiotic-resistant stains (*in vitro*)
- Provides proven, persistent antimicrobial protection for up to six hours without water or scrubbing (*in vitro*)
- Advanced, emollient-rich lotion base helps prevent dry, cracked hands, and maintains skin barrier integrity
- Compatible with latex and non-latex gloves. Latex gloves tested: Triflex, Encore. Non-latex gloves tested: Tactylon, Duraprene, Biogel

Suggested applications:

- Surgical hand antiseptic significantly reduces the number of microorganisms on the hands and forearms prior to surgery or patient care.
- Healthcare personnel hand antiseptic handwash to help reduce bacteria that potentially can cause disease - recommended for repeated use.

ORDERING INFORMATION:

Catalogue Number	Description	Items/Case
BELT CLIP	3M [™] Avagard [™] D Instant Hand Antiseptic, 88 mL	EACH
9221C	Belt Clip for $3M^{\mbox{\tiny TM}}$ Avagard $^{\mbox{\tiny M}}$ D Instant Hand Antiseptic, 88 mL	500
9222C	3M [™] Avagard [™] D Instant Hand Antiseptic, 500 mL	12
9226	3M [™] Avagard [™] D Molded Wall Bracket	

Catalogue Number	Description	Items/Case
9200C	Bottle, 16 fl. oz., 500 mL	1
9201C	$3M^{\mbox{\tiny TM}}$ Avagard $^{\mbox{\tiny TM}}$ Wall Bracket and Foot Pump	1
9203A	3M [™] Avagard [™] Replacement Foot Pump	1
9204C	3M [™] Avagard [™] Nail Cleaners	1



9

3M[™] Surgical Fixed Head Clipper

Compact, lightweight, rechargeable surgical clippers with fixed or pivoting blade assemblies for more efficient and safe hair removal.

- SHN 2007 recommends that patients should be advised to stop shaving, in the operative area, and one week before the surgery date.
- April 1999 CDC Guidelines recommend that if hair removal is required, the use of a surgical clipper immediately before an operation is one means to help reduce the risk of surgical site infections.
- Removes all types of head and body hair, wet or dry, to allow for standardization on one product for all hair removal needs.
- Unique design allows flexibility of clipping with either a pushing or a pulling motion.
- The cutting blade is on top of the stationary blade, away from the skin, so there's minimal risk of nicks or cuts.
- Single-use blades help reduce the risk of cross-contamination and are easy to snap on for use and pop off for disposal.
- Rechargeable clippers have a nickel metal hydride battery allowing multiple patients to be clipped from a single 12-14 hour recharge.
- Unique design is easy to clean and disinfect, wipe with a disinfectant or rinse under running water.

Suggested applications:

8

 Removing head and body hair from patients in preparation for procedures where hair removal is desired.

ORDERING INFORMATION:

Catalogue Number	Description	Items/Box	Items/Each
9670	3M [™] Single-use Blade Assembly (Black); for 3M [™] Surgical Clipper with Fixed Head - 9671	50	1
9671	3M [™] Surgical Clipper with Fixed Head - 9671; Battery or Cord Operated (charger/cord not included)	1	1
9672	3M [™] Clipper Charger/Cord; for 3M [™] Surgical Clipper with Fixed Head - 9671 (3.1 m cord)	1	1
9677	3M [™] Clipper Starter Kit; contains 3M [™] Surgical Clipper - 9671 and 3M [™] Clipper Charger/cord - 9672 (clipper blade assembly not included)	1	1
9669	3M TM Clipper Wall Bracket; for 3M TM Surgical Clipper with Fixed Head - 9671	1	1

3M[™] Surgical Pivoting Head Clipper

Now the 3MT[™] Surgical Clipper 9661 is available with a blade assembly to clip a full head of hair for neuro procedures. This combination of the Premium 40 assembly and pivoting head clipper allows for maximum blade contact at all times, no matter where you are clipping.

- SHN 2007 recommends that patients should be advised to stop shaving, in the operative area, and one week before the surgery date.
- Removes wet or dry hair
- Leaves minimal stubble length (0.5 0.7mm)
- Removes all types of body hair
- Individually sealed, single-use blades eliminate the need for cleaning the blade and reduce the risk of cross-contamination
- Blade assembly can be easily attached and removed with one hand for disposal.

Suggested applications:

 Removing head and body hair from patients in preparation for procedures where hair removal is desired.



Catalogue Number	Description	Items/Box	Items/Each
9600	3M [™] Single-Use Blade Assembly (purple); for 3M [™] Surgical Clipper with Pivoting Head - 9661	50	1
9661	3M [™] Surgical Clipper with Pivoting Head - 9661; Rechargeable	1	1
9662	3M [™] Clipper Charger; for 3M [™] Surgical Clipper with Pivoting Head - 9661	1	1
9667	3M [™] Clipper Starter Kit; contains 3M [™] Surgical Clipper - 9661 and 3M [™] Clipper Charger - 9662 (clipper blade assembly not included)	1	1

Solu I.V.[™] Maxi Swabsticks DIN 02240356

10 Meets the Safer Healthcare Now recommendation for central line insertion

- Contains 5 mls of 2% CHG and 70% alcohol
- Available in pink or clear solution



ORDERING INFORMATION:

Catalogue Number	Description	Items/Each
10208	Solu I.V. $^{\scriptscriptstyle\rm M}$ Maxi Swabstick 2% CHG and 70% alcohol - clear	30
10209	Solu I.V. $^{\scriptscriptstyle\rm M}$ Maxi Swabstick 2% CHG and 70% alcohol - tinted	30

Solu I.V.[™] Swabsticks DIN 02240356

IN 02240356

Designed for ease of use and specific to maintenance care of indwelling vascular access catheters

Solu I.V.[™] Swabsticks are available in three different concentrations:

- Solu I.V.[™] 0.5% CHG and 70% alcohol DIN 02060019
- Solu I.V.[™] 2% CHG and 70% alcohol DIN 02240356
- Solu I.V.[™] 2% CHG and alcohol free DIN 02242513

Catalogue Number	Description	Items/Each
10202	Solu I.V. $^{\mbox{\tiny TM}}$ Swabstick 0.5% CHG and 70% alcohol	50
10203	Solu I.V. [™] Swabstick 2% CHG and 70% alcohol	50
10207	Solu I.V. $^{\scriptscriptstyle\rm M}$ Swabstick 2% CHG and alcohol free	50





Solu I.V.™ Swabs DIN 02240356

12 Meets the Safer Healthcare Now recommendation for moving to CHG swabs for all applications including swabbing ports

Solu I.V. $\ensuremath{^{\scriptscriptstyle M}}$ Swabs are available in three different concentrations:

- Solu I.V.[™] 0.5% CHG and 70% alcohol DIN 02240356
- Solu I.V.[™] 2% CHG and 70% alcohol DIN 02240356
- Solu I.V.[™] 2% CHG and alcohol free DIN 02242513

ORDERING INFORMATION:

Catalogue Number	Description	Items/Each
10102	Solu I.V. $^{\scriptscriptstyle\rm TM}$ Swab 0.5% CHG and 70% alcohol - medium	200
10106	Solu I.V. $^{\scriptscriptstyle\rm TM}$ Swab 2% CHG and 70% alcohol - medium	100
10107	Solu I.V. $^{\scriptscriptstyle\rm TM}$ Swab 2% CHG and 70% alcohol - large	100
10108	Solu I.V. $^{\mbox{\tiny M}}$ Swab 2% CHG and alcohol free - large	100

3M[™] Steri-Drape[™] Angiography Drapes

Surgical drapes for venous access procedures, made of Absorbent Prevention Fabric*.

- Generous aperture area allows for proper positioning of the drape on the patient.
- Reliable, secure adhesive.
- 1192 drape has a large diamond shape aperture making landmarks easy to identify.

Suggested applications:

- Central and PICC Line Insertions.
- Brachial and Radial Artery Access.



ORDERING INFORMATION:

Catalogue Number			
1192	3M [™] Steri-Drape [™] Venous/Arterial Access Drape, absorbent prevention fabric, two adhesive patches with diamond shaped apertures	Overall size: 60 in. x 67 in. (152 cm x 170 cm) Aperture: 3 in x 8 in (7.6 cm x 20 cm)	12

*3M" Absorbent Fabric meets AAMI Standard PB70:2003 liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities and AAMI Technical Information Report TIR II:2005. Selection of surgical gowns and drapes in health care facilities.

$\mathbf{3}\mathbf{M}^{\text{\tiny TM}}$ Cavilon^{\text{\tiny TM}} No Sting Barrier Film

A versatile choice in your skin protection strategy.

- Forms a transparent protective film between the skin and potential irritants that cause breakdowns
- Sterile, non-cytotoxic, will not interfere with skin healing
- Alcohol-free, will not sting, even when applied to severely excoriated, reddened or denuded skin
- Compatible with CHG

Suggested applications:

14

- Prevents skin damage from tape and dressing adhesive
- Protects damaged skin from further breakdown

ORDERING INFORMATION:

Catalogue Number	Description	Size	ltems/ Each	Boxes/ Case
3343E	Cavilon [™] No Sting Barrier Film, Wand	1.0 mL wand (0.03 oz)	25	4
3344E	Cavilon [™] No Sting Barrier Film, Wipe	1.0 mL wipe (0.03 oz)	30	6
3345E	Cavilon [™] No Sting Barrier Film, Wand	3.0 mL wand (1 oz)	25	4



Combining the powerful antimicrobial activity of chlorhexidine gluconate with the simplicity of a Tegaderm[™] Dressing

- Proven to reduce skin flora under the dressing and inhibit re-growth at 7 days*
- Transparent window allows for continuous visualization of the insertion site
- Integrated design of dressing reduces application steps and minimizes potential for application errors

Suggested applications:

- Central Venous Catheters
- Peripherally Inserted Central Catheters
- Tunneled Central Venous Catheters
- Arterial Catheters
- Dialysis Catheters
- Midline Catheters
- Any percutaneous device

ORDERING INFORMATION:

Catalogue Number	Size	Dressings/ Box	Boxes/Case
1657R	3 1/2 in x 4 1/2 in (8.9 cm x 11.5 cm)	25	4
1659R	4 in x 6 1/8 in (10 cm x 15.5 cm)	25	4

*on healthy volunteers





3M[™] Tegaderm[™] I.V. Transparent Film Dressings

Reliable solutions for managing infusion therapy sites

- Notched design for a better seal and reduced edge-life around the catheter
- Sterile tape strips for anchoring hubs, lumens or tubing
- Soft cloth border for patient comfort and dressing durability
- Unique picture frame delivery

Suggested applications:

16

- Central Venous Catheters
- Peripherally Inserted Central Catheters
- Arterial Catheters
- Dialysis Catheters
- Midline Catheters
- Any percutaneous device

ORDERING INFORMATION:

Catalogue Number	Size	Dressings/Box	Boxes/Case
1610	2 in x 2 1/4 in (5 cm x 5.7 cm)	100	4
1633	2 3/4 in x 3 1/4 in (7 cm x 8.5 cm)	100	4
1635	3 1/2 in x 4 1/4 in (8.5 cm x 10.5 cm)	50	4
1650	4 in x 6 1/8 in (10 cm x 15.5 cm)	25	4
1655	3 1/2 in x 4 1/2 in (8.9 cm x 11.5 cm)	50	4



3M[™] Tegaderm[™] Transparent Film Dressing

A breathable dressing, impermeable to liquids, bacteria and viruses.*

- "Picture-frame" delivery—preferred by caregivers
- Secure adhesion, yet gentle to the skin
- Breathable film, allows oxygen in and moisture vapour out
- Barrier to external contaminants
- Conformable dressing flexes with skin for greater patient comfort



ORDERING INFORMATION:

Catalogue Number	Size	Dressings/Box	Boxes/Case
1622W	1 3/4 in x 1 3/4 in (4.4 cm x 4.4 cm)	100	4
1624W	2 3/8 in x 2 3/4 in (6 cm x 7 cm)	100	4
1626W	4 in x 4 3/4 in (10 cm x 12 cm)	50	4
1628	6 in x 8 in (15 cm x 20 cm)	10	8

*Laboratory testing has proven 3M[™] Tegaderm[™] and 3M[™] Tegaderm[™] HP Transparent Dressings provide a viral (HIV-1 and HBV) barrier while dressings remain intact without leakage.

3M[™] Tegaderm[™] HP Transparent Film Dressings

18 Dressings with quality, variety and versatility.

- Impermeable to liquids and bacteria
- Breathable, moisture vapour permeable
- Complete line of sizes and shapes with picture frame delivery

Suggested applications:

- Primary or secondary wounds
- Protection against skin breakdown
- Tegaderm[™] HP Film dressing has greater holding power even in the presence of moderate amounts of moisture
- Peripheral Catheters

ORDERING INFORMATION:

Catalogue Number	Size	Dressings/Box	Boxes/Case
9519HP	2 3/8 in x 2 3/8 in (6 cm x 6 cm)	100	4
9534HP	2 3/8 in x 2 3/4 in (6 cm x 7 cm)	100	4
9536HP	4 in x 4 3/4 in (10 cm x 12 cm)	50	4
9525HP	2 1/2in x 2 3/4in (6.5cm x 7cm)	100	4
9545HP (Oval)	1 7/8 in x 2 3/8 in (4.7 cm x 6 cm)	50	6
9546HP (Oval with Label)	4 in x 4 1/2 in (10 cm x 11.5 cm)	50	4
9548HP (Oval)	5 1/2 in x 6 1/2 in (14 cm x 16.5 cm)	10	8



$\mathbf{3M}^{\text{\tiny TM}}\,\mathbf{Tegaderm}^{\text{\tiny TM}}\,\mathbf{Transparent}\,\mathbf{Film}\,\mathbf{Roll}$

Tegaderm[™] Roll provides all the benefits of surgical tape, while providing a waterproof, bacterial and viral barrier.

- Waterproof backing provides a bacterial and viral barrier to outside contaminants and allows the patient to shower
- Can be cut and customized to desired lengthEasy-to-use delivery system
- Lasy-to-use derivery system
- Can be used for fixation of medical devices or infusion tubing
- Comformable fits body contours, permits free movement within stress on the skin

ORDERING INFORMATION:

Catalogue Number	Size	Rolls/Box	Boxes/Case
16002	2 in x 11 yd (5 cm x 10 m)	1	4
16004	4 in x 11 yd (10 cm x 10 m)	1	4
16006	6 in x 11 yd (15 cm x 10 m)	1	4



19

For Clinical and Sales information, contact your **3M Sales Representative. 1 800 410-6880** www.3M.com/ca/healthcare



Health Care 3M Realth Care 3M Canada Company 3M Centre, Building 275-4W-02 P.O. Box 5757 3M Centre, Building 275-4W-02 London, Ontario N6A 4T1 St. Paul, MN 55144-1000 Canada U.S. A. 3M Medica 1 800-364-3577 1-800-228-3957 D-41453 www.3m.ca www/M.com/healthcare Neuss, Germany

3M, Avagard, Cavilon, Tegaderm and Steri-Drape are trademarks are 3M. Used under license in Canada. Solu-I.V. is a trademark of Solumed, a 3M Canada Company. Please recycle. Printed in Canada. © 2010, 3M. All rights reserved. 0911-07919