



lively
wearable™

An introduction
to your new
Lively Wearable.



How it works

The Lively Wearable is a waterproof device you can wear on the GreatCall wristband, or around your neck when attached to our specially designed lanyard.

- Pair your Lively Wearable to your smartphone via Bluetooth.
- In any emergency, simply press the face of your Lively Wearable.
- The Wearable will connect you to 5Star service through the paired smartphone (within a 150 ft. range).
- 5Star Agents will then confirm your location, evaluate your situation and get you the help you need.



When to push the Lively Wearable Emergency Button

Potential emergencies:

- If you lose your balance and fall while you're out.
- If you mix up your medications and feel faint.
- If you slip in the shower.

Non-emergencies:

- If you lock yourself out of the house without your keys.
- If you're driving and get lost.
- If you feel nervous walking to your car.
- When in doubt, press the button and a highly trained 5Star Agent will be there to assist you right away.