

## ■ OmniPod™ System Federal Communications Commission (FCC) Notice Concerning Interference

---

The OmniPod System (both OmniPod and Personal Diabetes Manager or PDM) complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. These devices may not cause harmful interference.
2. These devices must accept any interference received, including interference that may cause undesirable operation.

Changes or modifications not expressly approved by Insulet could void the user's authority to operate the equipment.

Both OmniPod and PDM generate, use, and can radiate radio frequency energy, and may cause harmful interference to radio communications. There are no guarantees that interference will not occur in a particular installation. If the OmniPod System (both OmniPod and PDM) does cause harmful interference to radio and television reception, the interference may be corrected by one of the following measures:

- Move or relocate the OmniPod System (both OmniPod and PDM).
- Increase the distance between the OmniPod System (both OmniPod and PDM) and the other device that is emitting or receiving interference.

