



# Resuscitation supplies

**PHILIPS**

# Changes to this edition

April 2013

## New product added

### Multi Patient Disposable Cuffs – Multi Care

989803183311	Infant (orange)	41
989803183321	Pediatric (green)	41
989803183331	Small Adult (light blue)	41
989803183341	Adult (navy blue)	41
989803183351	Adult XL (navy blue with edge striping)	41
989803183361	Large Adult (burgundy)	41
989803183371	Thigh (grey)	41

### Single Care Disposable Cuffs

989803182281	Pediatric (green)	42
989803182291	Small Adult (light blue)	42
989803182301	Adult (navy blue)	42
989803182311	Adult XL (navy blue with edge striping)	42
989803182321	Large Adult (burgundy)	42

## Product revised

### Multifunction Electrode Pads

M3501A	Adult/Child with barrel connector	7
M3504A	Infant with barrel connector	7
M3713A	Adult/Child Plus with plug style connector	7
M3716A	Adult/Child Radiolucent with plug style connector	7
M3717A	Infant Plus with plug style connector	7
M3718A	Adult/Child Radiotransparent (Reduced Skin Irritation) with plug style connector	7
M3719A	Infant Radiotransparent (Reduced Skin Irritation) with plug style connector	7
989803166021	Adult/Child Preconnect	7

### Multi Patient Traditional Cuffs – Single Hose

40401A	Infant	36
40401B	Pediatric	36
40401C	Adult	36
40401D	Large Adult	36
40401E	Thigh	36

### Multi Patient Comfort Cuffs – Single Hose

M1572A	Pediatric, color coded: green tubing	37
M1573A	Small Adult, color coded: light blue tubing	37
M1574A	Adult, color coded: navy blue tubing	37
M1575A	Large Adult, color coded: burgundy tubing	37

**Multi Patient Easy Care Cuffs – Single Hose**

M4552B	Infant, orange	38
M4552B5	Infant, orange	38
M4553B	Pediatric, green	38
M4553B5	Pediatric, green	38
M4554B	Small Adult, light blue	39
M4554B5	Small Adult, light blue	39
M4555B	Adult, navy blue	39
M4555B5	Adult, navy blue	39
M4556B	Adult Extra-Long, navy blue	39
M4556B5	Adult Extra-Long, navy blue	39
M4557B	Large Adult, burgundy	39
M4557B5	Large Adult, burgundy	39
M4558B	Large Adult Extra-Long, burgundy	39
M4558B5	Large Adult Extra-Long, burgundy	40
M4559B	Thigh, grey	40
M4559B5	Thigh, grey	40

**Gentle Care Disposable Cuffs - Single Hose**

M4572B	Infant, orange	43
M4573B	Pediatric, green	43
M4574B	Small Adult, light blue	43
M4575B	Adult, navy blue	43
M4576B	Adult Extra-Long, navy blue	43
M4577B	Large Adult, burgundy	43
M4578B	Large Adult Extra-Long, burgundy	44
M4579B	Thigh, grey	44

**Blood Pressure Air Hoses**

M1598B	NIBP Air Hose	45
M1599B	NIBP Air Hose	45

**Product removed****Multifunction Electrode Pads**

M3502A	Adult Multifunction Pads (IEC – barrel connector)	8
M3503A	Infant/Child Multifunction Pads (IEC - barrel connector)	8

**Disposable Blood Pressure Kits**

M1567A	Single line blood pressure kit	26
M1568A	Dual line blood pressure kits	26
M1634A	Reusable adapter cable for M1567A and M1568A kits	26
M2271A	Holder for disposable pressure transducer	26
M2272A	IV Pole Mount	26

Dear Customer:

If you're looking for high-quality medical supplies to optimize the performance of your medical devices, you'll find them right here. Whether your focus is on quality patient care, reliability, or cost containment, you can turn to a single source – Philips Medical Supplies. Explore our complete family of supplies, from electrodes to SpO<sub>2</sub> sensors, for every area of care.

With Philips, you can buy your medical equipment and supplies from a single source, simply and easily. Use this guide to find the Philips supplies you need – or go to our website to purchase supplies: <http://www.philips.com/healthcarestore>.

## Using this guide

### Latex Cautions

Philips Medical Supplies currently distributes over 600 products, of which less than 9% contain natural rubber latex or dry natural rubber that could contact humans. Because these substances may cause allergic reactions, products in this guide that contain them are identified with the following symbol:

**LATEX**

**Caution:** Products that are identified with this symbol contain dry natural rubber or natural rubber latex that can come into contact with humans and may cause allergic reactions.

To the best of our knowledge, products not identified with this symbol are not known to contain natural rubber latex or dry natural rubber. The latex symbol indicates **only** which **medical supplies** contain natural rubber latex or dry natural rubber that can come into contact with humans and may cause allergic reactions. If you want to know whether a Philips **medical device** contains natural rubber latex or dry natural rubber, refer to the product labeling.

The FDA requires that products be labeled with human contacting latex content information when appropriate. Always refer to product labeling for accurate information on human contacting dry natural rubber or natural rubber latex.

### CE Marking

Products labeled with this symbol conform to the European Union (EU) directives as mandated by EU member state legislation and are available for sale in Europe to countries that accept CE marking with restrictions as may be noted. Medical devices labeled with this marking means the product is in conformity with the European

Medical Directive (MDD): 93/42/EEC Council Directive of 14 June 1993 on the approximation of the laws of the Member States concerning medical devices.

In order for a product to be able to be labeled with the CE marking, product quality and manufacturing quality systems must be documented to comply with certain standards. Compliance with the European directives and regulations ensures a standard level of performance for device function, safety, and health protection for users, patients, and others.

The CE symbol provided with some products in the Philips Medical Supplies Reference Guide and online store is intended only as an indication of availability in EU countries and is not intended as a legal product marking.

Countries that accept CE marked products include Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Russia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom. A CE mark does not necessarily indicate that a product is marketed or distributed in all of the countries listed above.

Non-medical products that comply with other harmonized directives may also bear the CE mark.

**Note:** Some products identified in this catalogue may not be licensed in accordance with Canadian law.

### Ordering Information

To order the medical supplies in this guide or to obtain more information or additional literature, contact your Philips Medical Supplies representative near you or visit our website: <http://www.philips.com/healthcarestore>

### Online Ordering

Online ordering is a simple process:

1. Log in using a secure password.
2. Enter the product number and quantity in the shopping cart.
3. Enter your PO or credit card number at check-out.
4. Your order is immediately transferred into our fulfillment system, where you can track it anytime.

If you are new to the web, we can guide you through the process on the phone in minutes. Many new users tell us that the ordering process is so intuitive that they do not need help, even the first time.

You can search for supplies in the online catalog three easy ways: by application category, Philips instrument, or keyword.

If you currently transmit via MMS/EDI or auto fax, the interactive system provides real-time information and services not available through one-way transmissions.

Online ordering is available to U.S. customers at the following url: <http://www.philips.com/newhealthcarestore>

# Resuscitation Supplies

## ALS Monitors-Defibrillators

ALS External Multifunction Electrode Pads and Cables	7
ALS External and Internal Paddle Sets	9
ALS Data Management Solutions	11
Q-CPR	12
ALS Accessories	13
ALS Batteries	14
Cases/Pouches	15
Capnography	16
ECG Cables and Leads	18
Patient Monitoring ECG Cables and Leads	20
ECG Electrodes	25
Paper	26
Invasive Blood Pressure	27
NIBP Cuffs and Air Hoses	36
SpO <sub>2</sub> Sensors and Cables	46
Temperature Probes	49

## AED Supplies

External AED Pads	52
AED Accessories and Kits	54
AED Data Management Solutions	55
AED Training Materials	56
AED Batteries	59
AED Cases/Mounts/Cabinets	60

# ALS Monitors – Defibrillators

## ALS External Multifunction Electrode Pads and Cables

### Multifunction Electrode Pads

Product #	Description	Packaging	Use With:
REVISÉ M3501A CE	Adult/Child with barrel connector	10 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ M3504A CE	Infant with barrel connector	5 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ M3713A CE	Adult/Child Plus with plug style connector	10 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ M3716A CE	Adult/Child Radiolucent with plug style connector	10 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ M3717A CE	Infant Plus with plug style connector	5 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ M3718A CE	Adult/Child Radiopaque (Reduced Skin Irritation) with plug style connector	10 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ M3719A CE	Infant Radiopaque (Reduced Skin Irritation) with plug style connector	5 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
REVISÉ 989803166021 CE	Adult/Child Preconnect	10 sets/case	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A



M3713A



M3717A



M3718A



989803166021

**External Multifunction Pads Cables, Adapters, Test Loads**


Product #	Description	Packaging	Use With:
989803171271 CE	Test Plug for M3508A Cable	1 pack = 1 test plug	HeartStart XL+ Defibrillator/Monitor 861290
M3507A CE	Hands-free Pads Cable (Barrel connector)	1 cable	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
M3508A CE	Hands-free Pads Cable (Plug connector)	1 cable	HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A
05-10000 CE	Medtronic™ Pads Adapter	1 adapter	Medtronic Physio-Control LifePak 500, 9, 10c, 11, 12, 20
05-10100 CE	Zoll™ Pads Adapter	1 adapter	Zoll 1200, 1400, 1600 and M-Series
05-10200 CE	Defibrillator Pads Adapter	1 adapter	Fits Philips hands-free cable M3507A to connect to: M3535A, M3536A, M4735A
M1781A CE	50-ohm test load for M3507A cable	1	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M3725A CE	HeartStart 50-ohm test load for M3507A. Can also be used in combination with 05-10200 adapter.	1	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4748A CE	HeartStart extension cable Length: 4" (10cm)	1 cable	HeartStart XL: M4735A; 989803171271 Test Plug



M3507A



M3508A



05-10000



05-10100



05-10200



## ALS External and Internal Paddle Sets

### External Anterior/Anterior Paddle Sets

Product #	Description	Packaging	Use With:
M3543A CE	HeartStart MRx Water Resistant External Paddles with PCI	1 set	HeartStart MRx: M3535A, M3536A
M4745A CE	Sterilizable External Paddles without PCI	1 set	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4746A CE	HeartStart XL External Paddles with PCI	1 set	HeartStart XL: M4735A

### Replacement Electrodes for External Paddle Sets

Product #	Description	Packaging	Use With:
M1789A CE	Rectangular Paddle Electrode Replacement	1 electrode	For use with M4746A External Paddle Set with PCI
M4759A	Rectangular Paddle Electrode Replacement	1 electrode	M3543A Paddle Set

### Switchless Reusable Internal Paddles (Discharge Switch on Defibrillator)

Product #	Description	Packaging	Use With:
M1741A CE	Adult Internal Paddles: 3", 7.5cm, switchless	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M1742A CE	Adult Internal Paddles: 2.4", 6.0cm, switchless	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M1743A CE	Pediatric Internal Paddles: 1.8", 4.5cm, switchless	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M1744A CE	Infant Internal Paddles: 1.1", 2.8cm, switchless	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A



M1741A-M1744A

**Switched Reusable Internal Paddles (Discharge Switch on Handle)**

Product #	Description	Packaging	Use With:
M4741A CE	Adult Internal Paddles: 3", 7.5cm, switched	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4742A CE	Adult Internal Paddles: 2.4", 6.0cm, switched	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4743A CE	Pediatric Internal Paddles: 1.8", 4.5cm, switched	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4744A CE	Infant Internal Paddles: 1.1", 2.8cm, switched	1 pair	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A

**Internal Paddle Adapters**

Product #	Description	Packaging	Use With:
M4740A CE	HeartStart Switchless Internal Paddle Adapter	1 adapter	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4748A CE	HeartStart Extension Cable	1 adapter	HeartStart XL: M4735A

**Sync Cables**

Product #	Description	Packaging	Use With:
M1783A CE	Sync Cable, 12-pin. Length: 8', 2.5m	1 cable	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M5526A CE	Sync Cable, 12-pin. Length: 25', 7.8m	1 cable	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A



M4741A-M4744A



M4740A

# ALS Data Management Solutions

## HeartStart Data Cards and Software

Product #	Description	Packaging	Use With:
989803147711 CE	HeartStart XL Data Card	1	HeartStart XL: M4735A
M3544A CE	HeartStart MRx Data Card Tray	1	HeartStart MRx: M3535A, M3536A
989803171261 CE	USB Data Drive	1 pack = 1 USB drive	HeartStart XL+ Defibrillator/Monitor 861290

# Q-CPR

## Q-CPR Compression Sensor and Accessories

Product #	Description	Packaging	Use With:
M4761A CE	Q-CPR Compression Sensor	1 sensor	HeartStart MRx: M3535A, M3536A
M4762A CE	Q-CPR Compression Sensor Adhesive Pads	10 pads	HeartStart MRx: M3535A, M3536A
(LATEX) M4763A CE	Pads/CPR Cable	1 cable	HeartStart MRx: M3535A, M3536A
989803162401 CE	CPR Meter – next generation Q-CPR	1 sensor	HeartStart MRx: M3535A, M3536A Software Rev. F or higher
989803163291 CE	CPR Meter Adhesive Pads – next generation Q-CPR	10 pads	HeartStart MRx: M3535A, M3536A Software Rev. F or higher
989803158661 CE	Pads/CPR Cable – next generation Q-CPR	1 cable	HeartStart MRx: M3535A, M3536A Software Rev. F or higher



M4761A and M4762A



989803162401

# ALS Accessories

## HeartStart MRx Accessories

Product #	Description	Packaging	Use With:
453564063841	MRx NIBP Calibration Kit	1 kit	HeartStart MRx: M3535A, M3536A
453564063851	MRx etCO <sub>2</sub> Calibration Kit	1 kit	HeartStart MRx: M3535A, M3536A
M3537A	Bed Rail/Roll Stand Hook for Defibrillator	1 hook	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M3549A	MRx Wide Bed Rail/Roll Stand Hook	1 hook	HeartStart MRx: M3535A, M3536A
M4737A	HeartStart MRx Monitor/Defibrillator Display Cover	1 cover	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5521A	Green Handle and label set	1 handle	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5522A	Blue Handle and label set	1 handle	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5523A	Yellow Handle and label set	1 handle	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5524A	Rose Handle and label set	1 handle	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5525A	Gray Handle and label set	1 handle	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5528A	Mounting Bracket	1 bracket	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
M5529A	DC Power Module	1	HeartStart MRx: M3535A, M3536A
<b>CE</b>			
989803171701	Bed Rail/Roll Stand Hook	1 pack = 1 bed rail/ roll stand mount	HeartStart XL+ Defibrillator/Monitor 861290
<b>CE</b>			
989803176541	HeartStart MRx Quick Disconnect DC Power Cable	1 pack = 1 complete set of cables	861593, 861494, 861495, M3536A



M5521A



M5528A



M5529A

# ALS Batteries

## HeartStart ALS Batteries

Product #	Description	Packaging	Use With:
M3516A CE	Sealed Lead Acid Battery	1 battery	HeartStart XL: M4735A
M3538A CE	Lithium Ion Battery with Capacity Gauge	1 battery	HeartStart MRx: M3535A, M3536A
M3506A CE	Battery Charger Adapter Tray	1 tray	HeartStart MRx: M3516A ALS
M4747A CE	HeartStart XL Battery Charger Kit Consists of: M3516A battery, M3506A charger tray and M3517A AC Power Module	1 kit	HeartStart XL: M4735A



M3516A



M3538A

## Cases/Pouches

### Cases and Pouches

Product #	Description	Packaging	Use With:
M3541A CE	MRx Complete Carrying Case (includes 3 accessory pouches, shoulder strap and sleeve)	1 case	HeartStart MRx: M3535A, M3536A
M4751A	HeartStart XL Accessory Pouch	1 pouch	HeartStart XL: M4735A
989803172501	Black MRx Carry Case	1 box = 1 case	HeartStart MRx M3535A, M3536A
989803171281 CE	Accessory Storage Bags	1 pack = 1 Therapy Accessory Pouch & 1 Monitoring Accessory Pouch	HeartStart XL+ Defibrillator/Monitor 861290
989803171291 CE	Cable Retention Straps	1 pack = 3 cable wraps	HeartStart XL+ Defibrillator/Monitor 861290
989803180871	MRx Black Complete Carrying Case	1 pack = 1 complete carrying case	M3536A



M3541A

# Capnography

## IntelliVue Microstream® CO<sub>2</sub> Solution

Product #	Description	Packaging	Use With:
989803177951	Smart CapnoLine H O <sub>2</sub> Adult	25 Filterlines	
CE			
989803177961	Smart CapnoLine H O <sub>2</sub> Adult Long	25 Filterlines	
CE			
989803177971	Smart CapnoLine H O <sub>2</sub> Pediatric	25 Filterlines	
CE			
989803177981	Smart CapnoLine H O <sub>2</sub> Pediatric Long	25 Filterlines	
CE			
989803178001	CapnoLine H O <sub>2</sub> Infant/Neonatal	25 Filterlines	
CE			
989803178011	CapnoLine H Infant/Neonatal Long	25 Filterlines	
CE			
989803178021	Nasal FilterLine Infant/Neonatal	25 Filterlines	
CE			
989803178031	Smart CapnoLine Guard	25 Filterlines	
CE			
989803178041	Smart CapnoLine Guard O <sub>2</sub>	25 Filterlines	
CE			
989803178051	Smart CapnoLine Guard O <sub>2</sub> Long	25 Filterlines	
CE			
989803178071	Hook and Loop Strap	1 box = 25	
CE			
989803179101	Nasal Filterline O <sub>2</sub> Adult	25 Filterlines	
CE			
989803179111	Nasal Filterline O <sub>2</sub> Adult Long	25 Filterlines	
CE			
989803179121	Nasal Filterline O <sub>2</sub> Pediatric	25 Filterlines	
CE			
M1920A	FilterLine Set, Adult/Pedi, Intubated	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06
CE			
M1921A	FilterLine H Set, Adult/Pediatric, Intubated (yellow)	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06
CE			



M1920A



M1921A



**IntelliVue Microstream® CO<sub>2</sub> Solution** *(continued)*

Product #	Description	Packaging	Use With:
M1923A CE	FilterLine H Set, Infant/Neonatal, Intubated (yellow)	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06
M2520A CE	Smart CapnoLine O <sub>2</sub> , Pediatric, Non-intubated dual purpose	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06
M2522A CE	Smart CapnoLine O <sub>2</sub> plus, Adult/intermed, Non-intubated dual purpose	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06
M2524A CE	Smart CapnoLine, Pediatric, Non-intubated single purpose	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06
M2526A CE	Smart CapnoLine plus, Adult/intermediate, Non-intubated single purpose	25 sets	HeartStart MRx: M3535A Opt A06, M3536A Opt A06



M1923A



M2520A



M2522A



M2524A



## ECG Cables and Leads

### 10-Lead ECG Cables and Leads (3+7)

Product #	Description	Packaging	Use With:
M3525A CE	10-lead ECG patient trunk cable for 3 and 7 electrode combination sets AAMI/IEC 12-pin connector Length: 2.7m	1 cable	HeartStart MRx: M3535A, M3536A
M3526A	3-electrode limb set, snap, AAMI	1 lead set	HeartStart MRx: M3535A, M3536A
M3527A	7-lead chest lead set, snap, AAMI	1 lead set	HeartStart MRx: M3535A, M3536A
M3528A CE	3-electrode limb lead set, snap, IEC	1 lead set	HeartStart MRx: M3535A, M3536A
M3529A CE	7-electrode chest lead set, snap, IEC	1 lead set	HeartStart MRx: M3535A, M3536A
M5530A CE	Combiner Plug for 3-Wire Lead Set	1 lead set	HeartStart MRx: M3535A, M3536A



M3525A

**10-Lead ECG Cables and Leads (5+5)**

Product #	Description	Packaging	Use With:
989803147691 CE	10-lead ECG patient trunk cable (5+5) AAMI/IEC 12-pin connector Length: 1.3m	1 cable	HeartStart MRx: M3535A, M3536A
989803158061	5-lead set, limb, AAMI	1 lead set	HeartStart MRx: M3535A, M3536A
989803158071	5-lead set, chest, AAMI	1 lead set	HeartStart MRx: M3535A, M3536A
989803158081	5-lead set, limb, IEC	1 lead set	HeartStart MRx: M3535A, M3536A
989803158091 CE	5-lead set, chest, IEC	1 lead set	HeartStart MRx: M3535A, M3536A



989803158061 and 989803158071



989803158081 and 989803158091

## Patient Monitoring ECG Cables and Leads

### IntelliVue 10-Lead (5+5) ECG Monitoring – AAMI/IEC

(for conventional diagnostic 12-lead ECG, use one limb 5-lead set and one chest 5-lead set)

Product #	Description	Packaging	Use With:
M1663A CE	10-Lead ECG Patient Trunk (5+5) – AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 6.5' (2m)	1 cable	HeartStart MRx: M3535A, M3536A
M1949A CE	10-Lead ECG Patient Trunk (5+5) – AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 8.9' (2.7m)	1 cable	HeartStart MRx: M3535A, M3536A, 860338, 860335
989803176161	Cbl 5-Lead, snap, shld, AAMI, limb rgd (5.2')	1 lead set	861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC
989803176171	Cbl 5-Lead, snap, shld, AAMI, chest, rgd (3.9')	1 lead set	861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC
989803176181	Cbl 5-Lead, snap, shld, IEC, limb, rgd (5.2')	1 lead set	861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC
989803176191	Cbl 5-Lead, snap, shld, IEC, chest, rgd	1 lead set	861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC



M1663A



M1949A

**5-Lead Sets for 10-Lead Monitoring – AAMI Color Coding**

Product #	Description	Packaging	Use With:
M1968A	ICU grabber, limb, AAMI	1 lead set	Use with M1663A or M1949A 10-lead patient trunk cable for HeartStart MRx: M3535A, M3536A, 860335, 860338
M1976A	ICU grabber, chest, AAMI	1 lead set	
M1644A	ICU snap, limb, AAMI	1 lead set	
M1602A	ICU snap, chest, AAMI	1 lead set	
M1973A	OR grabber, limb, AAMI	1 lead set	
M1979A	OR grabber, chest, AAMI	1 lead set	
989803173131 CE	5-Lead Set, Disposable, Bedside, AAMI, round	1 pack = 20 lead sets	Compatible Defibrillators M3535A (MRx monitor/ defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL)



M1968A

M1976A



M1973A



M1979A

**5-Lead Sets for 10-Lead Monitoring – IEC Color Coding**

Product #	Description	Packaging	Use With:
M1604A CE	ICU snap, chest, IEC	1 lead set	Use with M1663A or M1949A 10-lead patient trunk cable for HeartStart MRx: M3535A, M3536A, 860335, 860338
M1645A CE	ICU snap, limb, IEC	1 lead set	
M1971A CE	ICU grabber, limb, IEC	1 lead set	
M1974A CE	OR grabber, limb, IEC	1 lead set	
M1978A CE	ICU grabber, chest, IEC	1 lead set	
M1984A CE	OR grabber, chest, IEC	1 lead set	
989803174211 CE	5-Lead Set, Disposable, Bedside, IEC, round	1 pack = 20 lead sets	



M1971A



M1974A



M1978A



M1984A

**3-Lead Monitoring – AAMI/IEC**

Product #	Description	Packaging	Use With:
M1669A CE	3-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector. AAMI and IEC labels included. Length: 8.9' (2.7m)	1 cable	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A



M1669A

**3-Lead Sets for 3-Lead Monitoring – AAMI/IEC Color Coding**

Product #	Description	Packaging	Use With:
M1671A	ICU grabber, AAMI	1 lead set	Use with M1669A 3-Lead ECG Patient Trunk Cable, 860335, 860338
M1673A	ICU snap, AAMI	1 lead set	
M1675A	OR grabber, AAMI	1 lead set	
989803173121 CE	3-Lead Set, Disposable, Bedside, AAMI, round	1 pack = 20 lead sets	Compatible Defibrillators M3535A (MRx Monitor/ defib), M3536A (MRx Monitor/defib), M4735A (HeartStart XL)
M1672A CE	ICU grabber, IEC	1 lead set	Use with M1669A 3-Lead ECG Patient Trunk Cable, 860335, 860338
M1674A CE	ICU snap, IEC	1 lead set	
M1678A CE	OR grabber, IEC	1 lead set	
989803174201 CE	3-Lead Set, Disposable, Bedside, IEC, round	1 pack = 20 lead sets	Compatible Defibrillators M3535A (MRx Monitor/ defib), M3536A (MRx Monitor/defib), M4735A (HeartStart XL)



M1671A



M1672A



M1675A



M1678A



989803174201

**5-Lead Monitoring – AAMI/IEC**

Product #	Description	Packaging	Use With:
M1668A CE	5-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector. AAMI and IEC labels included. Length: 8.9' (2.7m)	1 cable	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A, 860335, 860338



M1668A

**5-Lead Sets for 5-Lead Monitoring – AAMI/IEC Color Coding**

Product #	Description	Packaging	Use With:
M1644A	ICU snap, limb, AAMI	1 lead set	Use with M1668A 5-Lead ECG Patient Trunk Cable for HeartStart MRx: M3535A, M3536A; 860335, 860338; HeartStart XL: M4735A
M1645A CE	ICU snap, limb, IEC	1 lead set	860338, 860335
M1968A	ICU grabber, limb, AAMI	1 lead set	Use with M1668A 5-Lead ECG Patient Trunk Cable for HeartStart MRx: M3535A, M3536A; 860335, 860338; HeartStart XL: M4735A
M1971A CE	ICU grabber, limb, IEC	1 lead set	860338, 860335
M1973A CE	OR grabber, limb, AAMI	1 lead set	Use with M1668A 5-Lead ECG Patient Trunk Cable for HeartStart MRx: M3535A, M3536A; 860335, 860338; HeartStart XL: M4735A
M1974A CE	OR grabber, limb, IEC	1 lead set	860338, 860335
989803173131 CE	5-Lead Set, Disposable, Bedside, AAMI, round	1 pack = 20 lead sets	Compatible Bedside Monitors M1001B, M1002A, M3000A, M3001A, M8105A (MP5), 863064 (VM6), 863065 (VM6), 863066 (VM8), 863068 (VM6), Compatible Defibrillators M3535A (MRx monitor/ defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL)
989803174211 CE	5-Lead Set, Disposable, Bedside, IEC, round	1 pack = 20 lead sets	Compatible Defibrillators M3535A (MRx monitor/ defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL)



M1644A



M1645A



M1968A



M1971A



M1973A



M1974A



989803174211



# ECG Electrodes

Product #	Description	Packaging	Use With:
989803148821 CE	Small Adult, radiolucent foam, solid gel ECG electrode Size: 1.4" x 1.6" (available in U.S. and Japan)	30/pouch 600/case	HeartStart MRx: M3535A, M3536A
M2202A CE	Adult radiotranslucent foam solid gel electrode Size: 2" (50mm) diameter	5/pouch 300/case	HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A
M4606A CE	ECG skin preparation paper	100 preps	All ECG Electrodes
M4612A CE	Adult solid gel ECG electrode Shape: teardrop	5/pouch 300/case	HeartStart MRx: M3535A, M3536A
M4613A CE	Adult solid gel ECG electrode Shape: teardrop	30/pouch 300/case	HeartStart MRx: M3535A, M3536A



989803148821



M2202A



M4606A



M4612A/M4613A

# Paper

## Recording Paper

Product #	Description	Packaging	Use With:
989803138171 CE	MRx Wide Printer Paper	10 rolls	M3535A, M3536A, M3536J Opt B04 (75mm printer)
989803138181 CE	MRx Wide Printer Paper	80 rolls	M3535A, M3536A, M3536J Opt B04 (75mm printer)
40457C CE	1-Channel Chem/Thermal Paper	10 rolls	M4735A, M3535A, M3536A, M3636J
40457D CE	1-Channel Chem/Thermal Paper	80 rolls	M4735A, M3535A, M3536A, M3536J



989803138181



40457C

## Invasive Blood Pressure

### TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits – Single Line

Product #	Description	Packaging	Use With:
989803177901 CE 0123	Single DPT TP4 60" (152cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Dash Tubing	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803181211 CE 0123	Single DPT TP4 60" (152cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Premium Stripe Tubing	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179771 CE 0123	Single DPT TP4 72" (183cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Dash Tubing	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803181221 CE 0123	Single DPT TP4 72" (183cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Premium Stripe Tubing	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179721 CE 0123	Single DPT TP4 84" (213cm), with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Dash Tubing	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803181231 CE 0123	Single DPT TP4 84" (213cm), with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Premium Stripe Tubing	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



98903177901



989803181211

**TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits – Single Line** *(continued)*

Product #	Description	Packaging	Use With:
989803179871 CE 0123	Single DPT TP4 9" (23cm) Patient Mount with Drip Chamber, One 3-Way Stopcock, 3ml/hr Flush Device	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803181141 CE 0123	Single DPT TP4 24" (61cm) Compartmental Pressure, One 4-Way Stopcock	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803179871



989803181141

**TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits with SAFESET™ Blood Sampling/Conservation System – Single Line**

Product #	Description	Packaging	Use With:
989803179741 CE 0123	Single DPT TP4 84" (213cm) w/SAFESET, 2 Cannula Sampling Ports, In-line Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179751 CE 0123	Single DPT TP4 84" (213cm) w/SAFESET, 2 Luer Sampling Ports, In-line Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179861 CE 0123	Single DPT TP4 24" (61cm) w/SAFESET, Patient Mount, 1 Cannula Sampling Port, In-Line Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179761 CE 0123	Single DPT TP4 60" (152cm) w/SAFESET, 1 Cannula Sampling Port, In-Line Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179781 CE 0123	Single DPT TP4 60" (152cm) w/SAFESET, 2 Cannula Sampling Ports, In-Line Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179791 CE 0123	Single DPT TP4 60" (152cm) w/SAFESET, 2 Luer Sampling Ports, In-Line Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803179751



98980179761



989803179781

**TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits – Multi-Lines**

Product #	Description	Packaging	Use With:
989803177911 CE 0123	Dual DPT TP4 72" (183cm) with Drip Chamber, Four 3-Way Stopcocks, Two 3ml/hr Flush Devices, Red/Blue Dash Tubing	10 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179731 CE 0123	Triple DPT TP4 60" (152cm) with Drip Chamber, Six 3-Way Stopcocks, Three 3ml/hr Flush Devices, Red/Blue/Yellow Dash Tubing	10 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A
989803181241 CE 0123	Dual DPT TP4 72" (183cm) with Drip Chamber, Four 3-Way Stopcocks, Two 3ml/hr Flush Devices, Red/Blue Premium Stripe Tubing	10 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803181251 CE 0123	Triple DPT TP4 60" (152cm) with Drip Chamber, Six 3-Way Stopcocks, Three 3ml/hr Flush Devices, Red/Blue/Yellow Premium Stripe Tubing	10 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A



989803177911



989803179731



989803181241



989803181251

**TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits with SAFESET™ Blood Sampling/Conservation System – Multi- Lines**

Product #	Description	Packaging	Use With:
989803179801 CE 0123	Dual DPT TP4 84" (213cm) w/SAFESET, (Red/Blue Dash Tubing), 2 Cannula Sampling Ports, In-Line Reservoir	10 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179811 CE 0123	Dual DPT TP4 84" (213cm) w/SAFESET, (Red/Blue Dash Tubing), 2 Luer Sampling Ports, In-Line Reservoir	10 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803179801



989803179811

**SAFESET™ Disposable Blood Sampling/Conservation System Only**

Product #	Description	Packaging	Use With:
989803179891 CE 0123	SAFESET kit w/ 84" (213cm) Tubing, Two In-Line Cannula Sampling Ports, 10ml In-Line Reservoir	20 kits	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803180851 CE 0123	SAFESET kit w/ 60" (152cm) Tubing, One In-line Cannula Sampling Ports, 10ml In-Line Reservoir	20 kits	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803180851



**TRANSPAC® 4 Reusable Transducer Cables**

Product #	Description	Packaging	Use With:
989803177921 CE 0123	Transpac 4 Cable for Philips Monitors 10' (3m)	1 cable	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179941 CE 0123	Transpac 4 Cable for Philips Monitors 15' (4.5m)	1 cable	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179951 CE 0123	Transpac 4 Cable for Philips Monitors 15' (4.5m), bulk package	20 cables	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803179941

**TRANSPAC® 4 Disposable Transducer Kit Accessories**

Product #	Description	Packaging	Use With:
989803177931 CE 0123	3-Position Transducer Mount, Reusable	28 mounts	TRANSPAC 4 Disposable Transducer Kits
989803177941 CE 0123	Transducer Holder for IV Pole, Reusable	5 holders	TRANSPAC 4 Disposable Transducer Kits
989803179901 CE 0123	SAFESET Mounting Bracket, Reusable	30 brackets	TRANSPAC 4 Disposable Transducer Kits and SAFESET System
989803179911 CE 0123	Single Transducer Mount, Reusable	30 mounts	TRANSPAC 4 Disposable Transducer Kits
989803179921 CE 0123	SAFESET Shielded Blunt Cannula	50 blunt cannulas	TRANSPAC 4 Disposable Transducer Kits and SAFESET System
989803179931 CE 0123	SAFESET Shielded Blunt Cannula w/ Blood Tube Holder	25 cannula with tube holders	TRANSPAC 4 Disposable Transducer Kits and SAFESET System
989803180861 CE 0123	Transducer Simulator, Reusable	1 simulator	TRANSPAC 4 Disposable Transducer Kits and Cables with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803177931



989803177941



989803179901



989803179921

**Memscap® Reusable Invasive Blood Pressure Solution**

Product #	Description	Packaging	Use With:
CPJ84046 CE	Pressure Transducer Holder	1 holder	CPJ840J6 transducer with HeartStart MRx: M3535A, M3536A, M3015A
CPJ84447 CE	Pressure Transducer Pole Mount	1 mount	CPJ840J6 transducer with HeartStart MRx: M3535A, M3536A, M3015A
CPJ840J6 CE	Reusable Pressure Transducer	1 transducer	HeartStart MRx: M3535A, M3536A, M3015A
CPJ84022 CE	Disposable Domes	50 domes	CPJ840J6 transducer with HeartStart MRx: M3535A, M3536A, M3015A
CPJ84026 CE	Monitoring Line Set	10 sets per bag	863065, 863066, 863068, 863087, 863088, M3536M, M3002A, M8102A, M3535A, M8105AT, M8105A, M3536A, M3015A, M3015B



CPJ840J6

## NIBP Cuffs and Air Hoses


### Multi-Patient Traditional Cuffs – Single Hose

Product #	Description	Packaging	Use With:
REVISED 40401A CE	Infant Bladder width = 2.4" (6cm) Bladder length = 4.7" (12cm) Limb circumference = 10 – 19cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED 40401B CE	Pediatric Bladder width = 3.5" (8.9cm) Bladder length = 7.1" (18cm) Limb circumference = 18 – 26cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED 40401C CE	Adult Bladder width = 4.8" (12.1cm) Bladder length = 9.1" (23cm) Limb circumference = 25 – 35cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED 40401D CE	Large Adult Bladder width = 6" (15.2cm) Bladder length = 13" (33cm) Limb circumference = 33 – 47cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED 40401E CE	Thigh Bladder width = 7" (17.8cm) Bladder length = 14.2" (36cm) Limb circumference = 46 – 66cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495



**Multi-Patient Comfort Cuffs – Single Hose**

Product #	Description	Packaging	Use With:
REVISED M1572A CE	Pediatric, color coded: green tubing Bladder width: 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED M1573A CE	Small Adult, color coded: light blue tubing Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference: 20.5 – 28cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED M1574A CE	Adult, color coded: navy blue tubing Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27 – 35cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED M1575A CE	Large Adult, color coded: burgundy tubing Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference: 34 – 43cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495



**Multi-Patient Easy Care Cuffs – Single Hose**

Product #	Description	Packaging	Use With:
REVISÉD M4552B CE	Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10 – 15cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4552B5 CE	Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10 – 15cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4553B CE	Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4553B5 CE	Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495



**Multi-Patient Easy Care Cuffs – Single Hose** (continued)

Product #	Description	Packaging	Use With:
REVISÉD M4554B CE	Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4554B5 CE	Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4555B CE	Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4555B5 CE	Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4556B CE	Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4556B5 CE	Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4557B CE	Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4557B5 CE	Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4558B CE	Large Adult Extra-Long, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495

**Multi-Patient Easy Care Cuffs – Single Hose** *(continued)*

Product #	Description	Packaging	Use With:
REVISED M4558B5 CE	Large Adult Extra-Long, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED M4559B CE	Thigh, grey Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference = 44 – 56cm	1 cuff	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED M4559B5 CE	Thigh, grey Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference = 44 – 56cm	5 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495





**Multi Patient Disposable Cuffs – Multi Care**

Product #	Description	Packaging	Use With:
<b>NEW</b> 989803183311 CE	Infant, orange Bladder width = 2.2" (5.5cm) Bladder length = 6.1" (15.4cm) Limb circumference = 10 – 15cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803183321 CE	Pediatric, green Bladder width = 3.1" (8cm) Bladder length = 8.5" (21.5cm) Limb circumference = 14 – 21.5cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803183331 CE	Small Adult, light blue Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference = 20.5 – 28cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803183341 CE	Adult, navy blue Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27 – 35cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803183351 CE	Adult XL, navy blue with edge striping Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27 – 35cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803183361 CE	Large Adult, burgundy Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference = 34 – 43cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803183371 CE	Thigh, grey Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm) Limb circumference: 42 – 54cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495

**Single Care Disposable Cuffs**

Product #	Description	Packaging	Use With:
<b>NEW</b> 989803182281 CE	Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm	20 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803182291 CE	Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm	20 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803182301 CE	Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm	20 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803182311 CE	Adult XL, navy blue with edge striping Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm	20 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
<b>NEW</b> 989803182321 CE	Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35 – 45cm	20 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495

**Gentle Care Disposable Cuffs – Single Hose**

Product #	Description	Packaging	Use With:
REVISÉD M4572B CE	Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10 – 15cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4573B CE	Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4574B CE	Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4575B CE	Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4576B CE	Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISÉD M4577B CE	Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35 – 45cm	10 cuffs	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495

**Gentle Care Disposable Cuffs – Single Hose** *(continued)*

<p>REVISED</p> <p>M4578B</p> <p>CE</p>	<p>Large Adult Extra-Long, burgundy</p> <p>Bladder width: 6.7" (17cm)</p> <p>Bladder length: 15" (38cm)</p> <p>Limb circumference = 35 – 45cm</p>	<p>10 cuffs</p>	<p>HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495</p>
<p>REVISED</p> <p>M4579B</p> <p>CE</p>	<p>Thigh, grey</p> <p>Bladder width: 8.3" (21cm)</p> <p>Bladder length 18.9" (48cm)</p> <p>Limb circumference = 44 – 56cm</p>	<p>10 cuffs</p>	<p>HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495</p>



**Blood Pressure Air Hoses**

Product #	Description	Packaging	Use With:
REVISED M1598B CE	NIBP Air Hose Length: 4.92' (1.5m)	1 cable	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495
REVISED M1599B CE	NIBP Air Hose Length: 9.84' (3.0m)	1 cable	HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495



M1599B

## SpO<sub>2</sub> Sensors and Cables

### Philips Disposable SpO<sub>2</sub> Sensors

Product #	Description	Packaging	Use With:
M1131A CE	Adult/Pediatric Finger Clip Disposable Sensor (single patient use) Length = 45cm Patient size: >44lb (>20kg)	20 sensors	HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278



M1131A

### Philips Reusable SpO<sub>2</sub> Sensors

Product #	Description	Packaging	Use With:
M1191B CE	Adult Finger Glove Sensor Length = 2m Patient size: >110lb (>50kg) 8-pin connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1191BL CE	Adult Finger Glove Sensor Length = 3m Patient size: >110lb (>50kg) 8-pin connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1191T CE	Adult Finger Glove Sensor Length = 0.45m Patient size: >110lb (>50kg) 9-pin (D-sub) connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1192A CE	Pediatric/Small-Adult Glove Sensor Length = 1.5m Patient size: 33lb to 110lb (15kg to 50kg) 8-pin connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278



M1191B



M1191T



M1192A

**Philips Reusable SpO<sub>2</sub> Sensors** (continued)

Product #	Description	Packaging	Use With:
M1192T CE	Pediatric/Small-Adult Glove Sensor Length = 0.45m Patient size: 33 – 110lb (15 – 50kg) 9 pin (D-sub) connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1194A CE	Pediatric/Adult Ear Clip Sensor Length = 1.5m Patient size: >88lb (>40kg) 8 pin connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1195A CE	Reusable Infant Finger Glove Sensor Length = 1.5m Patient size: 8.8 – 33lb (4 – 15kg) 8-pin connector	1 sensor	HeartStart MRx only: M3535A, M3536A, 863276, 863278
M1196A CE	Pediatric/Adult Finger Clip Sensor Length = 3m Patient size: >88lb (>40kg) 8 pin connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278
M1196S CE	Reusable Clip Adult SpO <sub>2</sub> Sensor Length = 2m Patient size: > 88lb (>40kg) 8-pin connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1196T CE	Pediatric/Adult Finger Clip Sensor Length 90cm Patient size: >88lb (>40kg) 9-pin (D-sub) connector	1 sensor	HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278
M1627A CE	Wristband Velcro straps for SpO <sub>2</sub> glove sensors	10 straps	M1191B/BL/T sensors, all ALS Defibrillators



M1194A



M1195A



M1196A

**SpO<sub>2</sub> Adapter Cables for Philips Sensors**

Product #	Description	Packaging	Use With:
M1941A CE	SpO <sub>2</sub> Extension Cable Adds additional length to 8-pin sensors with 8-pin sockets Length: 6.6' (2m)	1 cable	To add additional length to M1191B, M1192A, M1194A, M1195A sensors only, connecting to HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278
M1943A CE	SpO <sub>2</sub> Adapter Cable Adapts 9-pin sensors to 8-pin sockets Length: 3.3' (1m)	1 cable	Adapter cable for M1191T, M1192T and M1196T sensors only, connecting to HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278



M1941A



M1943A




# Temperature Probes

## Disposable Temperature Probes

Product #	Description	Packaging	Use With:
M1837A CE	Esophageal/Rectal Temperature Probe, 9 FR Length: 36" (91cm), sterilized, continuous monitoring Shelf life: at least 12 months	20 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21090A CE	Esophageal/Rectal Temperature Probe, 12 FR Length: 36" (91cm) sterilized, continuous monitoring Shelf life: at least 12 months	20 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21091A CE	Disposable sterile skin surface temperature probe Length: 30" (76cm) Shelf life: at least 12 months	1 box = 20 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21093A CE	Esophageal/Stethoscope Temperature Probe, 12 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months	20 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21094A CE	Esophageal/Stethoscope Temperature Probe, 18 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months	20 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21095A CE	Esophageal/Stethoscope Temperature Probe, 24 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months	20 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A



21090A



21091A



21094A

**Disposable Temperature Probes** (continued)

Product #	Description	Packaging	Use With:
21096A CE	Foley Catheter Temperature Probe, 16FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months	10 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21097A CE	Foley Catheter Temperature Probe, 18FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months	10 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
M2255A CE	Foley Catheter Temperature Probe, 14FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months	10 probes	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A



21096A



21097A



M2255A

**Disposable Temperature Probe Extension Cables**

Product #	Description	Packaging	Use With:
21082A CE	Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 10' (3m)	1 cable	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21082B CE	Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 5' (1.5m)	1 cable	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A



21082A

**Reusable Temperature Probes for Continuous Monitoring (meets EN12470 Standard)**

Product #	Description	Packaging	Use With:
21075A CE	Esophageal/Rectal Temperature Probe, 12FR (autoclavable) 2-prong plug, reusable, Length: 10' (3m)	1 probe	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21076A CE	Esophageal/Rectal Temperature Probe, 10FR (autoclavable) 2-prong plug, reusable, continuous monitoring Length: 10' (3m)	1 probe	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A
21078A CE	Skin Surface Temperature Probe (autoclavable) 0.37" (9.5mm) diameter 2-prong plug, reusable tape-on, continuous monitoring Length: 10' (3m)	1 probe	HeartStart MRx: M3535A Opt A06, M3536A Opt A06, 860335, 860338, M3015A



21075A



21076A



21078A

# AED Supplies

## External AED Pads

### External AED Pads and Adapters

Product #	Description	Packaging	Use With:
M5071A CE	Adult SMART Pads Cartridge	1 set	HeartStart HS1: M5066A, M5067A, M5068A
M5072A CE	Infant/Child SMART PADS Cartridge	1 set	HeartStart HS1: M5066A, M5067A, M5068A
989803139261 CE	SMART Pads II	1 set	HeartStart FRx: 861304, 861305
989803158211 CE	HeartStart Adult Defibrillator Pads	1 set	ForeRunner: E, S, EM FR series: 940010XX, 940020XX, 940030XX FR2 series: M3840A, M3841A, M3860A, M3861A
989803158221 CE	HeartStart Adult Defibrillator Pads	5 sets	ForeRunner: E, S, EM FR series: 940010XX, 940020XX, 940030XX FR2 series: M3840A, M3841A, M3860A, M3861A



M5071A



M5072A



989803139261

**External AED Pads and Adapters** *(continued)*

Product #	Description	Packaging	Use With:
M3870A CE	Infant/Child Reduced-Energy Defibrillator Pads	1 set	HeartStart FR2+: M3860A, M3861A
05-10000 CE	Medtronic™ Pads Adapter	1 adapter	989803158211, 989803158221
05-10100 CE	Zoll™ Pads Adapter	1 adapter	989803158211, 989803158221



M3870A



05-10000



05-10100

## AED Accessories and Kits

### AED Fast Response Kit

Product #	Description	Packaging	Use With:
68-PCHAT	Fast Response Kit	1 kit: 2 pairs/gloves, 1 CPR mask, 1 scissors, 1 prep razor, 1 paper towel, 1 zippered pouch	HeartStart HS1: M5066A, M5067A, M5068A ForeRunner AED: E, S, EM HeartStart FR: 940010XX, 940020XX, 940030XX HeartStart FR2: M3860A, M3861A, M3840A, M3841A HeartStart FRx: 861304, 861305

### AED Accessories

Product #	Description	Packaging	Use With:
989803139311 CE	Infant/Child Key. Under 8 years old (<55lb or <25kg)	1 key	HeartStart FRx: 861304, 861305
M3873A CE	FR2+ ECG Assessment Module, AAMI	1 module	HeartStart FR2+: M3860A, M3840A
M3874A CE	FR2+ ECG Assessment Module, IEC	1 set electrodes User guide	HeartStart FR2+: M3860A, M3840A
989803137771	Replacement ECG Monitoring Electrodes for the ECG Assessment Module	17 packs of 3 electrodes	M3873A FR2+ ECG Assessment Module – AAMI M3874A FR2+ ECG Assessment Module – IEC
861476	AED Awareness Posters 11" (28cm) w, 17" (43cm) h	1 package = 4 posters	861304 (FRx), M5066A (861282-OnSite), M3860A (861264-FR2), M3861A (861265-FR2)
861478	AED Signage Bundle. AED Wall Sign (989803170921), AED Awareness Placard (989803170901) and AED Awareness Poster 4-pack (861476 – English). Wall Sign: 9" (23cm) h, 6.1" (15cm) d, Awareness Placard: 10.25" (26cm) w, 4. Available in US only	1 bundle	FRx (861304), FR2 (M3860A/861264, M3861A/861265), OnSite (M5066A/861282)
861477	AED Wall Mount and Signage Bundle. Wall Mount: 10.5" (27cm) w, 8" (20cm) h, 6.9" (17cm) d; Wall Sign: 9" (23cm) h, 6.1" (15cm) d; Awareness Placard: 10.25" (26cm) w, 4.5" (11cm) h; Awareness Posters: 11" (28cm) w, 17" (43cm) h. Available in US only	1 bundle	FRx (861304), FR2 (M3860A/861264, M3861A/861265), OnSite (M5066A/861282)



989803139311



M3873A

# AED Data Management Solutions

## AED Data Cards and Data Management Accessories

Product #	Description	Packaging	Use With:
M3524A	Data Card Reader	1 card	CodeRunner SW: M3831A, M3833A; HeartStart FR2 series data card/tray: M3854A; HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A;
M3834A CE	HeartStart Event Review Software – Single License	1 CD with license	HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A HeartStart HS1: M5066A, M5067A HeartStart FRx: 861304, 861305
M3853A CE	AED Tray for Data Card	1 tray	HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A
M3854A CE	Data Card and Tray	1 card and tray	
SDCF05 CE	Data Card Adapter	1 adapter	
989803141811 CE	HeartStart Event Review Software – Organization-wide License PC	1 CD with license	HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A HeartStart HS1: M5066A, M5067A HeartStart FRx: 861304, 861305
ACT-IR CE	Infrared Data Cable	1 cable	HeartStart HS1: M5066A, M5067A, M5068A, HeartStart FRx: 861304, 861305 HeartStart Event Review Pro SW: 861276

# AED Training Materials

## AED Training Materials

Product #	Description	Packaging	Use With:
M5088A	Internal Manikin Adapter	1 set	HeartStart Trainer: M5085A and M5086A (for CPR/AED manikins) HeartStart HS1: M5073A Adult Training Pads Cartridge and M5074A Infant/Child Training Pads Cartridge
M5089A CE	External Manikin Adapter	1 set	HeartStart Trainer: M5085A and M5086A (for CPR/AED manikins) HeartStart HS1: M5073A Adult Training Pads Cartridge and M5074A Infant/Child Training Pads Cartridge
M5093A CE	Replacement Adult Training Pads with Replaceable Seal	1 set	HeartStart HS1: M5073A (Adult Training Pads Cartridge)
M5094A CE	Replacement Infant/Child Training Pads	1 set	HeartStart HS1: M5074A (Infant/Child Training Pads Cartridge)
M3855A CE	Charger for FR2 Training & Admin Pack	1 charger	M3864A (Training and Admin Pack for FR2 series)
M3864A CE	Training & Administration Pack for FR2 series Defibrillator	1 set	HeartStart FR2 and FR2+: M3860A, M3861A, M3840A, M3841A
861306 CE	HeartStart FRx Trainer	Trainer, case, 1 set pads, 1 external manikin adapter 1 IFU	989803139311 (Infant/child key) M5088A (Internal manikin adapter) M5089A (External Manikin Adapter) 989803139271 (HeartStart Training Pads II) 989803139291 (Replacement training pads II) M5090A (Adult pad placement guide) 989803139281 (Infant/child pad placement guide)
989803139271 CE	Training Pads II	1 cartridge	HeartStart FRx: 861304, 861305 HeartStart FRx Trainer: 861306



M3864A



861306



989803139271




**AED Training Materials** (continued)

Product #	Description	Packaging	Use With:
989803139281 CE	Infant/Child Pad Placement Guide	1 guide	HeartStart FRx Trainer: 861306 HeartStart Trainer: M5085A AED Trainer 2: M3752A
989803139321	HeartStart FRx Defibrillator Training Toolkit	1 kit (CD & DVD)	HeartStart FRx: 861304, 861305
989803150291	HeartStart FR2 Defibrillator Training Toolkit	1 kit (CD & DVD)	HeartStart FR2+: M3860A, M3861A, M3840A, M3841A
07-10900 CE	Defibrillator Training Pads	1 set	M3752A M3840A, M3841A, M3860A and M3861A when equipped with Training and Administration Pack (M3864A)
07-11000 CE	Defibrillator Trainer 2 Carrying Case	1 case	M3752A (AED Trainer 2)
M3752A CE	AED Trainer 2	FR2 AED trainer, 1-set of 07-10900, quick ref. card	07-10900, 07-11000, M3753A, M3754A, M3755A, M3756A, M3871A FR2 Infant/Child Training Pads
M3753A CE	Remote Control for AED Trainer 2	1 remote	M3752A (AED Trainer 2) M3756A (AED Little Anne Training System)
M3754A CE	Programming Kit for AED Trainer 2	CD-Rom training SW, serial port connector cable	M3752A (AED Trainer 2) M3756A (AED Little Anne Training System)
M3755A CE	AED Little Anne Training Pads	1 set	M3752A (AED Trainer 2) M3756A (AED Little Anne Training System)
M3756A CE (LATEX)	AED Little Anne Training System	Contains: AED Trainer 2, 1 remote, 1 AED Little Anne training manikin	M3755A Training pads for AED Little Anne
M3860-97800	Quick Reference Card, FR2 Defibrillator, English	1 card	FR2: M3860A, M3861A, M3868A, M3840A
M3871A	FR2 Infant/Child Training Pads	1 set	FR2: M3860A, M3861A, M3840A, M3841A when used with M3864A (Training and Admin Pack) M3752A (AED Trainer 2)



07-10900



M3871A

**AED Training Materials** (continued)

Product #	Description	Packaging	Use With:
M5066-89100	Training Toolkit	1 CD & DVD	HeartStart HS1: M5066A, M5067A
M5066-97800	Quick Reference Guide, OnSite	1 guide	HeartStart HS1: M5066A, M5067A, M5068A
M5073A CE	Adult Training Pads Cartridge	1 cartridge	HeartStart HS1: M5066A, M5067A, M5068A HeartStart Trainer: M5085A
M5074A CE	Infant/Child Training Pads Cartridge	1 cartridge	HeartStart HS1: M5066A, M5067A, M5068A HeartStart Trainer: M5085A
M5085A CE	HeartStart Trainer	1 case, one set of reusable adult training pads, one external manikin adapter	M5073A, M5074A, M5085-91900, M5087A, M5090A, 989803139281
989803139291	Replacement Training Pads II	1 set	HeartStart FRx Trainer: 861306
M5090A CE	Adult Pad Placement Guide	1 guide/per package	M5066A, M5068A
989803173221	HS1 Trainer, Cantonese	1 guide	M5085A
989803173231	HS1 Trainer, Traditional Chinese	1 guide	M5085A
989803173241	HS1 Trainer, Simplified Chinese	1 guide	M5085A
989803173341	HS1 Trainer, Korean	1 guide	M5085A
989803173461	FRx Trainer, Cantonese	1 guide	861306
989803173471	FRx Trainer, Traditional Chinese	1 guide	861306
989803173481	FRx Trainer, Simplified Chinese	1 guide	861306
989803173561	FRx Trainer, Korean	1 guide	861306



M5073A



M5074A



M5085A

# AED Batteries

## External AED Batteries

Product #	Description	Packaging	Use With:
989803136291 CE	FR2 Series Aviation Battery	1 battery	HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A
989803136301 CE	ForeRunner Aviation Battery	1 battery	ForeRunner AED: E, S, EM HeartStart FR AED: 940010XX, 940020XX, 94030XX
989803139301 CE	FRx Aviation Battery	1 battery	HeartStart FRx: 861304, 861305
BT1 CE	ForeRunner Disposable LitMnO <sub>2</sub> Battery Pack	1 battery	ForeRunner AED: E, S, EM HeartStart FR AED: 940010XX, 940020XX, 94030XX
M3848A CE	FR2+ Rechargeable Li Ion Battery	1 battery	HeartStart FR2+ Series: M3840A, M3841A, M3860A, M3861A, M3849A
M3849A CE	Charger for the Rechargeable Li Ion Battery	1 charger	HeartStart FR2+ Series: M3840A, M3841A, M3860A, M3861A, M3848A
M3863A CE	Long-Life LiMnO <sub>2</sub> Battery	1 battery	HeartStart FR2+ Series: M3840A, M3841A, M3860A, M3861A
M5070A CE	HeartStart Four-Year Battery	1 battery	HeartStart HS1: M5066A, M5067A, M5068A FRx: 861304, 861305
989803167281 CE	Lithium Ion Battery Module	1 Pack = 1 Battery	HeartStart XL+ Defibrillator/Monitor 861290



M3848A



M3863A



M5070A

## AED Cases/Mounts/Cabinets

### AED Wall Mounts/Accessories

Product #	Description	Packaging	Use With:
M3857A CE	Wall Mount Bracket	1 mount	HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM
M3859A CE	Secure-Pull Seal	1 box of 10	HeartStart FRx: 861304, 861305 HeartStart HS1: M5066A ForeRunner: E, S, EM HeartStart FR series: 940010XX, 940020XX, 94030XX HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A Wall Bracket: M3857A (Wall Mount Bracket) Cases: M3868A, M3869A, M5075A, M5076A, YC Cabinets: PFE7023D, PFE7024D
989803170891	AED Wall Mount Bracket	1 unit	FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282
989803170921	AED Wall Sign, red	1 sign	FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282
989803170931	AED Wall Sign, green	1 sign	FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282
989803170901	AED Awareness Placard, red	1 placard	FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282
989803170911	AED Awareness Placard, green	1 placard	FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282



M3857A



M3859A

**AED Cases/Wall Cabinets**

Product #	Description	Packaging	Use With:
989803133171	HeartStart Temperature Control Carrying Case	1 case	HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A
989803136531 CE	Defibrillator Cabinet (basic)	1 cabinet	HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM
989803139251 CE	FRx Carrying Case	1 case	HeartStart FRx: 861304, 861305
HC	ForeRunner Semi-Rigid Carrying Case	1 case	ForeRunner: E, S, EM HeartStart FR series: 940010XX, 940020XX, 94030XX
M3868A	Carrying Case for FR2 series	1 case	HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A
M3869A	Vinyl Carrying Case for FR2 Series	1 case	HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A
M5075A CE	Standard Carry Case for HeartStart HS1	1 case	HeartStart HS1: M5066A, M5067A, M5068A
M5076A CE	Slim carrying case for HeartStart HS1	1 case	HeartStart HS1: M5066A, M5067A, M5068A
M5087A	Replacement Carry Case, HeartStart Trainer	1 case	HeartStart Trainer: M5085A



M3868A



M5076A

**AED Cases/Wall Cabinets** (continued)

Product #	Description	Packaging	Use With:
PFE7023D CE	Defibrillator Cabinet (semi-recessed)	1 cabinet	HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM
PFE7024D CE	Defibrillator Cabinet (surface mount)	1 cabinet	HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM
YC	Carrying Case, Plastic Waterproof Shell	1 case	HeartStart FRx: 861304, 861305 ForeRunner: E, S, EM HeartStart FR series: 940010XX, 940020XX, 94030XX HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A, M3854A HeartStart HS1: M5066A, M5067A


# Index

05-10000.....	8	989803137771 .....	54	989803173341 .....	58
05-10000.....	53	989803138171 .....	26	989803173461 .....	58
05-10100.....	8	989803138181 .....	26	989803173471 .....	58
05-10100.....	53	989803139251 .....	61	989803173481 .....	58
05-10200.....	8	989803139261 .....	52	989803173561 .....	58
		989803139271 .....	56	989803174201 .....	23
07-10900.....	57	989803139281 .....	57	989803174211.....	22
07-11000.....	57	989803139291 .....	58	989803174211.....	24
		989803139301 .....	59	989803176161.....	20
68-PCHAT .....	54	989803139311 .....	54	989803176171.....	20
		989803139321 .....	57	989803176181.....	20
21075A .....	51	989803141811 .....	55	989803176191.....	20
21076A .....	51	989803147691 .....	19	989803176541 .....	13
21078A .....	51	989803147711 .....	11	989803177901 .....	27
21082A .....	50	989803148821 .....	25	989803177911 .....	30
21082B .....	50	989803150291 .....	57	989803177921 .....	33
21090A .....	49	989803158061 .....	19	989803177931 .....	34
21091A .....	49	989803158071 .....	19	989803177941 .....	34
21093A .....	49	989803158081 .....	19	989803177951 .....	16
21094A .....	49	989803158091 .....	19	989803177961 .....	16
21095A .....	49	989803158211 .....	52	989803177971 .....	16
21096A .....	50	989803158221 .....	52	989803177981 .....	16
21097A .....	50	989803158661 .....	12	989803178001 .....	16
		989803162401 .....	12	989803178011 .....	16
40401A .....	36	989803163291 .....	12	989803178021 .....	16
40401B .....	36	989803166021 .....	7	989803178031 .....	16
40401C .....	36	989803167281 .....	59	989803178041 .....	16
40401D .....	36	989803170891 .....	60	989803178051 .....	16
40401E .....	36	989803170901 .....	60	989803178071 .....	16
40457C .....	26	989803170911 .....	60	989803179101 .....	16
40457D .....	26	989803170921 .....	60	989803179111.....	16
		989803170931 .....	60	989803179121 .....	16
861306 .....	56	989803171261.....	11	989803179721 .....	27
861476.....	54	989803171271 .....	8	989803179731 .....	30
861477.....	54	989803171281 .....	15	989803179741 .....	29
861478.....	54	989803171291 .....	15	989803179751 .....	29
		989803171701 .....	13	989803179761 .....	29
453564063841 .....	13	989803172501 .....	15	989803179771 .....	27
453564063851 .....	13	989803173121 .....	23	989803179781 .....	29
		989803173131 .....	21	989803179791 .....	29
989803133171 .....	61	989803173131 .....	24	989803179801 .....	31
989803136291 .....	59	989803173221 .....	58	989803179811 .....	31
989803136301 .....	59	989803173231 .....	58	989803179861 .....	29
989803136531 .....	61	989803173241 .....	58	989803179871 .....	28

989803179891	.32	M1191T	.46	M1968A	.21
989803179901	.34	M1192A	.46	M1968A	.24
989803179911	.34	M1192T	.47	M1971A	.22
989803179921	.34	M1194A	.47	M1971A	.24
989803179931	.34	M1195A	.47	M1973A	.21
989803179941	.33	M1196A	.47	M1973A	.24
989803179951	.33	M1196S	.47	M1974A	.22
989803180851	.32	M1196T	.47	M1974A	.24
989803180861	.34	M1572A	.37	M1976A	.21
989803180871	.15	M1573A	.37	M1978A	.22
989803181141	.28	M1574A	.37	M1979A	.21
989803181211	.27	M1575A	.37	M1984A	.22
989803181221	.27	M1598B	.45	M2202A	.25
989803181231	.27	M1599B	.45	M2255A	.50
989803181241	.30	M1602A	.21	M2520A	.17
989803181251	.30	M1604A	.22	M2522A	.17
989803182281	.42	M1627A	.47	M2524A	.17
989803182291	.42	M1644A	.21	M2526A	.17
989803182301	.42	M1644A	.24	M3501A	.7
989803182311	.42	M1645A	.22	M3504A	.7
989803182321	.42	M1645A	.24	M3506A	.14
989803183311	.41	M1663A	.20	M3507A	.8
989803183321	.41	M1668A	.24	M3508A	.8
989803183331	.41	M1669A	.23	M3516A	.14
989803183341	.41	M1671A	.23	M3524A	.55
989803183351	.41	M1672A	.23	M3525A	.18
989803183361	.41	M1673A	.23	M3526A	.18
989803183371	.41	M1674A	.23	M3527A	.18
		M1675A	.23	M3528A	.18
ACT-IR	.55	M1678A	.23	M3529A	.18
		M1741A	.9	M3537A	.13
BT1	.59	M1742A	.9	M3538A	.14
		M1743A	.9	M3541A	.15
CPJ840J6	.35	M1744A	.9	M3543A	.9
CPJ84022	.35	M1781A	.8	M3544A	.11
CPJ84026	.35	M1783A	.10	M3549A	.13
CPJ84046	.35	M1789A	.9	M3713A	.7
CPJ84447	.35	M1837A	.49	M3716A	.7
		M1920A	.16	M3717A	.7
HC	.61	M1921A	.16	M3718A	.7
		M1923A	.17	M3719A	.7
M1131A	.46	M1941A	.48	M3725A	.8
M1191B	.46	M1943A	.48	M3752A	.57
M1191BL	.46	M1949A	.20	M3753A	.57



M3754A	.57	M4606A	.25	PFE7023D	.62
M3755A	.57	M4612A	.25	PFE7024D	.62
M3756A	.57	M4613A	.25		
M3834A	.55	M4737A	.13	SDCF05	.55
M3848A	.59	M4740A	.10		
M3849A	.59	M4741A	.10	YC	.62
M3853A	.55	M4742A	.10		
M3854A	.55	M4743A	.10		
M3855A	.56	M4744A	.10		
M3857A	.60	M4745A	.9		
M3859A	.60	M4746A	.9		
M3860-97800	.57	M4747A	.14		
M3863A	.59	M4748A	.8		
M3864A	.56	M4748A	.10		
M3868A	.61	M4751A	.15		
M3869A	.61	M4759A	.9		
M3870A	.53	M4761A	.12		
M3871A	.57	M4762A	.12		
M3873A	.54	M4763A	.12		
M3874A	.54	M5066-89100	.58		
M4552B	.38	M5066-97800	.58		
M4552B5	.38	M5070A	.59		
M4553B	.38	M5071A	.52		
M4553B5	.38	M5072A	.52		
M4554B	.39	M5073A	.58		
M4554B5	.39	M5074A	.58		
M4555B	.39	M5075A	.61		
M4555B5	.39	M5076A	.61		
M4556B	.39	M5085A	.58		
M4556B5	.39	M5087A	.61		
M4557B	.39	M5088A	.56		
M4557B5	.39	M5089A	.56		
M4558B	.39	M5090A	.58		
M4558B5	.40	M5093A	.56		
M4559B	.40	M5094A	.56		
M4559B5	.40	M5521A	.13		
M4572B	.43	M5522A	.13		
M4573B	.43	M5523A	.13		
M4574B	.43	M5524A	.13		
M4575B	.43	M5525A	.13		
M4576B	.43	M5526A	.10		
M4577B	.43	M5528A	.13		
M4578B	.44	M5529A	.13		
M4579B	.44	M5530A	.18		



© 2013 Koninklijke Philips N.V. All rights are reserved.  
Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands  
4522 962 95491 \* JUN 2013

**How to reach us**  
[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

**Product information**  
[www.philips.com/healthcarestore](http://www.philips.com/healthcarestore)