## Federal Communications Commission (FCC) Compliance Statement



Crossfire™ Integrated Arthroscopy System

**REF** 0475-000-000

## **FCC ID: SSH-XFC**

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

**Note:** FCC regulations provide that changes or modifications not expressly approved by Stryker Endoscopy could void your authority to operate this equipment.

Frequency of transmission: 13.56MHz

Type of frequency /

characteristics of 10% ASK

the modulation:

**Subcarrier:** 423.75kHz

Manchester coding

**Effective radiated power:** 50µW

