

SonoRemote Control User Guide Errata P12503-01

WARNING: To avoid the risk of electrical shock, do not use the SonoRemote control to control the ultrasound system when connected by USB cable. Use the USB cable only to charge the battery.

Please disregard references in the *SonoRemote Control User Guide* to controlling the ultrasound system when connected by USB cable. SonoSite does not currently support this configuration.

You can connect by USB cable to charge the battery.

For more information, please contact SonoSite or your SonoSite representative.

SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or
1-425-951-1200
F: 1-425-951-1201

SonoSite Ltd
Alexander House
40A Wilbury Way
Hitchin, Herts
SG4 OAP UK
T: +44-1462-444800
F: +44-1462-444801

03/2009
Copyright 2009 by
SonoSite, Inc.
All rights reserved.



P12503-01

