
# Respiratory Care Solutions across the continuum of care









# **Empower best practice**

Medline offers the depth and breadth of products to meet all of your needs with innovative solutions to enhance patient experience, improve outcomes and increase efficiency. Our solutions and products are clinically supported by our Respiratory Specialists and delivered by our unparalleled supply chain.

#### **Table of contents**

| lable of contents   |   |
|---|---|
| Airway management   |   |
| AG CUFFILL®<br>Endotracheal tubes, subglottic<br>Tracheostomy ties<br>Flex connector  | 5<br>6<br>7<br>7                                |
| Bronchial hygiene   |   |
| Volume incentive spirometers PEP therapy and inspiratory muscle   | 9   |
| Gas delivery & suction devices  |   |
| Air & oxygen flowmeters Medical gas hoses Oxygen regulators Suction regulators Suction aspirator  | 11<br>12<br>13<br>14–15<br>16                   |
| Humidification & filtration   |   |
| Bacterial/viral filters Bypass HME Heat & moisture exchangers Filtered heat & moisture exchangers Non-heated circuits Transport ventilator circuits Medline non-invasive ventilator circuits Pulmonary function testing filters & accessories | 18<br>19<br>20-21<br>21<br>22<br>22<br>22<br>23 |
| Medication delivery & asthma manage   | ment  |
| Medication nebulizers MDI spacers Peak flow meters Medication delivery accessories Masks, aerosol Nebulizer compressors   | 25-27<br>28-29<br>29<br>30<br>31<br>32-33       |
| Non-invasive interfaces   |   |
| Medline Flex™ non-invasive Ventilation masks  Medline Flex™ non-invasive ventilation masks accessories  DreamPort™ CPAP mask solution  CPAP masks AlRgel®  CPAP masks silicone  CPAP accessories  Go-PAP™ systems                             | 35<br>36<br>37<br>38–39<br>40<br>41<br>42       |
| O <sub>2</sub> -Max <sup>™</sup> systems  | 42-43   |

| O <sub>2</sub> & CO <sub>2</sub> disposables  |  |
|---|--|
| Southmedic OxyMask™ Cannulas, O₂ Cannulas, CO₂ Tubing, CO₂ Masks, CO₂ Tubing, O₂ Masks, O₂ Masks, o₂ Masks, aerosol Masks, tent & trach Bubble humidifiers O₂ accessories   | 45-46<br>47-48<br>49<br>50<br>51<br>52<br>53<br>54<br>54<br>55 |
| Pulse oximeters   |  |
| Fingertip pulse oximeters Handheld oximeters Pulse oximeter probes Pulse oximeter accessories   | 57<br>58<br>58<br>59   |
| Resuscitation   |  |
|   |  |
| Resuscitation bags, self inflating<br>Colorimetric CO <sub>2</sub> detectors<br>Resuscitation masks & accessories<br>Mouth-to-mask resuscitator<br>Hyperinflation bags  | 61–67<br>68<br>69<br>69<br>70                                  |
| Colorimetric CO <sub>2</sub> detectors<br>Resuscitation masks & accessories<br>Mouth-to-mask resuscitator   | 68<br>69<br>69   |
| Colorimetric CO <sub>2</sub> detectors Resuscitation masks & accessories Mouth-to-mask resuscitator Hyperinflation bags  Suction catheters & trach care  Closed suction Open suction catheters Open suction catheters Open suction kits Open suction kits with solutions Sterile solutions Trach care & cleaning trays Mucus/specimen traps | 68<br>69<br>69   |
| Colorimetric CO <sub>2</sub> detectors Resuscitation masks & accessories Mouth-to-mask resuscitator Hyperinflation bags  Suction catheters & trach care  Closed suction Open suction catheters Open suction catheters, sleeved Open suction kits Open suction kits with solutions Sterile solutions Trach care & cleaning trays             | 68<br>69<br>69<br>70<br>72–74<br>75<br>75<br>76<br>77<br>77    |

# Airway management

| AG CUFFILL®                    | اِ |
|--------------------------------|----|
| Endotracheal tubes, subglottic | (  |
| Tracheostomy ties              | -  |
| Flex connector                 | -  |

# AG CUFFILL® cuff pressure manometer

Now exclusive to Medline

**Precise airway cuff pressure is critical to minimizing complications associated with intubation.** Studies show that conventional methods caregivers use for cuff pressure readings are unreliable.¹ In fact, 52% of intubated patients have been found to be out of therapeutic cuff pressure range, putting them at risk for VAE or tracheal damage.²

**AG CUFFILL can help.** The pocket-sized and easy-to-use AG CUFFILL device is an accurate solution for both controlling cuff volume and measuring cuff pressure for cuffed ET tubes, trach tubes and LMAs.

**Simple:** Easily measure and adjust cuff pressure and volume of cuffed airways


**Accurate:** Digital volume and pressure readout, +/-2 (cmH20) accuracy with simple-one-handed calibration

Safe: Designed for single-patient use to help reduce the chances of cross-contamination risk

Accessible: Pocket-sized, light-weight, and stays at the patient's bedside

Cost effective: Priced to make regular cuff pressure checks and adjustments a new standard of care

| Item No.     | Description              | Pkg.  |
|--------------|--------------------------|-------|
| HCSCUFF0041  | AG Cuffill device, 10 cc | 10/bx |
| HCSCUFF0041H | AG Cuffill device. 10 cc | 1/ea  |






Our SubG subglottic endotracheal tubes are designed to help reduce the risk of nosocomial infections due to aspiration of oral and upper airway secretions.

- Polyurethane cuff helps to prevent channel formation
- High volume, low pressure barrel-shaped cuff for optimal sealing
- Kink-resistant tubing with easy-to-read insertion markings
- Subglottic suction port for removing secretions above the cuff
- Suction line cap keeps contaminants out of lumen

| Item No.  | Description              | Pkg.  |
|-----------|--------------------------|-------|
| DYNJ18860 | SubG ET tube, size # 6.0 | 10/bx |
| DYNJ18865 | SubG ET tube, size # 6.5 | 10/bx |
| DYNJ18870 | SubG ET tube, size # 7.0 | 10/bx |
| DYNJ18875 | SubG ET tube, size # 7.5 | 10/bx |
| DYNJ18880 | SubG ET tube, size # 8.0 | 10/bx |
| DYNJ18885 | SubG ET tube, size # 8.5 | 10/bx |

#### Advantages of the SubG ET tube

- **A.** Subglottic suction line and evacuation port for removal of pooled secretions above the cuff
- **B.** Barrel-shaped cuff made of polyurethane to aid in the reduction of microaspirations
- **C.** Thermosensitive material, a bullnose tip, and easy to read marking for ease in insertion





#### **Tracheostomy ties**

These soft, comfortable straps have adjustable hook-and-loop closures to keep tracheostomy tubes securely in place, ensuring uninterrupted airflow.

#### Premium tracheostomy ties

- Soft, no stretch cotton ensures patient safety and comfort while reducing skin irritation
- Offered in a wide range of sizes

#### 2-piece

| Item No. | Description                         | Pkg.  |
|----------|-------------------------------------|-------|
| HCS540   | Adult, fits up to 20" neck          | 10/bx |
| HCS541   | Pediatric, fits up to 18" neck      | 10/bx |
| HCS542   | Neonatal/infant, fits up to 9" neck | 10/bx |

#### 1-piece

| Item No. | Description                        | Pkg.  |
|----------|------------------------------------|-------|
| HCS597LA | Large adult, fits 13-19" neck      | 10/bx |
| HCS597M  | Small adult, fits 9–14" neck       | 10/bx |
| HCS597S  | Neonatal/pediatric, fits 7-9" neck | 10/bx |

#### Standard Medline tracheostomy ties

• Available in 2 sizes that fit a broad range of patients

#### 2-Piece

| Item No. | Description          | Pkg.  |
|----------|----------------------|-------|
| ZPP501   | Adult, 18¾" long     | 20/bx |
| ZPP501P  | Pediatric, 13½" long | 20/bx |

#### 1-piece

| Item No. | Description         | Pkg.  |
|----------|---------------------|-------|
| ZPP301   | Adult, 15" long     | 20/bx |
| ZPP301P  | Pediatric, 11" long | 20/bx |
| ZPP401   | Anti-disconnect     | 20/bx |

#### **Flex connector**

- Flexible connector is used on a ventilator
- · Connects to tubing

| Item No. | Description               | Pkg.  |
|----------|---------------------------|-------|
| HCS3222  | 15 x 22 mm flex connector | 50/cs |
| HCS3215  | 15 x 15 mm flex connector | 50/cs |



HCS540





# **Bronchial hygiene**

| 9 |
|---|
|   |

| PEP therapy an | d inspiratory musc | le 9 |
|----------------|--------------------|------|
|----------------|--------------------|------|

#### **Volume incentive spirometers**

Our volume incentive spirometer uses an intuitive design to help patients correctly perform and monitor their post-surgical breathing exercises.

- One way valve ensures patients inhale rather than exhale through the device
- Easy to use floating ball indicator and guidelines help patients inhale at the proper rate
- Sliding goal indicator can be easily adjusted allowing patients to monitor their progress
- Flexible tubing and dual sided calibrations allow use in various positions
- Available in 2,500 mL and 5,000 mL capacities

| Item No. | Description                   | Pkg.  |
|----------|-------------------------------|-------|
| HCS93250 | Incentive spirometer 2,500 mL | 12/cs |
| HCS93500 | Incentive spirometer 5,000 mL | 12/cs |

#### PEP therapy and inspiratory muscle

#### Threshold IMT

- · Provides consistent and specific pressure for inspiratory muscle strength and endurance training, regardless of how quickly or slowly patients breathe
- · Incorporates a flow-independent, one-way valve to provide consistent resistance and features an adjustable specific pressure setting (in cm H2O) to be set by a healthcare professional
- · When patients inhale through Threshold IMT, a spring-loaded valve provides resistance that exercises respiratory muscles through conditioning

| Item No. | Description                                | Pkg.  |
|----------|--|-------|
| HS730ML  | Threshold inspiratory muscle trainer (IMT) | 10/cs |

#### Threshold PFP

- Used for airway clearance, bronchial hygiene, or as an alternative to chest physical therapy
- Incorporates a flow-independent, one-way valve to ensure consistent resistance with adjustable specific pressure settings (in cm H2O) which are set by the healthcare professional
- When patients exhale through Threshold PEP, the resistive load creates positive pressure that helps open the airways and allows mucus to be expelled during "huff" coughing (forced expiratory technique)

| Item No. | Description   | Pkg.  |
|----------|---|-------|
| HS735ML  | Threshold positive expiratory pressure device (PEP) | 10/cs |



HCS93250



HS735ML

# Gas delivery & suction devices

| Air & oxygen flowmeters | 11    |
|-------------------------|-------|
| Medical gas hoses       | 12    |
| Oxygen regulators       | 13    |
| Suction regulators      | 14-15 |
| Suction aspirator       | 16    |

#### Air & oxygen flowmeters

#### Features:

- · Complete line of flowmeters from neonatal to adult
- Available in standard Thorpe style and Dial style
- MRI options are available
- One of the smallest and lightest brass bodied flowmeters available
- · Color coded inner tube and labeling



Dial/Click Style Flowmeter



Dual/Y-block Flowmeter Assembly



Thorpe Style



#### Accessories



Swivels/Christmas Trees HCS4700G (Green) HCS4700Y (Yellow)



Power Takeoff



#### Flowmeter ordering information RTFMXXXXX(YY)\*

#### Gas

O = Oxygen

#### LPM

A = Medical Air

1L = 1 LPM 3L = 3.5 LPM

2X = 200 CC

4L = 4 LPM (Dial only)

8L = 8 LPM

15 = 15 LPM

25 = 25 LPM (Dial only)

70 = 70 LPM

For dual Y-block flowmeter configurations: Replace RTFM prefix with RFMD

#### **Connection type**

AA = 1/8" FNPT

DH = DISS Handtight

DN = DISS Nut

OM = Ohmeda Male

CM = Chemetron Male

PB = Puritan Bennett

#### **Options**

PT = Power Takeoff

DI = Dial

M = MRI Compatible\*

\*Some restrictions apply

#### **Medical gas hose assemblies**

Medline offers custom-configured medical gas hose assemblies in four easy steps:

- 1/4" ID Non-Conductive
- 1. Select hose type
- 2. Select the inlet fitting (gas service and type)
- 3. Select the outlet fitting (gas service and type)
- 4. Specify the hose length

#### RT = Hose Assembly

#### 1. Hose type

OXY = Oxygen N2 = Nitrogen VAC = Vacuum CO<sub>2</sub> = Carbon Dioxide

N2O = Nitrous Oxide WAG = Waste Anesthetic Gas Disposal

AIR = Medical Air

#### 2. Inlet fitting (Connects to the Source)

OM: OHIO Male OF: OHIO Female Coupler DC: DISS Male with Check DN: DISS Nut/Nipple Female CM: Chemetron Male DH: DISS Handtight Female PB: Puritan-Bennett Male Coupler CF: Chemetron Female Coupler PF: Puritan-Bennett Female

#### 3. Outlet fitting (Connects to the Medical Device)


OM: OHIO Male OF: OHIO Female Coupler DC: DISS Male with Check DN: DISS Nut/Nipple Female CM: Chemetron Male DH: DISS Handtight Female PB: Puritan-Bennett Male Coupler CF: Chemetron Female Coupler PF: Puritan-Bennett Female

**4. Hose length,** in feet (e.g. 2 = 2-ft., 4 = 4-ft., 7 = 7-ft., etc.)

#### **EXAMPLE:**

RTOXYDHDC5 1/4" ID non-conductive oxygen hose assembly with a DISS handtight Female inlet and DISS Male outlet, 5' in length.

For custom configurations, contact your Medline Representative.



#### **Oxygen regulators**

For all medical and emergency applications.

- Fitted with brass sleeves around key parts
- Laser-etched knob and body

#### **Economy**

| Item No. | Description        | Flow range | Flow increments                          | Pkg. |
|----------|--------------------|------------|--|------|
| HCS8708M | 870 CGA connection | 0-8 LPM    | 0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8 | 1/ea |
| HCS8715M | 870 CGA connection | 0-15 LPM   | 0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15 | 1/ea |
| HCS8725M | 870 CGA connection | 0-25 LPM   | 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25     | 1/ea |
| HCS5408M | 540 CGA connection | 0-8 LPM    | 0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8 | 1/ea |
| HCS5415M | 540 CGA connection | 0-15 LPM   | 0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15 | 1/ea |



HCS8715M HCS5415M

#### **Super Mini**

| Item No.  | Description    | Flow range | Flow increments                          | Pkg. |
|-----------|----------------|------------|--|------|
| HCS8715MA | 870 connection | 0-15 LPM   | 0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15 | 1/ea |

#### **Accessories**

| Item No.  | Description   | Pkg.  |
|-----------|---|-------|
| HCS870W   | Oxygen cylinder wrench, small (CGA-870)                     | 10/bg |
| HCS870WC  | Oxygen cylinder wrench with security chain, small (CGA-870) | 10/bg |
| HCS870WLG | Oxygen cylinder wrench, large (CGA-870)                     | 25/bg |
| HCS166840 | Brass yoke washer seal for CGA-870 oxygen regulator         | 50/bg |



#### **EMS** regulators

Easily transport ventilator patients.

- Maximum inlet pressure: 3,000 PSI
- Chrome plated brass body
- Solid brass core

| Item No.     | 50 PSI power take-off | Hose barb | Liters per minute                       | Cylinder<br>style | Pin index | Pkg. |
|--------------|-----------------------|-----------|---|-------------------|-----------|------|
| CPRB34501C   | 1                     |           |   | CGA346            | Air       | 1/ea |
| CPRB87251B   | 0                     | •         | 0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25 | CGA870            | Oxygen    | 1/ea |
| CPRB87251B2C | 2                     | •         | 0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25 | CGA870            | Oxygen    | 1/ea |
| CPRB54251B2C | 2                     | •         | 0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25 | CGA540            | Oxygen    | 1/ea |
| CPRB54402C   | 2                     |           |   | CGA540            | Oxygen    | 1/ea |
| CPRB54501C   | 1                     |           |   | CGA540            | Oxygen    | 1/ea |



For specialty gas delivery products including blenders, analyzers, and more, please contact your  $\label{lem:mediane-sales-Representative} \ \text{or email } \textbf{ProductSupportRespiratory@Medline.com} \ \text{with your request.}$ 

## The choice is clear

With a unique transparent and disposable back, Medline regulators enhance efficiency, reduce costs and improve patient safety.

#### **Identify contamination instantly**

Up to 37% of suction regulators in hospitals are contaminated due to fluid breaches, and pathogens can spread back through the suction line within 48 hours. The Medline Suction Regulator has a clear back plate offering a view of the internal components of the regulator, so it's easy to see when it's contaminated.

1. Suction Regulators: A Potential Vector for Hospital Acquired Pathogens by Keith S. Kaye, MD, MPH; Dror Marchaim, MD; Chester Smialowicz, MD; Lauren Bentley, MSBME

#### There's no time for downtime

The clear back plate provides an internal view of the regulator, making it possible to see fluid breaches. Plus, cleaning is simple: the unit is easily taken apart, allowing you to discard the contaminated back plate and replace with a new component in minutes.



Step 1: Remove contaminated regulator



Step 2: Dispose



Step 3: Reassemble with only 1 NEW component in just minutes



#### Clear back plate

INFECTION CONTROL MEASURE

Remove from the outlet and identify contamination immediately

For specialty oxygen regulators (EMS, all-brass, etc.) please contact your Medline Sales Representative or email ProductSupportRespiratory@Medline.com with your request. **Suction regulators** 

#### **Enhance efficiency, reduce costs** and improve patient safety.

#### **Continuous models**

Continuous 3-mode (0-300mmHg) Pediatric 2-mode (0-160mmHg) Continuous high/surgical (0-760mmHg)

#### **Continuous/intermittent models**

Adult continous/intermittent (0-300mmHg) Pediatric continous/intermittent (0-160mmHg) Timing cycle: 16 seconds on/8 seconds off, can be adjusted on the back of the unit





Analog

**Inlet connection** 

DH = DISS Handtight

OM = Ohmeda Male

CM = Chemetron Male

PB = Puritan Bennett

AA = 1/8" FNPT

DN = DISS Nut

Digital

#### **Suction regulator ordering information**

#### RTXXXXX

#### **Display**

A = Analog D = Digital

#### Vacuum regulator type

- CI = Continuous/ Intermittent
- C3 = Continuous 3 Mode
- CH = Continuous High/ Surgical
- PI = Pediatric/Intermittent
- P2 = Pediatric 2 Mode

#### **Patient** connection

- A = 1/8" FNPT
- D = DISS Male
- B = Tubing Nipple
- V = 1/8" NPT Trap
- H = DISS Handtight Trap
- N = DISS Nut/ Nipple Trap

#### **Vacuum traps**

| Item No. | Description            |
|----------|------------------------|
| RTTRAPMT | 1/8" Male NPT Fitting  |
| RTTRAPDH | Female DISS Handtight  |
| RTTRAPDN | Female DISS Nut/Nipple |
|          | _                      |



Vacuum/Safety Traps

| Patient fittings    | Examples of wall ada | pters          |        |                 |
|---------------------|----------------------|----------------|--------|-----------------|
| DISS Male Hose Barb | Chemetron            | DISS Handtight | Ohmeda | Puritan-Bennett |

For MRI compatible please add "-M" to the end of the part number (available for adult/analog units only). For color coding, please specify in sales notes. Default color is white.

#### Suction aspirator and disposable canister kits

Our easy-to-use aspirators are ultra-quiet.

- Weight: 12.1-lb. (5.49 kg)
- Maximum suction: -21 inHg/520 mmHg
- Rugged construction is ideal for bedside carts
- Disposable 850 cc or 1,200 cc suction canister sets with float lids and 850 cc suction canister sets with filter lids.
- HCS7121, HCS7851, and HCS7852 are kits that include suction tubing, 1/each 6' and 18" length.

| Item No. | Description                | Pkg. |
|----------|----------------------------|------|
| HCS7000  | Portable suction aspirator | 1/ea |



HCS7000

#### Canister with float lid


| Item No. | Description                                | Pkg.  |
|----------|--|-------|
| HCS7850  | Canister with float lid: 850 cc            | 12/cs |
| HCS7851  | Canister with float lid & tubing: 850 cc   | 12/cs |
| HCS7120  | Canister with float lid: 1,200 cc          | 12/cs |
| HCS7121  | Canister with float lid & tubing: 1,200 cc | 12/cs |


#### Canister with filter lid

| Item No. | Description                               | Pkg.   |
|----------|---|--------|
| HCS7852  | Canister with filter lid & tubing: 850 cc | 50/cs  |
| HCS7853  | Canister with filter lid: 850 cc          | 100/cs |

#### Accessories

| Item No.  | Description   | Pkg. |
|-----------|---|------|
| HCS4053   | Hydrophobic filter set, barb both ends  | 3/pk |
| HCS701311 | Replacement gauge assembly: includes pressure gauge, pressure adjustment knob & barb tubing adapter | 1/ea |











Tubing included with HCS7851, HCS7852, and HCS7121

# Humidification & filtration

| Bacterial/viral filters                  | 18   |
|--|------|
| Bypass HME                               | 19   |
| Heat & moisture exchangers               | 20-2 |
| Filtered heat & moisture exchangers      | 21   |
| Non-heated circuits                      | 21   |
| Transport ventilator circuits            | 22   |
| Medline non-invasive ventilator circuits | 22   |
| Pulmonary function testing               | 23   |

#### **Bacterial/viral filters**

Medline's hydrophobic filters can be used for both respiratory and anesthesia applications. They have translucent housings to aid in viewing built up secretions or other obstructions.

#### Performance specifications

|                          | DYNJAABV1              | DYNJAABV100            | DYNJAAHEPA  |
|--------------------------|------------------------|------------------------|-------------|
| Description              | Bacterial/viral filter | Bacterial/viral filter | HEPA Filter |
| Tidal volume (mL)        | ≥250                   | ≥250                   | 200-1,500   |
| BFE (%)                  | 99.999                 | 99.999                 | 99.9999     |
| VFE (%)                  | 99.99                  | 99.99                  | 99.9999     |
| Volume (mL)              | 34                     | 34                     | 52          |
| Weight (g)               | 15.1                   | 15.1                   | 40          |
| Flow resistance @ 30 LPM | 0.71                   | 0.71                   | 1.48        |
| Flow resistance @ 60 LPM | 1.02                   | 1.02                   | 3.33        |
| CO₂ port                 | No                     | No                     | Yes         |
| Filter type              | Electrostatic          | Electrostatic          | HEPA        |
| Pkg.                     | 50/cs                  | 100/cs                 | 50/cs       |



#### Bypass heat and moisture exchanger

The AARC recommends that ventilator circuits should not be routinely changed and that passive humidifiers can be used for up to 48 hours. 1 It is widely understood that the opening of a respiratory circuit increases the risk of hospital acquired infection. Delivering aerosolized medication creates a challenge in maintaining a closed circuit.

The Medline Bypass Heat and Moisture Exchanger supports AARC clinical recommendations by eliminating the need to open the circuit during aerosol medication delivery for up to 72 hours.

#### Features:

- · Easy and intuitive to use
- Clear marking of "HME" and "Aerosol" modes
- Use for up to 72 hours

#### Performance specifications

- Tidal volume 250-1500 mL
- Moisture output 32.09 mg  $H_2O/L @ 500$  mL TV at 72 hours
- Resistance to flow 0.6 cm  $H_2O$  @ 30 L/min flow
- Dead space 29.2 mL
- Weight 39.3 g
- Duration of use up to 72 hours

| Item No.   | Description                  | Pkg.  |
|------------|------------------------------|-------|
| DYNJAABHME | 15 mm ID/22 mm OD x 22 mm ID | 25/cs |





1. AARC Evidence-Based Clinical Practice Guidelines, Care of the Ventilator Circuit And Its Relation to Ventilator-Associated Pneumonia. Respiratory Care, September, 2003, Volume 48, No 9.

# Heat and moisture exchangers and filtered heat and moisture exchangers

Our broad portfolio of HMEs and HMEFs meet your patient's humidification needs throughout the care continuum.

- Multiple tidal volume ranges and sizes to meet your clinical needs
- Non-filtered, electrostatic-filtered, and HEPA-filtered options available
- Angled options available to reduce the need for additional adapters and dead space
- Translucent housing for viewing of possible obstructions













|                                | DYNJAAHME36              | DYNJAAHME10              | DYNJAAHME12              | DYNJAAHME14              | DYNJAPHME20              |
|--------------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Description                    | HME, adult               | HME, compact             | HME, compact             | HME, compact             | HME, pediatric           |
| Туре                           | Straight                 | Straight                 | Straight                 | Angled                   | Straight                 |
| Connection                     | 220D/15ID x<br>22ID/150D |
| Vt (mL)                        | 200-1500                 | 120-1500                 | 120-1500                 | 120-1500                 | 90-250                   |
| Vd (mL)                        | 90                       | 32.1                     | 32.1                     | 34.5                     | 31.3                     |
| Weight (g)                     | 22                       | 14.7                     | 14.7                     | 16.8                     | 13.3                     |
| CO₂ port                       | Yes                      | Yes                      | No                       | Yes                      | Yes                      |
| Resistance (cm H₂O)            | 0.71 @ 30 LPM            | 0.6 @ 30 LPM             | 0.6 @ 30 LPM             | 0.5 @ 30 LPM             | 0.4 @ 30 LPM             |
| Moisture output<br>(Mg H₂O/L)* | 30 @ 500 mL              | 33.7 @ 500 mL            | 33.7 @ 500 mL            | 33.2 @ 500 mL            | 34.0 @ 500 mL            |
| Pkg.                           | 50/cs                    | 50/cs                    | 50/cs                    | 50/cs                    | 50/cs                    |
| HME media                      | Foam                     | Foam                     | Foam                     | Foam                     | Foam                     |









|                                | DYNJAPHME38   | DYNJAPHME34   | DYNJAAHME30   | DYNJAAHME40                              |
|--------------------------------|---------------|---------------|---------------|--|
| Description                    | HME, infant   | HME, neonatal | HME, trach    | HME, trach with suction and oxygen ports |
| Туре                           | Straight      | Straight      | Trach         | Trach                                    |
| Connection                     | N/A           | 15ID x 150D   | 15 mm ID      | 15 mm ID                                 |
| Vt (mL)                        | 75-500        | 15 – 50       | > 25          | >50                                      |
| Vd (mL)                        | 11.5          | 3.1           | 14            | 13.4                                     |
| Weight (g)                     | 13            | 4.1           | 3.5           | 7.1                                      |
| CO <sub>2</sub> port           | Yes           | No            | No            | No                                       |
| Resistance (cm H₂O)            | 2 @ 15 LPM    | 1.05 @ 10 LPM | 0.51@30LPM    | 0.42 @ 30 LPM                            |
| Moisture output<br>(Mg H₂O/L)* | 28.5 @ 250 mL | 24.8 @ 50 mL  | 25.9 @ 500 mL | 26 @ 500 mL                              |
| Pkg.                           | 50/cs         | 50/cs         | 50/cs         | 50/cs                                    |
| HME media                      | Paper         | Paper         | Paper         | Paper                                    |











|                             | DYNJAAHME1                | DYNJAAHME5                | DYNJAAHME9                | DYNJAAHME11               |
|-----------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
| Description                 | Filtered HME, large       | Filtered HME, large       | Filtered HME, compact     | Filtered HME, compact     |
| Туре                        | Straight                  | Straight                  | Straight                  | Straight                  |
| Connection                  | 22 M / 15 F x 22 F / 15 M | 22 M / 15 F x 22 F / 15 M | 22 M / 15 F x 22 F / 15 M | 22 M / 15 F x 22 F / 15 M |
| Vt (mL)                     | 150-1500                  | 150-1500                  | 120-1500                  | 120-1500                  |
| Vd (mL)                     | 40.7                      | 40.7                      | 31.3                      | 31.3                      |
| Weight (g)                  | 18.5                      | 18.5                      | 15                        | 15                        |
| CO₂ port                    | Yes                       | No                        | Yes                       | No                        |
| Resistance (cm H₂O)         | 0.8 @ 30 LPM              | 0.8 @ 30 LPM              | 2.1 @ 30 LPM              | 2.1 @ 30 LPM              |
| Moisture output (Mg H₂O/L)* | 34.1 @ 500 mL             | 34.1 @ 500 mL             | 33.6 @ 500 mL             | 33.6 @ 500 mL             |
| Filter media                | Electrostatic             | Electrostatic             | Electrostatic             | Electrostatic             |
| BFE(%)                      | 99.999                    | 99.999                    | 99.999                    | 99.999                    |
| VFE(%)                      | 99.99                     | 99.99                     | 99.99                     | 99.99                     |
| Pkg.                        | 50/cs                     | 50/cs                     | 50/cs                     | 50/cs                     |
| HME media                   | Foam                      | Foam                      | Foam                      | Foam                      |









|                             | DYNJAAHME13               | DYNJAPHME17               | DYNJAAHMEHEPA             | DYNJ603011                |
|-----------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
|                             |                           |                           |                           |                           |
| Description                 | Filtered HME, compact     | Filtered HME, pediatric   | Filtered HME, HEPA        | Filtered HME, flat, adult |
| Туре                        | Angled                    | Straight                  | Straight                  | Straight                  |
| Connection                  | 22 M / 15 F x 22 F / 15 M | 22 M / 15 F x 22 F / 15 M | 22 M / 15 F x 22 F / 15 M | 22 ID x 22 OD / 15 ID     |
| Vt (mL)                     | 120-1500                  | 90-1500                   | 200-1500                  | 200-1500                  |
| Vd (mL)                     | 33.7                      | 20.5                      | 52                        | 62                        |
| Weight (g)                  | 17.2                      | 13.6                      | 40                        | 37.2                      |
| CO <sub>2</sub> port        | Yes                       | Yes                       | Yes                       | Yes                       |
| Resistance (cm H₂O)         | 1.9 @ 30 LPM              | 2.4 @ 30 LPM              | 1.5 @ 30 LPM              | <1 @ 30 LPM               |
| Moisture output (Mg H₂O/L)* | 33.4 @ 500 mL             | 34.0 @ 250 mL             | 30 @ 500 mL               | 32.5 @ 500 mL             |
| Filter media                | Electrostatic             | Electrostatic             | HEPA                      | Electrostatic             |
| BFE(%)                      | 99.999                    | 99.999                    | 99.9999                   | >99.99%                   |
| VFE(%)                      | 99.99                     | 99.99                     | 99.9999                   | >99.99%                   |
| Pkg.                        | 50/cs                     | 50/cs                     | 50/cs                     | 50/cs                     |
| HME media                   | Foam                      | Foam                      | Foam                      | Foam                      |

#### Non-heated breathing circuits

Our family of non-heated breathing circuits support your mechanically ventilated patients with several configurations, which includes "Y" and parallel "Y" options.

#### Adult non-heated circuits

| Item No.   | Description  | Pkg.               |
|------------|--|--------------------|
| DYNJ300    | Adult dual limb circuit, wye without ports, uncuffed, 6-ft.  | 20/cs              |
| HCS1830    | Adult dual limb circuit, wye without ports, 6-ft.  | 20/cs              |
| HCS72R7002 | Adult dual limb circuit, uncuffed, wye without ports, 6-ft.  | 20/cs              |
| HCS1831    | Adult dual limb circuit, wye with temperature sensing port, 6-ft.                                      | 20/cs              |
| HCS1831P   | Adult dual limb circuit, parallel wye with temperature sensing port, 6-ft.                             | 20/cs              |
| HCS1571215 | Adult dual limb circuit, uncuffed, parallel wye with temperature sensing port, bacterial filter, 5-ft. | <sup>2</sup> 20/cs |
| HCS72R7401 | Adult dual limb circuit, uncuffed, parallel wye without ports, HMEF, and bacterial filter, 6-ft.       | 20/cs              |



| Item No. | Description  | Pkg.  |
|----------|--|-------|
| DYNJ100  | Infant dual limb circuit, uncuffed, parallel wye without ports, 5-ft.  | 20/cs |
| DYNJ200  | Pediatric dual limb circuit, uncuffed, parallel wye without ports 5-ft | 20/cs |

#### MRI-conditional transport ventilator circuit

We offer ready-to-use, transport ventilator circuits that include expandable tubing and a MRI-conditional PEEP valve

- Tubing is expandable up to 120" in length
- Circuit and PEEP valve are MRI-conditional
- Patients can be scanned with the device in operation under the following conditions: static magnetic field up to 3.0 T; maximum spatial gradient field less than or equal to 30 T/m

| Item No.    | Description  | Pkg.  |
|-------------|--|-------|
| DYNJAA300   | Transport ventilator circuit with bacterial filter | 20/cs |
| DYNJAA300NF | Transport ventilator circuit                       | 20/cs |

#### **Medline Non-Invasive Ventilation Circuits**

Medline offers flexible, smooth bore, non-invasive ventilation (NIV) circuits that are available in two configurations to meet your needs. These circuits are compatible with either vented or non-vented NIV masks with an anti-asphyxia valve and include a removable exhalation port and pressure line. The filter has a BFE and VFE of 99.99% (if included).

| Item No.   | Description   | Pkg.  |
|------------|---|-------|
| DYNJ500EP  | Single limb, smooth bore circuit with exhalation port, 6 ft.          | 10/ea |
| DYNJ500EPF | Single limb, smooth bore circuit with exhalation port and filter 6 ft | 10/ea |



DYNJ300

Medline can customize a ventilator circuit to your needs, contact your Medline Representative.



DYNJAA300



### Pulmonary function testing (PFT) filters and accessories

We offer multiple types of pulmonary function filters and accessories to meet your clinical needs. Our PFT filters have a BFE of 99.98% and VFE of 99.92% and mate with various PFT machines. They are packaged separately or kitted with a nose clip and mouthpiece.





DYNJ756SB DYNJ756ELB

#### **PFT filters**

| Item No.   | Description  | Fits the following devices   | Pkg.   |
|------------|--|--|--------|
| DYNJ756BMN | PFT filter kit, standard round patient end (includes nose clip and mouthpiece) | Koko, Kokolegend, Keystone, Koko-Trek, Cybermedic  | 100/cs |
| DYNJ756ELB | PFT filter, elliptical patient end   | Full Gas Systems,  | 100/cs |
| DYNJ756NB  | PFT filter kit, standard round patient end                                     | Collins 4g and 3g Full Gas Systems   | 100/cs |
| DYNJ756SB  | PFT filter kit, small round patient end  |  | 100/cs |
| DYNJ757BMN | PFT filter kit, standard round patient end (includes nose clip and mouthpiece) | Sensormedics, 2200, 2400, 2450, 2800 Vmax,   | 100/cs |
| DYNJ757ELB | PFT filter, elliptical patient end   | Keystone I & II, Pall, Cybermedics , Moose,<br>Gould 5004, New Vyaire Masterscreen                               | 100/cs |
| DYNJ757NB  | PFT filter kit, standard round patient end                                     | (Launched April 2019)  | 100/cs |
| DYNJ757SB  | PFT filter kit, small round patient end  |  | 100/cs |
| DYNJ758B   | PFT filter kit, standard round patient end                                     |  | 100/cs |
| DYNJ758BMN | PFT filter kit, standard round patient end (includes nose clip and mouthpiece) | Collins DC II and HDPFT 4000 Plethysmograph  | 100/cs |
| DYNJ758ELB | PFT filter, elliptical patient end   |  | 100/cs |
| DYNJ759B   | PFT filter kit, standard round patient end                                     |  | 100/cs |
| DYNJ759BMN | PFT filter kit, standard round patient end (includes nose clip and mouthpiece) | Incheck, Incheck Dial, Multispiro,<br>Creative Biomedics, Vitalograph,<br>Schiller, Microdirect, SDI (Pulmogard) | 100/cs |
| DYNJ759ELB | PFT filter, elliptical patient end   | Sermer, Price our ecc, SDF (Full togard)   | 100/cs |
| DYNJ760B   | PFT filter kit, standard round patient end                                     |  | 100/cs |
| DYNJ760BMN | PFT filter kit, standard round patient end (includes nose clip and mouthpiece) | Brentwood, Medical Graphics All  | 100/cs |
| DYNJ760ELB | PFT filter, elliptical patient end   |  | 100/cs |
| DYNJ780B   | PFT filter kit, standard round patient end                                     |  | 100/cs |
| DYNJ780BMN | PFT filter kit, standard round patient end (includes nose clip and mouthpiece) | Jaeger, Vyaire Masterscreen Series<br>(Previous Versions)  | 100/cs |
| DYNJ780ELB | PFT filter, elliptical patient end   |  | 100/cs |

#### **Nose clips**

| Item No. | Description                     | Pkg.   |
|----------|---------------------------------|--------|
| DYNJ750  | Nose clip                       | 100/bx |
| DYNJ751  | Nose clip with marshmallow pads | 100/bx |
| DYNJ752  | Large nose clip                 | 50/bx  |
| DYNJ755  | Nose clip and mouthpiece kit    | 100/bx |

#### Mouthpiece

| Item No. | Description                      | Pkg.  |
|----------|----------------------------------|-------|
| DYNJ754  | Mouthpiece, white, thermoplastic | 50/bx |



# Medication delivery & asthma management

| Medication nebulizers           | 25-27 |
|---------------------------------|-------|
| MDI spacers                     | 28-29 |
| Peak flow meters                | 29    |
| Medication delivery accessories | 30    |
| Masks, aerosol                  | 31    |
| Nebulizer compressors           | 32-33 |

## Medication nebulizers



#### Breath-actuated

#### **Aeromist BAN**

The Aeromist breath-actuated nebulizer delivers the performance you require from a high-efficiency nebulizer, with additional features to improve operational efficiency.

- Consistent drug delivery with minimal waste; demonstrated clinical efficacy in lab testing
- Fast treatment time improves patient satisfaction and saves caregiver time
- Limited drug loss to the environment reduces caregiver exposure to aerosol
- Connects directly with a face mask, no additional pieces needed
- Easily switch between breath-actuated and continuous mode
- Filtered BAN available: HCS4433FBAN



| Item No.    | Mouthpiece | Mask type | <b>Tubing length</b> | <b>Tubing connector</b> | Pkg.  |
|-------------|------------|-----------|----------------------|-------------------------|-------|
| HCS4433BAN  | •          |           | 7-ft.                | Standard                | 50/cs |
| HCSU4433BAN | •          |           | 7-ft.                | Universal               | 50/cs |
| HCS4435BAN  |            | Adult     | 7-ft.                | Standard                | 50/cs |
| HCSU4435BAN |            | Adult     | 7-ft.                | Universal               | 50/cs |
| HCS4436BAN  |            | Pediatric | 7-ft.                | Standard                | 50/cs |
| HCSU4436BAN |            | Pediatric | 7-ft.                | Universal               | 50/cs |

#### Filtered breath actuated nebulizer

| Item No.    | Mouthpiece | Mask type | <b>Tubing length</b> | <b>Tubing connector</b> | Pkg.  |
|-------------|------------|-----------|----------------------|-------------------------|-------|
| HCS4433FBAN | •          |           | 7-ft.                | Standard                | 50/cs |









HCS4436BAN

HCS4435BAN

## Medication nebulizers

#### Breath-enhanced

#### SideStream Plus® nebulizer

The SideStream Plus breath-enhanced nebulizer features an easy-action inspiratory valve that opens on inspiration to boost medication delivery and closes on exhalation, preserving medication and reducing waste. Less pressure is required to open the inspiratory valve, making treatments easier for pediatric and elderly patients.

- Proven SideStream five-jet venturi design provides a high respirable fraction and shortens treatment times
- Can be replaced every 6 months

| Item No.  | Description                                | Pkg.  |
|-----------|--|-------|
| HS870ML   | SideStream Plus nebulizer kit              | 10/cs |
| 1110824ML | SideStream Plus nebulizer kit (6' tubing)* | 50/cs |
| 1052521ML | SideStream Plus adult aerosol mask **      | 50/cs |

<sup>\*</sup> Acute-Care only tubing

#### **Entrainment**

#### SideStream® nebulizer

The SideStream® nebulizer offers treatment times up to twice as fast compared to conventional nebulizers with an improved respirable fraction.

- Diamond Jet provides a high respirable fraction and consistent delivery
- 4 minute average treatment time

#### **Disposable**

| Item No.  | Description   | Pkg.  |
|-----------|---|-------|
| HS800ML   | SideStream nebulizer kit  | 50/cs |
| 1028219ML | SideStream nebulizer (neb cup only, no additional components)     | 50/cs |
| 1028125ML | SideStream custom neb kit (2 nebulizer kits + 2 replacement cups) | 25/cs |
| 1036160ML | SideStream nebulizer kit w/ adult mask                            | 50/cs |
| 1041906ML | SideStream nebulizer kit w/ pediatric mask                        | 50/cs |
| 1025529ML | SideStream pediatric aerosol mask                                 | 50/cs |
| 1025531ML | SideStream adult aerosol mask                                     | 50/cs |

#### Reusable

| Item No.   | Description   | Pkg.  |
|------------|---|-------|
| HS860050ML | SideStream reusable nebulizer kit                           | 50/cs |
| 1072286ML  | SideStream reusable nebulizer kit w/ tucker the turtle mask | 10/cs |







<sup>\*\*</sup> Exclusive to SideStream Plus

## Medication nebulizers


#### Standard

**Upstream design** 

#### Handheld nebulizer kits

We offer a wide selection of reusable and disposable nebulizer kits for medication delivery.

- wasted medication





Mask

type

Disposable/

Pkg.

reusable\*

#### Reservoir Item# **Adapter** Valved tee Mouthpiece tube 11004405

| HCS4485    |             |             |   |   | 7-ft.  | Standard  | Adult     | Disposable | 50/cs |
|------------|-------------|-------------|---|---|--------|-----------|-----------|------------|-------|
| HCSU4485   |             |             |   |   | 7-ft.  | Universal | Adult     | Disposable | 50/cs |
| HCS4483    | 22 mm OD/ID |             | • | • | 7-ft.  | Standard  |           | Disposable | 50/cs |
| HCS448314  | 22 mm OD/ID |             | • | • | 14-ft. | Standard  |           | Disposable | 50/cs |
| HCSU4483   | 22 mm OD/ID |             | • | • | 7-ft.  | Universal |           | Disposable | 50/cs |
| HCSU448314 | 22 mm OD/ID |             | • | • | 14-ft. | Universal |           | Disposable | 50/cs |
| HCS4482    | 22 mm OD/ID |             | • |   | 7-ft.  | Standard  |           | Disposable | 50/cs |
| HCSU4482   | 22 mm OD/ID |             | • |   | 7-ft.  | Universal |           | Disposable | 50/cs |
| HCS4486    |             |             |   |   | 7-ft.  | Standard  | Pediatric | Disposable | 50/cs |
| HCSU4486   |             |             |   |   | 7-ft.  | Universal | Pediatric | Disposable | 50/cs |
| HCS4480*   |             |             |   |   | None   | None      |           | Disposable | 50/cs |
| HCS4483VT  |             | 22 mm OD/OD |   | • | 7-ft.  | Standard  |           | Disposable | 50/cs |
| HCS4481VT  |             | 22 mm OD/ID |   |   | 7-ft.  | Standard  |           | Disposable | 50/cs |
| HCS4482R   | 22 mm OD/ID |             | • |   | 7-ft.  | Standard  |           | Reusable   | 50/cs |
| HCS4483R   | 22 mm OD/ID |             | • | • | 7-ft.  | Standard  |           | Reusable   | 50/cs |

<sup>\*</sup>Neb cup only, no additional components

45° Angled

#### Side air intake design

HCS4482RR

| Item#   | Adapter     | Valved tee | Mouthpiece | Reservoir<br>tube | Tubing<br>length | Tubing connector | Mask<br>type | Disposable/<br>reusable* | Pkg.  |
|---------|-------------|------------|------------|-------------------|------------------|------------------|--------------|--------------------------|-------|
| HCS4493 | 22 mm OD/OD |            | •          | •                 | 7-ft.            | Standard         |              | Disposable               | 50/cs |
| HCS4492 | 22 mm 0D/0D |            | •          |                   | 7-ft.            | Standard         |              | Disposable               | 50/cs |
| HCS4495 | 22 mm 0D/0D |            |            |                   | 7-ft.            | Standard         | Adult        | Disposable               | 50/cs |

Disposable up to 2 weeks of use/reusable up to 6 months of use.



7-ft.

Standard

**Tubing** 

length

**Tubing** 

connector



Reusable

Side air intake design

50/cs

# **MDI spacers**

#### **OptiChamber® Diamond**

The OptiChamber® Diamond valved holding chamber's small and intuitive design enhances medication delivery and compliance for all patients.

- Detachable LiteTouch mask has a comfortable, soft-seal interface, designed to create ideal seal with minimal force
- · Anti-static valved holding chamber, smaller than most conventional chambers, intuitive design that enhances medication delivery
- Even under low pediatric pressures and flows, valves open freely
- Stepped mouthpiece helps patients transition from pediatric facemasks to mouthpiece; allows attachment to a 22 mm
- Mouthpiece and adapter can be removed from the chamber for easy cleaning

| Item No.   | Description                          | Pkg.  |
|------------|--------------------------------------|-------|
| 1098213ML  | OptiChamber Diamond only, 10 pack    | 10/cs |
| HCS1143035 | OptiChamber Diamond only, 50 pack    | 50/cs |
| 1098216ML  | OptiChamber Diamond with large mask  | 10/cs |
| 1098215ML  | OptiChamber Diamond with medium mask | 10/cs |
| 1098214ML  | OptiChamber Diamond with small mask  | 10/cs |



1098215ML



1098213ML

#### LiteTouch masks

The LiteTouch masks are designed to provide greater comfort and a more secure seal for increased treatment effectiveness.

- Mask guickly and comfortably conforms to facial contours, providing a soft, secure fit for increased compliance and treatment effectiveness
- Flexible seal makes treatment more comfortable and less technique dependent
- Ensures a proper fit with as little as one pound of force

| lt | em No.   | Description                           | Pkg.  |
|----|----------|---------------------------------------|-------|
| 10 | 083783ML | LiteTouch mask for MDI Spacer, large  | 10/cs |
| 10 | 083784ML | LiteTouch mask for MDI Spacer, medium | 10/cs |
| 10 | 083785ML | LiteTouch mask for MDI Spacer, small  | 10/cs |
|    |          |                                       |       |
|    |          |                                       | 1     |
|    |          |                                       |       |
|    |          | 2                                     |       |



#### **Medline MDI spacers**

Our spacers are designed for patients who have difficulty with metered dose inhalers (MDIs). With spacers, more medication is deposited in the lungs and less is deposited in the mouth and throat.

- A simple and economical way to improve effectiveness of inhaler
- Easy to use and helps avoid some of the unpleasant taste of inhaled medication
- Anti-static and compatible with any standard mask with a 22 mm fitting
- Masks are made of a soft silicone and fit snugly to ensure optimal drug delivery
- The collapsible model is portable and convenient

| Item No. | Description                    | Pkg.  |
|----------|--------------------------------|-------|
| HCSA100  | MDI spacer only                | 10/bx |
| HCSA200  | MDI spacer only, collapsible   | 10/bx |
| HCSA112  | MDI spacer with infant mask    | 10/bx |
| HCSA113  | MDI spacer with pediatric mask | 10/bx |
| HCSA114  | MDI spacer with adult mask     | 10/bx |





HCSA113 HCSA112



HCSA114

#### **Peak flow meters**

Peak flow meters help patients manage their asthma and gauge the effectiveness of their treatment program, which may help them predict their next attack.

#### Personal best

Lightweight, portable peak flow meter helps patients monitor their asthma. Handheld device measures Peak Expiratory Flow Rate so patients can initiate preemptive therapy in order to help prevent serious asthma attacks.

- Universal range allows the Personal Best peak flow meter to be used for all patients, from pediatric to adult
- Integrated case for storage
- Mouthpieces available separately

| Item No.   | Description                                    | Pkg.   |
|------------|--|--------|
| 1022973UML | Personal best peak flow meter, universal range | 10/cs  |
| HS714200ML | Plastic mouthpiece with one-way valve          | 200/cs |
| HS715100ML | Plastic pediatric mouthpiece                   | 100/cs |
| HS757100ML | Oval disposable plastic mouthpiece             | 100/cs |

#### Medline

- Individually calibrated to ensure accuracy
- Durable for years of reliable performance
- The pocket model offers superb convenience and portability
- Available in low range and full range versions—low range models are better for geriatric, pediatric or other patients with diminished flow rates

| Item No.  | Description  | Pkg.   |
|-----------|--|--------|
| HCSPF100  | Portable peak flow meter, full range,<br>60–720 LPM  | 12/bx  |
| HCSPF200  | Standard peak flow meter, full range,<br>60–850 LPM  | 12/bx  |
| HCSPFPEDI | Mini wright pediatric peak flow meter,<br>30–400 LPM | 100/cs |
| HCSPF220  | Disposable pediatric mouthpiece, adapter             | 20/bx  |
| HCSPF262  | Disposable one way valve mouthpiece, adult           | 100/cs |



1022973UML

HCSPF100



1-800-MEDLINE (633-5463) | medline.com 29

#### Valved tee adapters

- For inline medication delivery
- Available in 22 mm OD/OD and 22 mm OD/ID

| Item No. | Description | Pkg.  |
|----------|-------------|-------|
| HCSA1008 | 22 mm 0D/0D | 30/cs |
| HCSA1009 | 22 mm OD/ID | 30/cs |



#### **Southmedic inline MDI adapters**

- Single-patient-use adult metered dose inhaler (MDI) adapter ensures safe and complete delivery of prescribed aerosol medication
- Adapter remains inline to help reduce cross-contamination and improve patient comfort by reducing manipulation of circuit
- Arrows clearly mark the direction of flow; built-in tethered plug can be used to seal the breathing circuit when no MDI canister is in place
- Compatible with both plastic and metal canister stems
- Can be used with or without dosage counter

| Item No.    | Description                                | Pkg.  |
|-------------|--|-------|
| SMDIA1000ML | 22 mm ID $\rightarrow$ 22 mm OD / 15 mm ID | 25/cs |
| SMDIA1100ML | 22 mm OD / 15 mm ID $\rightarrow$ 22 mm ID | 25/cs |
| SMDIA1200ML | 22 mm OD $\rightarrow$ 22 mm OD / 15 mm ID | 25/cs |
| SMDIA1300ML | 15 mm OD $\rightarrow$ 15 mm ID            | 25/cs |





#### **Aerosol masks**

Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material. All aerosol masks feature an elongated fit for optimal patient fit.

- A secure fit with elastic straps
- Standard 22 mm connector
- Not made with natural rubber latex

#### Adult

| Item No. | Description                          | Pkg.  |
|----------|--------------------------------------|-------|
| HCS4630B | Aerosol mask with straight connector | 50/cs |
| HCS4630  | Aerosol mask with swivel connector   | 50/cs |

#### **Pediatric**

| Item No. Description |          | Description                        | Pkg.  |
|----------------------|----------|------------------------------------|-------|
|                      | HCS4631B | Aerosol mask with swivel connector | 50/cs |

#### SideStream

| Item No.  | Description  | Pkg.  |
|-----------|--|-------|
| 1044382ML | SideStream Tucker the Turtle pediatric aerosol mask* | 50/cs |
| HS880ML   | SideStream Sami the Seal pediatric aerosol mask**    | 50/cs |
|           | isposable and reusable                               |       |
| 1044382ML | HCS463   |       |

#### **Aeromist Compact**<sup>™</sup> **nebulizer compressors**

Our smallest nebulizer compressor. Choose from the standard white model or color options.

- Kits include a disposable nebulizer medicine cup, tee mouthpiece, 6" corrugate reservoir, 7-ft (2.13 m) of tubing, and 5 replacement filters
- Color choices: electric blue, neon green, neon pink, charcoal, or white
- Durable hard plastic shell
- · Warranty: 5 years
- Replacement Kits: HCS4483 (Disposable) and HCS4483R (Reusable)

#### **Specifications**

- Method of nebulization: dry piston compressor
- Electrical requirements: 120 VAC / 60 Hz
- Maximum jet nebulizer capacity: 10 mL
- Sound noise level: ≤55 dBA
- Operation air pressure: 9-16 PSI, Max air pressure: 35-50 PSI
- Operation air flow rate: 5-8 LPM

| Item No.     | Description   | Pkg. |
|--------------|---|------|
| HCS70004G    | Neon green compressor and nebulizer kit               | 8/cs |
| HCS70004BL   | Electric blue compressor and nebulizer kit            | 8/cs |
| HCS70004P    | Neon pink compressor and nebulizer kit                | 8/cs |
| HCS70004C    | Charcoal compressor and nebulizer kit                 | 8/cs |
| HCS70004     | White compressor and nebulizer kit                    | 8/cs |
| HCS70004ADLT | White compressor and nebulizer kit with adult mask    | 8/cs |
| HCS70004RD   | White compressor with reusable and disposable neb kit | 8/cs |
| HCS70004DP   | Sample pack, two of each color                        | 8/cs |
| HCS70004F    | Replacement filters                                   | 5/pk |
| HCS70004CAP  | Replacement cap for HCS70004                          | 8/pk |

#### **Aeromist Trek portable nebulizer compressor**

Better compliance with easy treatments, anytime and anywhere.

- Suitable for travel with a sleek design, extra small size, light weight and battery operation
- Durable hard plastic shell
- Compact and stable, with a small footprint
- Kit includes a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, mouthpiece, rechargeable battery, carry case, AC and DC adapters, adult and pediatric masks, and 5 replacement filters
- Warranty: 3 years for the compressor, 6 months for the battery

| Item No.   | Description  | Pkg. |
|------------|--|------|
| HCS2G0     | Portable battery-operated nebulizer compressor kit                                 | 1/ea |
| HCS2GONB   | Portable nebulizer compressor, battery not included                                | 1/ea |
| HCS2GOBATT | Replacement battery for portable nebulizer compressors (Items HCS2GO and HCS2GONB) | 1/ea |
| HCS2G0F    | Replacement filters  | 5/pk |
| HCS2GOAC   | Replacement AC adapter   | 1/ea |
| HCS2GODC   | Replacement car charger  | 1/ea |





Components included with each compressor



HCS2G0



#### **Aeromist Buddies**<sup>™</sup> **nebulizer compressors**

Kid-friendly nebulizers help put young patients at ease, while still providing safe and effective aerosol treatments.

- Kit includes disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, 6" reservoir tube, T-Mouthpiece, pediatric mask, and 5 replacement filters
- Durable hard plastic shell
- Warranty: 5 years

#### **Specifications**

- Method of nebulization: dry piston compressor
- Electrical requirements: 120 VAC / 60 Hz
- Maximum nebulizer capacity: 10 mL
- Sound noise level: ≤60 dBA
- Operation air pressure: 10-17 PSI, Max air pressure: 35-50 PSI
- Operation air flow rate: 4-7.5 LPM

| Item No.     | Description  | Pkg. |
|--------------|--|------|
| HCSMONKEY    | Monkey compressor and nebulizer kit                          | 8/cs |
| HCSMONKEYCAP | Replacement filter cap for HCSMONKEY                         | 8/pk |
| HCSFISHNEB*  | Fish compressor and nebulizer kit                            | 8/cs |
| HCSPUPPY*    | Puppy compressor and nebulizer kit                           | 8/cs |
| HCSPUPPYCAP  | Replacement filter cap for HCSPUPPY                          | 8/pk |
| HCSPUPPYRD*  | Puppy compressor with disposable and reusable nebulizer kits | 8/pk |
| HCSTURTLE    | Turtle compressor and nebulizer kit                          | 8/cs |
| HCSTURTLECAP | Replacement filter cap for HCSTURTLE                         | 8/pk |
| HCS70004F    | Replacement filters  | 5/pk |



All compressors include a complementary nebulizer kit.



# Non-invasive interfaces

| Medline Flex® non-invasive Ventilation masks                            | 35    |
|---|-------|
| Medline Flex <sup>™</sup> non-invasive<br>ventilation masks accessories | 36    |
| DreamPort™ CPAP mask solution   | 37    |
| CPAP masks AIRgel®  | 38-39 |
| CPAP masks silicone   | 40    |
| CPAP accessories  | 41    |
| Go-PAP <sup>™</sup> systems   | 42    |
| 02-Max <sup>™</sup> systems   | 42-43 |

#### Medline Flex<sup>™</sup> non-invasive ventilation (NIV) adult full face masks with FlexiPlane™ technology

Mask fit and minimizing leaks is critical to success with non-invasive ventilation. Our NIV adult full face masks are designed to offer patient comfort and reduce the risk of skin breakdown. They are available in small, medium, and large sizes as well as vented and non-vented options to meet your needs. Medline NIV masks include the following features:

- · Dual swivel, flexible elbow with FlexiPlane technology which allows it to flex in the 3rd plane which is designed to reduce pressure points on the patient's face
- · Large, adjustable silicone forehead pad designed to reduce pressure points and aid in stability
- · Soft and flexible cushion to reduce the need to overtighten the mask.
- · Comfortable headgear with clips for quick release

| Item No.  | Description   | Pkg.  |
|-----------|---|-------|
| HCSFFVAS  | Small, vented, anti-asphyxia valve(AAV) and headgear      | 10/cs |
| HCSFFVAM  | Medium, vented, anti-asphyxia valve(AAV) and headgear     | 10/cs |
| HCSFFVAL  | Large, vented, anti-asphyxia valve(AAV) and headgear      | 10/cs |
| HCSFFNVAS | Small, non-vented, anti-asphyxia valve(AAV) and headgear  | 10/cs |
| HCSFFNVAM | Medium, non-vented, anti-asphyxia valve(AAV) and headgear | 10/cs |
| HCSFFNVAL | Large, non-vented, anti-asphyxia valve(AAV) and headgear  | 10/cs |
| HCSFFNVS  | Small, non-vented, and headgear                           | 10/cs |
| HCSFFNVM  | Medium, non-vented, and headgear                          | 10/cs |
| HCSFFNVL  | Large, non-vented, and headgear                           | 10/cs |



**HCSFFVAM** 

#### Medline non-invasive ventilation (NIV) adult nasal masks

NIV adult nasal masks are available in small, medium, and large sizes and designed to improve patient comfort.

- · Offered in vented and non-vented options
- · Silicone forehead pad designed to reduce pressure points and aids in stability
- · Comfortable headgear with clips for quick release

| Item No. | Description   | Pkg.  |
|----------|---|-------|
| HCSNVAS  | Small, vented, anti-asphyxia valve(AAV) and headgear      | 10/cs |
| HCSNVAM  | Medium, vented, anti-asphyxia valve(AAV) and headgear     | 10/cs |
| HCSNVAL  | Large, vented, anti-asphyxia valve(AAV) and headgear      | 10/cs |
| HCSNNVAS | Small, non-vented, anti-asphyxia valve(AAV) and headgear  | 10/cs |
| HCSNNVAM | Medium, non-vented, anti-asphyxia valve(AAV) and headgear | 10/cs |
| HCSNNVAL | Large, non-vented, anti-asphyxia valve(AAV) and headgear  | 10/cs |



HCSNVAL

#### Medline non-invasive ventilation (NIV) circuits

Medline offers flexible, smooth bore, NIV circuits that are available in two configurations to meet your needs. These circuits are compatible with either vented or non-vented NIV masks with an anti-asphyxia valve and include a removable exhalation port and pressure line. The filter has a BFE of 99.99% and VFE of 99.9% (if included).

| Item No.   | Description   | Pkg.  |
|------------|---|-------|
| DYNJ500EP  | Single limb, smooth bore circuit with exhalation port, 6 ft.            | 10/ea |
| DYNJ500EPF | Single limb, smooth bore circuit with exhalation port and filter, 6 ft. | 10/ea |



DYNJ500EP

#### Medline Flex<sup>™</sup> non-invasive ventilation (NIV) adult mask accessories

We offer replacement 4-point, breathable headgear for both full face and nasal NIV masks.

- · Breathable material that stretches to offer comfort and aid in mask cushion seal
- · Includes clips for quick release

| Item No. | Description                                 | Pkg.   |
|----------|---|--------|
| HCSFFHG  | Headgear with clips for full face NIV masks | 100/cs |
| HCSNHG   | Headgear with clips for nasal NIV masks     | 100/cs |



DYNJ500EPF

#### **DreamPort**<sup>™</sup> **CPAP mask solution**

The DreamPort CPAP mask solution from Bleep LLC, is an innovative CPAP device designed to eliminate leaks, and improve patient comfort and compliance.



- It is one-size-fits-all and is customizable to each patient's needs
- Adhesive, hypoallergenic surgical tape on the DreamPorts create an airtight seal
- The DreamPorts pair with the DreamWay that ensures a strapless connection to standard-sized CPAP tubing
- · Lightweight and easy to use

| Item No.   | Description                        | Pkg.  |
|------------|------------------------------------|-------|
| CPAP100383 | Bleep DreamPort sleep solution kit | 24/cs |
| CPAP100382 | Bleep DreamPorts, 32 ports per box | 54/cs |
| CPAP100387 | Bleep DreamWay CPAP interface      | 24/cs |
| CPAP100434 | Bleep DreamPreps, 100 wipes        | 8/cs  |



watch a video about the DreamPort, visit



CPAP100383 (includes 1 DreamWay Interface and 1 box of DreamPorts)



CPAP100382





CPAP100434

#### AIRgel® CPAP masks

Amazing comfort and customized fit technology. The patented AIRgel replaceable cushioning minimizes discomfort and skin irritation while effectively sealing airflow. It is easily detached and reattached for easy cleaning, saving time and money.

- · Proprietary AlRgel cushioning minimizes discomfort and skin irritation while effectively sealing airflow
- The flexible and formable mask provides a custom fit

#### Aura nasal mask

- The shapeable shell of each mask can be easily molded to fit individual needs
- Flexible spacebar is designed to provide superior stability without touching the skin
- Active headgear connector moves with the patient to help maintain an effective seal

| Item No.  | Description            | Mask Size | Pkg. |
|-----------|------------------------|-----------|------|
| CPAP50176 | Adjustable gel cushion | Small     | 1/ea |
| CPAP50177 | Adjustable gel cushion | Medium    | 1/ea |
| CPAP50178 | Adjustable gel cushion | Large     | 1/ea |

#### Mojo<sup>®</sup> full-face gel mask

- Full-face design in four sizes
- Adjustable tension wheel minimizes discomfort while ergonomic cushions adapt to unique facial features for a comfortable fit
- Magnetic headgear clips hold firmly and disconnects with one hand easy on-and-off

| Item No.  | Description      | Headgear style | Mask size | Pkg. |
|-----------|------------------|----------------|-----------|------|
| CPAP50835 | Vented           | One size       | Small     | 1/ea |
| CPAP50836 | Vented           | One size       | Medium    | 1/ea |
| CPAP50837 | Vented           | One size       | Large     | 1/ea |
| CPAP50838 | Vented           | One size       | XL        | 1/ea |
| CPAP55045 | Headgear<br>Only | One size       | One size  | 1/ea |



Headgear included

CPAP50177





#### **AIRgel® CPAP masks**

#### iQ Blue® nasal gel mask

- One-size nasal mask design fits a wide variety of patients
- Clearest field of vision possible with no forehead pad required patients can keep their glasses on while reading or watching TV
- Easy on-and-off thanks to quick release, breathable headgear

| Item No.  | Mask style    | Headgear style | Mask size | Pkg. |
|-----------|---------------|----------------|-----------|------|
| CPAP50655 | Vented        | Stable fit     | One size  | 1/ea |
| CPAP55044 | Headgear only | Stable fit     | One size  | 1/ea |

#### MiniMe® 2 pediatric gel mask

- Nasal mask in four sizes, designed for children ages 2 to 12
- Headgear with 5 adjustment points
- Soft, breathable headgear for all-night comfort
- Available in vented and non-vented designs
- Adjustable cushion with memory wire for a true custom fit

| Item No.   | Description | Headgear size | Mask size | Pkg.  |
|------------|-------------|---------------|-----------|-------|
| CPAP60216  | Vented      | X-Small/small | X-Small   | 10/cs |
| CPAP602172 | Vented      | X-Small/small | Small     | 10/cs |
| CPAP602182 | Vented      | Medium/large  | Medium    | 10/cs |
| CPAP60221  | Vented      | Medium/large  | Large     | 10/cs |
| CPAP60256  | Non-vented  | X-Small/small | X-Small   | 10/cs |
| CPAP602572 | Non-vented  | X-Small/small | Small     | 10/cs |
| CPAP602582 | Non-vented  | Medium/large  | Medium    | 10/cs |
| CPAP60261  | Non-vented  | Medium/large  | Large     | 10/cs |
| CPAP55055  | Headgear    | X-Small/small | N/A       | 10/cs |
| CPAP55056  | Headgear    | Medium/large  | N/A       | 10/cs |

#### Ascend® masks

Patented AIRgel replaceable cushion is easily detached and reattached for easy cleaning and replacement. Leak ports are built-in, so a separate exhalation port is not needed.

- Minimizes pressure points and skin irritation
- Ergonomic cushion technology delivers an effective seal for the most comfortable fit
- Flexible spacebar with a unique bendable design provides superior stability, holding a seal without contacting the skin
- · Lightweight, breathable headgear slides off easily with one hand
- 360° rotating swivel for optimal tube positioning and easier mobility



| Item No.  | Description   | Pkg. |
|-----------|---|------|
| CPAP50825 | Starter kit with small, medium, large cushions and headgear | 1/ea |
| CPAP50813 | Full-face mask, small                                       | 1/ea |
| CPAP50814 | Full-face mask, medium                                      | 1/ea |
| CPAP50815 | Full-face mask, large                                       | 1/ea |
| CPAP50207 | Replacement cushion, small                                  | 1/ea |
| CPAP50208 | Replacement cushion, medium                                 | 1/ea |
| CPAP50209 | Replacement cushion, large                                  | 1/ea |
| CPAP55051 | CPAP headgear replacement, one size                         | 1/ea |



CPAP50655





CPAP50815

#### **CPAP** nasal mask

| Item No.  | Description   | Pkg. |
|-----------|---|------|
| CPAP50174 | Starter kit with small, medium, large cushions and headgear | 1/ea |
| CPAP50187 | Nasal mask, small   | 1/ea |
| CPAP50188 | Nasal mask, medium  | 1/ea |
| CPAP50189 | Nasal mask, large   | 1/ea |
| CPAP50201 | Replacement cushion, small                                  | 1/ea |
| CPAP50202 | Replacement cushion, medium                                 | 1/ea |
| CPAP50203 | Replacement cushion, large                                  | 1/ea |
| CPAP55051 | CPAP headgear replacement, one size                         | 1/ea |

#### **Silicone CPAP masks**

A consistent seal without excessive tightness, delivering the prescribed therapy with enough comfort to help improve compliance.

- OmniClip and silicone forehead pad for a custom fit at the bridge of the nose for patients of all head shapes and sizes
- SilentVent elbow flushes out exhaled air quietly, without an additional inline exhalation port
- Headgear attaches with convenient clips, includes a handy wash bag to keep it clean and fresh
- Cushions, forehead pads, elbows, cones and clips are individually replaceable
- · Available in four sizes, with and without headgear

| Item No.                                     | Size                                | Description  | Pkg. |
|--|-------------------------------------|--|------|
| CPAP9400<br>CPAP9401<br>CPAP9402<br>CPAP9403 | Small<br>Medium<br>Large<br>X-Large | <ul> <li>Full face mask with head strap</li> <li>SilentVent anti-asphyxia elbow, 22 mm male</li> <li>OmniClip with silicone forehead pad</li> </ul>                  | 1/ea |
| CPAP9404<br>CPAP9405<br>CPAP9406<br>CPAP9407 | Small<br>Medium<br>Large<br>X-Large | <ul> <li>Full face mask, no head strap</li> <li>SilentVent anti-asphyxia elbow, 22 mm male</li> <li>OmniClip with silicone forehead pad</li> </ul>                   | 1/ea |
| CPAP9408<br>CPAP9409<br>CPAP9410<br>CPAP9411 | Small<br>Medium<br>Large<br>X-Large | <ul> <li>Nasal mask with head strap</li> <li>SilentVent elbow, 22 mm male</li> <li>Silicone forehead pad</li> <li>OmniClip with silicone forehead pad</li> </ul>     | 1/ea |
| CPAP9412<br>CPAP9413<br>CPAP9414<br>CPAP9415 | Small<br>Medium<br>Large<br>X-Large | <ul> <li>Nasal mask, no head gear strap</li> <li>SilentVent elbow, 22 mm male</li> <li>Silicone forehead pad</li> <li>Omniclip with silicone forehead pad</li> </ul> | 1/ea |



CPAP9403





#### Aloha nasal pillow systems

Medline offers the Aloha nasal pillow system for CPAP therapy patients. This lightweight, nasal pillow system is quiet and offers a stable and comfortable seal.

| Item No.    | Description   | Pkg. |
|-------------|---|------|
| SUHCM014    | Includes: one small, medium and large pillow headgear | 1/ea |
| SUHCMP014SP | Replacement nasal pillow for Aloha, small             | 1/ea |
| SUHCMP014MP | Replacement nasal pillow for Aloha, medium            | 1/ea |
| SUHCMP014LP | Replacement nasal pillow for Aloha, large             | 1/ea |
| SUHHG014    | Replacement headgear for Aloha                        | 1/ea |



SUHCM014

#### **CPAP** accessories

#### **Tubing**

| Item No.   | Description                        | Pkg.  |
|------------|------------------------------------|-------|
| SUH0610    | 6-ft. tubing with flex connectors  | 10/bx |
| SUHTUB008  | 8-ft. tubing with flex connectors  | 1/ea  |
| SUHTUB010  | 10-ft. tubing with flex connectors | 1/ea  |
| SUHTUB0018 | 18-ft. tubing with flex connectors | 1/ea  |

#### **Chin straps**

| Item No. | Description | Pkg. |
|----------|-------------|------|
| SUHCS001 | Standard    | 1/ea |
| SUHCS003 | Adjustable  | 1/ea |
| SUHCS004 | Deluxe      | 1/ea |
| SUHCS006 | Premium     | 1/ea |



SUH0610



SUHCS006

#### **Go-PAP**<sup>™</sup> **CPAP** systems

The Go-PAP CPAP System offers emergency CPAP with integrated nebulization using the flow from an oxygen flowmeter or cylinder.

- Connects to barbed connector and can generate up to 90 LPM flow
- Includes 3 options for PEEP of 5.0, 7.5 and 10.0 cm  $H_{\rm 2}{\rm 0}$
- Integrated nebulization without interrupting the delivered flow to the patient



| Item No.      | Description  | Pkg.  |
|---------------|--|-------|
| CPAP3134602   | With BiTrac™ ED mask and standard headgear, adult medium                                 | 10/cs |
| CPAP3134602N  | With BiTrac ED mask, standard headgear and nebulizer canister, adult medium              | 10/cs |
| CPAP3134602NC | With BiTrac ED mask, standard headgear, nebulizer canister and Neb-Connect, adult medium | 10/cs |
| CPAP3134603   | With BiTrac ED mask and standard headgear, adult large                                   | 10/cs |
| CPAP3134603N  | With BiTrac ED mask, standard headgear and nebulizer canister, adult large               | 10/cs |
| CPAP3134603NC | With BiTrac ED mask, standard headgear, nebulizer canister and Neb-Connect, adult large  | 10/cs |

#### O<sub>2</sub>-Max<sup>™</sup> CPAP systems

The O<sub>2</sub>-Max CPAP System offers emergency CPAP with integrated nebulization using the flow from an oxygen flowmeter or cylinder.

- Connects to DISS port and can generate up to 140 LPM flow
- Includes 3 options for PEEP of 5.0, 7.5 and 10.0 cm  $H_2O$
- Integrated nebulization without interrupting the delivered flow to the patient



#### O<sub>2</sub>-MAX CPAP system with 3-set CPAP (5.0, 7.5 and 10.0 cm H<sub>2</sub>0)

| Item No.      | Description   | Pkg.  |
|---------------|---|-------|
| CPAP3137555X  | With BiTrac ED mask, adult large  | 10/cs |
| CPAP3137555XN | With BiTrac ED mask and nebulizer canister, adult large                       | 10/cs |
| CPAP3137555X1 | With BiTrac ED mask and Ohmeda Quik-Connect™, adult large                     | 10/cs |
| CPAP3137555N1 | With BiTrac ED mask, Ohmeda Quik-Connect and nebulizer canister, adult large  | 10/cs |
| CPAP3137555X2 | With BiTrac ED mask and Chemtron Quik-Connect, adult large                    | 10/cs |
| CPAP3137556X  | With BiTrac ED mask, adult medium   | 10/cs |
| CPAP3137556XN | With BiTrac ED mask and nebulizer canister, adult medium                      | 10/cs |
| CPAP3137556X1 | With BiTrac ED mask, Ohmeda Quik-Connect and nebulizer canister, adult medium | 10/cs |
| CPAP31375561  | With BiTrac ED mask and Ohmeda Quik-Connect, adult medium                     | 10/cs |
| CPAP3137556X2 | With BiTrac ED mask and Chemtron Quik-Connect, adult medium                   | 10/cs |

#### O<sub>2</sub>-MAX CPAP system with 5-set CPAP (5.0, 7.5, 10.0, 12.5 and 15 cm H<sub>2</sub>0)

| Item No.      | Description   | Pkg.  |
|---------------|---|-------|
| CPAP3137554X  | With BiTrac ED mask, adult large  | 10/cs |
| CPAP3137554X1 | With BiTrac ED mask, Ohmeda Quik-Connect and nebulizer canister, adult large  | 10/cs |
| CPAP3137554XN | With BiTrac ED mask and nebulizer canister, adult large                       | 10/cs |
| CPAP3137557X  | With BiTrac ED mask, adult medium   | 10/cs |
| CPAP3137557X1 | With BiTrac ED mask, Ohmeda Quik-Connect and nebulizer canister, adult medium | 10/cs |
| CPAP3137557XN | With BiTrac ED mask and nebulizer canister, adult medium                      | 10/cs |

O<sub>2</sub>-MaxCPAP systems continues on the next page

#### **O<sub>2</sub>-Max**<sup>™</sup> **CPAP systems** (continued)

#### O<sub>2</sub>-MAX CPAP system (with fixed 10 cm H<sub>2</sub>0 CPAP)

| Item No.      | Description   | Pkg.  |
|---------------|---|-------|
| CPAP3137535X  | With BiTrac ED mask, adult large  | 10/cs |
| CPAP3137538X  | With BiTrac ED mask, adult medium   | 10/cs |
| CPAP3137538XN | With BiTrac ED mask and nebulizer canister, adult medium                      | 10/cs |
| CPAP31375381  | With BiTrac ED mask and Ohmeda Quik-Connect, adult medium                     | 10/cs |
| CPAP3137538X1 | With BiTrac ED mask, Ohmeda Quik-Connect and nebulizer canister, adult medium | 10/cs |

#### O<sub>2</sub>-MAX Trio CPAP system (Can be set to 30%, 60% or 90% Fio<sub>2</sub>)

| Item No.      | Description   | Pkg.  |
|---------------|---|-------|
| CPAP3138007X  | 3-SET PEEP system with BiTrac ED mask with head strap, adult medium   | 10/cs |
| CPAP3138017X  | 3-SET PEEP system with BiTrac ED mask with head strap, adult large  | 10/cs |
| CPAP3138007X1 | 3-SET PEEP system with BiTrac ED mask with head strap, Ohmeda Quik-Connect and nebulizer canister, adult medium | 10/cs |
| CPAP3138107X  | 3-SET PEEP system with BiTrac ED mask with full cap headgear, adult medium                                      | 10/cs |
| CPAP3138017X1 | 3-SET PEEP system with BiTrac ED mask with head strap, Ohmeda Quik-Connect and nebulizer canister, adult large  | 10/cs |
| CPAP3138117X  | 3-SET PEEP system with BiTrac ED mask with full cap headgear, adult large                                       | 10/cs |
| CPAP3138117XN | 3-SET PEEP system with BiTrac ED mask with full cap headgear, small volume nebulizer, adult large               | 10/cs |
| CPAP31381171  | 3-SET PEEP system with BiTrac ED mask with full cap headgear, adult large                                       | 10/cs |
| CPAP3138008X  | 5-SET PEEP system with BiTrac ED mask with head strap, adult medium   | 10/cs |
| CPAP3138108X  | 5-SET PEEP system with BiTrac ED mask with full cap headgear, adult medium                                      | 10/cs |
| CPAP3138018X  | 5-SET PEEP system with BiTrac ED mask with head strap, adult large  | 10/cs |
| CPAP3138118X  | 5-SET PEEP system with BiTrac ED mask with full cap headgear, adult large                                       | 10/cs |
| CPAP3138003X  | 10 cm H <sub>2</sub> O fixed PEEP system with BiTrac ED mask with head strap, adult medium                      | 10/cs |
| CPAP3138013X  | 10 cm H <sub>2</sub> O fixed PEEP system with BiTrac ED mask with head strap, adult large                       | 10/cs |

#### O<sub>2</sub>-Max accessories

| Item No.      | Description   | Pkg.  |
|---------------|---|-------|
| CPAP3137028   | Small mask for Go-PAP™ and O <sub>2</sub> -Max CPAP systems | 10/cs |
| CPAP3137029   | Medium mask for Go-PAP and O <sub>2</sub> -Max CPAP systems | 10/cs |
| CPAP3137030   | Large mask for Go-PAP and O <sub>2</sub> -Max CPAP systems  | 10/cs |
| CPAP313111111 | Quik-Connect™ Ohmeda DISS                                   | 10/cs |
| CPAP3138000   | O <sub>2</sub> -Max Trio with 6" tubing and DISS            | 10/cs |

# O<sub>2</sub> & CO<sub>2</sub> disposables

| Southmedic OxyMask <sup>®</sup> | 45-4  |
|---------------------------------|-------|
| Cannulas, O <sub>2</sub>        | 47-48 |
| Cannulas, CO <sub>2</sub>       | 49    |
| Tubing, CO <sub>2</sub>         | 50    |
| Masks, CO <sub>2</sub>          | 51    |
| Tubing, O <sub>2</sub>          | 52    |
| Masks, O <sub>2</sub>           | 53    |
| Masks, aerosol                  | 54    |
| Masks, tent & trach             | 54    |
| Bubble humidifiers              | 55    |
| O <sub>2</sub> accessories      | 55    |

### **OxyMask**<sup>™</sup>

#### Everyone breathes easier with OxyMask.

Introducing OxyMask—a unique open-style oxygen mask by Southmedic. Its proprietary design allows one mask to provide a broad range of oxygen concentrations while increasing patient comfort, reducing safety risks for patients, and helping to improve overall clinical and operational efficiencies



#### **Patient experience**



Reduce risk of CO<sub>2</sub> re-breathing

Easier communication

Ability to sip fluids to take oral medications

Fewer interruptions in care

#### **Caregiver efficiency**



Reduced risk of errors Only one product to learn Standardization of care

Time saving

#### **Operational efficiency**



Reduce SKUs Reduce O2 consumption Save storage space Reduce waste

#### **OxyMask**

Open mask design that delivers oxygen flows from 1–15 liters per minute at a range of 24–90% FiO<sub>2</sub>

| Item No.      | Size                 | Universal<br>tubing length | Pkg.  |
|---------------|----------------------|----------------------------|-------|
| OM11258ML     | Adult                | 7-ft.                      | 25/bx |
| OP11258ML     | Adult Plus           | 7-ft.                      | 25/bx |
| OK11258ML     | Pediatric            | 7-ft.                      | 25/bx |
| OT10258ML*    | Infant               | 7-ft.                      | 25/bx |
| OM112514ML    | Adult                | 14-ft.                     | 25/bx |
| OP112514ML    | Adult Plus           | 14-ft.                     | 25/bx |
| OAPLUS11258** | Adult Plus<br>OxyArm | 7-ft.                      | 25/bx |

<sup>\*</sup>Infant Mask delivers from ¼ L/min- 6+ at 22%-65% FiO2

#### OxyMask<sup>™</sup> Non-PVC

| Item No.   | Size       | Universal<br>Tubing length | Pkg.  |
|------------|------------|----------------------------|-------|
| OM70258ML  | Adult      | 7-ft.                      | 25/bx |
| 0P70258ML  | Adult Plus | 7-ft.                      | 25/bx |
| OK70258ML  | Pediatric  | 7-ft.                      | 25/bx |
| OT70258ML* | Infant     | 7-ft.                      | 25/bx |

<sup>\*</sup>Infant Mask delivers from ¼ L/min- 6+ at 22%-65% FiO2

<sup>\*\*</sup>OxyArm™ delivers 24–50% FiO2

#### OxyMulti-Mask™

- Open mask design that delivers oxygen flows from 1–15 liters per minute at a range of 23-83% FiO<sub>2</sub>
- Delivers oxygen through standard small bore tube using 22mm adapter or can be used with aerosol nebulizers

| Item No.     | Description                              | Universal<br>tubing length | Pkg.  |
|--------------|--|----------------------------|-------|
| OHH10258ML   | Adult mask with 22 mm adapter            | 7-ft.                      | 25/bx |
| OPHH10258ML  | Adult Plus mask with 22 mm adapter       | 7-ft.                      | 25/bx |
| OKHH10258ML* | Pediatric mask with 22 mm adapter        | 7-ft.                      | 25/bx |
| OAT1025ML    | Adult trach mask                         | 7-ft.                      | 25/bx |
| OAT1025CSML  | Adult trach mask with comfort neck strap | 7-ft.                      | 25/bx |

<sup>\*</sup>Pediatric mask delivers 22-93% FiO2



OHH10258ML

#### OxyMask aerosol

Closed mask with open bottom design and right angle swivel adapter for optimal aerosol delivery, even when patient is laying down or sitting.

| Item No.  | Description                       | Pkg.  |
|-----------|-----------------------------------|-------|
| OMN1025ML | Adult mask with 22 mm adapter     | 25/bx |
| OTN1025ML | Infant mask with 22 mm adapter    | 25/bx |
| OKN1025ML | Pediatric mask with 22 mm adapter | 25/bx |



#### OxyMask ETCO<sub>2</sub>

Open design mask that delivers oxygen flows from 1–15 liters per minute at a range of 24-65% FiO<sub>2</sub> with end-tidal CO<sub>2</sub> monitoring.

#### With female pigtail adapter

| Item No.   | Size       | Universal tubing length | Female pigtail adapter length | Pkg.  |
|------------|------------|-------------------------|-------------------------------|-------|
| OM21258ML  | Adult      | 7-ft.                   | 8-in.                         | 25/bx |
| OP21258ML  | Adult Plus | 7-ft.                   | 8-in.                         | 25/bx |
| OK21258ML  | Pediatric  | 7-ft.                   | 8-in.                         | 25/bx |
| OM212514ML | Adult      | 14-ft.                  | 15-in.                        | 25/bx |

#### With male gas sampling line

| Item No.     | Size       | Universal tubing<br>length | Male gas sampling<br>line length | Pkg.  |
|--------------|------------|----------------------------|----------------------------------|-------|
| OM21258SLMML | Adult      | 7-ft.                      | 8-ft.                            | 25/bx |
| OP21258SLMML | Adult Plus | 7-ft.                      | 8-ft.                            | 25/bx |
| OK21258SLMML | Pediatric  | 7-ft.                      | 8-ft.                            | 25/bx |
| OM212514SLM  | Adult      | 14-ft.                     | 15-ft.                           | 25/bx |



OM21258SLMML



#### SuperSoft oxygen cannulas

Even short term use of an oxygen cannula can lead to the development of pressure ulcers on the ears and cheeks. The Medline SuperSoft nasal cannula is made of an advanced material that is soft, flexible and non-tacky. Also features curved, non-flared nasal prongs. You will notice a difference right away:

- Promote patient compliance with added comfort patients are less apprehensive about wearing a cannula that feels soft and comfortable
- Not made with natural rubber latex

#### Adult

| Item No.  | <b>Tubing length</b>             | Connector | Pkg.  |
|-----------|----------------------------------|-----------|-------|
| HCS4510S  | No tubing (cannula harness only) | Male      | 50/cs |
| HCSS4504S | 4-ft.                            | Standard  | 50/cs |
| HCSU4504S | 4-ft.                            | Universal | 50/cs |
| HCSS4514S | 7-ft.                            | Standard  | 50/cs |
| HCSS4516S | 14-ft.                           | Standard  | 50/cs |
| HCSU4516S | 14-ft.                           | Universal | 50/cs |
| HCSS4515S | 25-ft.                           | Standard  | 25/cs |
| HCSU4515S | 25-ft.                           | Universal | 25/cs |

#### **Pediatric**

| Item No.  | <b>Tubing length</b> | Connector | Pkg.  |
|-----------|----------------------|-----------|-------|
| HCSS4518S | 7-ft.                | Standard  | 50/cs |
| HCSU4518S | 7-ft.                | Universal | 50/cs |

# HCSU4517S

Crush-resistant tubing This multi-channel interior feature helps to maintain a strong airflow

#### Infant

| Item No.  | <b>Tubing length</b> | Connector | Pkg.  |
|-----------|----------------------|-----------|-------|
| HCSS4517S | 7-ft.                | Standard  | 50/cs |
| HCSU4517S | 7-ft.                | Universal | 50/cs |

#### Neonatal

| Item No.  | <b>Tubing length</b> | Connector | Pkg.  |
|-----------|----------------------|-----------|-------|
| HCSS4520S | 7-ft.                | Standard  | 50/cs |
| HCSU4520S | 7-ft.                | Universal | 50/cs |

#### Preemie

| Item No.  | <b>Tubing length</b> | Connector | Pkg.  |
|-----------|----------------------|-----------|-------|
| HCSS4521S | 7-ft.                | Standard  | 50/cs |
| HCSU4521S | 7-ft.                | Universal | 50/cs |

#### **Soft-Touch oxygen cannulas**

- · Curved, non-flared nasal prongs
- · Available in a wide variety of sizes and lengths
- · Kink-resistant tubing
- Not made with natural rubber latex

#### Adult

| Item No.   | <b>Tubing length</b>                | Connector | Pkg.  |
|------------|-------------------------------------|-----------|-------|
| HCS4510    | No tubing<br>(cannula harness only) | Standard  | 50/cs |
| HCS4504B   | 4-ft.                               | Standard  | 50/cs |
| HCSU4504B  | 4-ft.                               | Universal | 50/cs |
| HCS4514    | 7-ft.                               | Standard  | 50/cs |
| HCS4516B   | 14-ft.                              | Standard  | 50/cs |
| HCSU4516B  | 14-ft.                              | Universal | 50/cs |
| HCS4515B   | 25-ft.                              | Standard  | 25/cs |
| HCSU4515B  | 25-ft.                              | Universal | 25/cs |
| HCS02KIT1* | Oxygen set up kit                   | N/A       | 10/cs |

<sup>\*</sup>Contains two adult soft-touch cannulas, Christmas tree adater, empty bubble humidifier, oxygen tubing connector and 25-ft. green oxygen tubing



#### **Pediatric**

| Item No. | Tubing length | Connector | Pkg.  |  |
|----------|---------------|-----------|-------|--|
| HCS4518  | 7-ft.         | Standard  | 50/cs |  |
| HCSU4518 | 7-ft.         | Universal | 50/cs |  |

#### Infant

| Item No. | Tubing length | Connector | Pkg.  |
|----------|---------------|-----------|-------|
| HCS4517  | 7-ft.         | Standard  | 50/cs |
| HCSU4517 | 7-ft.         | Universal | 50/cs |



- Kink-resistant tubing with standard and universal oxygen connectors
- Straight non-flared nasal prongs
- Not made with natural rubber latex

| Item No.  | Tubing length | Connector | Pkg.  |
|-----------|---------------|-----------|-------|
| HCS4511B  | 7-ft.         | Standard  | 50/cs |
| HCSU4511B | 7-ft.         | Universal | 50/cs |
| HCS4512B  | 14-ft.        | Standard  | 50/cs |
| HCSU4512B | 14-ft.        | Universal | 50/cs |
| HCS4513B  | 25-ft.        | Standard  | 25/cs |
| HCSU4513B | 25-ft.        | Universal | 25/cs |



#### **High-Flow SuperSoft oxygen cannulas**

- Adult liter flow from 1-15 LPM
- For use with high-flow bubble humidifiers
- Not made with natural rubber latex

#### Green tubing

| Item No.    | <b>Tubing length</b> | Connector | Pkg.  |
|-------------|----------------------|-----------|-------|
| HCS4514HFG  | 7-ft.                | Standard  | 25/cs |
| HCS4516HFG  | 14-ft.               | Standard  | 25/cs |
| HCS4515HFG  | 25-ft.               | Standard  | 10/cs |
| HCS4514HFKG | High flow oxygen kit | Standard  | 7/cs  |



#### **Soft-Touch oxygen cannulas** with foam ear covers

- Pre-attached foam ear covers to help reduce skin irritation and increase patient comfort
- Curved, non-flared nasal prongs
- · Kink-resistant tubing
- Not made with natural rubber latex

#### Clear tubing

| Item No.  | <b>Tubing length</b> | Connector | Pkg.  |
|-----------|----------------------|-----------|-------|
| HCS4514FM | 7-ft.                | Standard  | 25/cs |
| HCS4516FM | 14-ft.               | Standard  | 25/cs |
| HCS4515FM | 25-ft.               | Standard  | 25/cs |



HCS4514FM

#### **SuperSoft CO<sub>2</sub> cannulas**

Our comprehensive line of CO<sub>2</sub> nasal sampling cannulas features a "Pull N' Peel" design to resist tangling, while also providing accurate, uninterrupted CO₂ sampling and oxygen delivery. Universal O₂ connecters save you from searching for adaptors.

- Universal O<sub>2</sub> connecter is designed to fit any flowmeter without the need for adaptors
- Male and female luer connecters are available to meet caregiver needs

#### CO<sub>2</sub> Nasal sampling cannulas with male luer

| Item No.  | Description      | Gas sampling<br>line length | UOM    |
|-----------|------------------|-----------------------------|--------|
| HCS4560S  | Adult, 7-ft.     | 4" Pigtail                  | 25/cs  |
| HCS4562S  | Pediatric, 7-ft. | 7-ft.                       | 25/cs  |
| HCS4568BS | Adult, 7-ft.     | 7-ft.                       | 100/cs |
| HCS4568S  | Adult, 7-ft.     | 7-ft.                       | 25/cs  |
| HCS4571S  | Adult, 7-ft.     | 7-ft.                       | 25/cs  |
| HCS4565S  | Adult, 10-ft.    | 10-ft.                      | 25/cs  |
| HCS4565BS | Adult, 10-ft.    | 10-ft.                      | 100/cs |
| HCS4566BS | Adult, 14-ft.    | 14-ft.                      | 100/cs |
| HCS4566S  | Adult, 14-ft.    | 14-ft.                      | 25/cs  |

#### CO<sub>2</sub> Nasal sampling cannulas with female luer

| Item No.   | Description      | Gas sampling<br>line length | UOM    |
|------------|------------------|-----------------------------|--------|
| HCS4561S   | Adult, 7-ft.     | 4" Pigtail                  | 25/cs  |
| HCS4563S   | Pediatric, 7-ft. | 7-ft.                       | 25/cs  |
| HCS4568FBS | Adult, 7-ft.     | 7-ft.                       | 100/cs |
| HCS4568FS  | Adult, 7-ft.     | 7-ft.                       | 25/cs  |
| HCS4572S   | Adult, 7-ft.     | 7-ft.                       | 25/cs  |
| HCS4565FS  | Adult, 10-ft.    | 10-ft.                      | 25/cs  |

#### CO<sub>2</sub> Oral sampling cannulas with male luer

| Item No. | Description  | Gas sampling<br>line length | UOM   |
|----------|--------------|-----------------------------|-------|
| HCS4573S | Adult, 7-ft. | 7-ft.                       | 25/cs |


#### CO<sub>2</sub> Oral sampling cannulas with female luer

| Item No. | Description  | Gas sampling<br>line length | UOM   |
|----------|--------------|-----------------------------|-------|
| HCS4574S | Adult, 7-ft. | 7-ft.                       | 25/cs |



Oxygen delivery through the bottom ports









Oral scoop sampling for mouth-breathing patients.



4-port nasal sampling allows for dual nare CO, sampling and 0, delivery.

#### MicroFilter CO<sub>2</sub> sampling

- Intended for short- term applications
- Proven to measure end-tidal CO<sub>2</sub> within 2 mmHg of CapnoLine® Plus
- 4 port design available with or without oral sampling scoop

## Accurate CO<sub>2</sub> sampling you can trust

#### Cannulas

| Item No    | Sampling<br>type | Size      | Universal<br>tubing<br>length | MicroFilter<br>line | Pkg.   |
|------------|------------------|-----------|-------------------------------|---------------------|--------|
| HCS4607M0  | Oral/ nasal      | Adult     | 7-ft                          | 7-ft                | 25/cs  |
| HCS4607MN  | Nasal            | Adult     | 7-ft                          | 7-ft                | 25/cs  |
| HCS4707MN  | Nasal            | Pediatric | 7-ft                          | 7-ft                | 25/cs  |
| HCS4607MOB | Oral/ nasal      | Adult     | 7-ft                          | 7-ft                | 100/cs |
| HCS4607MNB | Nasal            | Adult     | 7-ft                          | 7-ft                | 100/cs |
| HCS4614M0  | Oral/ nasal      | Adult     | 14-ft                         | 14-ft               | 25/cs  |
| HCS4614MN  | Nasal            | Adult     | 14-ft                         | 14-ft               | 25/cs  |
| HCS4714MN  | Nasal            | Pediatric | 14-ft                         | 14-ft               | 25/cs  |



#### **ET Tube adapters**

• In line sampling line with 22mm ID x 15mm ID adapter intended for short- term applications for adult and pediatrics

| Item No.  | MicroFilter line | Pkg.  |
|-----------|------------------|-------|
| HCS4607ET | 7-ft             | 25/cs |

#### **Accessories**

• In line sampling line intended for short-term application for adult and pediatric

| Item No.  | MicroFilter line | Pkg.   |
|-----------|------------------|--------|
| HCS4607GS | 7-ft             | 25/ cs |



#### CO<sub>2</sub> masks

#### Capnoxygen® mask by Southmedic

• Closed mask allows for oxygen delivery and end-tidal CO₂ monitoring

#### With female pigtail adapter

| Item No. | Size       | Universal<br>tubing length | Female pigtail adapter length | Pkg.  |
|----------|------------|----------------------------|-------------------------------|-------|
| C0201ML  | Adult      | 7-ft.                      | 8-in.                         | 25/bx |
| C0203ML  | Adult plus | 7-ft.                      | 8-in.                         | 25/bx |
| C0204ML  | Pediatric  | 7-ft.                      | 8-in.                         | 25/bx |

#### With male gas sampling line

| Item No.   | Size       | Universal tubing length | Male gas sampling<br>line length | Pkg.  |
|------------|------------|-------------------------|----------------------------------|-------|
| CO201SLMML | Adult      | 7-ft.                   | 8-ft.                            | 25/bx |
| CO203SLMML | Adult plus | 7-ft.                   | 8-ft.                            | 25/bx |



#### **Oxygen tubing**

- Unique, relaxed, kink-resistant, low-memory material helps prevent kinking
- Selected products in green color to provide greater visual contrast to assist in patient safety
- For standard and high flow applications
- Not made with natural rubber latex

#### Clear

| Item No.   | Tubing length   | Connector | Pkg   |
|------------|-----------------|-----------|-------|
| HCS4507    | 7-ft. (2.1 m)   | Standard  | 50/cs |
| HCSU4507   | 7-ft. (2.1 m)   | Universal | 50/cs |
| HCS4524    | 14-ft. (4.3 m)  | Standard  | 50/cs |
| HCSU4524   | 14-ft. (4.3 m)  | Universal | 50/cs |
| HCS4525    | 25-ft. (7.6 m)  | Standard  | 25/cs |
| HCSU4525   | 25-ft. (7.6 m)  | Universal | 25/cs |
| HCS4540FT  | 40-ft. (1.0 m)  | Standard  | 10/cs |
| HCS455010  | 50-ft. (15.2 m) | Standard  | 10/cs |
| HCSU455010 | 50-ft. (15.2 m) | Universal | 10/cs |

#### Green

| Item No.   | Tubing length   | Connector | Pkg   |
|------------|-----------------|-----------|-------|
| HCS4507G   | 7-ft. (2.1 m)   | Standard  | 50/cs |
| HCSU4507G  | 7-ft. (2.1 m)   | Universal | 50/cs |
| HCS4524G   | 14-ft. (4.3 m)  | Standard  | 50/cs |
| HCSU4524G  | 14-ft. (4.3 m)  | Universal | 50/cs |
| HCS4525G   | 25-ft. (7.6 m)  | Standard  | 10/cs |
| HCSU4525G  | 25-ft. (7.6 m)  | Universal | 25/cs |
| HCS4540FTG | 40-ft. (1.0 m)  | Standard  | 25/cs |
| HCS4550GB  | 50-ft. (15.2 m) | Standard  | 10/cs |
| HCSU4550GB | 50-ft. (15.2 m) | Universal | 10/cs |










#### **Corrugated aerosol tubing**

| Item No. | Description   | Pkg   |
|----------|---|-------|
| HSK600   | Corrugated, segmented tubing, 100-ft. (30.5 m), opaque with cuttable cuff every 6-in. (15.2 cm) | 1/cs  |
| HSK626   | Corrugated, segmented tubing, 6-ft. (1.8 m), opaque   | 50/cs |
| HCS4541  | Corrugated tubing, 6-in. (15.2 cm), blue  | 50/cs |



HCS4541



HSK600

#### Oxygen masks

Choose from a large selection of masks for all your oxygen-delivery needs. All oxygen masks feature an elongated fit for optimal patient fit.

- Soft, anatomical form for patient comfort
- A secure fit with elastic straps
- Convenient choice of universal or standard connectors
- Not made with natural rubber latex
- All oxygen masks comes with 7' of 3 channel kink-resistant clear tubing

#### Adult

| Item No.   | Description  | Connector | Reservoir<br>bag | Safety<br>vent | UOM   |
|------------|--|-----------|------------------|----------------|-------|
| HCS4600B   | Medium concentration                                       | Standard  |                  |                | 50/cs |
| HCSU4600B  | Medium concentration                                       | Universal |                  |                | 50/cs |
| HCS4640B   | Non-rebreather   | Standard  | •                | •              | 50/cs |
| HCSU4640B  | Non-rebreather   | Universal | •                | •              | 50/cs |
| HCS120E    | Venturi mask with 6 jet dilluters                          | Standard  |                  |                | 50/cs |
| HCSU120E   | Venturi mask with 6 jet dilluters                          | Universal |                  |                | 50/cs |
| HCS121E    | Venturi mask with single jet dilluter                      | Standard  |                  |                | 50/cs |
| HCSU121E   | Venturi mask with single jet dilluter                      | Universal |                  |                | 50/cs |
| HCS123E    | Venturi mask with single jet dilluter and corrugate tubing | Standard  |                  |                | 50/cs |
| HCS4670B   | Non-rebreather   | Standard  | •                |                | 50/cs |
| HCSU4670B  | Non-rebreather   | Universal | •                |                | 50/cs |
| HCS4680B*  | 3-In-1 mask  | Standard  | •                | •              | 50/cs |
| HCSU4680B* | 3-In-1 mask  | Universal | •                | •              | 50/cs |
|            |  |           |                  |                |       |



HCS121E

HCSU4640B Non-rebreather with safety vent (Reservoir bag and tubing not pictured)



HCS4670B Non-rebreather (Reservoir bag and tubing not pictured)

#### **Pediatric**

| Item No.   | Description          | Connector | Reservoir<br>bag | Safety<br>vent | UOM   |
|------------|----------------------|-----------|------------------|----------------|-------|
| HCS4601B   | Medium concentration | Standard  |                  | •              | 50/cs |
| HCSU4601B  | Medium concentration | Universal |                  | •              | 50/cs |
| HCS4642B   | Non-rebreather       | Standard  | •                | •              | 50/cs |
| HCSU4642B  | Non-rebreather       | Universal | •                | •              | 50/cs |
| HCS4675B   | Non-rebreather       | Standard  | •                |                | 50/cs |
| HCSU4675B  | Non-rebreather       | Universal | •                |                | 50/cs |
| HCS4685B*  | 3-ln-1 mask          | Standard  | •                | •              | 50/cs |
| HCSU4685B* | 3-In-1 mask          | Universal | •                | •              | 50/cs |



Comes with bag and tubing

<sup>\*3-</sup>in-1 masks come packaged with two 1-way valves, reservoir bag, 7' of oxygen tubing,oxygen connector, and oxygen mask

#### **Aerosol masks**

 $\label{thm:extra} \mbox{Extra patient comfort comes from an anatomical form}$ and a soft, light and comfortable vinyl material. All aerosol masks feature an elongated fit for optimal patient fit.

- A secure fit with elastic straps
- Standard 22 mm connector
- Not made with natural rubber latex

#### Adult

| Item No. | Description                          | Pkg.  |
|----------|--------------------------------------|-------|
| HCS4630B | Aerosol mask with straight connector | 50/cs |
| HCS4630  | Aerosol mask with swivel connector   | 50/cs |

#### **Pediatric**

| Item No. | Description                        | Pkg.  |
|----------|------------------------------------|-------|
| HCS4631B | Aerosol mask with swivel connector | 50/cs |

#### **Tent and trach masks**

#### Adult

| Item No.   | Description                          | Pkg.  |
|------------|--------------------------------------|-------|
| HCS4632B   | Face tent with 22 mm connector       | 50/cs |
| HCS4632BCT | Face tent with 6 ft corrugate tubing | 25/cs |
| HCS245E    | Trach mask                           | 50/cs |

#### Pediatric

| Item No. | Description  | Pkg.  |
|----------|--|-------|
| HCS4621B | Trach mask   | 50/cs |
| HCS4761  | 22 mm aerosol adapter,<br>works with HCS4632B and HCS4621B | 50/cs |

#### SideStream

| Item No.  | Description  | Pkg.  |
|-----------|--|-------|
| 1044328ML | SideStream Tucker the Turtle pediatric aerosol mask* | 50/cs |
| HS880ML   | SideStream Sami the Seal pediatric aerosol mask**    | 50/cs |

<sup>\*</sup>only SideStream disposable and reusable

<sup>\*\*</sup>SideStream and SideStream Plus







HCS4630B



HCS4631B



HCS4630





HCS245E

#### **Bubble humidifiers**

Choose from a wide variety of sizes and configurations.

#### Prefilled

| Item No. | Description          | Pkg.  |
|----------|----------------------|-------|
| HCS00350 | 350 mL, with adaptor | 20/cs |
| HCS00550 | 550 mL, with adaptor | 12/cs |
| HCS00640 | 650 mL, with adaptor | 10/cs |

#### Empty bubble humidifier

| Item No.    | Description                             | Pkg.  |
|-------------|---|-------|
| HCS64377    | 300 mL, 6 PSI                           | 50/cs |
| VTB6900     | 350 mL, hi-flow 6 PSI                   | 10/cs |
| HCSO2KIT1*  | Oxygen setup kit                        | 10/cs |
| HCS4514HFG* | High flow kit: 1 HCS4514HFG*, 1 VTB6900 | 7/cs  |

 $<sup>^{\</sup>star}$ Contains two adult soft-touch cannulas, christmas tree adapter, empty bubble humidifier, oxygen tubing connector, and 25' green oxygen tubing.





HCS00350

HCS00550





VTB6900

Medline can customize an oxygen set up kit to your needs, contact your Medline

Representative.

HCS64377

#### **Accessories**

Our complete selection includes a wide variety of everyday connectors and adaptors.

| Item No. | Description   | Pkg.   |
|----------|---|--------|
| HCS4700G | Green nut nipple adapter                              | 50/cs  |
| HCS4700Y | Yellow nut nipple adapter                             | 50/cs  |
| HCS4700C | Clear nut nipple adapter                              | 50/cs  |
| HCS4701  | Rigid connector for oxygen tubing                     | 200/cs |
| HCS4702  | Swivel connector for oxygen tubing                    | 25/cs  |
| HCS0440  | EZ wrap cushions for cannulas                         | 50/cs  |
| AGISH749 | Y splitter connector                                  | 5/cs   |
| DYNJ750  | Oxygen nose clip                                      | 100/cs |
| HCS1642  | Pressure line adapter, 5–7 mm ID                      | 50/cs  |
| HCS2503  | Disposable water trap                                 | 25/cs  |
| HCS974   | Aerosol drainage bag with "T" adapter hanger          | 50/cs  |
| HCS974Y  | Aerosol drainage bag with "Y" adapter hanger          | 50/cs  |
| HCS676   | 15" humidifier connector                              | 50/cs  |
| HCS1790  | 21" humidifier connector                              | 50/cs  |
| HCS8986  | Humidifier connector with rigid black elbow connector | 1/ea   |
| HSK962A  | Circuit connector, 22 OD / 15 ID 22 OD                | 50/cs  |
| HCS4672T | Aerosol tee adapter                                   | 50/cs  |





## Pulse oximeters

| Fingertip pulse oximeters  | 57 |
|----------------------------|----|
| Handheld oximeters         | 58 |
| Pulse oximeter probes      | 58 |
| Pulse oximeter accessories | 59 |

#### **Fingertip pulse oximeters**

Get quick and precise SpO<sub>2</sub> and pulse readings in seconds with our high-quality fingertip pulse oximeters. Small and lightweight, these clip sensors automatically take readings when a finger is inserted.

- Fit a wide range of finger sizes
- · Long battery life and are solidly constructed for worry-free use
- All models include AAA batteries, lanyard and warranty





#### Basic finger pulse oximeter

- Small and lightweight, this clip sensor automatically takes readings when a finger is inserted and gives a precise reading in seconds
- Features a digital LCD screen that displays SpO2, pulse rate and pulse bar and lights up for easy reading in all lighting conditions
- Includes AAA batteries, lanyard and warranty
- Rx only
- · 1 year warranty

HCSM70B, 1/ea



#### Soft-touch finger pulse oximeter

- Soft padding enhances patient comfort
- SpO<sub>2</sub>, pulse rate, perfusion index, pulse bar, and waveforms
- Backlit OLED display with 7 viewing modes
- · Adjustable brightness
- RX Only
- 1 year warranty

HCSM70S, 1/ea



#### Pediatric finger pulse oximeter

- · Child-friendly frog design makes kids smile
- SpO<sub>2</sub>, pulse rate, pulse bar, plethysmograph
- Patented multi-directional display has 6 viewing modes
- Bright, digital OLED with adjustable readout sizes
- Soft padding for patient comfort
- · Hinged battery cover
- 2 year warranty

HCSM70P, 1/ea

#### **Handheld oximeters**

Dependable and easy-to-use pulse oximeters with alarms for spot checks and continuous monitoring.

#### Medline MD300

- SpO₂ and pulse rate
- Adjustable display: 2 readout sizes and 10 brightness levels
- Includes 1 adult finger probe
- Automatic power off
- Includes protective silicone cover

| Item No. | Description  | Pkg  |
|----------|--|------|
| HCSMD300 | Handheld spot check oximeter with 30" (70 cm) adult finger probe | 1/ea |



#### **Medline MD400**

- SpO<sub>2</sub>, pulse rate, pulse bar, waveform
- Adjustable audible and visual alarms
- Two display modes
- Comfortable finger clip probe
- Data transfer to computer via USB cable
- Low battery indicator

| ı | tem No.  | Description  | Pkg  |
|---|----------|--|------|
| ŀ | HCSMD400 | Handheld continuous economy pulse oximeter with adult finger probe | 1/ea |



#### **Replacement probes**

| Item No.  | Description  | Pkg  |
|-----------|--|------|
| HCSMD400A | Adult fingertip style, for HCSMD400 and HCSMD300     | 1/ea |
| HCSMD400P | Pediatric fingertip style, for HCSMD400 and HCSMD300 | 1/ea |
| HCSMD400I | Infant wrap style, for<br>HCSMD400 and HCSMD300      | 1/ea |







## Resuscitation

| Resuscitation bags, self inflating | 61-67 |
|------------------------------------|-------|
| Colorimetric CO₂ detectors         | 68    |
| Resuscitation masks & accessories  | 69    |
| Mouth-to-mask resuscitator         | 69    |
| Hyperinflation bags                | 70    |

#### **Manual resuscitators**

Innovations provide clinicians with improved ease of use. They are available in infant, pediatric and adult sizes, and come in drawstring bag color-coded to indicate size.

#### A better bag

- Delivers up to 99% FiO2
- Maneuverable the body of the bag swivels 360° for multiple clinician hand positions, helping to reduce clinician fatigue
- Translucent the translucent, lightweight bag body and components allow for visualization of vomitus or occlusions in the system
- Textured provides an enhanced, non-slip grip
- Repositionable the patient port swivels 360°, allowing the clinician freedom to move around the patient without placing excessive torque or stress on the tracheal tube or mask
- Manometer allows for a clear visual reading of the centimeters of water pressure being delivered to the patient

#### An improved mask

• Firmly sealing-premium material provides a soft cushion and a firm seal

|                              | Adult<br>(mask size #5)     | Small adult<br>(mask size #5) | Pediatric<br>(mask size #3) | Infant<br>(mask size #2)           |
|------------------------------|-----------------------------|-------------------------------|-----------------------------|------------------------------------|
| Bag volume                   | 1,900 mL                    | 1,500 mL                      | 1,000 mL                    | 300 mL                             |
| Stroke volume                | 950-960 mL                  | 640-660 mL                    | 600-620 mL                  | 140-170 mL                         |
| Dead space with mask         | 10 mL                       | 10 mL                         | 10 mL                       | 10 mL                              |
| Inspiratory resistance       | <5 cm H₂O                   | <5 cm H₂O                     | <5 cm H <sub>2</sub> O      | <5 cm H₂O                          |
| Exhalation resistance        | $<$ 5 cm $H_2O$ 50 LPM flow | $<$ 5 cm $H_2O$ 50 LPM flow   | $<$ 5 cm $H_2O$ 50 LPM flow | <5 cm H <sub>2</sub> O 50 LPM flow |
| Attainable delivery pressure | ≥80 cm H <sub>2</sub> 0     | ≥80 cm H <sub>2</sub> 0       | ≥80 cm H <sub>2</sub> 0     | ≥80 cm H <sub>2</sub> 0            |
| Ventilatory frequency        | ≥60 cm bpm (no load)        | ≥60 cm bpm (no load)          | ≥60 cm bpm (no load)        | ≥90 cm bpm (no load)               |
| Body mass use range          | ≥40 kg (88 lbs)             | ≥40 kg (88 lbs)               | ≤40 kg (88 lbs)             | ≤10 kg (22 lbs)                    |
| Minimal tidal volume         | 600 mL                      | 600 mL                        | 500 mL                      | 150 mL                             |



#### **Resuscitation bags**

#### **Adult resuscitation bags**



| Item No.     | Description  | Mask Size |
|--------------|--|-----------|
| CPRM1116     | Bag reservoir, medium adult mask, and 7-ft. oxygen tubing  | 5         |
| CPRM11163    | Bag reservoir, small, medium and large adult masks, and 7-ft. oxygen tubing                              | 4,5,6     |
| CPRM1116B    | Bag reservoir, pre-inflated medium adult mask, and 7-ft. oxygen tubing                                   | 5         |
| CPRM1116C    | Bag reservoir, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                                  | 5         |
| CPRM1116C2   | Bag reservoir, CO₂ detector, medium and large adult masks, and 7-ft. oxygen tubing                       | 5, 6      |
| CPRM1116F    | Bag reservoir, filter, medium adult mask, and 7-ft. oxygen tubing  | 5         |
| CPRM1116FC   | Bag reservoir, filter, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                          | 5         |
| CPRM1116FM   | Bag reservoir, manometer, filter, medium adult mask, and 7-ft. oxygen tubing                             | 5         |
| CPRM1116FMT  | Bag reservoir, manometer, filter, medium adult mask, and 14-ft. oxygen tubing                            | 5         |
| CPRM1116FP   | Bag reservoir, PEEP valve, filter, medium adult mask, and 7-ft. oxygen tubing                            | 5         |
| CPRM1116FPC  | Bag reservoir, PEEP valve, filter, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing              | 5         |
| CPRM1116FPM  | Bag reservoir, manometer, PEEP valve, filter, medium adult mask, and 7-ft. oxygen tubing                 | 5         |
| CPRM1116L    | Bag reservoir, large adult mask, and 7-ft. oxygen tubing   | 5         |
| CPRM1116M    | Bag reservoir, manometer, pre-inflated medium adult mask, and 7-ft. oxygen tubing                        | 5         |
| CPRM1116MC   | Bag reservoir, manometer, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                       | 5         |
| CPRM1116P    | Bag reservoir, PEEP valve, medium adult mask, and 7-ft. oxygen tubing                                    | 5         |
| CPRM1116PC   | Bag reservoir, PEEP valve, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                      | 5         |
| CPRM1116PCM  | Bag reservoir, manometer, PEEP valve, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing           | 5         |
| CPRM1116PD5  | Bag reservoir, PEEP valve, medium adult mask, and 7-ft. oxygen tubing                                    | 5         |
| CPRM1116PM   | Bag reservoir, manometer, PEEP valve, medium adult mask, and 7-ft. oxygen tubing                         | 5         |
| CPRM1116PMT  | Bag reservoir, manometer, PEEP valve, medium adult mask, and 14-ft. oxygen tubing                        | 5         |
| CPRM1116PW   | Bag reservoir, pre-attached PEEP valve, medium adult mask, and 7-ft. oxygen tubing                       | 5         |
| CPRM1116T    | Bag reservoir, medium adult mask, and 14-ft. oxygen tubing   | 5         |
| CPRM1116WM   | Bag reservoir, pre-attached medium adult mask, and 7-ft. oxygen tubing                                   | 5         |
| CPRM1126     | Tubing reservoir, medium adult mask, and 7-ft. oxygen tubing   | 5         |
| CPRM1126C    | Tubing reservoir, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                               | 5         |
| CPRM1126F    | Tubing reservoir, filter, medium adult mask, and 7-ft. oxygen tubing                                     | 5         |
| CPRM1126FC   | Tubing reservoir, filter, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                       | 5         |
| CPRM1126FL   | Tubing reservoir, filter, large adult mask, and 7-ft. oxygen tubing                                      | 6         |
| CPRM1126FM   | Tubing reservoir, manometer, filter, medium adult mask, and 7-ft. oxygen tubing                          | 5         |
| CPRM1126FP   | Tubing reservoir, PEEP valve, filter, medium adult mask, and 7-ft. oxygen tubing                         | 5         |
| CPRM1126FPC  | Tubing reservoir, PEEP valve, filter, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing           | 5         |
| CPRM1126FPM  | Tubing reservoir, manometer, PEEP valve, filter, medium adult mask, and 7-ft. oxygen tubing              | 5         |
| CPRM1126FPWM | Tubing reservoir, manometer, pre-attached PEEP valve, filter, medium adult mask, and 7-ft. oxygen tubing | 5         |
| CPRM1126FPWT | Tubing reservoir, pre-attached PEEP valve, filter, medium adult mask, and 14-ft. oxygen tubing           | 5         |
| CPRM1126M    | Tubing reservoir, manometer, medium adult mask, and 7-ft. oxygen tubing                                  | 5         |
| CPRM1126P    | Tubing reservoir, PEEP valve, medium adult mask, and 7-ft. oxygen tubing                                 | 5         |
| CPRM1126PC   | Tubing reservoir, PEEP valve, CO₂ detector, medium adult mask, and 7-ft. oxygen tubing                   | 5         |

#### Small adult resuscitation bags

| Item No.  | Description   | Mask Size |
|-----------|---|-----------|
| CPRM4416  | Bag reservoir, medium adult mask, and 7-ft. oxygen tubing               | 5         |
| CPRM4416M | Bag reservoir, manometer, medium adult mask, and 7-ft. oxygen tubing    | 5         |
| CPRM4416P | Bag reservoir, PEEP valve, medium adult mask, and 7-ft. oxygen tubing   | 5         |
| CPRM4426  | Tubing reservoir, medium adult mask, and 7-ft. oxygen tubing            | 5         |
| CPRM4426M | Tubing reservoir, manometer, medium adult mask, and 7-ft. oxygen tubing | 5         |

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.

#### Manual resuscitator ordering guide

#### Adult resuscitation bags

|              | Rese | rvoir style | Tubing | glength | Pressure     | CO₂      | CO <sub>2</sub> |            |                   |           |      |
|--------------|------|-------------|--------|---------|--------------|----------|-----------------|------------|-------------------|-----------|------|
| Item No.     | Bag  | 40" tubing  | 7-ft.  | 14-ft.  | relief valve | detector | FLOCAP          | Peep valve | Exhalation filter | Manometer | Pkg. |
| CPRM1116     | •    |             | •      |         |              |          |                 |            |                   |           | 6/cs |
| CPRM11163    | •    |             | •      |         |              |          |                 |            |                   |           | 6/cs |
| CPRM1116B    | •    |             | •      |         |              |          |                 |            |                   |           | 6/cs |
| CPRM1116C    | •    |             | •      |         |              | •        |                 |            |                   |           | 6/cs |
| CPRM1116C2   | •    |             | •      |         |              | •        |                 |            |                   |           | 6/cs |
| CPRM1116F    | •    |             | •      |         |              |          |                 |            | •                 |           | 6/cs |
| CPRM1116FC   | •    |             | •      |         |              | •        |                 |            | •                 |           | 6/cs |
| CPRM1116FM   | •    |             | •      |         |              |          |                 |            | •                 | •         | 6/cs |
| CPRM1116FMT  | •    |             |        | •       |              |          |                 |            | •                 | •         | 6/cs |
| CPRM1116FP   | •    |             | •      |         |              |          |                 | •          | •                 |           | 6/cs |
| CPRM1116FPC  | •    |             | •      |         |              | •        |                 | •          | •                 |           | 6/cs |
| CPRM1116FPM  | •    |             | •      |         |              |          |                 | •          | •                 | •         | 6/cs |
| CPRM1116L    | •    |             | •      |         |              |          |                 |            |                   |           | 6/cs |
| CPRM1116M    | •    |             | •      |         |              |          |                 |            |                   | •         | 6/cs |
| CPRM1116MC   | •    |             | •      |         |              | •        |                 |            |                   | •         | 6/cs |
| CPRM1116P    | •    |             | •      |         |              |          |                 | •          |                   |           | 6/cs |
| CPRM1116PC   | •    |             | •      |         |              | •        |                 | •          |                   |           | 6/cs |
| CPRM1116PCM  | •    |             | •      |         |              | •        |                 | •          |                   | •         | 6/cs |
| CPRM1116PD5  | •    |             | •      |         |              |          |                 | •          |                   |           | 6/cs |
| CPRM1116PM   | •    |             | •      |         |              |          |                 | •          |                   | •         | 6/cs |
| CPRM1116PMT  | •    |             |        | •       |              |          |                 | •          |                   | •         | 6/cs |
| CPRM1116PW   | •    |             | •      |         |              |          |                 | •          |                   |           | 6/cs |
| CPRM1116T    | •    |             |        | •       |              |          |                 |            |                   |           | 6/cs |
| CPRM1116WM   | •    |             | •      |         |              |          |                 |            |                   | •         | 6/cs |
| CPRM1126     |      | •           | •      |         |              |          |                 |            |                   |           | 6/cs |
| CPRM1126C    |      | •           | •      |         |              | •        |                 |            |                   |           | 6/cs |
| CPRM1126F    |      | •           | •      |         |              |          |                 |            | •                 |           | 6/cs |
| CPRM1126FC   |      | •           | •      |         |              | •        |                 |            | •                 |           | 6/cs |
| CPRM1126FL   |      | •           | •      |         |              |          |                 |            | •                 |           | 6/cs |
| CPRM1126FM   |      | •           | •      |         |              |          |                 |            | •                 | •         | 6/cs |
| CPRM1126FP   |      | •           | •      |         |              |          |                 | •          | •                 |           | 6/cs |
| CPRM1126FPC  |      | •           | •      |         |              | •        |                 | •          | •                 |           | 6/cs |
| CPRM1126FPM  |      | •           | •      |         |              |          |                 | •          | •                 | •         | 6/cs |
| CPRM1126FPWM |      | •           | •      |         |              |          |                 | •          | •                 | •         | 6/cs |
| CPRM1126FPWT |      | •           |        | •       |              |          |                 |            | •                 |           | 6/cs |
| CPRM1126M    |      | •           | •      |         |              |          |                 |            |                   | •         | 6/cs |
| CPRM1126P    |      | •           | •      |         |              |          |                 | •          |                   |           | 6/cs |
| CPRM1126PC   |      | •           | •      |         |              | •        |                 | •          |                   |           | 6/cs |

#### Small adult resuscitation bags

|           | Rese | rvoir style | Tubing | g length | Pressure | CO <sub>2</sub> | CO <sub>2</sub> |            |                   |           |      |
|-----------|------|-------------|--------|----------|----------|-----------------|-----------------|------------|-------------------|-----------|------|
| Item No.  | Bag  | 40" tubing  | 7-ft.  | 14-ft.   |          |                 |                 | Peep valve | Exhalation filter | Manometer | Pkg. |
| CPRM4416  | •    |             | •      |          |          |                 |                 |            |                   |           | 6/cs |
| CPRM4416M | •    |             | •      |          |          |                 |                 |            |                   | •         | 6/cs |
| CPRM4416P | •    |             | •      |          |          |                 |                 | •          |                   |           | 6/cs |
| CPRM4426  |      | •           | •      |          |          |                 |                 |            |                   |           | 6/cs |
| CPRM4426M |      | •           | •      |          |          |                 |                 |            |                   | •         | 6/cs |

#### **Resuscitation bags**



#### **Pediatric resuscitation bags**

| Item No.     | Description  | Mask Size |
|--------------|--|-----------|
| CPRM2216     | Bag reservoir, pediatric mask, and 7-ft. oxygen tubing   | 3         |
| CPRM2216C    | Bag reservoir, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing   | 3         |
| CPRM2216F    | Bag reservoir, filter, pediatric mask, and 7-ft. oxygen tubing   | 3         |
| CPRM2216FC   | Bag reservoir, filter, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing                                 | 3         |
| CPRM2216FM   | Bag reservoir, manometer, filter, pediatric mask, and 7-ft. oxygen tubing                                    | 3         |
| CPRM2216FM3  | Bag reservoir, manometer, filter, infant, child and small adult masks, and 7-ft. oxygen tubing               | 2, 3, 4   |
| CPRM2216FP   | Bag reservoir, PEEP valve, filter, pediatric mask, and 7-ft. oxygen tubing                                   | 3         |
| CPRM2216FPC  | Bag reservoir, PEEP valve, filter, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing                     | 3         |
| CPRM2216FPM  | Bag reservoir, manometer, PEEP valve, filter, pediatric mask, and 7-ft. oxygen tubing                        | 3         |
| CPRM2216FPM3 | Bag reservoir, manometer, PEEP valve, filter, infant, child and small adult masks, and 7-ft. oxygen tubing   | 2, 3, 4   |
| CPRM2216M    | Bag reservoir, manometer, pediatric mask, and 7-ft. oxygen tubing  | 3         |
| CPRM2216M3   | Bag reservoir, manometer, infant, child and small adult masks, and 7-ft. oxygen tubing                       | 2, 3, 4   |
| CPRM2216P    | Bag reservoir, PEEP valve, pediatric mask, and 7-ft. oxygen tubing   | 3         |
| CPRM2216P3M3 | Bag reservoir, manometer, 30mm PEEP valve, infant, pediatric, and small adult masks, and 7-ft. oxygen tubing | 2, 3, 4   |
| CPRM2216PM6  | Bag reservoir, PEEP valve, manometer, pre-inflated large adult mask, and 7-ft. oxygen tubing                 | 6         |
| CPRM2216PC   | Bag reservoir, PEEP valve, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing                             | 3         |
| CPRM2216PCM  | Bag reservoir, manometer, PEEP valve, CO <sub>2</sub> detector, pediatric mask, and 7-ft. oxygen tubing      | 3         |
| CPRM2216PM   | Bag reservoir, manometer, PEEP valve, pediatric mask, and 7-ft. oxygen tubing                                | 3         |
| CPRM2216PMT  | Bag reservoir, manometer, PEEP valve, pediatric mask, and 14-ft. oxygen tubing                               | 3         |
| CPRM22162    | Bag reservoir, infant and pediatric mask, and 7-ft. oxygen tubing  | 2, 3      |
| CPRM22163    | Bag reservoir, neonatal, infant and pediatric mask, and 7-ft. oxygen tubing                                  | 1, 2, 3   |
| CPRM2226     | Tubing reservoir, pediatric mask, and 7-ft. oxygen tubing  | 3         |
| CPRM2226B    | Tubing reservoir, pre-inflated pediatric mask, and 7-ft. oxygen tubing                                       | 3         |
| CPRM2226C    | Tubing reservoir, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing                                      | 3         |
| CPRM2226F    | Tubing reservoir, filter, pediatric mask, and 7-ft. oxygen tubing  | 3         |
| CPRM2226F2   | Tubing reservoir, filter, pediatric and small adult masks, and 7-ft. oxygen tubing                           | 3, 4      |
| CPRM2226FC   | Tubing reservoir, filter, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing                              | 3, 4      |
| CPRM2226FM   | Tubing reservoir, manometer, filter, pediatric mask, and 7-ft. oxygen tubing                                 | 3         |
| CPRM2226FM3  | Tubing reservoir, manometer, filter, infant, child and small adult masks, and 7-ft. oxygen tubing            | 2, 3, 4   |
| CPRM2226FP   | Tubing reservoir, PEEP valve, filter, pediatric mask, and 7-ft. oxygen tubing                                | 3         |
| CPRM2226FPC  | Tubing reservoir, PEEP valve, filter, CO₂ detector, pediatric mask, and 7-ft. oxygen tubing                  | 3         |
| CPRM2226FPM  | Tubing reservoir, manometer, PEEP valve, filter, pediatric mask, and 7-ft. oxygen tubing                     | 3         |
| CPRM2226FPWM | Tubing reservoir, manometer, pre-attached PEEP valve, filter, pediatric mask, and 7-ft. oxygen tubing        | 3         |
| CPRM2226M    | Tubing reservoir, manometer, pediatric mask, and 7 ft. oxygen tubing   | 3         |
| CPRM2226P    | Tubing reservoir, PEEP valve, pediatric mask, and 7 ft. oxygen tubing  | 3         |
| CPRM2226PC   | Tubing reservoir, PEEP valve, CO₂ detector, pediatric mask, and 7 ft. oxygen tubing                          | 3         |
| CPRM2226PW   | Tubing reservoir, pre-attached PEEP valve, pediatric mask, and 7 ft. oxygen tubing                           | 3         |

If you do not see your preferred configuration, **Medline can customize a resuscitator bag to your preference.** 

#### Manual resuscitator ordering guide

#### **Pediatric resuscitation bags**

|              | Rese | rvoir style | Tubing | length |                        | 60                          | 60            |            |                   |           |      |
|--------------|------|-------------|--------|--------|------------------------|-----------------------------|---------------|------------|-------------------|-----------|------|
| Item No.     | Bag  | 40" tubing  | 7-ft.  | 14-ft. | Pressure relief valve  | CO <sub>2</sub><br>detector | CO₂<br>FLOCAP | Peep valve | Exhalation filter | Manometer | Pkg. |
| CPRM2216     | •    |             | •      |        | 40 cm H₂0              |                             |               |            |                   |           | 6/cs |
| CPRM2216C    | •    |             | •      |        | 40 cm H <sub>2</sub> 0 | •                           |               |            |                   |           | 6/cs |
| CPRM2216F    | •    |             | •      |        | 40 cm H₂0              |                             |               |            | •                 |           | 6/cs |
| CPRM2216FC   | •    |             | •      |        | 40 cm H <sub>2</sub> 0 | •                           |               |            | •                 |           | 6/cs |
| CPRM2216FM   | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            | •                 | •         | 6/cs |
| CPRM2216FM3  | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            | •                 | •         | 6/cs |
| CPRM2216FP   | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          | •                 |           | 6/cs |
| CPRM2216FPC  | •    |             | •      |        | 40 cm H <sub>2</sub> 0 | •                           |               | •          | •                 |           | 6/cs |
| CPRM2216FPM  | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          | •                 | •         | 6/cs |
| CPRM2216FPM3 | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          | •                 | •         | 6/cs |
| CPRM2216M    | •    |             | •      |        | 40 cm H₂0              |                             |               |            |                   | •         | 6/cs |
| CPRM2216M3   | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            |                   | •         | 6/cs |
| CPRM2216P    | •    |             | •      |        | 40 cm H₂0              |                             |               | •          |                   |           | 6/cs |
| CPRM2216P3M3 | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          |                   | •         | 6/cs |
| CPRM2216PC   | •    |             | •      |        | 40 cm H₂0              | •                           |               | •          |                   |           | 6/cs |
| CPRM2216PCM  | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          |                   | •         | 6/cs |
| CPRM2216PM   | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          |                   | •         | 6/cs |
| CPRM2216PM6  | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          |                   | •         | 6/cs |
| CPRM2216PMT  | •    |             |        | •      | 40 cm H₂0              |                             |               | •          |                   | •         | 6/cs |
| CPRM22162    | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            |                   |           | 6/cs |
| CPRM22163    | •    |             | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            |                   |           | 6/cs |
| CPRM2226     |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            |                   |           | 6/cs |
| CPRM2226B    |      | •           | •      |        | 40 cm H₂0              |                             |               |            |                   |           | 6/cs |
| CPRM2226C    |      | •           | •      |        | 40 cm H <sub>2</sub> 0 | •                           |               |            |                   |           | 6/cs |
| CPRM2226F    |      | •           | •      |        | 40 cm H₂0              |                             |               |            | •                 |           | 6/cs |
| CPRM2226F2   |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            | •                 |           | 6/cs |
| CPRM2226FC   |      | •           | •      |        | 40 cm H₂0              | •                           |               |            | •                 |           | 6/cs |
| CPRM2226FM   |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            | •                 | •         | 6/cs |
| CPRM2226FM3  |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            | •                 | •         | 6/cs |
| CPRM2226FP   |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          | •                 |           | 6/cs |
| CPRM2226FPC  |      | •           | •      |        | 40 cm H <sub>2</sub> 0 | •                           |               | •          | •                 | •         | 6/cs |
| CPRM2226FPM  |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          | •                 | •         | 6/cs |
| CPRM2226FPWM |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          | •                 | •         | 6/cs |
| CPRM2226M    |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               |            |                   | •         | 6/cs |
| CPRM2226P    |      | •           | •      |        | 40 cm H <sub>2</sub> 0 |                             |               | •          |                   |           | 6/cs |
| CPRM2226PC   |      | •           | •      |        | 40 cm H <sub>2</sub> 0 | •                           |               | •          |                   |           | 6/cs |
| CPRM2226PW   |      | •           | •      |        | 40 cm H₂0              |                             |               | •          |                   |           | 6/cs |

#### **Resuscitation bags**



#### Infant resuscitation bags

| Item No.      | Description   | Mask Size |
|---------------|---|-----------|
| CPRM3312      | Bag reservoir, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM33122     | Bag reservoir, neonatal and infant masks, and 7-ft. oxygen tubing   | 1, 2      |
| CPRM33123     | Bag reservoir, neonatal, infant and pediatric mask, and 7-ft. oxygen tubing   | 1, 2, 3   |
| CPRM3312C     | Bag reservoir, CO₂ detector, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312F     | Bag reservoir, filter, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312FC    | Bag reservoir, filter, CO₂ detector, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312FM    | Bag reservoir, manometer, filter, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3312FM123 | Bag reservoir, manometer, filter, neonatal, infant and pediatric masks, and 7-ft. oxygen tubing                                   | 1, 2, 3   |
| CPRM3312FP    | Bag reservoir, PEEP valve, filter, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312FPC   | Bag reservoir, PEEP valve, filter, CO₂ detector, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312FPM   | Bag reservoir, manometer, PEEP valve, filter, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3312FPM3  | Bag reservoir, manometer, PEEP valve, filter, pre-inflated neonatal and infant masks, and 7-ft. oxygen tubing                     | 0, 1, 2   |
| CPRM3312M     | Bag reservoir, manometer, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3312M3    | Bag reservoir, manometer, neonatal, infant and pediatric masks, and 7-ft. oxygen tubing   | 1, 2, 3   |
| CPRM3312P     | Bag reservoir, PEEP valve, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312P3M3  | Bag reservoir, manometer, 30mm PEEP valve, premie pre-inflated, neonatal valved, and infant valved masks, and 7-ft. oxygen tubing | 0, 1, 2   |
| CPRM3312PC    | Bag reservoir, PEEP valve, CO₂ detector, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312PM    | Bag reservoir, manometer, PEEP valve, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3312PM1   | Bag reservoir, manometer, PEEP valve, neonatal mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3312PM2   | Bag reservoir, manometer, PEEP valve, pre-inflated infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3312PM3   | Bag reservoir, manometer, PEEP valve, neonatal, infant and pediatric masks, and 7-ft. oxygen tubing                               | 1, 2, 3   |
| CPRM3312PMC2  | Bag reservoir, manometer, PEEP valve, CO₂ detector, neonatal and infant masks, and 7-ft. oxygen tubing                            | 1, 2      |
| CPRM3312PMTD1 | Bag reservoir, manometer, PEEP valve, neonatal mask, and 14-ft. oxygen tubing   | 1         |
| CPRM3312PMX   | Bag reservoir, manometer, PEEP valve, without mask, and 7-ft. oxygen tubing   |           |
| CPRM3312PWM   | Bag reservoir, manometer, pre-attached PEEP valve, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3322      | Tubing reservoir, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322C     | Tubing reservoir, CO₂ detector, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322F     | Tubing reservoir, filter, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322FC    | Tubing reservoir, filter, CO₂ detector, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322FM    | Tubing reservoir, manometer, filter, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3322FP    | Tubing reservoir, PEEP valve, filter, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322FPC   | Tubing reservoir, PEEP valve, filter, CO₂ detector, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322M     | Tubing reservoir, manometer, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3322MC    | Tubing reservoir, manometer, CO₂ detector, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3322MG    | Tubing reservoir, manometer, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3322P     | Tubing reservoir, PEEP valve, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322PC    | Tubing reservoir, PEEP valve, CO₂ detector, infant mask, and 7-ft. oxygen tubing  | 2         |
| CPRM3322PM    | Tubing reservoir, manometer, PEEP valve, infant mask, and 7-ft. oxygen tubing   | 2         |
| CPRM3412P     | Bag reservoir, PEEP valve, preemie and infant masks, and 7-ft. oxygen tubing  | 0, 2      |
|               |   |           |

If you do not see your preferred configuration, **Medline can customize a resuscitator bag to your preference.** 

#### Manual resuscitator ordering guide

#### Infant resuscitation bags

|               | Rese | rvoir style | Tubing | glength | _                                   |                          |               |            |                   |           |       |
|---------------|------|-------------|--------|---------|-------------------------------------|--------------------------|---------------|------------|-------------------|-----------|-------|
| Item No.      | Bag  | 40" tubing  | 7-ft.  | 14-ft.  | Pressure relief valve               | CO <sub>2</sub> detector | CO₂<br>FLOCAP | Peep valve | Exhalation filter | Manometer | Pkg.  |
| CPRM3312      | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            |                   |           | 12/cs |
| CPRM3312C     | •    |             | •      |         | 40 cm H <sub>2</sub> 0              | •                        |               |            |                   |           | 12/cs |
| CPRM3312F     | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            | •                 |           | 12/cs |
| CPRM3312FC    | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            | •                 |           | 12/cs |
| CPRM3312FM    | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            | •                 | •         | 12/cs |
| CPRM3312FM123 | •    |             | •      |         | 40 cm H₂0                           |                          |               |            | •                 | •         | 12/cs |
| CPRM3312FP    | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          | •                 |           | 12/cs |
| CPRM3312FPC   | •    |             | •      |         | 40 cm H <sub>2</sub> 0              | •                        |               | •          | •                 |           | 12/cs |
| CPRM3312FPM   | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          | •                 | •         | 12/cs |
| CPRM3312FPM3  | •    |             | •      |         | $40 \text{ cm H}_2\text{O}$         |                          |               | •          | •                 | •         | 12/cs |
| CPRM3312M     | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            |                   | •         | 12/cs |
| CPRM3312M3    | •    |             | •      |         | $40 \text{ cm H}_2\text{O}$         |                          |               |            |                   | •         | 12/cs |
| CPRM3312P     | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   |           | 12/cs |
| CPRM3312P3M3  | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PC    | •    |             | •      |         | 40 cm H <sub>2</sub> 0              | •                        |               | •          |                   |           | 12/cs |
| CPRM3312PM    | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PM1   | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PM2   | •    |             | •      |         | $40\;cm\;H_{\scriptscriptstyle 2}0$ |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PM3   | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PMC2  | •    |             | •      |         | $40\;cm\;H_{\scriptscriptstyle 2}0$ | •                        |               | •          |                   | •         | 12/cs |
| CPRM3312PMTD1 | •    |             |        | •       | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PMX   | •    |             | •      |         | $40\;cm\;H_{\scriptscriptstyle 2}0$ |                          |               | •          |                   | •         | 12/cs |
| CPRM3312PWM   | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   | •         | 12/cs |
| CPRM33122     | •    |             | •      |         | $40\;cm\;H_{\scriptscriptstyle 2}0$ |                          |               |            |                   |           | 12/cs |
| CPRM33123     | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            |                   |           | 12/cs |
| CPRM3322      |      | •           | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            |                   |           | 12/cs |
| CPRM3322C     |      | •           | •      |         | 40 cm H <sub>2</sub> 0              | •                        |               |            |                   |           | 12/cs |
| CPRM3322F     |      | •           | •      |         | $40\;cm\;H_2O$                      |                          |               |            | •                 |           | 12/cs |
| CPRM3322FC    |      | •           | •      |         | 40 cm H <sub>2</sub> 0              | •                        |               |            | •                 | •         | 12/cs |
| CPRM3322FM    |      | •           | •      |         | $40\;cm\;H_{\scriptscriptstyle 2}0$ |                          |               |            | •                 |           | 12/cs |
| CPRM3322FP    |      | •           | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          | •                 |           | 12/cs |
| CPRM3322FPC   |      | •           | •      |         | $40 \text{ cm H}_2\text{O}$         | •                        |               | •          | •                 |           | 12/cs |
| CPRM3322M     |      | •           | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            |                   | •         | 12/cs |
| CPRM3322MC    |      | •           | •      |         | 40 cm H₂0                           | •                        |               |            |                   | •         | 12/cs |
| CPRM3322MG    |      | •           | •      |         | 40 cm H <sub>2</sub> 0              |                          |               |            |                   | •         | 12/cs |
| CPRM3322P     |      | •           | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   |           | 12/cs |
| CPRM3322PC    |      | •           | •      |         | 40 cm H <sub>2</sub> 0              | •                        |               | •          |                   |           | 12/cs |
| CPRM3322PM    |      | •           | •      |         | 40 cm H₂0                           |                          |               | •          |                   | •         | 12/cs |
| CPRM3412P     | •    |             | •      |         | 40 cm H <sub>2</sub> 0              |                          |               | •          |                   |           | 12/cs |

#### **Colorimetric CO<sub>2</sub> detectors**

Medline's colorimetric CO<sub>2</sub> detectors offer a simple way to indicate the presence of carbon dioxide during a patient's breathing cycle.

It's quick response of changes in color between blue and yellow make it easy to visualize and determine whether exhaled CO<sub>2</sub> is present.

- Our CO₂ detectors are long lasting and can be used for up to 24 hours
- Both CO₂ detectors have a small internal dead space. Item No. CPRM2400 is 16 mL and CPRM2401 is 3 mL
- Offered in two sizes to meet your needs in supporting patients weighing 1 kg or greater

| Item No. | Description   | Pkg.  |
|----------|---|-------|
| CPRM2400 | Colorimetric CO₂ detector. For patients who weigh ≥10 kg (≥22-lbs.)               | 10/cs |
| CPRM2401 | Colorimetric Pediatric CO₂ detector. For patients who weigh 1–15 kg (2.2–33-lbs.) | 10/cs |





#### **Resuscitation masks and accessories**

#### **Resuscitation masks**

Our resuscitation masks provide an ergonomic cushion design for an optimal seal and clear cone to visualize patient breathing.

- Removable hook ring for hand-held or head-strap application (hook rings included on sizes 3-6)
- No valve prevents mask deflation
- · Available in all sizes, including preemie

| Item No. | I      | Description          | Pkg.  |
|----------|--------|----------------------|-------|
| CPRMAS   | K11R 1 | Neonatal, size 1     | 40/cs |
| CPRMAS   | K12R I | nfant, size 2        | 40/cs |
| CPRMAS   | K13R ( | Child, size 3        | 40/cs |
| CPRMAS   | K14R 9 | 5mall adult, size 4  | 40/cs |
| CPRMAS   | K15R 1 | Medium adult, size 5 | 20/cs |
| CPRMAS   | K16R l | _arge adult, size 6  | 20/cs |

Medline offers a wide variety of masks which are not shown here. Please contact us for more information.



#### **Resuscitation accessories**

Our high quality resuscitation accessories are designed to work with all Medline resuscitation devices and allow you to customize to patient need.

| Item No.   | Description  | Pkg.  |
|------------|--|-------|
| CPRMP19    | 5-20 cm H2O adjustable PEEP valve,<br>19 mm ID and 22 mm OD connections        | 20/cs |
| CPRMP30    | 5-20 cm H2O adjustable PEEP valve,<br>19 mm ID, 22 mm OD, and 30 mm ID adapter | 20/cs |
| CPRMFILTER | Expiratory filter  | 50/cs |



#### Mouth-to-mask resuscitator mask

This ready-to-use CPR mask is a highly effective one-way filtering system that shields the rescuer and victim from the risk of cross-contamination.

- A clear plastic design allows for easy placement, and the soft inflated cuff improves the fit and patient comfort
- Standard 15 mm fittings are compatible with any resuscitator
- Standard mask size accommodates children and adults.
- Easy-to-understand illustrations printed on insert
- · Includes hard plastic transport case

| Item No. | Description        | Pkg.  |
|----------|--------------------|-------|
| HCS64184 | With one-way valve | 12/cs |


Medline offers a wide variety of masks which are not shown here. Please contact us for more information.



#### **Hyperinflation systems**

Help prevent or treat lung collapse with a hyperinflation system.

- Bag design uses a "dipped" process to enhance tactile sensitivity and improve control
- Uniform inflation and the highest level of performance
- 100% compliant with ISO 5362 standards
- Variety of sizes to more accurately fit every patient



#### **Hyperinflation systems**

7-ft. (2.13 m) of O₂ tubing included\*

| Item No.      | Breathing bag                      | O <sub>2</sub> connector | Mask                       | Manometer | Pkg.  |
|---------------|------------------------------------|--------------------------|----------------------------|-----------|-------|
| DYNJAAHY005   | 0.5 L                              | Standard                 | N/A                        |           | 20/cs |
| DYNJAAHY600   | 0.5 L with tail-end<br>bleed valve | Universal                | N/A                        |           | 50/cs |
| DYNJAAHY30    | 0.5 L                              | Standard                 | 1 Neonatal                 |           | 20/cs |
| DYNJAAHY29    | 0.5 L                              | Universal                | 2 Infant<br>(pre-inflated) |           | 20/cs |
| DYNJAAHY01    | 1 L                                | Standard                 | N/A                        |           | 20/cs |
| DYNJAAHY601   | 1 L with tail-end<br>bleed valve   | Universal                | N/A                        |           | 50/cs |
| DYNJAAHY02    | 2 L                                | Standard                 | N/A                        |           | 20/cs |
| DYNJAAHY03    | 3 L                                | Standard                 | N/A                        |           | 20/cs |
| DYNJAAHY03M5  | 3 L                                | Universal                | 5 Adult                    |           | 20/cs |
| DYNJAAHY5005  | 0.5 L                              | Standard                 | N/A                        | •         | 50/cs |
| DYNJAAHY5053P | 0.5 L                              | Universal                | N/A                        | •         | 50/cs |
| DYNJAAHY32    | 0.5 L                              | Standard                 | 1 Neonatal                 | •         | 50/cs |
| DYNJAAHY33    | 0.5 L                              | Standard                 | 2 Infant                   | •         | 50/cs |
| DYNJAAHY35    | 0.5 L                              | Universal                | 2 Infant                   | •         | 50/cs |
| DYNJAAHY34    | 0.5 L                              | Universal                | 2 Infant<br>(pre-inflated) | •         | 50/cs |
| DYNJAAHY501   | 1 L                                | Universal                | N/A                        | •         | 50/cs |
| DYNJAAHY502   | 2 L                                | Universal                | N/A                        | •         | 50/cs |
| DYNJAAHY503   | 3 L                                | Universal                | N/A                        | •         | 50/cs |

<sup>\*</sup>Custom options are available





Available in all sizes, both inflatable and pre-inflated.



DYNJAAHY5005

## Suction catheters and trach care

| Closed suction                   | 72-74 |
|----------------------------------|-------|
| Open suction catheters           | 75    |
| Open suction catheters, sleeved  | 75    |
| Open suction kits                | 76    |
| Open suction kits with solutions | 77    |
| Sterile solutions                | 77    |
| Trach care & cleaning trays      | 78    |
| Mucus/specimen traps             | 79    |



#### **ClearPro**<sup>™</sup> **closed suction catheter**

The power to set—and follow—best practice.

ClearPro provides the performance you expect without a time-based change-out recommendation on the label. This gives you the freedom to set your own change-out policy to meet both clinical and regulatory recommendations—and drive down costs.

Plus, it includes advanced features to help provide better care:

- Isolated cleaning chamber for rapid irrigation
- Pressure wiper seal to help maintain PEEP

#### Additional details:


- Kitted with all necessary adapters and connectors
- Does not contain DEHP
- Available in French sizes 10–16 with catheter depth markings and color-coded by ISO standards

#### Endotracheal


| Item No.  | Description                | Pkg.  | Thumb valve color |
|-----------|----------------------------|-------|-------------------|
| DYNCPDS10 | Double swivel elbow, 10 Fr | 20/cs | Black             |
| DYNCPDS12 | Double swivel elbow, 12 Fr | 20/cs | White             |
| DYNCPDS14 | Double swivel elbow, 14 Fr | 20/cs | Green             |
| DYNCPDS16 | Double swivel elbow, 16 Fr | 20/cs | Orange            |
| DYNCPTP10 | T-Piece adapter, 10 Fr     | 20/cs | Black             |
| DYNCPTP12 | T-Piece adapter, 12 Fr     | 20/cs | White             |
| DYNCPTP14 | T-Piece adapter, 14 Fr     | 20/cs | Green             |
| DYNCPTP16 | T-Piece adapter, 16 Fr     | 20/cs | Orange            |

#### Tracheostomy

| Item No.   | Description                | Pkg.  | Thumb valve color |
|------------|----------------------------|-------|-------------------|
| DYNCPDS10T | Double swivel elbow, 10 Fr | 20/cs | Black             |
| DYNCPDS12T | Double swivel elbow, 12 Fr | 20/cs | White             |
| DYNCPDS14T | Double swivel elbow, 14 Fr | 20/cs | Green             |
| DYNCPDS16T | Double swivel elbow, 16 Fr | 20/cs | Orange            |
| DYNCPTP10T | T-Piece adapter, 10 Fr     | 20/cs | Black             |
| DYNCPTP12T | T-Piece adapter, 12 Fr     | 20/cs | White             |
| DYNCPTP14T | T-Piece adapter, 14 Fr     | 20/cs | Green             |
| DYNCPTP16T | T-Piece adapter, 16 Fr     | 20/cs | Orange            |



#### **Colored thumb valves**



#### **Accessories**







## **Medline closed suction catheter systems**

#### **Adult**

- Labeled for facilities to determine change out policy
- Soft tip helps reduce damage to mucous membranes
- Transparent connector for easy observation
- High-quality, translucent EVA sleeve helps reduce cross-contamination and facilitates operation when advancing the catheter

#### **Endotracheal**

| Item No.  | Description                | Pkg.  | Thumb valve color |
|-----------|----------------------------|-------|-------------------|
| DYNCSDS12 | Double swivel elbow, 12 Fr | 10/cs | White             |
| DYNCSDS14 | Double swivel elbow, 14 Fr | 10/cs | Green             |
| DYNCSTP12 | T-Piece adapter, 12 Fr     | 10/cs | White             |
| DYNCSTP14 | T-Piece adapter, 14 Fr     | 10/cs | Green             |

#### Tracheostomy

| Item No.   | Description                | Pkg.  | Thumb valve color |
|------------|----------------------------|-------|-------------------|
| DYNCSDS12T | Double swivel elbow, 12 Fr | 10/cs | White             |
| DYNCSDS14T | Double swivel elbow, 14 Fr | 10/cs | Green             |
| DYNCSTP12T | T-Piece adapter, 12 Fr     | 10/cs | White             |
| DYNCSTP14T | T-Piece adapter, 14 Fr     | 10/cs | Green             |

# DYNCSDS12T

Also available with T-Piece adapter

#### **Pediatric**

- Colored depth markings for clear viewing in low-light settings
- Available with an elbow or Y-style modular patient adapter
- Three adapter sizes to ensure compatibility
- Labeled for facilities to determine change out policy
- Compatible with all modes of ventilation

#### Elbow adapter

| Item No. | Description         | Pkg.  | Thumb valve color |
|----------|---------------------|-------|-------------------|
| DYNCPE5  | Elbow adapter, 5 Fr | 10/cs | Black             |
| DYNCPE6  | Elbow adapter, 6 Fr | 10/cs | Light green       |
| DYNCPE7  | Elbow adapter, 7 Fr | 10/cs | Purple            |
| DYNCPE8  | Elbow adapter, 8 Fr | 10/cs | Light blue        |

#### **Y-Connector**

| Item No. | Description       | Pkg.  | Thumb valve color |
|----------|-------------------|-------|-------------------|
| DYNCPY5  | Y-Connector, 5 Fr | 10/cs | Black             |
| DYNCPY6  | Y-Connector, 6 Fr | 10/cs | Light green       |
| DYNCPY7  | Y-Connector, 7 Fr | 10/cs | Purple            |
| DYNCPY8  | Y-Connector, 8 Fr | 10/cs | Light blue        |



Elbow adapter & Y-Connector products come with three adapter sizes in each kit

## **Open suction catheters**

Our D419 series catheters offers a full range of straight packaging.

- Available in a complete range of continuous adult and pediatric sizes
- All Medline suction catheters feature easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while minimizing trauma potential

#### Single sterile-straight pack

| Item No.  | French size | Pkg.   |
|-----------|-------------|--------|
| DYND41906 | 6 Fr        | 100/cs |
| DYND41908 | 8 Fr        | 100/cs |
| DYND41900 | 10 Fr       | 100/cs |
| DYND41901 | 12 Fr       | 100/cs |
| DYND41902 | 14 Fr       | 100/cs |
| DYND41903 | 16 Fr       | 100/cs |
| DYND41904 | 18 Fr       | 100/cs |

#### Coiled

| Item No.   | French size | Pkg.  |
|------------|-------------|-------|
| DYND41906C | 6 Fr        | 50/cs |
| DYND41908C | 8 Fr        | 50/cs |
| DYND41900C | 10 Fr       | 50/cs |
| DYND41901C | 12 Fr       | 50/cs |
| DYND41902C | 14 Fr       | 50/cs |
| DYND41903C | 16 Fr       | 50/cs |
| DYND41904C | 18 Fr       | 50/cs |

#### Poly-cath

| Item No.  | French size | Pkg.   |
|-----------|-------------|--------|
| DYND40990 | 10 Fr       | 100/cs |
| DYND40992 | 14 Fr       | 100/cs |

#### Continuous valve

| Item No.  | French size | Pkg.   |
|-----------|-------------|--------|
| DYND41960 | 10 Fr       | 100/cs |
| DYND41962 | 14 Fr       | 100/cs |
| DYND41964 | 18 Fr       | 100/cs |

## **Sleeved suction catheters**

Our D407 series open-suction catheters are contained in a sleeve to help avoid cross-contamination.

| Item No.   | French Size | Pkg.  |
|------------|-------------|-------|
| DYND40708F | 8 Fr        | 50/cs |
| DYND40700F | 10 Fr       | 50/cs |
| DYND40701F | 12 Fr       | 50/cs |
| DYND40702F | 14 Fr       | 50/cs |
|            |             |       |



## **Kits and trays**

Our D409 series kits and trays provide all the necessities for airway clearance.

- Irrigation solutions can be poured into the tray or the pop-up cup
- All kits and trays contain Medline Aloetouch® 3G gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to help minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while helping minimizing trauma potential

### Suction catheter mini tray with gloves

| Item No.  | French Size | Pkg.   |
|-----------|-------------|--------|
| DYND40986 | 6 Fr        | 100/cs |
| DYND40988 | 8 Fr        | 100/cs |
| DYND40980 | 10 Fr       | 100/cs |
| DYND40981 | 12 Fr       | 100/cs |
| DYND40982 | 14 Fr       | 100/cs |
| DYND40983 | 16 Fr       | 100/cs |
| DYND40984 | 18 Fr       | 100/cs |
|           |             |        |

#### Suction catheter pack with gloves and pop-up cup

| Item No.  | French Size | Pkg.  |
|-----------|-------------|-------|
| DYND40976 | 6 Fr        | 50/cs |
| DYND40978 | 8 Fr        | 50/cs |
| DYND40970 | 10 Fr       | 50/cs |
| DYND40971 | 12 Fr       | 50/cs |
| DYND40972 | 14 Fr       | 50/cs |
| DYND40973 | 16 Fr       | 50/cs |
| DYND40974 | 18 Fr       | 50/cs |

#### Suction catheter mini tray with gloves, poly-cath

| Item No.   | French Size | Pkg.   |
|------------|-------------|--------|
| DYND48980A | 10 Fr       | 100/cs |
| DYND48982  | 14 Fr       | 100/cs |

Medline manufactures over 125,000 sterile procedure kits every day.














DYND40973









## **Suction catheter trays with sterile solutions**

Our D414 series suction catheters come with everything needed for airway clearance, including a choice of irrigation fluids.

- All kits and trays contain Medline Aloetouch® 3G gloves
- · Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while helping to minimize trauma potential

#### Water

| Item No.  | Description   | Pkg.  |
|-----------|---|-------|
| DYND41476 | Suction catheter tray with sterile water, 100 mL; 6 Fr  | 36/cs |
| DYND41478 | Suction catheter tray with sterile water, 100 mL; 8 Fr  | 36/cs |
| DYND41470 | Suction catheter tray with sterile water, 100 mL, 10 Fr | 36/cs |
| DYND41471 | Suction catheter tray with sterile water, 100 mL; 12 Fr | 36/cs |
| DYND41472 | Suction catheter tray with sterile water, 100 mL; 14 Fr | 36/cs |
| DYND41474 | Suction catheter tray with sterile water, 100 mL; 18 Fr | 36/cs |

#### Saline

| Item No.    | Description  | Pkg.  |
|-------------|--|-------|
| DYND41478SA | Suction catheter tray with sterile saline, 100 mL; 8 Fr  | 36/cs |
| DYND41470SA | Suction catheter tray with sterile saline, 100 mL; 10 Fr | 36/cs |
| DYND41442   | Suction catheter tray with sterile saline, 100 mL; 14 Fr | 36/cs |
| DYND40442   | Suction catheter tray with saline and drape; 14 Fr       | 36/cs |

#### **Sterile solutions**

These convenient, ready-to-use containers reduce the contamination risks and potential waste associated with bulk containers of sterile water and saline.

- · Ideal for device irrigation, suction, dressing moistening and wound debridement
- Cup size provides the proper amount for suctioning

#### Water

| Item No.    | Description             | Pkg.  |
|-------------|-------------------------|-------|
| RDI30295    | Sterile water, 100 mL   | 48/cs |
| DYND40570   | Sterile water, 110 mL   | 48/cs |
| PCS1550     | Sterile water, 250 mL   | 24/cs |
| DYND500MLW  | Sterile water, 500 mL   | 18/cs |
| DYND1000MLW | Sterile water, 1,000 mL | 6/cs  |

#### Saline

| Item No.    | Description                    | Pkg.  |
|-------------|--------------------------------|-------|
| RDI30296    | Sterile saline, 0.9%, 100 mL   | 48/cs |
| DYND40540   | Sterile saline, 0.9%, 110 mL   | 48/cs |
| PCS1650     | Sterile saline, 0.9%, 250 mL   | 24/cs |
| DYND500MLS  | Sterile saline, 0.9%, 500 mL   | 18/cs |
| DYND1000MLS | Sterile saline, 0.9%, 1,000 mL | 6/cs  |













DYND500MLS RD130296

# **Tracheostomy care and cleaning trays**

We offer a wide variety of tracheostomy care and cleaning trays in standard and customized configurations to meet diverse clinical needs.

- Available in economy, basic and full-compartment kits
- May be customized

Did you know?
We can build a custom kit to your exact specifications.
Medline is a leading manufacturer of custom procedure kits.

|                        |       | 2850 AN | DA SS ON THE | 0200400 | 000 000 000 000 000 000 000 000 000 00 | 880 A | SZSOO ON THE O | 000000000000000000000000000000000000000 | OF O | 04/10 OF 105/1 |       |       | K ON CO |       |
|------------------------|-------|---------|--------------|---------|--|-------|----------------|---|--|----------------|-------|-------|---------|-------|
| Tray compartments      | 3     | 1       | 1            | 1       | 3                                      | 3     | 3              | 3                                       | 3  | 3              | 1     | 3     | 3       | 3     |
| AloeTouch gloves       | 2     |         | 2            | 2       | 2                                      | 2     | 2              | 2                                       | 2  | 2              |       | 2     |         | 2     |
| Drape                  | 1     |         | 1            | 1       | 1                                      | 1     | 1              |   | 1  | 1              |       | 1     | 1       | 1     |
| Pipe cleaners          |       | 2       | 2            |         |  |       |                |   |  | 3              | 3     | 3     | 3       | 3     |
| Non-woven gauze 4 x 4  | 4     | 2       | 2            | 2       | 2                                      | 2     | 2              | 2                                       | 1  | 3              | 4     | 1     | 3       | 3     |
| Cotton tip applicator  | 2     | 2       | 5            | 5       | 2                                      | 5     | 2              | 2                                       | 4  | 2              | 2     | 2     | 3       | 3     |
| Pliable cleaning brush | 1     | 1       | 1            |         | 1                                      | 1     | 1              |   | 1  | 1              | 1     | 1     | 1       | 1     |
| Twill tape 34"         | 1     | 1       | 1            | 1       | 1                                      | 1     | 1              | 1                                       | 1  |                | 1     | 1     | 1       | 1     |
| 6-ply trach dressing   | 1     |         | 1            | 1       | 1                                      | 1     | 1              | 1                                       | 1  | 1              | 1     | 1     | 1       | 1     |
| Hydrogen peroxide      |       |         |              |         | 2                                      | 2     |                | 1                                       |  | 1              |       |       |         |       |
| Saline 100 mL          |       |         |              |         | 1                                      |       |                |   | 1  | 1              |       |       |         |       |
| Water 100 mL           |       |         |              |         |  |       |                | 1                                       |  |                |       |       |         |       |
| Suction catheter 14 fr |       |         |              |         |  |       | 1              |   |  | 1              |       |       |         |       |
| Forceps                |       |         |              |         |  |       |                |   |  |                |       |       | 1       | 1     |
| Pkg                    | 20/cs | 30/cs   | 20/cs        | 20/cs   | 20/cs                                  | 20/cs | 20/cs          | 20/cs                                   | 20/cs                                    | 20/cs          | 30/cs | 20/cs | 20/cs   | 20/cs |





DYND40582

## **Mucus/specimen traps**

Our mucus and specimen traps have transport caps for more secure specimen collection.


- Sterile plastic vials, graduated in 1cc increments
- Includes threaded cap, tubing not made from natural rubber latex and patient ID label

#### Sterile DeLee mucus trap

| Item No.  | Description                                       | Pkg.  |
|-----------|---|-------|
| DYND44108 | Trap with catheter and cap, capacity 20 cc; 8 Fr  | 50/cs |
| DYND44110 | Trap with catheter and cap, capacity 20 cc; 10 Fr | 50/cs |

#### Sterile mucus specimen trap

| Item No.  | Description                              | Pkg.  |
|-----------|--|-------|
| DYND44140 | Trap with transport cap, capacity; 40 cc | 50/cs |
| DYND44180 | Trap with transport cap, capacity; 80 cc | 50/cs |



# **Rest assured**

We are stocking up on top respiratory products to meet your seasonal needs



# **Medication delivery**

## **MDI spacers and adapters**

| Item No.             | Description                                       |  |  |  |
|----------------------|---|--|--|--|
| MDI spacer           |   |  |  |  |
| HCS1143035           | OptiChamber Diamond MDI Spacer, 50 pack           |  |  |  |
| 1098213ML            | OptiChamber Diamond MDI Spacer, 10 pack           |  |  |  |
| HCSA100              | MDI Spacer  |  |  |  |
| MDI spacer with mask |   |  |  |  |
| 1098214ML            | OptiChamber Diamond MDI Spacer with small mask    |  |  |  |
| 1098215ML            | OptiChamber Diamond MDI Spacer with medium mask   |  |  |  |
| 1098216ML            | OptiChamber Diamond MDI Spacer with large mask    |  |  |  |
| MDI adapter          |   |  |  |  |
| SMDIA1000ML          | MDI inline adapter, 22 mm ID -> 22mm OD/15 mm ID  |  |  |  |
| SMDIA1100ML          | MDI inline adapter, 22 mm OD/15 mm ID -> 22mm ID  |  |  |  |
| SMDIA1200ML          | MDI inline adapter, 22 mm OD -> 22 mm OD/15 mm ID |  |  |  |
| SMDIA1300ML          | MDI inline adapter, 15 mm OD -> 15 mm ID          |  |  |  |





SMDIA1000ML

## **Medication nebulizers**

| Item No.         | Description  |
|------------------|--|
| Nebulizer with r | mouthpiece   |
| HCS4483          | Nebulizer kit with tee mouthpiece, 6" reservoir, 7' tubing and standard connector          |
| HCS4483R         | Reusable nebulizer kit with tee mouthpiece, 6" reservoir, 7' tubing and standard connector |
| Nebulizer with r | mask   |
| HCS4485          | Nebulizer kit with mask, 7' tubing and standard connector, adult                           |
| HCS4486          | Nebulizer kit with mask, 7' tubing and standard connector, pediatric                       |
| Nebulizer BAN    |  |
| HCS4433BAN       | Aeromist BAN nebulizer with tee mouthpiece, 7' tubing                                      |
| HCS4433FBAN      | Aeromist BAN nebulizer with tee mouthpiece, filter, 7'tubing                               |
| Nebulizer comp   | ressor   |
| HCS70004         | Nebulizer compressor with disposable nebulizer medicine cup, 7' tubing and mouthpiece      |
| Valved tee adap  | oter   |
| HCSA1009         | Valved tee adapter, 22mm OD/ID   |



HCS4483

 $A vailable\ in\ Universal\ tubing\ options.\ For\ additional\ options\ please\ visit\ medline.com.$ 

# O<sub>2</sub> disposables

# Oxygen cannulas

| Item No.                      | Description  |
|-------------------------------|--|
| Standard O <sub>2</sub> cann  | nulas  |
| HCS4511B                      | Standard nasal cannula, 7' tubing and standard connector, adult                    |
| HCS4512B                      | Standard nasal cannula, 14' tubing and standard connector, adult                   |
| HCS4513B                      | Standard nasal cannula, 25' tubing and standard connector, adult                   |
| Soft-touch O <sub>2</sub> car | nnulas   |
| HCS4514                       | Soft-touch nasal cannula with 7' tubing and standard connector, adult              |
| HCS4516B                      | Soft-touch nasal cannula with 14' tubing and standard connector, adult             |
| HCS4515B                      | Soft-touch nasal cannula with 25' tubing and standard connector, adult             |
| SuperSoft O <sub>2</sub> can  | nnulas   |
| HCSS4514S                     | SuperSoft oxygen cannula, 7' tubing and standard connector, adult                  |
| HCSS4516S                     | SuperSoft oxygen cannula, 14' tubing and standard connector, adult                 |
| HCSS4515S                     | SuperSoft oxygen cannula, 25' tubing and standard connector, adult                 |
| High-flow O <sub>2</sub> can  | nulas  |
| HCS4514HFG                    | SuperSoft high-flow oxygen cannula, standard connector and 7' green tubing, adult  |
| HCS4516HFG                    | SuperSoft high-flow oxygen cannula, standard connector and 14' green tubing, adult |
| HCS4515HFG                    | SuperSoft high-flow oxygen cannula, standard connector and 25' green tubing, adult |



HCSS4514S



HCS4515HFG

# Oxygen and aerosol masks

| Item No.             | Description   |
|----------------------|---|
| O <sub>2</sub> masks |   |
| HCS120E              | Venti oxygen mask with nose clip, color diluters, 7' tubing and standard connector, adult                                 |
| HCS121E              | Venti oxygen mask with single dial, 7' tubing and standard connector, adult   |
| HCS4600B             | Medium concentration oxygen mask, 7' tubing and standard connector, adult   |
| HCS4640B             | Partial non-rebreathing oxygen mask with reservoir bag, safety vent, check valve, 7' tubing and standard connector, adult |
| HCS4680B             | Oxygen mask 3-in-1 with 7' tubing and standard connector, adult   |
| OxyMask™             |   |
| OM11258ML            | OxyMask with 7' universal tubing, adult   |
| OP11258ML            | OxyMask with 7' universal tubing, adult plus  |
| Aerosol mask         |   |
| HCS4630B             | Aerosol mask with straight connector, adult   |



## **Prefilled bubble humidifier**

| Item No.         | Description                         |
|------------------|-------------------------------------|
| Prefilled bubble | humidifier                          |
| HCS00550         | Prefilled bubble humidifier, 550 ml |

 $A vailable in Universal \ tubing \ options. For \ additional \ options \ please \ visit \ medline.com.$ 

# O<sub>2</sub> disposables, continued

# Oxygen & aerosol tubing

| Item No.              | Description   |
|-----------------------|---|
| Aerosol tubing        |   |
| HSK600                | 100' L aerosol tubing with 22 mm cuffs                      |
| HSK626                | 6' L aerosol tubing with 22 mm cuffs                        |
| O <sub>2</sub> tubing |   |
| HCS4507               | Clear kink-resistant oxygen tubing, 7', standard connector  |
| HCS4507G              | Green kink-resistant oxygen tubing, 7', standard connector  |
| HCS4524               | Clear kink-resistant oxygen tubing, 14', standard connector |
| HCS4524G              | Green kink-resistant oxygen tubing, 14', standard connector |



# **Pulse oximeter**

| Item No.                 | Description                                       |  |
|--------------------------|---|--|
| Fingertip pulse oximeter |   |  |
| HCSM70B                  | Basic fingertip pulse oximeter                    |  |
| HCSM70S                  | Soft-touch fingertip pulse oximeter, OLED display |  |
| Handheld pulse oximeter  |   |  |
| HCSMD300                 | Handheld pulse oximeter with 1 adult finger probe |  |



HCSM70S

# **Ventilation products**

## **Closed suction**

| Item No.     | Description   |
|--------------|---|
| ETT length   |   |
| DYNCPDS14    | ClearPro double swivel closed suction catheter with endotracheal tube connection, 14 fr |
| DYNCPTP14    | ClearPro t-piece closed suction catheter with endotracheal tube connection, 14 fr       |
| DYNCSDS14    | Double-swivel closed suction catheter with endotracheal tube, 14 fr                     |
| DYNCSTP14    | Closed suction catheter with t-piece and endotracheal tube, 14 fr                       |
| Trach length |   |
| DYNCPDS14T   | ClearPro double swivel closed suction catheter with tracheostomy connection, 14 fr      |
| DYNCPTP14T   | ClearPro t-piece closed suction catheter with tracheostomy connection, 14 fr            |
| DYNCSDS14T   | Double-swivel closed suction catheter for tracheostomy care, 14 fr                      |
| DYNCSTP14T   | Closed suction catheter with t-piece for tracheostomy care, 14 fr                       |



DYNCPDS14

Available in Universal tubing options. For additional options please visit medline.com.

# **Ventilation products, continued**

# Heat and moisture exchanger & filtration

| Item No.                | Description  |
|-------------------------|--|
| Bacterial/viral filters |  |
| DYNJAABV1               | Bacterial/viral filter   |
| DYNJAAHEPA              | HEPA filter  |
| Heat and moisture ex    | cchangers  |
| DYNJAABHME              | Bypass HME   |
| DYNJAAHME12             | HME, compact   |
| Filtered heat and mo    | isture exchangers  |
| DYNJAAHME1              | Filtered HME, with CO <sub>2</sub> port  |
| DYNJAAHME11             | Filtered HME, compact, with CO <sub>2</sub> port                                       |
| DYNJAAHME13             | Compact, angled, with CO <sub>2</sub> port   |
| DYNJAAHME5              | Filtered HME   |
| DYNJAAHME9              | Filtered HME, compact, with CO <sub>2</sub> port                                       |
| DYNJAAHMEHEPA           | HEPA filtered HME, with CO <sub>2</sub> port   |
| Cuffill®                |  |
| HCSCUFF0041             | AG Cuffill Device, 10 cc   |
| HCSCUFF0041H            | AG Cuffill Device, 10 cc   |
| Non-invasive masks      |  |
| HCSFFNVAM               | Flex™ full face NIV mask, adult medium, non-vented, AAV, and headgear                  |
| HCSFFNVM                | Flex™ full face NIV mask, adult medium, non-vented, and headgear                       |
| HCSFFVAM                | Flex™ full face NIV mask, adult medium, vented, AAV, and headgear                      |
| HCSNVAM                 | Nasal NIV mask, adult medium, vented, AAV, and headgear                                |
| Resuscitation bags      |  |
| CPRM1116                | Adult manual resuscitation bag with medium adult mask                                  |
| CPRM1116F               | Adult manual resuscitation bag with filter and medium adult mask with 7' oxygen tubing |
| Resuscitation access    | ories  |
| CPRM2400                | Colorimetric CO <sub>2</sub> detector  |
| CPRMFILTER              | Expiratory filter for Medline resuscitation bags                                       |
| CPRMP30                 | PEEP valve, adjustable, 19mm ID with 22mm OD, 30mm ID connections                      |
| SubG endotracheal t     | ube  |
| DYNJ18870               | Subg ET tube, size # 7.0   |
| DYNJ18875               | Subg ET tube, size # 7.5   |
| DYNJ18880               | Subg ET tube, size # 8.0   |
| DYNJ18885               | Subg ET tube, size # 8.5   |
| Suction canister with   | float lid  |
| HCS7850                 | Canister with float lid: 850 cc  |
| HCS7851                 | Canister with float lid & tubing: 850 cc   |
| Non-heated circuits     |  |
| DYNJ300                 | Adult dual limb circuit, wye without ports, uncuffed, 6'                               |



DYNJAABV1





CPRM1116F



CPRM2400

RTAC3AAA



# **Minor capital**


| Item No.                   | Description                                     |
|----------------------------|---|
| HCS7000                    | Portable suction aspirator                      |
| Suction regulators         | Various configurations available—call for quote |
| 02 & air flowmeters        | Various configurations available—call for quote |
| Oxygen hoses & accessories | Various configurations available—call for quote |
| O2 regulators              | Various configurations available—call for quote |

| A  | D                                 | M   | R                              |
|--|-----------------------------------|---|--------------------------------|
| Adapters                                       | Detectors, CO <sub>2</sub>        | Manual Resuscitators, 61-67                             | Regulators                     |
| Aerosol, 54                                    | Colorimetric, 68, 83              | Masks   | Oxygen, 13, 83                 |
| MDI, 30, 80                                    | Disposable Canister Kits, 16, 83  | Aerosol, 31, 54, 81                                     | Suction, 14–15, 83             |
| Valved Tee, 30, 80                             | DreamPort CPAP Mask Solution, 37  | CO₂, Capnoxygen, 51                                     | Resuscitation, 61–70           |
| Aeromist BAN Nebulizer, 25                     | bream ore erra mask solution, sr  | CPAP, 38–41   | Accessories, 69                |
| Aerosol  | E                                 | Solution, 37  | Bags, 62–67, 83                |
|  |                                   | Face Tent, 54   |                                |
| Masks, 31, 46, 54, 81                          | Endotracheal                      | Gel, Full-Face, 39<br>LifeTouch, 28                     | Masks, 61                      |
| Tubing, 52, 82                                 | Holder, 6                         | Mouth-to-Mask Resuscitator, 69                          | Resuscitators                  |
| AG Cuffill®, 5, 83                             | Suction Catheters, 72, 74, 82     | Nasal, 38–40  | Manual, 61–67                  |
| Air Flowmeters, 11                             | Tubes, Subglottic, 5, 83          | Non-Invasive  | Mouth-To-Mask, 69              |
| Airway Management, 5–7                         | Exchangers, Heat and              | Ventilation, 35–36, 83                                  |                                |
| Aspirator, Suction, 16                         | Moisture, 20-21, 83               | Oxygen, 45-46, 53, 81                                   | S                              |
| Asthma Management, 25–33                       | Bypass, 19, 83                    | Resuscitation, 69                                       | SideStream                     |
|  | Filtered, 21, 83                  | Southmedic OxyMask, 45-46                               | Disposable, 26                 |
| В  | Trach, 20                         | Trach, 54   | Masks, 31, 54                  |
| Bacterial Filters, 18, 83                      | Hach, 20                          | MDI   | Nebulizer, 26                  |
| Bags, Resuscitation, 62–67, 83                 | г                                 | Adapters, 30  | Reusable, 26                   |
| Breathing Circuits,                            | F                                 | Spacers and Masks, 29, 80                               | SideStream Plus, 26            |
| Non-Heated, 22, 83                             | Filters                           | Medication Delivery, 25–33                              | Sleeved Suction Catheters, 75  |
|  | Bacterial/Viral, 18, 83           | MicroFilter CO <sub>2</sub> Cannulas, 50                |                                |
| Bronchial Hygiene, 8                           | Heat and Moisture                 | Mouthpiece, Disposable, 23<br>MRI-Conditional Transport | Southmedic OxyMask, 45–46, 8   |
| Bubble Humidifiers, 55                         | Exchangers, 20-21, 83             | Ventilator Circuit, 22                                  | Spacers, MDI, 28–29, 80        |
| Bypass Heat and Moisture                       | Pulmonary Function                | Mucus Traps, 79   | Specimen Traps, 79             |
| Exchanger, 19, 83                              | Testing (PFT), 23                 | Mucus Haps, 19  | Spirometers, Incentive, 9      |
|  | Filtration, 19–23                 | M   | Sterile Solutions, 77          |
| C  | Finger Pulse Oximeters, 57, 82    | N   | SubG Endotracheal Tubes, 6, 83 |
| Canister Kits, Suction, 16, 83                 | Flex Connector, 7                 | Nasal Pillow Systems, 41                                | Suction                        |
| Cannulas                                       |                                   | Nebulizers, 25–27, 80                                   | Aspirator, 16, 83              |
| CO <sub>2</sub> , 49–50                        | Flowmeters, Air & Oxygen 11, 83   | Compressors, 32–33, 80<br>Handheld, 25-27, 80           | Canister Kits, 16              |
| Oxygen, 47–48, 81                              | Flowmeters, Peak 29               | Non-Heated Breathing                                    | Catheters, 72–77               |
| , ,  | Foam Ear Cover Cannulas, 48       | Circuits, 28, 83  | Devices, 11–16                 |
| Capnoxygen Masks, 51                           |                                   | Non-Invasive  | Regulators, 14–15, 83          |
| Catheters, Suction                             | G                                 | Interfaces, 35–43                                       | 11050101013, 11 13, 03         |
| ClearPro, 72–73, 82                            | Gas                               | Ventilation Masks, 35–36, 83                            | -                              |
| Closed, 72–74                                  | Delivery, 11–16                   | Nose Clips, 23  | Т                              |
| Open, 75–77                                    | Hose Assemblies, 12               |   | Tent Masks, 54                 |
| Kits, 76-78                                    | Go-PAP CPAP Systems, 42           | 0   | Ties, Tracheostomy, 7          |
| Chin Straps, CPAP, 41                          | do 1711 ci711 Systems, 12         | O₂ Disposables, 45–55                                   | Tracheostomy, 72, 74           |
| Clips, Nose, 23                                | Н                                 | O₂-Max CPAP Systems, 42–43                              | Care & Cleaning Trays, 78      |
| CO <sub>2</sub>                                |                                   | Open Suction Catheters, 75–77                           | Trach HME, 20                  |
| Cannulas, 49-50                                | Handheld Oximeters, 58            | OptiChamber Diamond, 28, 80                             | Masks, 54                      |
| Detectors                                      | Heat and Moisture                 | Oximeters   | Ties, 7                        |
| Colorimetric, 68, 83                           | Exchangers, 20–21, 82             | Finger Pulse, 57, 82                                    | Traps                          |
| Disposables, 45–55                             | Bypass, 19, 83                    | Handheld, 58  | Mucus, 79                      |
| Masks, 46, 51                                  | Filtered, 21, 83                  | Pulse, 57–58  | Specimen, 79                   |
| Colorimetric CO <sub>2</sub> Detectors, 68, 83 | Trach, 20                         | Oxygen  | Vacuum Trap, 15                |
| Compressors, Nebulizer, 32–33, 80              | High Flow Oxygen Cannulas, 48, 81 | Accessories, 55   | Tube                           |
| Portable, 32                                   | Hose Assemblies, Gas, 12          | Cannulas, 47–48, 81<br>Cylinder Wrench, 13              | SubGlottic, 6, 83              |
|  | Humidification, 18–23             | Flowmeters, 11, 83                                      |                                |
| Corrugated Aerosol Tubing, 52                  | Humidifiers, Bubble, 55           | Masks, 45–46, 53, 81                                    | Tubing                         |
| COVID-19 Key Respiratory                       | Prefilled, 55, 81                 | Regulators, 13, 83                                      | Corrugated Aerosol, 52         |
| Codes, 80-83                                   | Hyperinflation Systems, 70        | Tubing, 52, 82  | CPAP, 41                       |
| CPAP   | Tryperiillation Systems, 70       |   | Oxygen, 52, 82                 |
| Accessories, 41                                |                                   | P   |                                |
| DreamPort Mask Solution, 37                    | 1                                 | Peak Flow Meter, 29                                     | V                              |
| Masks, 38-41                                   | Incentive Spirometers, 9          | PEEP Valves, 69, 83                                     | Vacuum Traps, 15               |
| Full-Face, 38-40                               | Interfaces, Non-Invasive, 35–43   | PEP therapy and   | Valved Tee Adapters, 30        |
| Nasal, 38-41                                   |                                   | inspiratory muscle , 9                                  | Ventilation Masks, 35–36, 83   |
| Go-PAP, O2-Max Nasal Pillow                    |                                   | Portable Nebulizer                                      | Ventilator Circuit, 22         |
| Systems, 42–43                                 |                                   | Compressor, 32  | Flex Connector, 7              |
| CPR Mask, 69                                   |                                   | Pressure Monitoring                                     | Viral Filters, 18, 83          |
| Cuff Pressure Controller                       |                                   | AG CUFFILL, 5, 83                                       | vii ai i iicci 3, 10, 03       |
| CUFFILL AG, 5, 83                              |                                   | Probes  | 10/                            |
| Cylinder Wrench, Oxygen, 13                    |                                   | Oximeter, 58  | W                              |
| Symmetry Tremen, Oxygen, 15                    |                                   | Pulmonary Function Testing (PFT)<br>Filters, 23         | Washer, Brass Yoke, 13         |
|  |                                   | Pulse Oximeters, 57–58                                  | Wrench, Oxygen Cylinder, 13    |
|  |                                   | Finger, 57, 82  |                                |
|  |                                   | Probes, 58  |                                |
|  |                                   |   |                                |

| Notes |  |  |
|-------|--|--|
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |

| Notes |  |  |
|-------|--|--|
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |

| Notes |  |
|-------|--|
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |





Medline Industries, Inc.

Three Lakes Drive, Northfield, IL 60093 Medline United States | 1-800-MEDLINE (633-5463) medline.com | info@medline.com



**Medline Canada** 

1-800-268-2848 | medline.ca | canada@medline.com Medline México

01-800-831-0898 | medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2021 Medline Industries, Inc. IntelliCuff is the trademark of Hamilton Medical. Respironics and Sidestream are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. LiteTouch Respironics and LiteTouch are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. OptiChamber Respironics and OptiChamber are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. Respironics and Personal Best are trademarks of Koninklijke Philips N.V. and its affiliates, all rights reserved. Bleep, Dreamports, DreamWay, and DreamPreps are trademarks of Bleep LLC. AlRgel, Ascend, VeraSeal 2, Mojo, MiniMe, MiniMe 2, and Innova are trademarks of Sleepnet. AlRgel, Ascend, VeraSeal 2, Mojo, MiniMe, MiniMe 2, and Innova are trademarks of Sleepnet. Citrus II CPAP Mask Wipes is a trademark of Beaumont Products Inc. Go-PAP and O2-Max are trademarks of Pulmodyne. FloCap is a trademark of Sunmed Medical. The Oxymask trademark is the property of Southmedic Inc., and is used under license by Medline Industries, Inc. MKT19W1759222 / LIT085 / 15