

Gel and Foam Positioners

Support pressure relief, comfort and stability



Take the pressure off.

To safely position patients during surgery and rehab therapy, you need a reliable support system.

Designed in a variety of shapes, sizes and levels of resiliency, our gel and foam positioners help take the pressure off. They elevate and stabilize patients to reduce the risk of pressure injuries and other complications—so you can handle patients with care.

Under pressure.

\$26.8 billion

U.S. health system spends on HAPIs¹

1 in 5

patients affected by pressure injuries each year²

\$20,900 - \$151,700

Cost of individual patient care per pressure injury³

Head and Neck Positioners 4

Upper Extremity Positioners 8

Lower Extremity Positioners 10

Multipurpose Positioners 13

Gel Specialty Positioners 16

Bean Bag Positioners 17

Foam Positioning Wedge 18

Support Surfaces 20

Foam Surgical Positioning Kits 22

ECO Foam Positioners 23

Take a strong position against pressure injuries.

Recent guidelines from the National Pressure Ulcer Advisory Panel reinforce the importance of reducing pressure, friction and shear in preventing pressure injuries.⁴ Our positioners were designed to help meet these goals.



Gel positioners

- Impervious to body fluids
- Designed to last
- Reusable and easy to clean
- Perfect balance of quality and versatility

Foam positioners

- Accessible and economical
- Fulfill procedural requirements without compromising care



Head and Neck Positioners

Provide support and stability where patients need it most during surgery.

Gel Positioners

Donut Head Pad

- Helps reduce pressure on back of the head during supine procedures
- Supports head and neck while protecting ears during lateral procedures

Item No.	Description	Pkg.
PCD082000	Adult 8" x 3" x 1.75" (20.3 x 7.6 x 5.1 cm)	1/ea
PCD082001	Adult Bariatric, 8" x 3" x 3" (20.3 x 7.6 x 7.6 cm)	1/ea
PCD082004	Pediatric, 5.25" x 2.25" x 1.25" (14 x 5.7 x 3.2 cm)	1/ea
PCD082005	Neonatal, 3.5" x 1.5" x 0.75" (8.9 x 3.8 x 1.5 cm)	1/ea



PCD082004

Horseshoe Donut Head Pad

- Helps protect the patient's face during surgery
- Enables easy access to anesthesia tubes during prone-position surgeries

Item No.	Description	Pkg.
PCD082006	Adult, 8" x 3" x 3" (20.3 x 7.6 x 7.6 cm)	1/ea
PCD082007	Adult, 8" x 3" x 1.75" (20.3 x 7.6 x 5.1 cm)	1/ea
PCD082009	Pediatric, 5.25" x 2.25" x 1.25" (13.3 x 5.7 x 3.2 cm)	1/ea
PCD082008	Neonatal, 3.5" x 1.5" x 0.75" (8.9 x 3.8 x 1.9 cm)	1/ea



PCD082007

Lateral Head Pad with Center Dish

- Centering dish lets pad conform to head for support and protection
- Enhances patient comfort and helps protect against added pressure

Item No.	Description	Pkg.
PCD082023	Adult, 1" x 9" x 10" (2.5 x 22.9 x 25.4 cm)	1/ea
PCD082025	Pediatric, 1" x 6" x 7" (2.5 x 15.2 x 17.8 cm)	1/ea
PCD082027	Neonatal, 5" x 4" x 1" (12.7 x 10.2 x 2.5 cm)	1/ea



PCD082025

Prone Position Head Rest

- Comfortably rests patient's head face-down with respiratory tubes in place
- Center opening provides protects patient's nose and mouth

Item No.	Description	Pkg.
PCD082010	11" x 9" x 6" (27.9 x 22.9 x 15.2 cm)	1/ea



PCD082010

Richard Type Head Rest

- Ideal for side or supine-positioned patients
- T-shaped cutout helps protect ears when used laterally, or back of head in supine position

Item No.	Description	Pkg.
PCD082021	9" x 8" x 3" (22.9 x 20.3 x 7.6 cm)	1/ea



PCD082021

Supine Head Pad with Center Dish

- Best for positioning patient during supine procedures, such as eye surgery
- Center dish enhances comfort and provides protection from added pressure

Item No.	Description	Pkg.
PCD082022	Adult, 1.5" x 9" x 10" (3.8 x 22.9 x 25.4 cm)	1/ea
PCD082024	Pediatric, 1.5" x 6" x 7" (3.8 x 15.2 x 17.8 cm)	1/ea



PCD082022

Supine Head Rest

- Contoured design cushions the patient's head while providing stability
- Allows for easy surgical access during face and eye procedures

Item No.	Description	Pkg.
PCD082020	Contoured 7" x 7" x 3.25" (17.8 x 17.8 x 8.3 cm)	1/ea



PCD082022

Head and Neck Positioners

Foam Positioners

Foam Head Positioner

- Resilient, single-use foam helps reduce pressure while maximizing support
- Protects facial features and provides comfort during surgery
- Accommodates endotracheal tubes on both sides for easier positioning
- Comfort Touch Head Rest has a top layer of soft foam for increased comfort

Item No.	Description	Pkg.
NON081143	Slotted, 8.5" x 8" x 4" (21.6 x 20.3 x 10.2 cm)	24/cs
NON081143E	Slotted, Economy, 8.5" x 8" x 4" (21.6 x 20.3 x 10.2 cm)	24/cs
NON081144	Adult, Compressed, 9" x 8" x 4.5" (22.9 x 20.3 x 11.4 cm)	24/cs
NON081144E	Adult, Economy, Compressed, 9" x 8" x 4.5" (22.9 x 20.3 x 11.4 cm)	24/cs
NON081144NC	Adult, Non-Compressed, 9" x 8" x 4.5" (22.9 x 20.3 x 11.4 cm)	24/cs
NON081145	Child, 7" x 6.5" x 3.5" (17.8 x 16.5 x 8.9 cm)	24/cs
NON081146	Frame Head Rest, 10.75" x 9.75" x 6" (27.3 x 24.8 x 15.2 cm)	28/cs
NON081147	Comfort Touch Frame Head Rest, Two-Layer Design, 10.75" x 9.75" x 6" (27.3 x 24.8 x 15.2 cm)	28/cs
NON081147P	Pediatric, Comfort Touch Frame Head Rest, Two-Layer Design, 10.75" x 9.75" x 6" (27.3 x 24.8 x 15.2 cm)	28/cs
NON081148NC	Adult, Extra-Firm Density, Non-Compressed, 9" x 8" x 4.5" (22.9 x 20.3 x 11.4 cm)	24/cs



NON081143



NON081144



NON081147

Reusable Neck Support

- Helps provide anatomically correct support to the neck and head
- Relieves tension in shoulders; helps reduce side-to-side movement
- Cradle molds itself to the head for personalized support
- Nylex cover is easy to clean; can be used long-term

Item No.	Description	Pkg.
MSC04160	Neck Positioner Cradle, 15" x 7" x 3.5" (38.1 x 17.8 x 8.9 cm)	1/ea
MSC04165	Head and Neck Positioner, 9" x 8" x 4.5" (22.9 x 20.3 x 11.4 cm)	1/ea
MSC04175	Neck Stabilizer, 19" x 14" x 6" (48.3 x 35.6 x 15.2 cm)	1/ea

Memory Foam Pillow

- Made of visco elastic memory foam designed for extra comfort
- 100% cotton slip cover provides patients with support and comfort

Item No.	Description	Pkg.
NON084241	Foam Pillow, Temp	6/cs



NON084241

Ring Positioner

- Multipurpose cushions in a variety of diameters, including multi-ring style
- Positions head, ear or iliac crest; also supports elbows, heels and knees

Item No.	Description	Pkg.
NON081139	Ring Head Positioner, 7" x 1.5" (17.8 x 3.8 cm)	24/cs
NON081140	Multi-Ring	36/cs
NON081140E	Multi-Ring, Economy	24/cs
NON081141	9" OD, 4" ID, x 2" (22.9/10.2 x 5 cm)	32/cs
NON081141E	Economy, 9" OD, 4" ID, x 2" (22.9/10.2 x 5 cm)	1/ea
NON081142	7" OD, 2" ID, x 2" (17.8/5 x 5 cm)	72/cs
NON0811421	7" OD, 3.25" ID x 1.5" (17.8/8.3 x 3.81 cm)	72/cs
NON081149	Contoured, 9" x 3.25" (22.9 x 8.3 cm)	32/cs



NON081140



NON081141



NON081149

Upper Extremity Positioners

Select sizes and styles to support arms, elbows, hands and shoulders.

Gel Positioners

Armboard Pad

- Effective in relieving pressure on the hand, forearm, elbow and bicep
- Provide additional stability to ensure proper positioning
- Flat and contoured pads fit most surgical table armboards

Item No.	Description	Pkg.
PCD085000	Flat, 15" x 4" x 0.25" (38.1 x 10.1 x 0.6 cm)	1/ea
PCD085001	Flat, 25" x 7" x 1" (63.5 x 17.8 x 2.5 cm)	1/ea
PCD085002	Flat, 30" x 10" x 0.5" (76.2 x 25.4 x 1.3 cm)	1/ea
PCD085005	Flat, 25" x 5" x 0.375" (63.2 x 12.7 x 0.95 cm)	1/ea
PCD085008	Flat, 26" x 6" x 0.37" (66 x 15.2 x 0.9 cm)	1/ea
PCD085009	Flat, 15" x 6" x 0.37" (38.1 x 15.2 x 0.9 cm)	1/ea
PCD085010	Contoured, 20" x 5.75" x 2.25" (51 x 14.6 x 5.7 cm)	1/ea



Ulnar and Brachial Nerve Protectors

- Ulnar nerve cushions cradle the elbow's bony prominences
- Designed to cradle ulnar nerve area; protects against skin breakdown

Item No.	Description	Pkg.
PCD085011	Ulnar/Brachial Nerve Pad, 17" x 11" x 0.375" (43.18 x 27.94 x 0.95cm)	1/ea
PCD085012	Diamond Ulnar Nerve Pad, 15.75" x 6" x 0.375" (40 x 15.2 x 0.95 cm)	1/ea
PCD085013	Ulnar/Brachial Nerve Pad, 21" x 24" x 0.25" (53.3 x 61 x 0.6 cm)	1/ea
PCD085015	Ulnar/Brachial Nerve Pad, 24" x 19" x 0.25" (61 x 48.3 x 0.6 cm)	1/ea
PCD085016	Ulnar/Brachial Nerve Pad, 27" x 19" x 0.25" (61 x 68.6 x 0.6 cm)	1/ea
PCD085018	Ulnar/Brachial Nerve Pad, 15" x 21" x 0.357" (38.1 x 53.3 x 0.95 cm)	1/ea



Foam Positioners

Armboard Pads and Arm Cradles

- High-quality convoluted foam offers support and comfort
- Designed to fit most surgical table armboards; can be cut for a better fit
- Cradles flex to accommodate a variety of angles for best positioning practice

Item No.	Description	Pkg.
NON081343	Armboard Pads, Pair, 2" x 8" x 20" (5 x 20.3 x 50.8 cm)	24/cs
NON081343E	Armboard Pads, Pair, Economy 2" x 8" x 20" (5 x 20.3 x 50.8 cm)	24/cs
NON081343L	Armboard Pads, Pair, Long 2" x 8" x 26" (5 x 20.3 x 66 cm)	12/cs
NON081346	Armboard Pads, Convoluted, 16" x 6" x 2" (40.6 x 15.2 x 5 cm)	24/cs
NON081344	Arm Cradles, Pair, 3" x 5" x 24" (7.6 x 12.7 x 61 cm)	24/cs
NON081344E	Arm Cradles, Pair, Economy 3" x 5" x 24" (7.6 x 12.7 x 61 cm)	24/cs
NON081344NC	Arm Cradles, Pair, Non-Comp 3" x 5" x 24" (7.6 x 12.7 x 61 cm)	12/cs
NON081451	Arm Cradle, 15.25" x 12.5" x 8" (38.7 x 31.7 x 20.3 cm)	12/cs
NON081657	Arm Cradles, Arilla Cuff	5/cs



NON081343



NON081344




NON081451

Ulnar Nerve Cushions

- Designed to cradle sensitive ulnar nerve area; protects against skin breakdown
- Convoluted foam offers pressure relief
- Hook-and-loop strap provides secure placement
- Great for use during rehabilitation

Item No.	Description	Pkg.
NON081340	Ulnar Nerve Cushions, Pair	72/cs
NON081340E	Ulnar Nerve Cushions, Pair, Economy	72/cs
NON081341	Ulnar Nerve Cushions, Firm, Pair	36/pr



NON081340



NON081340E

Lower Extremity Positioners

Select sizes and styles to support heels, ankles, legs, knees and feet.

Gel Positioners

Heel and Ankle Positioner

- Bony prominences are commonly affected by pressure injuries
- Helps protect heels and ankles during exams and surgery

Item No.	Description	Pkg.
PCD084002	Heel Protector, Clam Shell	1/ea
PCD084004	Heel Support, Cup Style	1/ea
PCD085026	Heel Elevator	1/ea



PCD084002



PCD084004

Heel and Elbow Dish

Item No.	Description	Pkg.
PCD085025	Heel and Elbow Dish, 3.25" x 1.5" x 0.75" (8.3 x 3.8 x 2 cm)	1/ea



PCD085025

Leg and Knee Positioners

Item No.	Description	Pkg.
PCD084009	Candy Cane Stirrup Pad, 2" x 28" x 0.25" (5 x 66.04 x 0.63 cm)	1/ea
PCD084010	Candy Cane Stirrup Pad, Pair, 2" x 26" x 0.25" (5 x 66.04 x 0.63 cm)	1/pr
PCD084011	Ankle Stirrup Strap, Small, Pair, 12.5" x 1.87" x 0.375" (32 x 4.8 x 1 cm)	1/pr
PCD084012	Ankle Stirrup Strap, Large, Pair, 17.37" x 1.87" x 0.37" (44 x 4.8 x 1 cm)	1/pr
PCD087000	Knee Crutch Pad, Small, 11" x 9" x 0.37" (27.9 x 22.9 x 0.9 cm)	1/ea
PCD087004	Knee Crutch Pad, Large, 20" x 20" x 0.5" (50.8 x 50.8 x 1.3 cm)	1/ea
PCD087005	Knee Crutch Pad, Extra-Large, 27" x 11" x 0.5" (68.6 x 28 x 1.3 cm)	1/ea
PCD088006	Knee Crutch, Large, Bariatric, 20" x 20" x 1" (50.8 x 50.8 x 2.5 cm)	1/ea



PCD084011

Disposable Stirrup Straps

- Budget-friendly disposable foam stirrup straps
- Aid in positioning while maintaining patient comfort

Item No.	Description	Pkg.
NON52710	Disposable Stirrup Straps without Slip Ring, 19" x 3.5" (48.2 x 8.9 cm)	240/cs
NON52712	Disposable Stirrup Straps with Slip Ring, 19" x 3.5" (48.2 x 8.9 cm)	240/cs
NON081458	Disposable Stirrup Pad, 9" x 23.5" x 2" (22.9 x 59.7 x 5 cm)	10/cs



NON52712

Lithotomy Boot Pads

Item No.	Description	Pkg.
PCD084006	Small (Allen Boots)	1/ea
PCD084007	Large (Allen Boots)	1/ea



PCD084006

Gel Positioners Care and Cleaning:

- Use standard operating room cleaning agents. Bleach and long soaks not recommended. Avoid sharp instruments or dragging pads and positioning devices across sharp edges.
- Gel pads may be warmed in warming cabinets, heating blankets or by immersing in water to make the pads approach normal human body temperature. Do not heat over 140°F.
- Do not attempt to sterilize or autoclave. Store flat at room temperature. Do not move patients while on the pads.

Lower Extremity Positioners

Foam Positioners

Convuluted Foam Heel Protectors

- Anatomically-designed pads protect against increased pressure
- Made with convuluted foam for air circulation and comfort
- Hook-and-loop strap ensures a secure fit around the ankle

Item No.	Description	Pkg.
NON081440	Heel Protectors, Pair	36/pr



NON081440

Lower Extremity Surgical Positioners

- Provides specialized support for legs, knees, heels, ankles and feet
- Offers stability, protection and comfort for a variety of procedures

Item No.	Description	Pkg.
NON081442	Foot Drop Stop, Strap, 7" x 7" x 6.25" (17.8 x 17.8 x 15.9 cm)	24/cs
NON081443	Leg/Arm Positioner, 17" x 11" x 8" (43.1 x 28 x 20.3 cm)	6/cs
NON081445	Elevating Foot Cradle Positioner, 9" x 9" x 14" (22.9 x 22.9 x 35.6 cm)	4/cs
NON081446	Arthroscopic Knee Holder	10/cs
NON081447	Well Leg Holder	18/cs
NON081457	Wedge, Leg, with Strap, 29" x 11" x 11" (73.7 x 28 x 28 cm)	1/ea
NON081459	Foam Leg Elevator, 12" x 33" x 12" (30.5 x 83.8 x 30.5 cm)	1/ea
MSCLEGWELL	Nylex Leg Elevator, 12" x 33" x 12" (30.5 x 83.8 x 30.5 cm)	1/ea



NON081446



NON081447

Total Knee Elevator

- Provides specialized support for legs, knees, heels, ankles and feet
- Offers stability, protection and comfort for a variety of procedures

Item No.	Description	Pkg.
NON081455	Total Knee Elevator, 27" x 15.5" x 12.5" (68.5 x 39.3 x 31.75 cm)	1/ea



NON081445

Multipurpose Positioners

Get the support you need for a variety of surgical procedures.

Gel Positioners

- Used to support the torso during prone position surgeries
- Can also be used under knees and ankles, or along the abdomen

Chest Rolls

Item No.	Description	Pkg.
PCD08363	Flat-Bottom 6" x 3" x 3" (15.2 x 7.6 x 7.6 cm)	1/ea
PCD08383	Flat-Bottom 8" x 3" x 3" (20.3 x 7.6 x 7.6 cm)	1/ea
PCD083004	Flat-Bottom 12" x 4" x 4" (30.4 x 10.2 x 10.2 cm)	1/ea
PCD083006	Flat-Bottom 12" x 6" x 6" (30.4 x 15.2 x 15.2 cm)	1/ea
PCD083008	Flat-Bottom 16" x 4" x 4" (40.6 x 10.2 x 10.2 cm)	1/ea
PCD083009	Flat-Bottom 16" x 6" x 6" (40.6 x 15.2 x 15.2 cm)	1/ea
PCD083010	Flat-Bottom 20" x 4" x 4" (50.8 x 10.2 x 10.2 cm)	1/ea
PCD083011	Flat-Bottom 20" x 6" x 6" (50.8 x 15.2 x 15.2 cm)	1/ea
PCD083103	Flat-Bottom 10" x 3" x 3" (25.4 x 7.6 x 7.6 cm)	1/ea
PCD083133	Flat-Bottom 12" x 3" x 3" (30.4 x 7.6 x 7.6 cm)	1/ea
PCD083134	Dome, 20" x 6" x 3" (50.8 x 15.2 x 7.6 cm)	1/ea
PCD083135	Dome, 16" x 6" x 3" (40.6 x 15.2 x 7.6 cm)	1/ea
PCD083136	Dome, 14" x 6" x 3" (35.6 x 15.2 x 7.6 cm)	1/ea
PCD083137	Dome, 12" x 6" x 3" (30.5 x 15.2 x 7.6 cm)	1/ea
PCD083312	Contoured 12" x 6" x 5.5" (30.4 x 16.5 x 12.7 cm)	1/ea
PCD083316	Contoured 16" x 6" x 5.5" (40.6 x 15.9 x 12.7 cm)	1/ea
PCD083320	Contoured 20" x 6" x 5.5" (50.8 x 15.9 x 12.7 cm)	1/ea



PCD083006



PCD083316

Pediatric Chest Rolls

Item No.	Description	Pkg.
PCD08342	Flat-Bottom 4" x 2" x 2" (10.2 x 5.1 x 5.1 cm)	1/ea
PCD08362	Flat-Bottom 6" x 2" x 2" (15.2 x 5.1 x 5.1 cm)	1/ea
PCD08382	Flat-Bottom 8" x 2" x 2" (20.3 x 5.1 x 5.1 cm)	1/ea
PCD083102	Flat-Bottom 10" x 2" x 2" (25.4 x 5.1 x 5.1 cm)	1/ea
PCD083122	Flat-Bottom 12" x 2" x 2" (30.5 x 5.1 x 5.1 cm)	1/ea



PCD08342

Multipurpose Positioners

Gel Positioners

Gel Sandbag Positioners

- Used to aid in laminectomies and ankle support operations
- Good for neck extensions during surgeries and ophthalmic procedures
- Available in a wide selection of sizes

Item No.	Description	Pkg.
PCD086000	0.5 lb., 1" x 6" x 1.5" (2.5 x 15.2 x 3.8 cm)	1/ea
PCD086001	1 lb., 1.5" x 6" x 2" (3.8 x 15.2 x 5 cm)	1/ea
PCD086002	1 lb., 1" x 8" x 1.5" (2.5 x 20.3 x 3.8 cm)	1/ea
PCD086003	1 lb., 10" x 4.5" x 1.5" (25.4 x 11.4 x 3.8 cm)	1/ea
PCD086004	1.5 lb., 12" x 4.5" x 1.5" (30.5 x 11.4 x 3.8 cm)	1/ea
PCD083015	2.5 lb., 10" x 4.5" x 1.5" (25.4 x 11.4 x 3.8 cm)	1/ea
PCD086008	3 lb., 10" x 4.5" x 2" (25.4 x 11.4 x 5 cm)	1/ea
PCD086010	4 lb., 12" x 4.5" x 2" (30.4 x 11.4 x 5 cm)	1/ea
PCD083017	6.5 lb., 14" x 5" x 2" (35.6 x 12.7 x 5.1 cm)	1/ea
PCD086011	7 lb., 14" x 5" x 2.5" (35.5 x 12.7 x 6.35 cm)	1/ea
PCD086005	8 lb., 2" x 18" x 6" (5.1 x 45.7 x 15.2 cm)	1/ea
PCD086012	10 lb., 16" x 12" x 1.5" (40.6 x 30.4 x 3.8 cm)	1/ea



PCD086008

Toboggan Positioner

- Makes specific positions easier while maintaining comfort and safety
- Can be used with toboggan gel pads to protect arms and legs

Item No.	Description	Pkg.
PCD087020	Gel Pad for Toboggan, 17" x 8" x 0.375" (43.2 x 20.3 x 1 cm)	1/ea
PCD087021	Toboggan Shield with Gel Pad	1/ea
PCD087022	Toboggan Shield, 15.75" x 9.5" x 0.375" (40 x 24.1 x 1 cm)	1/ea
PCD087023	Arm/Leg Gel Pad Only, 20" x 8" x 0.375" (50.8 x 20.3 x 1 cm)	1/ea
PCD087024	Gel Pad & Shield Set, Curved	1/ea
PCD087025	Arm/Leg Shield Only, Curved, 19.75" x 5.5" x 0.125" (50.2 x 14 x 0.3 cm)	1/ea



PCD087021

Foam Positioners

Roll Positioners and Bolsters

- Supports the head, neck, arm, knee, leg, foot and other parts of the body
- Also available with Nylex covers that can be easily wiped clean

Item No.	Description	Pkg.
MSC04155	Smooth Roll, 15" x 5" Diameter, Covered (38.1 x 12.7 cm)	1/ea
MSC04158	Smooth Roll, 36" x 7" Diameter, Covered (91.4 x 17.8 cm)	1/ea
MSC04150	Convuluted Roll, 15" x 5" Diameter, Covered (38.1 x 12.7 cm)	1/ea
NON081243	Smooth Roll, 16" x 5" Diameter (40.6 x 12.7 cm)	12/cs
NON081246	Bolster, Pair, 7" x 18" (17.8 x 45.7 cm)	6/pr
NON081264	Rectangular Foam Block, Uncovered, 16" x 6" x 3" (40.6 x 15.2 x 7.6 cm)	12/cs



MSC04150



NON081246

Disposable Foam Armboard Straps

- Strong yet soft foam straps with hook-and-loop closures
- A cost-effective choice

Item No.	Description	Pkg.
NON081540	Armboard Strap, 1.5" x 32" (3.8 x 81.3 cm)	48/cs
NON55299	Armboard Strap, 3" x 30" (7.6 x 76.2 cm)	24/cs



NON081540

Multipurpose Positioning Straps

- Ideal for a wide range of uses
- Soft, yet extremely strong; ensures proper positioning for security and comfort
- Adjusts quickly and easily to fit most patients
- Features durable plastic rings for quick, easy adjustment
- Straps provide secure attachment to the OR table or gurney

Disposable

Item No.	Description	Pkg.
NON081541	Multipurpose Positioning Strap, 4" x 30" (10.1 x 76.2 cm)	12/cs
DYNDSTRAP1	Positioning Strap, 3.25" x 43.5" (8.3 x 110.5 cm)	50/cs
DYNDSTRAP1P	Positioning Strap, 1PC, 5" x 72" (12.7 x 182.9 cm)	20/cs
NON548001	Positioning Strap, Metal Rings, 3" x 30" (7.6 x 76.2 cm)	24/cs



NON081541



DYNDSTRAP1

Reusable

Item No.	Description	Pkg.
NON081544	Positioning Strap, Metal Rings 4" x 60" (10.2 x 152.4 cm)	6/cs
NON081545	Patient Positioning Strap, Reusable	4/cs



DYNDSTRAP1P

Gel Specialty Positioners

Position patients comfortably and safely.

Gel Trapezoid Positioner

- Used to separate the patient's legs
- Can also be used to elevate knees, feet and more

Item No.	Description	Pkg.
PCD083027	Gel Trapezoid Positioner, 10" x 7" x 3.5" (25.4 x 17.8 x 8.9 cm)	1/ea



PCD083027

Gel Lateral Positioner

- Two-piece positioner creates a stable, safe cradle of upper body and arm
- Base is constructed of a lightweight, high-density foam
- Easy to move and put in place

Item No.	Description	Pkg.
PCD083026	Gel Lateral Positioner	1/ea



PCD083026

Cover for Gel Positioner Chest Roll

- Can be wrapped around open rolls for superb cushioning during surgeries
- Hook-and-loop closure holds securely yet easy to adjust

Item No.	Description	Pkg.
PCD083032	Cover for Gel Positioner Axillary Roll	1/ea

Bean Bag Positioners

Customize support to improve surgical positioning.

Traditional Bean Bag Positioners

- Constructed with durable polyurethane film
- Can be hooked up to a suction line to remove air and customize the shape
- Clean with standard OR cleaners; do not sterilize

Item No.	Description	Pkg.
PCDBB1020	20" x 10" (51 x 25 cm)	1/ea
PCDBB2020	20" x 20" (51 x 51 cm)	1/ea
PCDBB3040	Shoulder Cutout, 40" x 30" (102 x 76 cm)	1/ea
PCDBB4036	Shoulder Cutout, 40" x 36" (102 x 91 cm)	1/ea
PCDBB4646	Bariatric with Shoulder Cutout, 46" x 46" (117 x 117 cm)	1/ea
PCDBB3076	Full-Length, 76" x 30" (193 x 76 cm)	1/ea
PCDBBVALVE	Valve Replacement Kit	1/ea



PCDBB2020



PCDBB4036



PCDBB1020G



PCDBB3040G



PCDBB4646G

Gel Bean Bag Positioners

- Reusable gel bean bags made with a viscoelastic polymer gel layer and durable black polyurethane film
- Easy to clean
- Easily attaches to wall suction, aspirator pumps or any OR vacuum source

Item No.	Description	Pkg.
PCDBB1020G	20" x 10" (51 x 25 cm)	1/ea
PCDBB2020G	20" x 20" (51 x 51 cm)	1/ea
PCDBB3040G	Shoulder Cutout, 40" x 30" (102 x 76 cm)	1/ea
PCDBB4036G	Shoulder Cutout, 40" x 36" (102 x 91 cm)	1/ea
PCDBB4646G	Bariatric with Shoulder Cutout, 46" x 46" (117 x 117 cm)	1/ea

Foam Positioning Wedges

Versatile stability for surgeries, imaging and rehabilitation.

Disposable Positioning Wedges

- Versatile, disposable and available in multiple sizes
- Works for a variety of patient and procedural needs
- Radiolucent wedges won't interfere with crucial testing

Item No.	Description	Pkg.
NON081244	7.5" x 7.5" x 17" (19 x 19 x 43.2 cm)	8/cs
NON081230	8" x 11" x 18" (20.3 x 28 x 71.1 cm)	6/cs
NON081245	8" x 11" x 22" (20.3 x 28 x 55.9 cm)	8/cs
NON081254	10" x 7" x 2.5" (25.4 x 17.8 x 6.4 cm)	16/cs
NON081233	16" x 11" x 5" (40.6 x 27.9 x 12.7 cm)	8/cs
NON081242	18" x 18" x 26" (45.7 x 45.7 x 66 cm)	4/cs
NON081255	Spinal, 21.25" x 7" x 7" (54 x 17.8 x 17.8 cm)	16/cs
NON081265	23" x 11" x 11" (58.4 x 28 x 28 cm)	12/cs



NON081265

Abduction Pillows

- Feature straps to hold legs in proper position for optimum healing
- Variety of sizes to accommodate different patient frames
- Nylex cover wipes down easily; treated with antimicrobial*

Item No.	Description	Pkg.
MSC019014	Reusable T-Shaped with Straps, 5" x 6.5" x 8.5" (12.7 x 16.5 x 21.6 cm)	1/ea
MSC04143	Reusable Small, 18" x 12" x 6" (45.7 x 30.5 x 15.2 cm)	1/ea
MSC04145	Reusable Medium, 22" x 15" x 6" (55.9 x 38 x 15.2 cm)	1/ea
MSC04147	Reusable Large, 25" x 18" x 6" (63.5 x 45.7 x 15.2 cm)	1/ea
NON081448	Disposable, Small, 6" x 12" x 18" (15.2 x 30.5 x 45.7 cm)	6/cs
NON081449	Disposable, Medium, 6" x 15" x 22" (30.5 x 38 x 55.9 cm)	6/cs
NON081450	Disposable, Large, 6" x 18" x 25" (30.5 x 45.7 x 63.5 cm)	6/cs



MSC019014



NON081449

* This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

Nylex®-Covered Foam Positioners Care and Cleaning:

- Do not use abrasive and/or phenolic-based cleaning agents.
- **Routine care:** Wipe surface with neutral detergent and lukewarm water; rinse with water and dry thoroughly.
- **Washing:** Wash at no higher than 120° F. Tumble dry at low heat. Do not use bleach.
- **Disinfection:** Wipe surface with a quaternary solution at recommended dilution ratios. Neutral solutions less than 10pH are preferred.



MSC04105



MSCWGRNDTOP



MSC2451508

Reusable Nylex®-Covered Positioning Wedges

- Sewn-on, easy-to-clean Nylex cover is treated with antimicrobial*
- Flame-retardant and fluid-resistant
- Radiolucent, so they won't interfere with crucial testing
- Available in a variety of wedges and angles

Item	Description	Pkg.
MSC018312	10" x 7" x 2.5" (25.4 x 17.8 x 6.4 cm)	2/cs
MSC018317	10" x 10" x 10" (25.4 x 25.4 x 25.4 cm)	1/cs
MSC012812	12" x 8" x 12" (30.5 x 20.3 x 30.5 cm)	4/cs
MSC0161565	15.5" x 10.5" x 5.5" (39.4 x 26.7 x 14 cm)	4/cs
MSC04110	16" x 11" x 5" (66 x 27.9 x 39.4 cm)	1/ea
MSC0161210	16" x 12" x 10" (40.6 x 30.5 x 25.4 cm)	2/cs
MSCENDOWDG	18" x 10" x 10" (45.7 x 25.4 x 25.4 cm)	1/ea
MSC04100	19" x 7.5" x 7.5" (48.3 x 19 x 19 cm)	1/ea
CBC11X15X25	19" x 15" x 11" (48.3 x 38.1 x 27.9 cm)	2/cs
MSC04115	21" x 21" x 7" (53.3 x 53.3 x 17.8 cm)	1/ea
MSC04116	21" x 21" x 10" (53.3 x 53.3 x 24.4 cm)	2/cs
MSC04120	21" x 21" x 12.5" (53.3 x 53.3 x 31.8 cm)	1/ea
MSC04105	22" x 11" x 8" (55.9 x 28 x 20.3 cm)	1/ea
MSC04105FNS	22" x 11" x 8" (55.9 x 28 x 20.3 cm)	2/cs
MSCM1004280A	22" x 11" x 11" (55.9 x 28 x 28 cm)	1/ea
CBC8X22X22	22" x 22" x 8" (55.9 x 55.9 x 20.3 cm)	2/cs

Item	Description	Pkg.
MSC2451508	Platform, 24.5" x 15" x 8" (62.2 x 38.1 x 20.3 cm)	1/cs
CBC6X18X24	24" x 18" x 6" (61 x 45.7 x 15.2 cm)	2/cs
MSC04121	24" x 24" x 7" (61 x 61 x 17.8 cm)	2/cs
CBC04010	24" x 24" x 10" (61 x 61 x 25.4 cm)	2/cs
MSC04122	24" x 24" x 10" (61 x 61 x 25.4 cm)	2/cs
MSC04123	24" x 24" x 12" (61 x 61 x 30.5 cm)	2/cs
MSC04114	26" x 18" x 18" (66 x 45.7 x 45.7 cm)	2/cs
MSC0271715	27" x 17" x 15" (68.6 x 43.2 x 38.1 cm)	1/ea
MSC04103	28" x 17" x 7" (71.1 x 43.2 x 17.8 cm)	2/cs
MSC301009	Platform, 30" x 10" x 9" (76.2 x 25.4 x 22.9 cm)	1/cs
MSC04106	30" x 11" x 7" (76.2 x 27.9 x 17.8 cm)	2/cs
MSC018316	30" x 9" x 9" (76.2 x 22.9 x 22.9 cm)	1/cs
CBC102MP113	34" x 12" x 10" (86.3 x 30.5 x 25.4 cm)	2/cs
MSC04101	36" x 12" x 7" (91.4 x 30.5 x 17.8 cm)	2/cs
MSC04113	39" x 12" x 6" (99 x 30.5 x 15.2 cm)	2/cs
MSCBARIWEDGE	48" x 12.5" x 12" (121.9 x 31.8 x 30.5 cm)	1/ea
MSC04127	48" x 20" x 20" (121.9 x 50.8 x 50.8 cm)	1/ea
MSC04141	Popliteal Support with Nylex Cover	2/ea
MSCWDGRNDTOP	Round Top Wedge, 12" x 21" x 7.375" (30.5 x 53.3 x 18.8 cm)	1/cs
MSCWDGARCH	Foam Wedge, Nylex Cover with Gel Layer, 23" x 17" x 17" (58.4 x 43.2 x 43.2 cm)	1/ea

* This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

Sterile Foam Bumps

- Bumps help prevent intraoperative injury from pressure, obstruction and stretching of an extremity
- Rectangular design allows for the extremity to be situated at various heights
- Can be used during shoulder abduction, lateral shoulder arthroscopy, incision closure of total hip or total knee, vascular surgery, extremity trauma or carpal tunnel

Item No.	Description	Pkg.
DYNJBUMP01	Flat, 10" x 6" x 4" (25.4 x 15.2 x 10.2)	28/cs
DYNJBUMP02	Flat, 10" x 4" x 4.5" (25.4 x 10.2 x 11.4 cm)	24/cs
DYNJBUMP03	Round Top, 10" x 6" x 6" (25.4 x 15.2 x 15.2 cm)	12/cs
DYNJBUMP04	Round Top, 10" x 4" x 4.5" (25.4 x 10.2 x 11.4 cm)	24/cs
DYNJBUMP06	Round Top, 12" x 7" x 7" (30.5 x 17.8 x 17.8 cm)	6/cs
DYNJBUMP07	Bolster, 18" x 7" x 3.5" (45.7 x 17.8 x 8.9 cm)	6/cs



DYNJBUMP01



DYNJBUMP04

Support Surfaces

Create comfortable surfaces for patient procedures.

Gel Pads

- Optimum pads for procedures in the exam room and OR, or for rehabilitation
- Ideal for most patient handling applications
- Pads ensure patient comfort and insulate from cold, hard table surfaces

Item No.	Description	Pkg.
PCD081000	Amsco Table Pad Set	3/cs
PCD081001	Head Table Pad, 20" x 9" x 0.5" (51 x 22.8 x 1.3 cm)	1/ea
PCD081002	Torso Table Pad, 39.5" x 20" x 0.5" (100.3 x 51 x 1.3 cm)	1/ea
PCD081003	Foot Table Pad with Cutout, 20" x 20" x 0.5" (51 x 51 x 1.3 cm)	1/ea
PCD081004	Gel Table Pad, 51.5" x 20.5" x 0.5" (130.8 x 52.1 x 1.3 cm)	1/ea
PCD081005	Foot Table Pad, 20" x 20" x 0.5" (51 x 51 x 1.3 cm)	1/ea
PCD081006	Full Table Pad Wide, 72" x 24" x 0.5" (61 x 183 x 1.3 cm)	1/ea
PCD081007	Full Table Pad, 72" x 20" x 0.5" (183 x 51 x 1.3 cm)	1/ea
PCD081008	3/4 Table Pad with Cutout, 46" x 20" x 0.5" (117 x 51 x 1.3 cm)	1/ea
PCD081009	3/4 Table Pad, 46" x 20" x 0.5" (117 x 51 x 1.3 cm)	1/ea
PCD081010	OR Gel Table Pad, 20" x 48" x 1" (50.8 x 122 x 2.5 cm)	1/ea
PCD081011	OR Gel Table Pad, 12" x 37" x 1" (30.5 x 94 x 2.5 cm)	1/ea
PCD083000	Hip Table Pad, 20" x 20" x 0.5" (51 x 51 x 1.3 cm)	1/ea
PCD084008	Foot Fracture Table Pad, 15" x 4" x 0.375" (38.1 x 10.16 x 0.95 cm)	1/ea
PCD088015	Cysto Table Pad, 24" x 30" x 0.5" (61 x 76.2 x 1.3 cm)	1/ea
PCD088017	Stirrup Pad, 12" x 12" x 0.25" (30.5 x 30.5 x 0.6 cm)	1/ea
PCD089002	Pediatric Pad, Small, 11" x 9" x 0.37" (28 x 23 x 0.95 cm)	1/ea
PCD089001	Pediatric Pad, Medium, 17" x 11" x 0.25" (43 x 28 x 0.64 cm)	1/ea



PCD081003

Specialty Gel Pads

Item No.	Description	Pkg.
PCD088007	Wilson Frame Pad, Set, 29" x 10" x 0.5" (74 x 25 x 1.3 cm)	1/ea
PCD088008	Andrews Frame Knee Pad, Set, 22" x 16" x 1" (56 x 40.6 x 2.5 cm)	1/ea
PCD088010	Mayo Stand Pad Cover, 21" x 18.5" x 0.62" (53 x 47 x 1.5 cm)	1/ea
PCD088011	Kit, Patch, Gel Repair	4/cs
PCD088013	Wilson Frame Pad, 29" x 10" x 0.5" (74 x 25 x 1.3 cm)	1/ea



PCD088009

Gel Peg Board Sets

- Gel overlays provide comfort and support when positioned on peg boards
- Can be easily cleaned and reused

Item No.	Description	Pkg.
PCD088020	Peg Board Pad, Full Set	1/ea
PCD088021	Peg Board Pad, Short Post	1/ea
PCD088022	Peg Board Pad, Tall Post	1/ea
PCD088023	Peg Board Pad, Overlay	1/ea



PCD088021



PCD088023

Foam Convuluted OR Pads

- High-quality convoluted foam offers support and comfort
- Cradles flex to accommodate a variety of positioning angles

Item No.	Description	Pkg.
NON082440	OR Utility Pad, 2" x 12" x 20" (5 x 30 x 50 cm)	12/cs
NON082442	OR Table Pad, 2" x 20" x 72" (5 x 50 x 189 cm)	12/cs
NON082444	OR Table Pad, 2" x 16" x 20" (5 x 40.6 x 50.8 cm)	12/cs
NON082446	OR Table Pad, 5" x 20" x 70" (12.7 x 50.8 x 177.8 cm)	6/cs
NON08071Z	OR Flat Pad, 12" x 15" x 0.5" (30 x 38.1 x 1.3 cm)	10/cs



NON082440

Foam Rubber Pads

- Foam rubber pads provide maximum safety and comfort
- Comfort Touch foam pad features soft foam that contours to patient for increased comfort

Item No.	Description	Pkg.
NON61054	Foam Rubber Pad, 36" x 21" x 1" (91.4 x 53.3 x 2.5 cm)	1/ea
NON58351	Foam Rubber Pad, 36" x 21" x 0.25" (91.4 x 53.3 x 0.6 cm)	1/ea
NON2029	Comfort Touch OR Pad, 29" x 20" x 1.5" (73.6 x 50.8 x 3.8 cm)	1/ea



NON2029

Foam Surgical Positioning Kits

Save resources, storage space and set-up time.

Customizable Kits

Item No.	Description	Pkg.
NON081641	Prone Positioning Kit	8/cs
NON081646	Knee Positioning Kit	10/cs
NON081647	Schlein Kit	6/cs
NON081647BP	Schlein, Back Pad	6/cs
NON081652	Jackson Frame Kit	12/cs
NON081652XC	Jackson Frame Kit, No Cushion	12/cs
NON081653	Schlein Frame Kit, Chin Strap	10/cs
NON081654	Wilson Frame Kit	6/cs
NON081660	Hana Table Frame Kit	6/cs
NON081660XC	Hana Table Frame Kit, No Headrest/Straps	6/cs
NON1094387	Robotic Kit	6/cs

Jackson Frame Kit

Kit includes:

- 4 Bouffant Covers
- Bouffant Pillow
- Comfort Touch Frame Head Rest
- Arm Cradles, 1 Pair



Wilson Frame Positioning Kit

Kit includes:

- 1 slotted head rest with memory foam layer
- 1 pair of arm cradles
- 2 pad sleeves
- 1 transverse bar cover
- 1 base cover drape



Schlein Kit

Kit includes:

- Schlein Pad
- Large Pad
- Small Pad
- Strap Harness
- Chin Strap
- Belt Strap, Small
- Belt Strap, Large

ECO Foam Positioners

Eco-friendly, body-friendly positioning.



Eco-friendly Foam Positioners

- Made with 10% renewable resources to reduce non-renewable energy consumption
- Consist of soy-based polyols vs. petroleum-based
- Physically equivalent to standard polyether and polyester foams

Item No.	Description	Pkg.
NON081340ECO	Ulnar Nerve Protector	1/pr
NON081141ECO	Ring Head Positioner, 9" (22.9 cm) diameter	1/ea
NON081343ECO	Foam Arm Board Pads, 2" x 8" x 20" (5.1 x 20.3 x 50.8 cm)	1/pr
NON081344ECO	Arm Cradle Positioner, 3" x 5" x 24" (7.6 x 12.7 x 61 cm)	1/pr
NON081245ECO	Disposable Foam Wedge, 8" x 11" x 22" (20.3 x 27.9 x 55.9 cm)	1/ea



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093
Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996 | medline.ca | canada@medline.com
Medline México
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us



References: **1.** Padula W, Delarmente B. The national cost of hospital-acquired pressure injuries in the United States. *Int Wound J.* 2019;1-7 doi:10.1111/iwj.13071. [Epub ahead of print] **2.** Padula WV, Mishra MK, Makic MB, Sullivan PW. Improving the Quality of Pressure Ulcer Care with Prevention: a cost-effectiveness analysis. *Med Care.* 2011; 49(4):385-92. doi: 10.1097/MLR.0b013e31820292b3. **3.** Agency for Healthcare Research and Quality, Rockville, MD. Available at: <https://www.ahrq.gov/professionals/systems/hospital/pressureulcertoolkit/putool1.html>. Accessed March 14, 2019. **4.** Wounds International. Pressure Ulcer Prevention: pressure, shear, friction and microclimate in context. Available at: <https://www.woundsinternational.com/resources/details/international-review-pressure-ulcer-prevention-pressure-shear-friction-and-microclimate-context>. Accessed January 7, 2020.