














Phonak Earmould & Ear Plug Order Form Account No: 200 _____

| | |
|-------------------|--|
| Hearing Clinic: | Client Name: |
| Delivery Address: | Date Ordered: Date Required: |
| Contact Name: | VIP: 24hr in house \$50 per order (+GST) |
| Contact Number: | SEND IMPRESSIONS WITH ORDER FORM TO: |
| Email: | Phonak: 12 Inglewood Place, Baulkham Hills, NSW 2153 Ph: 1800 809 321 F: 02 9634 8373 E: orders@phonak.com.au |

Step 1: Enter the Audiogram (for AOV calculation)

| | |
|---|--|
| Left <input type="radio"/> 250 <input type="radio"/> 500 <input type="radio"/> 1k <input type="radio"/> 2k <input type="radio"/> 3k <input type="radio"/> 4k dB Hz | Right <input type="radio"/> 250 <input type="radio"/> 500 <input type="radio"/> 1k <input type="radio"/> 2k <input type="radio"/> 3k <input type="radio"/> 4k dB Hz |
|---|--|

Step 2: Choose your Earmould or Ear Plug option

| | | | | | | | | |
|--|--|---|--|--|---|---|--|--------------------------|
| SlimTip – Hollow (HA) L R  with Std xR with Power xR with Micro/SlimTube OK to change to cShell if SlimTip not possible due to ear geometry | SlimTip – Solid (HA) L R  with Std xR with Micro/SlimTube OK to change to cShell if SlimTip not possible due to ear geometry | SlimTip – Solid / Soft L R  with Std xR with Power xR with Micro/SlimTube OK to change to cShell if SlimTip not possible due to ear geometry | Canal Mould L R  | Semi Skeleton (HA only) L R  | Canal Lock (HA only) L R  | Skeleton (HA only) L R  | | |
| Half Shell L R  | Carved Shell L R  | Uncarved Shell L R  | Ultra Carved Shell L R  | cShell (RIC only) L R  | CROS B* *Compatible with CROS II-13 Only L R  | Custom Ear Plugs – Noise Protection Silicon Plug: L R Handling Options: No Handle Fixed carry cord | Custom Ear Plugs – Swim Protection Silicon Plug: L R Handling Options: No Handle Fixed carry cord | Special Requests: |

Step 3: Choose your Wax Management System (cShell only)

NOTE: Cerustop will default unless alternative is selected.

| | | | | |
|----------|---|---|---------------|---|
| Cerustop | Smart Guard* <small>*not for SP & UP power levels.</small> | Waxbuster** <small>**size permitting</small> | Extended Tube | OK to change to Cerustop if SmartGuard or WaxBuster does not fit. |
|----------|---|---|---------------|---|

Step 4: Choose your Earmould material

| | | | | | | | | | | | | | | |
|--|--|--|------------------------------------|------------------------------|---------------------------------|---------------------------------|--------------------------------|---------------------------------|------------------------------|-------------------------------|-------------------------------|--------------------------------|--------------------------------|--------------------------------|
| Hard Acrylic (HA) L R Colour: Clear Rosa Tint Beige | Silicone L R S70 (Firm; available in clear & rosa tint only) S40 (Soft; available in all colours. Refer to colour chart) | Colour Chart: Please choose Earmould Colour <table border="0" style="width: 100%; font-size: small;"> <tr> <td><input type="checkbox"/> Rosa Tint</td> <td><input type="checkbox"/> Red</td> <td><input type="checkbox"/> Orange</td> <td><input type="checkbox"/> Purple</td> </tr> <tr> <td><input type="checkbox"/> Clear</td> <td><input type="checkbox"/> Yellow</td> <td><input type="checkbox"/> Tan</td> <td><input type="checkbox"/> Pink</td> </tr> <tr> <td><input type="checkbox"/> Blue</td> <td><input type="checkbox"/> Green</td> <td><input type="checkbox"/> Black</td> <td><input type="checkbox"/> White</td> </tr> </table> | <input type="checkbox"/> Rosa Tint | <input type="checkbox"/> Red | <input type="checkbox"/> Orange | <input type="checkbox"/> Purple | <input type="checkbox"/> Clear | <input type="checkbox"/> Yellow | <input type="checkbox"/> Tan | <input type="checkbox"/> Pink | <input type="checkbox"/> Blue | <input type="checkbox"/> Green | <input type="checkbox"/> Black | <input type="checkbox"/> White |
| <input type="checkbox"/> Rosa Tint | <input type="checkbox"/> Red | <input type="checkbox"/> Orange | <input type="checkbox"/> Purple | | | | | | | | | | | |
| <input type="checkbox"/> Clear | <input type="checkbox"/> Yellow | <input type="checkbox"/> Tan | <input type="checkbox"/> Pink | | | | | | | | | | | |
| <input type="checkbox"/> Blue | <input type="checkbox"/> Green | <input type="checkbox"/> Black | <input type="checkbox"/> White | | | | | | | | | | | |

Step 5: Choose your Venting (mm) (AOV will default unless alternative is chosen)

| | | | | | | | | |
|-----|----------|-----|-----|-----|-----|-----|-----|------|
| AOV | Occluded | 0.8 | 1.0 | 1.5 | 2.0 | 2.5 | 3.0 | Open |
| L | L | L | L | L | L | L | L | L |
| R | R | R | R | R | R | R | R | R |

Step 6: Choose your Tubing for Earmoulds & SlimTip options for BTE Devices (not all options are compatible with earmoulds/shells)

| Std Earmould Tubing | microTube <small>micro CORE BTEs</small> | SlimTube <small>CORE BTE, Milo Plus, Baseo Q</small> | SlimTube HEs <small>Spice, Quest, & Venture BTE's</small> | Power SlimTubes <small>Naida S & Q Only</small> |
|---------------------|---|---|--|--|
| L R | OL OR | OL OR | OL OR | OL OR |
| Thick 13 | 1L 1R | 1L 1R | 1L 1R | 1L 1R |
| Thick Dry 13 | 2L 2R | 2L 2R | 2L 2R | 2L 2R |
| Double Wall 13 | 3L 3R | 3L 3R | 3L 3R | 3L 3R |
| 3mm Libby Dry | | | | |

OR Step 7: Choose your xReceiver for RIC Devices (not all options are compatible with earmoulds/shells. xR in cShells cannot be changed. cShell must be returned if new xR is needed)

| Std Receiver <small>For Slim Tip or cShell</small> | Power Receiver <small>For Slim Tip or cShell</small> | Super Power Receiver <small>For cShell only</small> | Super Power+ Receiver <small>For Naida S & Q and Audéo V</small> | Ultra Power Receiver <small>For Audéo B</small> |
|---|---|--|---|--|
| OL OR | OL OR | OL OR | OL OR | OL OR |
| 1L 1R | 1L 1R | 1L 1R | 1L 1R | 1L 1R |
| 2L 2R | 2L 2R | 2L 2R | 2L 2R | 2L 2R |
| 3L 3R | 3L 3R | 3L 3R | 3L 3R | 3L 3R |

Step 8: Choose Additional options (not all options are available for all models)

| | | | |
|----------|-------------|--------------|-------------------|
| L R | L R | L R | L R |
| Elbow | Removal Pin | Top Lock ON | Tapered Canal Tip |
| Tubelock | Belled Bore | Top Lock OFF | Concha Lock |

