

Greetings from the Rocky Mountains!

My name is Zubin Meshginpoosh and I am the new President & Chief Commercial Officer at Westone Laboratories. Although I have officially only been a member of the distinguished Westone team since November, my top priority from day one has always been to learn all the ins and outs of Westone's business as quickly as possible. While I may have over twenty years of experience working within the healthcare industry, I am still eager to learn from you as to how we can further improve our partnerships to better serve you and your patients. I, along with the rest of the remarkable team here at Westone, am nothing short of ecstatic to get to work with you as we bring in 2019!

It is hard to believe just how far the company has come since its humble beginnings in 1959 when Ron Morgan founded Westone in his family cabin in the mountains of Divide, Colorado. While 2018 was another great year for Westone, I look forward to helping improve all aspects of Westone's business to ensure that our 60th year of operation is our best yet!

Since November, I have had the opportunity to meet with an abundance of our customers, and it is easy to see why Westone has been so successful for the past 60 years. We have the most dedicated customers in the world and we will never be able to thank you enough for all of the years of support you have shown us. Several customers have been with us since 1959, and it is your commitment and dedication to Westone that has kept us thriving for so long.

I look forward to meeting you in the coming year at conferences and in the field to further grow our business and truly earn the title of "The In-Ear Experts" in 2019 and beyond! In the meanwhile, if you would like to contact me directly I would love to hear from you! My contact information is listed below.

Thank you again for all of your support,



Zubin Meshginpoosh
President & Chief Commercial Officer
Westone Laboratories
zubinm@westone.com
Cell: 248.231.6340

The *Finest* impression materials are available here.


Our impression materials are formulated to pick up every detail and feature in the ear, cure quickly and provide the best possible blueprint from which any custom-fit product can be produced.





Westone EZ Mold offers the most accurate and stable instant earmold material on the market! EZ Mold earpiece material was developed to fit any need, whether it's a temporary hearing instrument earpiece, or industrial and recreational hearing protection. In times of immediate need, look no further than Westone's EZ Mold!

- Product Features:
- Vinyl polysiloxane material
 - High viscosity
 - 7-10 minute set-up time
 - NRR 26



EZ Mold Singles

- Set includes:
- Color Choice
 - 24pk EZ Mold Singles
 - Oto-Ease .5oz Bottle
 - Guide
 - EZ Mold Patient Care Guides



EZ Mold Singles 24pk / Approx. 12 pairs

10331	24pk - Beige	\$72.00
10332	24pk - Pink	\$72.00
10333	24pk - Purple	\$72.00
10334	24pk - Orange	\$72.00
10335	24pk - Yellow	\$72.00
10336	24pk - Green	\$72.00
10337	24pk - Blue	\$72.00
10338	24pk - Brown	\$72.00
10339	24pk - Red	\$72.00
10341	24pk - Forest Green	\$72.00

EZ Mold Singles 22pk / Assorted Colors / Approx. 11 pairs

10342	22pk - Assorted Colors	\$72.00
-------	------------------------	---------

EZ Mold Tub Sets

Set includes:

- 1 EZ Mold 8oz Tub – Base (Part A)
- 1 EZ Mold 8oz Tub – Catalyst (Part B)
- Measuring spoon set (A & B)
- 20 Oto-Ease, single use packets
- Process Guide
- 30 Patient care guides



10389

EZ Mold 8oz Tub Set / Approx. 20 pairs

10381	8oz Tub Set - Beige	\$137.50
10382	8oz Tub Set - Pink	\$137.50
10383	8oz Tub Set - Purple	\$137.50
10384	8oz Tub Set - Orange	\$137.50
10385	8oz Tub Set - Yellow	\$137.50
10386	8oz Tub Set - Green	\$137.50
10387	8oz Tub Set - Blue	\$137.50
10388	8oz Tub Set - Brown	\$137.50
10389	8oz Tub Set - Red	\$137.50
10391	8oz Tub Set - Forest Green	\$137.50

EZ Mold 4oz Tub Set / Approx. 15 pairs

10349	4oz Tub Set - Beige	\$79.25
10350	4oz Tub Set - Pink	\$79.25
10351	4oz Tub Set - Purple	\$79.25
10352	4oz Tub Set - Orange	\$79.25
10353	4oz Tub Set - Yellow	\$79.25
10354	4oz Tub Set - Green	\$79.25
10355	4oz Tub Set - Blue	\$79.25
10356	4oz Tub Set - Brown	\$79.25
10357	4oz Tub Set - Red	\$79.25
10359	4oz Tub Set - Forest Green	\$79.25

EZ Mold Accessories and Tools

20180	Syringe 1pc, 10 Tips	\$8.95
42300	Pouch - 1pc	\$3.50
42346	Black Cord - 1pc	\$2.35
42350	Screw Cord Installation Tool	\$5.50

10375	Egger Lacquer L Nano 100ml	\$65.50
20312	Foam 9 mm - Black, 50pk	\$4.95
20313	Foam 11 mm - White, 50pk	\$4.95
20314	Foam 13 mm - Blue, 50pk	\$4.95



20180



42300



10375



20312



42346



42350



20313



20314

Siliclone®

Westone's SilliClone® impression material is the quintessential choice for any hearing healthcare professional in search of the best material available for taking ear impressions. Made from a vinyl polysiloxane material, Westone's SiliClone cartridges come in both DS-50 and DM-50 cartridge configurations, ensuring the material is compatible with any impression gun.

Available in 3 distinct materials: Standard SiliClone • SiliClone Firm • SiliClone X-Act

SiliClone®

SiliClone is manufactured to our own strict specifications, is highly accurate, and dimensionally stable.


- Soft vinyl polysiloxane material
- Low viscosity
- Shore A hardness of 25
- 4-7 minute set-up time
- Approximately 4-6 impressions per cartridge

Use model DMA 50 Gun

10101	SiliClone 48ml - 1 pc, No Tips	\$8.15
10102	SiliClone 48ml - 1 pc, 4 Tips	\$9.45
10103	SiliClone 48ml - 8pk, No Tips	\$62.00
10104	SiliClone 48ml - 8pk, 32 Tips	\$72.00

Use model DS 50 Gun

10151	SiliClone S-50 - 1 pc, No Tips	\$8.40
10152	SiliClone S-50 - 1 pc, 4 Tips	\$9.70
10153	SiliClone S-50 - 8pk, No Tips	\$64.00
10154	SiliClone S-50 - 8pk, 32 Tips	\$74.00



SiliClone® X-ACT

X-act combines the accuracy and dimensional stability of Westone's S-50 cartridge material with the expansion properties of our Pink Silicast syringe material. It is easily expressed from a manual impression gun. X-act's flow characteristics are more like those seen in a higher viscosity syringe material than the typical cartridge material.

- Vinyl polysiloxane material
- Medium to High viscosity*
- Shore A hardness of 20-22
- 3-4 minute set-up time
- Approximately 4-6 impressions per cartridge

*Used with DS 50 Impression Gun

Use model DS 50 Gun

10111	SiliClone X-ACT - 1pc, No Tips	\$9.45
10112	SiliClone X-ACT - 1pc, 4 Tips	\$11.00
10113	SiliClone X-ACT - 8pk, No Tips	\$74.00
10114	SiliClone X-ACT - 8pk, 32 Tips	\$84.00



SiliClone® FIRM

SiliClone FIRM is designed for those who want the ear canal expansion properties of our Pink Silicast but want to take advantage of the benefits of a cartridge impression material delivery system.

- Vinyl polysiloxane material
- High viscosity
- Shore A hardness of 25-30
- 6-10 minute set-up time
- Approximately 4-6 impressions per cartridge

Use model DS 50 Gun

10181	SiliClone FIRM - 1pc, No Tips	\$9.25
10182	SiliClone FIRM - 1pc, 4 Tips	\$10.75
10183	SiliClone FIRM - 8pk, No Tips	\$71.00
10184	SiliClone FIRM - 8pk, 32 Tips	\$81.00



SILICONE® Singles

Silicone Singles®

Silicone Singles are the most convenient, accurate, and unique silicone impression system on the market! One Silicone Single package contains enough impression material to make one impression.

- Vinyl polysiloxane material
- High viscosity
- Shore A hardness of 24-30
- 3-5 minute set-up time

10200	Silicone Single	\$1.75
10201	Silicone Singles - 10pk	\$13.20
10203	Silicone Singles - 24pk	\$31.00
10204	Silicone Singles - 50pk	\$63.00
10205	Silicone Singles - 100pk	\$120.00



Blend® Powder and Liquid Impression Material

This powder and liquid impression material is the only option we carry for patients with silicone allergies.

10500	24ct	\$35.75
10503	100ct	\$145.00



Silicast® PINK

Pink Silicast® Silicone Impression Material

Pink Silicast is for those who prefer to measure out their own impression material. This economical product is the same great material as our Silicone Singles, but in bulk form. Tubs come with color-coded dosing spoons.

- Vinyl polysiloxane material
- High viscosity
- Shore A hardness of 32
- 3-5 minute set-up time

10301	335g Tub Set, Approx. 45 Impressions	\$52.25
10302	525g Tub Set, Approx. 70 Impressions	\$81.00



10301

Silicast® BLUE

Blue Silicast® Silicone Impression Material

This highly accurate and dimensionally stable economical impression material is easy to use at a low per-impression cost. Blue Silicast will not shrink, melt, or deform. Each tub includes a tube of hardener and a dosing spoon.

- Addition-Cure Silicone material
- High viscosity
- Shore A hardness of 30
- 5-8 minute set-up time
- Approximately 50 impressions per tub

10401	Elastomer, 700g Tub	\$60.00
10402	Hardener, 1 Tube	\$15.00



10402



10401

Impression Kits

SiliClone® S-50 Starter Pack

The Starter Pack contains all the essential materials you will need to take S-50 impressions.

Includes:

- S-50 Cartridge Impression Gun (Model DS 50)
- SiliClone S-50 Impression Material Cartridges, 4pk, 16 mixing tips
- PROS Oto-Dams System, 4 sizes, 6 each (7 mm, 9 mm, 11 mm, 13 mm)

10132	Starter Pack	\$82.50
-------	--------------	---------



Impression Taking Traveler Kit

The ultimate impression kit for the mobile hearing healthcare professional.

Includes:

- Heine Mini 3000 Otoscope, hard-shell case, specula, 50pk (4 mm)
- Wahl Hair Trimmer
- S-50 Cartridge Impression Gun
- SiliClone S-50 Impression Material, 8pk, 32 mixing tips
- PROS Oto-Dams, 4 sizes, 24 each (7 mm, 9 mm, 11 mm and 13 mm)
- LED Earlite
- LED Earlite PROS Tips, 4pk
- Oto-Ease Lubricant, 100 individual doses
- Purell Hand Sanitizer
- Alcohol Prep Pads
- 25 zip-lock bags, 4" x 6" (10.2 cm x 15.2 cm)
- COS (Colorado) Shipping boxes
- Crush-proof carrying case, waterproof/airtight

10166	Traveler Kit	\$485.00
-------	--------------	----------

10168	Traveler Kit - International	\$485.00
-------	------------------------------	----------



Impression Taking Starter Kit

All the materials any beginner needs to start taking impressions.

Includes:

- Heine Mini 3000 Otoscope, hard-shell case, replacement specula, 50pk (4 mm)
- Wahl Trimmer
- S-50 Cartridge Impression Gun
- SiliClone S-50 Impression Material, 8pk, 32 mixing tips
- PROS Oto-Dams System, 4 sizes, 24 ea (7 mm, 9 mm, 11 mm and 13 mm)
- LED Earlite
- LED Earlite PROS Tips, 4pk
- Oto-Ease Lubricant, 100 individual doses
- Purell Hand Sanitizer
- Alcohol Prep Pads, 200ct
- Zip-lock bags, 25ct, 4" x 6" (10.2 cm x 15.2 cm)
- COS (Colorado Springs) Shipping boxes



10161	Starter Kit	\$365.00
10163	Starter Kit - International	\$365.00

Impression Taking Refill Kit

Everything needed to restock the Impression Taking Traveler and Starter Kit.

Includes:

- Heine Otoscope Specula, 50pk (4 mm)
- SiliClone S-50 Impression Material, 8pk, 32 mixing tips
- PROS Oto-Dams, 4 sizes, 24 each (7 mm, 9 mm, 11 mm and 13 mm)
- LED Earlite PROS Tips, 4pk
- Oto-Ease lubricant, 100 individual doses
- Purell Hand Sanitizer
- Alcohol Prep Pads
- Zip-lock bags, 25ct, 4" x 6" (10.2 cm x 15.2 cm)
- COS (Colorado Springs) Shipping boxes




10162	Refill Kit	\$145.00
10164	Refill Kit - International	\$145.00

Impression Guns

Westone® Practicon BFC3 Powered Impression Gun

Making a great impression is now easier and faster than ever with our Powered Impression Gun. (Better Faster Comfort, Control, Compatibility) dispenses any type of 1:1 or 2:1 automix material at the press of a button. No repetitive trigger squeezing to wear out your hand and throw off your aim. Works with SB-50 cartridges. Rechargeable battery dispenses up to sixty 50ml cartridges on a full charge. Convenient, user-friendly, handheld design fits smoothly into normal practice workflow.


20115 Impression Gun \$299.99



48ml Cartridge Impression Gun

Manual 48ml Cartridge Impression Gun is lightweight and easy to use. Model DMA 50 is designed for all 48ml impression material cartridges.


20111 Impression Gun \$40.00



S-50 Cartridge Impression Gun

Manual S-50 Cartridge Impression Gun is lightweight and easy-to-use. Model DS 50 is designed for use with all S-50 cartridges.

20112 Impression Gun \$49.00



Trigger Insert for Manual Impression Gun

Trigger Insert shortens the throw of the trigger, making manual impression guns easier to use.

20113 Trigger Insert \$2.50



Mixing Tips

SiliClone Cartridge Mixing Tips

Designed for 48 ml cartridges, this tip completely mixes the impression material as it is expressed into the ear. Discount pricing available on quantities of three or more.

20141	Adult, 5.4 mm, 32pk	\$19.00
-------	---------------------	---------



SiliClone FIRM Cartridge Mixing Tips

Designed for SiliClone FIRM or any other S-50 impression material cartridges, these tips completely mix the impression material as it is expressed into the ear. Discount pricing available on quantities of three or more.


20148	6.5 mm, 32pk	\$20.50
-------	--------------	---------



SiliClone S-50 Cartridge Mixing Tips

Designed for S-50 cartridges, these tips completely mix the impression material as it is expressed into the ear. Discount pricing available on quantities of three or more.

20143	Pediatric, 4.2 mm, 32pk	\$20.50
-------	-------------------------	---------



20146	Adult, 5.4 mm, 32pk	\$20.50
-------	---------------------	---------




Syringes

Silicast Syringe

Ideal for use with Silicone Singles and Pink Silicast impression materials. The unique design of the Silicast Syringe allows impression material to flow with less pressure.

Use silicone impression material only. Powder and liquid impression material will damage the Silicast Syringe.

20131	Silicast Syringe	\$18.25
20151	Replacement O-rings, 5pk	\$2.75



20131



20151

Polyethylene Syringe

The economical, high-quality polyethylene syringe can be used with any silicone or powder and liquid impression material.

20121	Polyethylene Syringe	\$4.50
-------	----------------------	--------



Disposable Extension Tips

Use with both the Silicast Syringe and Polyethylene Syringe to take impressions of smaller ears. Disposable extension tips help guard against infection.

20152 Tips, 25pk \$7.75



The BEST Syringe

The compact design of The BEST Syringe is suitable for any size hand. The syringe holds enough impression material for one complete ear impression. To help guard against infection, the specula are considered disposable tips.


Use silicone impression material only.

Pediatric

20182	Syringe 1pc, 10 tips	\$8.95
20184	Syringe 4pk, 40 Tips	\$35.00
20198	Retainer Ring, Red, 1pc	\$1.25

Standard

20181	Syringe 1pc, 10 Tips	\$8.95
20183	Syringe 4pk, 40 Tips	\$35.00
20197	Retainer Ring, Gray, 1pc	\$1.25



Clear

20180	Syringe 1pc, 10 Tips	\$8.95
20186	Syringe 4pk, 40 Tips	\$35.00
20196	Retainer Ring, Translucent Red, 1pc	\$1.25

Assortments

20185	Syringe 4pk, 20 Std Tips / 20 Peds Tips	\$35.00
20195	Assorted O-rings, 2pk	\$0.35
23121	Pediatric Specula, 2.5 mm - 50pk	\$3.50
23131	Specula, 4 mm - 50pk	\$4.00



Accessories

Metal Spatula

Metal spatula with wooden handle is ideal for mixing BLEND or silicone impression materials.

20171	Spatula, 1pc	\$5.50
-------	--------------	--------



SPLEAD Pad

SPLEAD pads can be used as a mixing surface for silicone impression materials and other lab supplies, such as adhesives, and earpiece or shell build-up materials.

- 35 Wax-impregnated sheets per pad
- Features a foam, non-slip surface on the bottom

20161	Pad, 2pk	\$13.50
-------	----------	---------



Bite Blocks

These one-time-use styrofoam spacers hold a patient's mouth open while taking open jaw impressions. Westone recommends open jaw/bite block impressions for hearing protection devices and musicians' products.

20251	Bite Block, 100pk	\$33.95
20170	Bite Block, 25pk	\$11.50



Demonstration Ears

These anatomically correct models of the external ear are ideal for learning how to take impressions, or to show patients how earpieces and hearing instruments interact with the ear.

Available as a pair (right and left) only.

20221	Demonstration Ears - One pair	\$25.00
-------	-------------------------------	---------



Wahl® Ear, Nose, Brow Trimmer

Battery operated trimmer is perfect for trimming hair away from the ear canal. Steel blades rinse clean under running water.

Includes:

- Rotary head
- Reciprocating head
- Detailing head
- 1 AA battery

20241 Ear, Nose, Brow Trimmer \$12.50



Duco® Cement

Duco Cement is the ideal adhesive for gluing powder and liquid impressions into shipping boxes. Duco Cement dries quickly and helps to ensure that powder and liquid impressions are not damaged during shipment.

20261 1oz Tube \$2.95



LevelEar Impression Marker

Accurately marks ear impressions for hearing instruments that include directional microphones. The built-in bubble level and corresponding straightedge help to easily create impressions before they cure.

Laboratories can use these marks as a perfect reference for placing directional microphones.

20175 Impression Marker, 1pc \$8.25



Oto-Dams

PROS® Pressure Relief Oto-Dam System

These revolutionary oto-dams are manufactured with a ventilation tube that relieves air pressure that can build up while impression material is injected into the ear and when the cured impression is removed. Each pack includes a specially designed slotted earlite tip.



24 Pack

20301	Assortment, 6 each	\$17.50
20302	7 mm - Green	\$17.50
20303	9 mm - Black	\$17.50
20304	11 mm - Beige	\$17.50
20305	13 mm - Blue	\$17.50

Slim PROS® Pressure Relief Oto-Dam System

Vented just like our standard PROS oto-dams but half as thick, Slim PROS are perfect for deep canal impressions. Each package includes one PROS slotted earlite tip.




24 Pack

20306	Assortment, 6 each	\$17.50
20307	9 mm - Black	\$17.50
20308	11 mm - Beige	\$17.50
20309	13 mm - Blue	\$17.50

Foam Oto-Dams

Inexpensive and designed to fit any patient, these foam ear dams feature a string removal filament.



50 Pack

20311	7 mm - White	\$5.35
20312	Foam 9 mm - Black, 50pk	\$4.95
20313	Foam 11 mm - White, 50pk	\$4.95
20314	Foam 13 mm - Blue, 50pk	\$4.95

Slim Stops

Slim Stops are approximately half the thickness of foam ear dams, and are perfect for deep canal impressions.

50 Pack




20321	9 mm - Black	\$5.35
20322	11 mm - White	\$5.35
20323	13 mm - Blue	\$5.35



50 Pack Cotton Oto-Dams

Cotton oto-dams are ideal for deep canal impressions, and are pre-tied.

20331	Size 4 - Mini	\$6.20
20332	Size 2 - Small	\$6.20
20333	Size 1 - Medium	\$6.20
20334	Size 0 - Large	\$6.20
20335	Size 3 - Small/Mini	\$6.20



Your *Premier* selection for all clinical supplies.

At Westone, we pride ourselves on being a one-stop-shop to meet any and all of your audiological supply needs. From the highest quality Otoscopes and Stethoscopes to the latest test equipment and everything in-between, all of our clinical products are guaranteed to leave not only you, but your patients 100% satisfied.



Otoscopies

Heine® Mini 3000 Otoscope LED Illumination

Compact modern pocket Otoscope with maintenance-free LED illumination

Features:

- Maintenance-free illumination: LED cannot be changed
- LED Thermal management for constant light output
- Illumination with Excellent color rendering
- Twice as bright as Xenon Halogen instruments
- Exclusive Battery Performance Indicator: know when your AA batteries need changing
- Exclusive fade-out function reduces illumination intensity
- Operation time up to 10 hours, battery powered

Includes:

- Unique attachment clip with integrated on/off switch
- High-quality handle is shock proof and features a sturdy no-slip grip
- (2) AA batteries
- (5) 2.5 mm and (5) 4 mm diameter AllSpec disposable tips
- 5 year warranty (excludes bulb)

23068	Head Only	\$195.00
23069	Otoscope	\$245.00



Heine® Beta 400 LED Fiber Optic Otoscope

This otoscope has been designed with an all metal housing and a scratch resistant glass window which guarantee a lifetime of maintenance-free use. LED bulb with virtually unlimited hours of working life (up to 100,000 hours). Gives the necessary overview and perfect orientation during examinations.

Features:

- 4.2 x Magnification
- All-metal housing
- LED HQ with continuous brightness control

Bundle Includes:

- BETA 400 LED Fiber Optic Otoscope with Heine AV Connector
- 1 set (4 pieces) of Reusable Tips: 2.4/3/4/5 mm diameter
- 5 each of AllSpec disposable tips 2.5 and 4 mm
- Hard Case (Note: Power Supply will not fit inside the Case)
- LED Bulb
- 5 Year warranty

23085	Otoscope	\$913.00
-------	----------	----------



Heine® Mini 3000 Otoscope

Compact modern pocket Otoscope features a unique design with direct illumination.

Features:

- Enhanced Heine XHL® Xenon Halogen Technology produces a very bright concentrated white light for perfect illumination with 100% more light than conventional bulbs
- Swiveling view window with 3x magnification and optimized casing surface for razor-sharp images and minimal reflection

Includes:

- Unique attachment clip with integrated on/off switch
- High-quality handle is shock proof and features a sturdy no-slip grip
- (2) AA batteries
- 1 set (4 pieces) reusable tips
- (5) 2.5 mm and (5) 4 mm diameter AllSpec disposable tips
- 5 year warranty (excludes bulb)



23070	Otoscope - Black	\$130.00
23071	Otoscope - Blue	\$130.00
23079	Replacement Bulb, X-01.88.110, 2.5 V	\$30.00



Heine® Mini 3000 Otoscope Fiber-Optic Illumination

Compact modern pocket Otoscope comes with Fiber-Optic (F.O.) illumination.

Features:

- Enhanced Heine XHL® Xenon Halogen Technology produces a very bright concentrated white light for perfect illumination with 100% more light than conventional bulbs
- Swiveling view window with 3x magnification and optimized casing surface for razor-sharp images and minimal reflection

Includes:

- Unique attachment clip with integrated on/off switch. Turns off automatically when returned to the pocket; switch is guaranteed for up to 20,000 cycles
- High-quality handle is shock proof and features a sturdy no-slip grip
- (2) AA batteries
- 4 reusable tips, 2 pediatric, 2 large
- (5) 2.5 mm and (5) 4 mm diameter AllSpec disposable tips
- 5 year warranty (excludes bulb)



23080	Otoscope - Black	\$220.00
23081	Otoscope - Blue	\$220.00
23083	Head Only - Black	\$150.00
23089	Replacement Bulb Fiber-Optic, X-01.88.105, 2.5 V	\$30.00

Heine® AllSpec® Tip Dispenser

Convenient wall-mounted unit holds four sleeves of either 2.5 mm or 4 mm specula. Mount on the examination room wall for quick access. Includes 50 each of 2.5 mm and 4 mm specula.

23141	Tip Dispenser	\$42.00
-------	---------------	---------



Heine® Otoscope Hard Case

Protect your otoscope from damage with a durable, hard shell case. Fits any Heine Mini 2000 or 3000 Otoscope, or the Heine Mini 2000 or 3000 (FO) Fiber Optic Otoscope.

23100	Otoscope Hard Case	\$31.50
-------	--------------------	---------



Replacement Halogen Bulbs for Heine® Otoscopes

Replacement bulb for the Mini 2000 Heine Otoscope, Fibr lux, Mini 2000 Fiber optic, and Beta 200 Heine Otopscopes.

23107	X0188056 Mini 2000 Bulb	\$26.00
23108	X0188041 Fibr lux Bulb	\$26.00
23109	X0188037 Mini 2000 FO Bulb	\$26.00
23111	X0288078 3.5V Bulb	\$34.00



23107



23108



23109



23111

Replacement Specula for Heine® Otoscopes

Replacement specula for the Heine Otoscope.

Discount pricing:

- 2 or more = \$3.13 each
- 20 or more = \$1.82 each

23121	Pediatric Specula, 2.5 mm - 50pk	\$3.50
23131	Specula, 4 mm - 50pk	\$4.00



23121



23131

Earnimals

Animal shaped graphics that fit on most otoscopes to help relax younger patients. Assorted packages only.

Available in two sizes:

- Standard - fits most otoscopes with disposable specula
- Large - fits operating-type otoscopes.

23191	Standard, Assorted 6pk	\$3.85
23192	Standard, Assorted 24pk	\$12.50
23193	Large, 5pk	\$4.00



KleenSpec® Replacement Specula

23112	Replacement Specula, 5 mm	\$4.00
23113	Replacement Specula, 3 mm	\$4.00
23114	Replacement Specula, 4 mm	\$4.00



Welch Allyn® Otoscope

The 2.5V Halogen PocketScope™ Otoscope with Throat Illuminator is small, convenient, and easy to use.

- Halogen light for true tissue color and consistent, long-lasting illumination
- Detachable otoscope head creates throat illuminator or all-purpose penlight
- Fiber optics: cool light, no reflections, no obstructions
- Wide-angle viewing lens
- Sealed system for pneumatic otoscope
- Two AA batteries (included)
- 1 year warranty

23049	Replacement Otoscope Bulb	\$25.00
23050	2.5V Pocket Otoscope	\$229.00




Welch Allyn® Replacement Specula

Replacement specula for standard Welch Allyn otoscopes.

- Discount pricing:
- 2 or more = \$3.50 each
 - 25 or more = \$2.40 each

23115	Specula, 2.75 mm - 34pk	\$4.00
23119	Specula, 4.25 mm - 34pk	\$4.00



Welch Allyn® Ni-CAD Batteries

Rechargeable 3.5 V Ni-CAD replacement batteries for Welch Allyn Otoscopes.

When ordering replacement batteries please be sure to check the color code on the otoscope handle to ensure that the correct replacement battery is ordered.

23117	Battery 72300	\$45.00
23118	Battery 72200	\$45.00



Earlites

LED Earlite Batteries

Replacement 1.5 V batteries for LED Earlite

24109	Battery, 3pk	\$7.00
-------	--------------	--------



Westone® Earlite

The 18 Lumens light is concentrated in a small tip that fits most ears. (2) AAA batteries and pocket clip included. Works with tips 24141, 24142, 24143, 24144 and 24145.

24102	Earlite with 2 Tips	\$7.55
24151	Earlite, Westone Logo	\$5.95
24152	Earlite, Unprinted/No Logo	



24102



24151



24152



24151

E.N.T. Pocket Light

Powerful xenon blub is 400% brighter than ordinary blubs and is perfect for oto-dam placement.

- Includes:
- (2) AAA batteries
 - One straight Earlite tip

23153	Pocket Light, Black	\$22.95
23154	Pocket Light, Blue	\$22.95



E.N.T. Pocket Light - Replacement Xenon Light Bulb

Replacement Xenon bulb for ENT Pocket Light

23155	Replacement Xenon Light Bulb	\$7.95
-------	------------------------------	--------



E.N.T. Pocket Light, Cerumen Management Set

Using a powerful Xenon bulb, this handy light is 400% brighter than those using ordinary bulbs.

Compatible with all Westone earlite tips 24121, 24122, 24123, 24124 and 24125.

Includes:

- E.N.T. Pocket Light
- (2) AAA batteries
- Xenon bulb
- (2) Lighted SofTouch cerumen loops
- (2) Lighted stainless steel cerumen loops

23156	Replacement Lens Cap	\$3.95
23161	ENT Pocket Earlite	\$32.95



E.N.T. Pocket Light, Deluxe Set

The small, lightweight ENT Pocket Light features a bright xenon light for effective examination of the ear. Its versatility, durability, and ease of use make it ideal for both clinical and field settings.

Includes:

- E.N.T. Pocket Light
- (2) AAA batteries
- Xenon bulb
- (2) Lighted SofTouch cerumen loops
- (2) Lighted stainless steel cerumen loops
- (3) Disposable adult specula (4 mm)
- Spare lamp module

23151	Pocket Light, Black	\$42.95
23152	Pocket Light, Blue	\$42.95



Dino-Lite Digital Earscope

A digital microscope solutions ideal for capturing and saving images of the ear. This easy-to-use, compact video otoscope connects to a computer or laptop via USB.

Features:

- Resolution 640 x 480 pixels
- Connects to USB 2.0
- Easy to load and use Software
- 4 white LEDs built in
- Video rate: up to 30 frames-per-second
- Body length: 4" (10.2 cm)
- Body diameter: 1.25" (3.2 cm)
- Cable length: 6' (180 cm)

Includes:

- Dino-Lite Digital Earscope
- Dino Capture application software
- User's manual
- (5) Welch Allyn KleenSpec disposable specula 5 mm
- Deluxe storage bag
- 1 year warranty

23201	Digital Earscope	\$249.00
-------	------------------	----------



Krypton Replacement Bulb

Krypton replacement bulb for the Oto-Probe II and the Home Care Ear Examination Kit

24105	Replacement Bulb	\$1.75
-------	------------------	--------



Earlite Tips

Replacement Earlite Tips

Replacement tips that fit the Disposable Earlite, the Oto-Probe II and the ENT Pocket Light.

24121	Earlite Tip, 4pk	\$1.95
24122	PROS Earlite Tip, 4pk	\$5.15
24139	Earlite Tip, 100pk	\$25.00



24121



24122

Cerumen Management Tips

Three different styles of cerumen management tips designed to fit the Disposable Earlite, Oto-Probe II, and ENT Pocket Light.

24123	Scoop Tips, 4pk	\$2.75
24124	Stainless Steel Loop Tip, Small/Large - 2pk	\$10.50
24125	Nylon Loop, Small/Large - 2pk	\$10.50



24123



24124



24125

Welch Allyn® Earlite Tips

Replacement tips for use with Welch Allyn Earlites.

24131	Straight Tip - 4pk	\$2.25
24133	Plastic Scoop - 4pk	\$2.75



24131



24133

LED Earlite Cerumen Management Tips

Cerumen Management tips can be used to facilitate cerumen removal. Sized to fit the Coast LED Earlite. Will work with earlite 24102, 24151 and 24152.

24143	Cerumen Management Tip, 4pk	\$4.50
24144	Stainless Steel Loop, Small/Large - 2pk	\$12.00
24145	Nylon Loop, Small/Large - 2pk	\$12.00



24143



24144



24145

LED Earlite Replacement Tips

Replacement tips to fit Rayovac LED Earlite, Coast LED earlite and LED lenser.

24138	PROS Tip, 100pk	\$150.00
24140	Straight Tip, 100pk	\$55.00
24141	Straight Tip, 4pk	\$3.00
24142	PROS Tip, 4pk	\$6.00



Conductivity Enhancers

NuPrep® Skin Prep Gel

Nuprep Skin Prep effectively lowers impedance to improve your tracings. Its mild abrasive formula improves conductivity and helps achieve maximum efficiency with your equipment.

30681	4oz Tube, 3pk	\$21.95
30683	25g Tube, 6pk	\$12.90



Ten20® Conductive Paste

Ten20 contains the right balance of adhesiveness and conductivity, enabling the electrodes to remain in place while allowing the transmittance of electrical signals.

30685	4oz Jar	\$18.50
30687	4oz Tube	\$18.80




Equipment

Wells Dental Gold Chrome Finishing Machine

Modify earmolds of any material with this high-speed bench grinder. Includes a Wells 3/32" high-speed spindle that allows quick tool changes while the lathe is running. Low speed = 12,000 RPM. High speed = 24,000 RPM.

Note: Ships direct from the manufacturer. Please allow 2-4 weeks for product delivery. Additional shipping charges will apply.



41101	Gold Chrome Finishing Machine	\$1,475.00
-------	-------------------------------	------------

Buffalo Dental V35 Hand piece System

This 120 V electric system features a no-frills, multi-speed control console teamed with the Buffalo Magic Wand™ gray hand piece. This combination provides a heavy-duty system that offers durability, high torque, and economical pricing. Maximum speed is 35,000 RPM for modifying earpieces of any material.

Note: Ships direct from the manufacturer. Please allow 2-4 weeks for product delivery. Additional shipping charges will apply.



41112	Replacement Hand piece	\$385.00
41113	V35 Hand piece System	\$425.00

Handler™ Lathe Shield and Light

A lighted aluminum splash shield that mounts alongside the motor for increased protection from flying debris. The shield measures 12" h x 6" w x 10" d. Includes one insert pan.

Please specify if the shield is for the right or left side.

Note: Ships direct from the manufacturer. Please allow 2-4 weeks for product delivery. Additional shipping charges will apply.




41122	Replacement Safety Shield	\$28.00
41123	Lathe Shield and Light	\$95.00

Handler™ 26A Red Wing Lathe

Durable lathe is designed to provide years of service. Features a 1/4 hp, two-speed grinder and polisher with standard dental, tapered shafts. Low speed = 1725 RPM. High speed = 3450 RPM.

Note: Ships direct from the manufacturer. Please allow 2-4 weeks for product delivery. Additional shipping charges will apply.



41121	26A Red Wing Lathe	\$350.00
-------	--------------------	----------

Handler™ Repair and Modification Kit

The perfect solution for hearing instrument repair and earpiece modification.

Includes:

- Handler 26A Red Wing Lathe
- 5/32 inch brass and steel chuck, right side
- Circular saw
- Aluminum oxide points- No:s 75, 89 and 91
- Pinvise and 4 chucks
- 5 inch 45 ply muslin buff
- 5 SA, SS anti-magnetic tweezer
- Pair of safety goggles
- 6 Ball Burs one of each of the following sizes: .055 inch, .083 inch, .091 inch, .098 inch, .114 inch, .122 inch, Tapered spindle, left side



Note: This product ships directly from the manufacturer. Please allow 2-4 weeks for product delivery. Additional shipping charges will apply.

41131	Kit	\$435.00
41132	Kit without lathe/bench grinder	\$135.00
41133	Ball Burs only - 6pk	\$12.00

Handler™ Porta-Vac 550 Dust Collector

Compact solution for the dust and debris created by in-office modifications. With 225 CFM of powerful suction, the Porta-Vac collects dust particles 0.5 microns and larger. Compatible with either hand-held or bench grinders, the lighted unit also features a removable clean-out tray, replaceable acrylic shield, and reusable/replaceable filters.

Note: This product ships directly from the manufacturer. Please allow 2-4 weeks for product delivery. Additional shipping charges will apply.



41124	Porta-Vac 550 Dust Collector	\$365.00
-------	------------------------------	----------

Grinding Tools & Accessories

Handler™ 8T Chuck

Replacement chuck for the Red Wing Lathe. Accepts bits and burrs with a 1/4" shaft or less. Right side only.

41138 Handler 8T Chuck \$27.50



Modification Tool Combo Kit for Hard and Soft Earpieces

Kit contains all the tools required for hard and soft earpiece modification.

Includes:

- Tube puller - Monofilament
- Tube puller - Wire Loop
- Blue grinding stone
- Reamer for #13 tube
- #33 Drill Bit - 1/8" shank
- #52 Drill Bit - 3/32" shank
- White grinding stones #4 and #6
- Burr .037" (.94 mm)
- Burr.033"(.84 mm)
- Core Drill for tubes
- Core Drill for large bore
- Core Drill for venting
- Mandrel and screw
- Small polishing buffs - 3pk

41141 Combo Kit \$62.50



Modification Tool Kit for Hard Earpieces

Contains everything necessary to modify hard earpieces.

Includes:

- Tube puller - Monofilament
- Reamer for #13 tube
- #33 Drill Bit - 1/8" shank
- #52 Drill Bit - 3/32" shank
- Burr .037" (.94 mm)
- Burr.033"(.84 mm)
- Mandrel and screw
- Small polishing buffs - 3pk
- White grinding stones #4 and #6

41142 Kit for Hard Earpieces \$35.00




Modification Tool Kit for Soft Earpieces

Contains everything necessary to modify soft earpieces.

Includes:

- Tube Puller - Wire Loop
- Core Drill for tubes
- Core Drill for large bores
- Core Drill for venting
- White grinding stones #6
- Blue grinding stone
- Burr .033" (.84 mm)



41143	\$42.50
-------	---------

Tapered Spindle for Bench Mounted Grinders

Tapered spindle press fits on bench grinder shafts. The threaded design allows quick and easy replacement of polishing wheels.


41136	Tapered Spindle - Left side	\$10.50
41137	Tapered Spindle - Right side	\$10.50



Blue Grinding Stone

Coarse blue grinding stone with a 3/32" shank. Perfect for modifying all silicone earmolds.

41201	Blue Grinding Stone	\$8.50
-------	---------------------	--------



Fast-Kut Grinding Stone

Medium coarse grinding stone with a 3/32" shank. Perfect for modifying vinyl and acrylic earmolds.

42110	#4, Small	\$9.50
42111	#6, Large	\$9.50



White Grinding Stone

Medium coarse grinding stone with 3/32" shank for modifying vinyl and acrylic earpieces.

42112	#4, Small	\$5.00
42113	#5, Medium	\$5.00
42114	#6, Medium Coarse - Large	\$5.00



Super Fine Carbide Cutter

Super fine carbide cutter perfect for precision shell modification.

42115	Carbide Cutter	\$27.95
-------	----------------	---------



Excel® Assorted Drill Bits and Index

Includes 20 bits ranging from #61 to #80. Features a clear plastic dome over the indexed base.

42131	Assorted Drill Bits and Index	\$18.95
-------	-------------------------------	---------



High-Speed Drill Bit

High-speed steel drill bit.


42141	#33 Steel Bit	\$2.00
42144	#41 Steel Bit	\$4.00
42147	#52 Steel Bit, 3/32" shank	\$3.00



Eyelit Kit with Drill

Small eyelets for installation in the faceplate of an ITE hearing instrument. The eyelet kit contains five eyelets and a drilling tool. Perfect for use with all ITE Oto-Clips

42121	Kit with drill, 5pk	\$10.50
42125	Eyelit, 1pc	\$0.85



Excel® Swivel Head Pin Vise

This Excel® Pinvise with 2 Collets is perfect for hearing instrument repair and other detail work. Features two double-end collets. Accepts a wide range of drill bits from #30 to #80.

42132	Swivel Head Pin Vise	\$5.50
-------	----------------------	--------



Precision 6 Piece Screwdriver Set

Includes four flat-tip screwdrivers and two Phillips tips.

41186	6pc Set	\$4.00
-------	---------	--------



Precision 11 Piece Screwdriver Set

Six flat-tip screwdrivers, three Phillips tips, one needle-point tip, and one magnet tip.

41187	11pc Set	\$5.50
-------	----------	--------



OptiVisor® HB Magnifier

Magnifying visor is available in a choice of three magnifying powers, with an adjustable headband and a flip-up visor.

41151	1.75x Power Lens, Working distance: 14" (35.5 cm)	\$35.00
41152	2x Power Lens, Working distance: 10" (25.4 cm)	\$38.50
41153	2.5x Power Lens, Working distance: 8" (20.3 cm)	\$38.50



Truline Safety Glasses

A close, stylish fit is available with this fashionable, lightweight wrap-around lens. Features a soft, flexible nose piece and fold-up temples with soft, rubber sleeves for comfort. Complies with ANSI.

41192	Safety Glasses, 1 pair	\$6.50
-------	------------------------	--------



Excel® End Nipper Pliers

The perfect tool for trimming tubing from the end of a canal after the earmold has been retubed.

42130	Pliers	\$6.50
-------	--------	--------



Needle Nose Pliers

Needle nose pliers for a multitude of uses in the hearing healthcare clinic.

42133	5" Needle Nose Pliers	\$6.50
41180	4.5" Needle Nose Pliers	\$12.75



42133



41180

Excel® Sprue Cutter

A unique design easily trims tubing flush to the earpiece.

42134	4.5" Sprue Cutter	\$6.50
41177	Euro Tool 5" Flush Cutter	\$15.50



42134



41177

Alligator Forceps

Standard 5.5" stainless steel alligator forceps with serrated jaw. Not intended for invasive medical procedures.

41179 Stainless Steel, 5.5" \$29.00



Excel 6" Curved point Tweezers

Bent tip, precision stainless steel tweezers for working on hearing instruments and other small parts.

41181 6" Tweezers \$5.25



Excel® Tweezers

Fine-tipped tweezers for a variety of clinical uses.

42135 4.75" Tweezers \$4.00



Stainless Fine Point Tweezers

The fine point stainless steel tweezers are ideal for a variety of clinical uses.

41178 Anti-Magnetic 4.25" \$8.25



Probe Set - Stainless Steel

An incredibly handy tool set that belongs in every home and professional toolbox. Will get into tiny cracks and crevices other tools can't for probing, material removal.

41176 Probe Set, 5.5"- 6.5" - 4pk \$19.75



Excel® Soft Grip Razor Knife

Soft grip razor knife is great for trimming a variety of materials and features a rear locking device and safety cap. Shipped with six blades.

42136	Razor Knife	\$4.75
42153	Blades, 5pk	\$1.35




Excel® Precision Cutting Kit

Protect work surfaces with this 3-ply, self-healing cutting mat. Includes a sharp Rite-Kut knife with safety cap.

Cutting mat dimensions are 5 1/2" x 9".

42137	Kit	\$8.25
-------	-----	--------



Excel® Scissors

Sharp scissors suitable for general use in the clinic.

42139	3.5" Scissors	\$3.75
-------	---------------	--------



Excel® Extra Hand with Magnifier

Two spring-loaded alligator clips mounted on ball joints will hold parts at any angle. The 2X magnifier enables the user to view even the smallest projects.

42138	Magnifier	\$16.00
-------	-----------	---------



Excel® Sanding Stick with Extra Belt

Ideal for minor modification and smoothing of any earmold. Includes an extra sanding belt.

42151	Sanding Stick with Extra Belt	\$5.00
-------	-------------------------------	--------



Buffing Wheels

Ideal for polishing modified earmolds and hearing instrument shells. The 5 inch diameter muslin buffing wheels attach to spindle-equipped bench grinders. Available in three different thicknesses.

42171	12 ply	\$3.50
42175	20 ply	\$4.15
42177	50 Ply	\$4.65



Small Polishing Buffs

The 1 1/4" polishing buffs are sized for use with hand-held grinders that will accept a 3/32" shaft.

42181	Small, 3pk	\$4.75
-------	------------	--------



Mandrel and Screw

Designed for 1 1/4" polishing buffs, the mandrel and screw fit any hand-held grinding tool that accepts a 3/32" shaft tool.


42191	3/32"	\$2.00
-------	-------	--------



Nelson Tools™ Replacement Brushes

Replacement brushes for The Cleaning Station & Maintenance Caddy. Each brush is coded with a colored dot on the end of the handle.

42245	Black dot, best for vents on ITE instruments	\$3.50
42246	Blue dot, best for battery compartments	\$3.50
42247	Brown dot, best for vents and earpiece tubing	\$3.50
42248	Red dot, best for mic openings and hearing instruments	\$3.50
42249	Tan dot, best for 3mm and 4mm tubing and delicate parts	\$3.50
42250	White dot, best for hearing instruments and battery contacts	\$3.50
42251	Yellow dot, best for small vents	\$3.50



Buffing Compound

A non-abrasive jeweler's polish used for hand or wheel buffing of earpieces and hearing instruments.

42161 .5lb bar \$8.50



Nelson Tools™ Tool Kit

The maintenance caddy contains useful tools for cleaning and maintaining hearing instruments and earpieces. Features an impact-resistant, foam-filled carrying case.

Includes:

- 1 Large black bristle brush
- 6 White bristle brushes, including vent sizes
- 2 Wire loops
- 3 Wire picks

42231 Tool Kit \$30.00



Nelson Tools™ Wire Picks

Replacement picks for the Nelson Maintenance Caddy or Cleaning Station. Each pick is coded with a colored dot on the end of the handle.

42252	Gold dot, .012 gauge, Removes cerumen baskets	\$3.50
42253	Pink dot, .009 gauge, Removes debris from receiver ports and microphones	\$3.50
42254	Silver dot, .028 gauge, Removes batteries from instruments	\$3.50



Nelson Tools™ Wire Loops

Replacement loops for the Nelson Maintenance Caddy or Cleaning Station. Each loop is coded with a colored dot on the end of the handle.

42255	Green dot, Removes cerumen from ITE and BTE instruments	\$3.50
42256	Orange dot, Heavy gauge, Removes cerumen from ITE and BTE	\$3.50



Stethoscopes

Eko Core Digital Stethoscope

53115	Eko Core Digital Stethoscope	\$299.99
53116	Attachment Only	\$199.00

In the Box:

- Eko Core Attachment
- Cardiology-grade analog stethoscope
- Pediatric & adult diaphragms
- Wall Charger and Charging cable
- 1 Year Eko Warranty




Eko Core Digital Stethoscope

The Eko Digital Stethoscope is a complete electronic stethoscope with digital & analog capabilities that wirelessly pairs to the free Eko Mobile App for recording, visualization, sharing, and annotation of sounds.



53116



View it. Hear it.
Record it. Share it.
Stream it. Analyze it.
Save it to an EHR.

Connect to your Eko Core Digital Stethoscope wirelessly with mobile and web apps.

Westone® Listening Stethoscope

Popular with hearing healthcare professionals and the parents of hearing impaired children, this inexpensive and easy-to-use stethoscope comes with two listening tips. One tip is a probe tip for use with ITE hearing instruments and the other a listening cup for use with BTE or body worn hearing instruments. Available in Black, Blue or Red.

53111	Black, 1pc	\$18.00
53112	Blue, 1pc	\$18.00
53113	Red, 1pc	\$18.00



ADC Adscope Vistascope

Unique chest piece design features eight interchangeable diaphragms that provide up to twice the amplification of ordinary acoustic stethoscopes.

- Includes:
- Large and small PVC eartips
 - Instructions
 - Lifetime warranty

53121	Stethoscope, 1pc	\$57.50
53123	Replacement Eartips - Large, Black	\$3.75
53124	Eartip Adapter	\$3.75



Steth-o-Mate

Replaces the hard plastic tips found on most stethoscopes that make them nearly impossible for use by medical professionals who wear hearing instruments. This soft silicone eartip allows a person wearing a CIC or ITC hearing instrument to use a stethoscope without removing the hearing instrument.

53311	.450" (11.43 mm) dia. for CICs, 2pk	\$23.00
53312	.700" (17.78 mm) dia. for ITCs, 2pk	\$23.00



Steth-o-Mate for 3M® Littman® Stethoscopes

Same great product as the original Steth-o-Mate but sized for Littman stethoscopes. Available in black and as a pair only.

53313	.450" (11.43 mm) dia. for CICs, 2pk	\$23.00
53314	.700" (17.78 mm) dia. for ITCs, 2pk	\$23.00



Lightweight Listening Stethoscope

Check the performance of hearing instruments with this lightweight and inexpensive all plastic stethoscope. Features a listening cup for use with BTE earpieces, ITE hearing instruments and body worn hearing devices.

53114 Stethoscope, 1pc \$13.00



Adapt-a-Scope

Large metal male adapter designed for use with stethoscope tubing. The adapter fits into a standard snap ring like the one used in Westone's Style No.1 Earmold.

53351 Metal Adapter \$4.00



Y Adapter for Duoscope

Y Adapter designed for use with stethoscope tubing. Splits a monaural signal into a binaural signal.

53361 Y Adapter \$2.85



ITE Probe Listening Tip

Slips into the sound bore of a hearing instrument and allows the hearing healthcare professional to listen to the instruments output.

53211 Listening Tip \$1.90




Listening Tube

Economical hearing instrument listening tube is perfect for mobile professionals and the parents of hearing impaired children.

Features an 18" (45.7 cm) long tube with a soft medium size eartip and a listening cup.

53225 Listening Tube 4.75



E-Scope®

Cardionics E-Scope® II, Electronic Stethoscope

E-Scope® II Electronic Stethoscopes use state of the art technology to amplify heart, breath, and Korotkoff sounds. Murmurs and faint blows are much more distinguishable when compared to acoustic stethoscopes. All controls are located on the chest piece.

Includes:

- (2) Types of eartips (Clinical Model)
- Specialist adult diaphragm, stainless steel
- Soft pack carrying case
- Operator's manual
- Accessory Kit
- (1) AAA Battery
- Lifetime Technical Support
- 1 Year warranty: parts and labor



53142 Clinical Stethoscope \$325.00

E-Scope® II, Electronic Stethoscope Hearing Impaired Model

E-Scope® II Electronic Stethoscopes use state of the art technology to amplify heart, breath, and Korotkoff sounds. Murmurs and faint blows are much more distinguishable when compared to acoustic stethoscopes. All controls are located on the chest piece.

Includes:

- Specialist adult diaphragm, stainless steel
- Specialist adult bell, stainless steel
- Soft-pack carrying case
- Belt clip, swivels for ease-of-use
- Operator's manual
- Accessory Kit
- (1) AAA battery
- Lifetime Technical Support
- 1 Year warranty: parts and labor



53143 Hearing Impaired Stethoscope \$389.00

E-Scope® Dual Audio Cable

This cable connects the E-Scope to a hearing instrument audio shoe. The audio shoe is a plastic snap-on fitting, designed to direct the sound through the hearing instrument, and is available from the hearing instrument manufacturer or hearing healthcare professionals.



53147 Audio Cable \$65.00

E-Scope® Patch Cord

The Patch cord connects the Cardionics E-Scope to the personal audio cable of the cochlear implant or streaming device.

53150	Patch Cord	\$30.00
-------	------------	---------



E-Scope® Large Headphones

Large headphones with 3.5mm connector are intended for use with BTE-type hearing instruments, and can also be worn with ITE or "open fit" type BTE hearing instruments.

53152	Large Headphones	\$90.00
-------	------------------	---------



E-Scope® Traditional Headphones

The Traditional Headphones have a separate volume control and are worn over-the-head. These headphones can be worn with ITE or "open fit" type BTE hearing aids.

53155	Traditional Headphones	\$53.00
-------	------------------------	---------



E-Scope® Convertible Headphones

The Convertible Headphones will convert from over-the-head to behind-the-head by moving the headphone earpieces. These headphones can be worn with ITE or "open fit" type BTE hearing aids.

53156	Convertible Headphones	\$60.00
-------	------------------------	---------



Ear Tips

Grason™ OAE-Tymp Eartip

The unique design of the Grason OAE-Tymp Eartip allows a single eartip to fit a range of ears. Latex free.

30500A	7-11 mm, 100pk	\$22.50
30500B	12-16 mm, 100pk	\$22.50
30502	3-6 mm, 100pk	\$22.50



30500A



30500B



30502

Grason™ Single-Use Eartips Starter Pack

Economical and easy-to-use eartips eliminate the need for cleaning and sterilization of testing equipment and provide a positive seal. Contains 300 tips.

30501	300 Tips - 12 sizes, 25 each	\$80.50
-------	------------------------------	---------



Grason™ Single-Use Eartip Refill Packs

30504	4 mm, 100pk	\$19.75
30505	5 mm, 100pk	\$19.75
30506	6 mm, 100pk	\$19.75
30507	7 mm, 100pk	\$19.75
30508	8 mm, 100pk	\$19.75
30509	9 mm, 100pk	\$19.75
30510	10 mm, 100pk	\$19.75
30511	11 mm, 100pk	\$19.75
30512	12 mm, 100pk	\$19.75
30513	13 mm, 100pk	\$19.75
30514	14 mm, 100pk	\$19.75
30516	16 mm, 100pk	\$19.75



30504



30505



30506



30507



30508



30509



30510



30511



30512



30513



30514



30516

Grason™ Single-Use Pediatric Eartip Set

Disposable pediatric eartips designed to fit ECP, SNS and SGS probes. Latex free.

Includes:
50 of each of the following sizes: 3.5 mm, 4.5 mm, 5.5 mm, and 6.5 mm.

30503	300 Tips - 4 sizes, 50 each	\$45.00
-------	-----------------------------	---------



Grason™ Single-Use Pediatric Eartip Refills

Disposable pediatric eartips designed to fit ECP, SNS and SGS probes.

30517	3.5 mm, 100pk	\$19.75
30518	4.5 mm, 100pk	\$19.75
30519	5.5 mm, 100pk	\$19.75
30520	6.5 mm, 100pk	\$19.75



Test Tips

Disposable test tip used for hearing instrument evaluation. Allows patients to hear the benefits of signal amplification.

30601	Size 1, 1pc	\$1.75
30602	Size 2, 1pc	\$1.75
30603	Size 3, 1pc	\$1.75
30604	Size 4, 1pc	\$1.75
30605	Size 5, 1pc	\$1.75
30606	Test Tips, 5pk	\$7.00



Generic BTE Temporary Mold Single-Use Replacement Eartips

Refill eartip for the ER-13E Temporary BTE Kit - ER 13-14, 3 mm horn, yellow tip / green stem.


30625	Single-Use Eartips - 36pk	\$22.50
-------	---------------------------	---------



E-A-RLink® Eartips

Refill eartips for the ER-13E Temporary BTE Kit.

30622	3A, Adult - 50pk	\$22.50
30623	3B, Small - 50pk	\$22.50
30624	3C, Jumbo - 24pk	\$18.50
30626	14D2, Infant - 20pk	\$22.50
30627	3A, Adult - 750pk	\$220.00
30628	3B, Small - 750pk	\$220.00



Comply™ Pro Fitting System

Contains everything needed to start using Comply Canal Tips in your practice.

Includes:

- 10 Standard canal tips, 5 sizes (two of each size)
- 10 Slim canal tips, (two of each size)
- 10 Sound tube adapters

30521	Fitting System	\$62.00
-------	----------------	---------



Comply™ Canal Tubes

Tubes for use with any size Comply Canal Tips.

30561	Canal Tubes, 15pk	\$34.00
30562	Sound Tubes, 30pk	\$67.00



Comply™ Canal Tips

Twelve-packs of Standard, Short, and Slim Comply Canal Tips

30531	Standard, 0 Vents - 12pk	\$19.25
30532	Standard, 1 Vent - 12pk	\$19.25
30533	Standard, 2 Vents - 12pk	\$19.25
30534	Standard, 3 Vents - 12pk	\$19.25
30535	Standard, 4 Vents - 12pk	\$19.25



Comply™ Canal Tips, Slim

Twelve-packs of Slim Comply Canal Tips.

30541	Slim, 0 Vents, Black - 12pk	\$19.25
30542	Slim, 1 Vent - 12pk	\$19.25
30543	Slim, 2 Vents - 12pk	\$19.25
30544	Slim, 3 Vents - 12pk	\$19.25
30545	Slim, 4 Vents - 12pk	\$19.25



Vented Snap Tip

Disposable test tip used for hearing instrument evaluation. Allows patients to hear the benefits of signal amplification.

30596	9 mm, Medium - 10pk	\$14.00
-------	---------------------	---------



E-A-RTemp™ Temporary Eartip

Soft foam eartips can be used for hearing instrument evaluation, as temporary earmolds for patients, or as inexpensive disposable demonstration tips.

30609	#13A standard with tube, 25pk	\$41.50
-------	-------------------------------	---------



Supra-Aural Earphone Cushions

These audiometric headset cushions meet the requirements specified in ANSI S3.6, ICE 645, and ISO 389-1.

30642	Blue, 1pc	\$7.50
30643	Red, 1pc	\$7.50
30646	Red/Blue Set, 2pk	\$14.50



Pediatric Supra-Aural Earphone Cushion

Pediatric Supra-Aural Earphone Cushion provides compatibility between smaller ears and the audiometric headset. The reduced diameter helps eliminate contact with the skin covering the cranial bone and helps give better test results, as outlined in ANSI S3.6 Section 9.1. Calibration remains the same as a standard earphone cushion.

30651	Black, 1pc	\$8.75
30652	Blue, 1pc	\$8.75
30653	Red, 1pc	\$8.75
30656	Red/Blue Set, 2pk	\$17.50



Bone Conduction Headband

Holds bone oscillator in place during bone conduction testing. Adjustable Velcro fasteners.

54204	Headband, 24in	\$36.00
-------	----------------	---------



Universal Hearing Screener

Designed to maintain consistent output in all frequencies, with 40 dB presented at 500, 1000, 2000 and 4000 Hz. The Universal Hearing Screener is ideal for quick screenings in a variety of satellite locations such as primary care doctor offices, senior centers and nursing homes. Operates on two AAA batteries (included), with automatic shutoff when batteries are depleted. Perfect for your field bag or order them with custom labels for distribution to referring physicians.

30679	Universal Hearing Screener	\$19.95
-------	----------------------------	---------



TRUE-FIT™ FOAM TIPS

Replaceable True-Foam ear-tips are so comfortable that they can be worn for hours and hours. Westone's ear-tips enhance the performance in-ear earphones by ensuring superior noise isolation, and a superb fit.

5 pair / 10 pack \$7.50

62800	Orange, 15 mm
62801	Black, 12.6 mm
62802	Red, 12.6 mm
62803	Blue, 11 mm
62804	Green, 11 mm
62821	Combo (1 pair each size)



62800

Length
15.5 mm

Diameter
15 mm



62802

Length
15.4 mm

Diameter
12.6 mm



62803

Length
14.9 mm

Diameter
11 mm



62801

Length
11 mm

Diameter
12.6 mm



62804

Length
11 mm

Diameter
11 mm



62821

1 TRUE-FIT FOAM


Smooth eartip surface increases comfort and longevity. Plush medium recovery foam ensures easy insertion and maximum seal.

2 ACOUSTIC ISOLATION

Color coded smooth bore stem with direct-coupled exit port delivers accurate sound while flexible design contours to your ear for superior noise isolation.

3 EARVAULT™ PROFILE

Extended five-tip size range was developed by profiling over 15 million ears and our EarVault database that contains anatomical 3D scans of the human ear. True-fit foam tips allow an individual fit for excellent retention, comfort, and sound isolation.



TRUE-FIT™ + STAR™ FIT KIT

10 pair / 20 pack

\$7.50


62810	Fit Kit (1 pair each size)
-------	----------------------------

STAR™ SILICONE TIPS


Westone's replaceable Star Silicone eartips are available in 5 different sizes ensuring superior noise isolation and a superb fit. Westone ear-tips fit a variety of popular in-ear monitors, which allows you to enjoy the best ear-tips on the market with your current IEMs.

5 pair / 10 pack \$7.50


62805	Orange, 13.5 mm
62806	Red, 12 mm
62807	Black, 12 mm
62808	Blue, 10.5 mm
62809	Green, 10.5 mm
62822	Combo (1 pair each size)




62805




62806



62808



62807



62809

1 PATENTED STAR DESIGN



Multiple STAR flex zones allow the tip to conform to your individual ear anatomy. High-purity silicone construction delivers durability and biocompatibility.

2 ACOUSTIC ISOLATION

Color coded smooth bore stem with direct-coupled exit port delivers accurate sound while flexible design contours to your ear for superior noise isolation.

3 EARVAULT™ PROFILE

Extended five-tip size range was developed by profiling over 15 million ears and our EarVault database that contains anatomical 3D scans of the human ear. Star silicone tips allow an individual fit for excellent retention, comfort, and sound isolation.



62810

Screw End Cords

Screw-in Cord with Handles

Handy cords with screw-in anchors attach to earmolds and earplugs. The 14" long cord is useful for hanging around the neck when earmolds/earplugs are not in use.

The unique screw-in anchor also serves as a removal/insertion handle.



42291	Cord with Red & Blue Handles - 1pc	\$2.95
42293	Handles only, 2pk	\$2.10

Nylon Cord with Red and Blue Screw Ends

Nylon cord with red and blue anchors will screw into soft vinyl and silicone earpieces. Length: 26" (66 cm).

Installation tool sold separately.



42346	Black Cord - 1pc	\$2.35
42346-BL	Blue Cord - 1pc	\$2.35
42346-GR	Green Cord - 1pc	\$2.35
42346-PK	Pink Cord - 1pc	\$2.35
42346-PU	Purple Cord - 1pc	\$2.35
42346-RD	Red Cord - 1pc	\$2.35



Screw Cord Installation Tool

Installation tool for mounting Westone screw cords and anchors in earpieces and earplugs.

42350	Tool - 1pc	\$5.50
-------	------------	--------



Nylon Cord with Clear Screw Ends

Nylon cord with clear anchors, screw into silicone and vinyl earpiece. Length: 26" (66 cm).

Installation tool sold separately.

42349	Black Cord - 1pc	\$2.35
42349-BL	Blue Cord - 1pc	\$2.35
42349-O	Orange Cord - 1pc	\$2.35
42349-PK	Pink Cord - 1pc	\$2.35
42349-PU	Purple Cord - 1pc	\$2.35
42349-R	Red Cord - 1pc	\$2.35
42349-YL	Yellow Cord - 1pc	\$2.35



Clear Plastic Cord with Clear Screw Ends

Clear plastic cord helps keep earmolds and earplugs from being lost. Clear anchors screw into soft vinyl and silicone earmolds and earplugs. Length: 26" (66 cm).

Installation tool sold separately.

42348	Clear Cord - 1pc	\$2.35
-------	------------------	--------



Clear Plastic Cord with Red & Blue Screw Ends

This clear, plastic 26" cord helps keep earmolds from being lost. The cord is designed to screw into both soft vinyl and silicone earmolds. Unique screw-in anchors are colored red to indicate the right side and blue to indicate the left.

42345	Red & Blue - 1pc	\$2.35
-------	------------------	--------



Replacement Cord Screw Ends

Replacement red (right) and blue (left) ends for Westone screw cords. Two per package.

42347	Replacement Screw Ends - 2pk	\$1.25
-------	------------------------------	--------



Pop End Cords

Replacement Nylon Cord with Pop Ends

Nylon "pop" cord replacement for earpieces manufactured with preset anchors. Length: 26" (66 cm) Features right angle mushroom anchors.

42342	Black Cord - 1pc	\$2.35
42342-BL	Blue Cord - 1pc	\$2.35
42342-GR	Green Cord - 1pc	\$2.35
42342-OR	Orange Cord - 1pc	\$2.35
42342-PK	Pink Cord - 1pc	\$2.35
42342-PU	Purple Cord - 1pc	\$2.35
42342-R	Red Cord - 1pc	\$2.35
42342-YL	Yellow Cord - 1pc	\$2.35
42343	Black Cord - 10pk	\$15.25
42330	Black Cord - 100pk	\$150.00



Clear Plastic Cord with Pop Ends

Plastic "pop" cord replaces lost or damaged cord used on AquaNot swim plugs or other earmolds manufactured with preset anchors.

Length: 26" (66 cm)
Features clear mushroom anchors.

42340	Clear Cord - 10pk	\$15.25
42341	Clear Cord - 1pc	\$2.35



Black Nylon Bolo Cord with Pop Ends

Nylon bolo cords feature a bolo cinch and right-angle mushroom anchors. Length: 24" (61 cm)


42344	Black Cord - 1pc	\$2.75
-------	------------------	--------



Metal Detectable Plastic Cord

Replacement blue metal plastic "pop" cord for earpieces manufactured with preset anchors. Used in the food processing industry for its detection capability. Length: 29" (74 cm)

20150	Blue, 1pc	\$2.35
20154	Blue, 10pk	\$21.50



Replacement Angled Pop Ends

Clear angled end for Westone pop cord.

42357	Angled Pop End - 1pc	\$0.50
-------	----------------------	--------



ADCO PRODUCTS

ADCO ITE and CIC Modication Kit

Everything you need to adjust the fit of ITE and CIC hearing instruments.

- Includes:
- Soft ADCO Powder
 - ADCO Liquid
 - ADCO Sheen
 - ADCO Sheen Thinner
 - Hearing Instrument Clamp
 - Mixing supplies
 - Detailed instructions.

42495	Modification Kit	\$58.50
-------	------------------	---------



Hard ADCO Addon

Brush on the clear liquid Hard ADCO Addon for minor buildup of hard earmolds and hearing instrument shells. Use with Hard ADCO Addon Thinner.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

42465	.5oz Bottle	\$6.75
42467	Thinner, .5oz Bottle	\$3.00



42465



42467

Soft ADCO Addon

Brush on Soft ADCO Addon for minor buildup of soft earpiece surfaces. The clear liquid leaves a flexible finish. Use with Soft ADCO Addon Thinner.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

42461	.5oz Bottle	\$6.75
42463	Thinner, .5oz Bottle	\$3.00



42461



42463

ADCO Build

ADCO Build is used for major buildup of hard earpieces and hearing instruments. We recommend applying ADCO Sheen Sealer to the modified area of the earpiece/hearing instrument.

42471	.5oz Liquid / 1oz Powder	\$24.00
42473	.5oz Liquid only	\$11.75
42475	1oz Liquid / 2oz Powder	\$40.50
42478	2oz Powder only	\$18.50



42471



42473



42475



42478

NOTE: Eartips and other products used in the ear are non-refundable due to hygiene/health concerns.

ADCO Sheen

Brush on the clear liquid Hard ADCO Addon for minor buildup of hard earmolds and hearing instrument shells. Use with Hard ADCO Addon Thinner. THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

42481	.5oz Bottle	\$5.75
42483	Thinner - .5oz Bottle	\$3.00



42481



42483

ADCO Economy Ki

Includes all the components needed to build up soft and hard earmolds or hearing instruments.

- Includes:
- ADCO Build Powder
 - ADCO Build Liquid
 - Hard ADCO Addon
 - Hard ADCO Addon Thinner

- Soft ADCO Addon
- Soft ADCO Addon Thinner
- ADCO Sheen
- Mixing supplies and instructions

42491	Economy Kit	\$43.50
-------	-------------	---------



Bionox®

Bionix® Lighted Placement Tool

With its bright light, depth guide, and ergonomically designed handle, the Lighted Placement Tool allows audiologists to safely and accurately place Oto-Dams in the ear.

- Includes:
- 15 lighted placement tips
 - 1 light source
 - 1 magnification lens

24110	Tool, 15pk	\$49.50
-------	------------	---------



Curettes

Bionix® Lighted Ear Curette with Magnification

The only lighted curette that articulates up to 90 degrees to assist in the removal of cerumen and foreign bodies from the ear and the ear canal. 3 mm wide tip.



Articulating

42518	Articulating - 50pk	\$95.50
-------	---------------------	---------

CeraPik

42519	CeraPik - 50pk	\$79.85
42532	CeraPik - 200pk	\$310.00

InfantScoop

42520	InfantScoop - 50pk	\$79.85
42530	InfantScoop - 200pk	\$300.00

AngleLoop

42521	AngleLoop - 50pk	\$79.85
42526	AngleLoop - 200pk	\$300.00

MicroLoop

42522	MicroLoop - 50pk	\$79.85
42527	MicroLoop - 200pk	\$300.00

WaveCurette

42523	WaveCurette - 50pk	\$79.95
42528	WaveCurette - 200pk	\$300.00

FlexLoop

42524	FlexLoop - 50pk	\$79.85
42529	FlexLoop - 200pk	\$300.00

Variety Packs

42525	Variety - 50pk	\$79.85
42531	Variety - 200pk	\$300.00

50 Variety Pack Includes:
1 Light Source and 1 Magnification Lens

200 Variety Pack Includes:
4 Light Sources and 4 Magnification Lens

Accessories

42534	Light Source Handle - 1pc	\$41.50
42535	Magnification Lens - 5pk	\$17.25
42505	Lighted Forceps 10pk	\$97.00



Bionix® Safe Ear Curettes™

The only curette that articulates up to 90 degrees to assist in the removal of cerumen and other foreign bodies residing within the ear. 3 mm wide tip.

50 Count Box

42510	Articulating, Purple	\$54.50
42511	AngleLoop, Red	\$46.00
42512	CeraSpoon, Yellow	\$46.00
42513	FlexLoop, White	\$46.00
42514	VersaLoop, Purple	\$46.00
42515	InfantScoop, Blue	\$46.00
42516	MicroLoop, Green	\$46.00
42517	ControLoop, Orange	\$46.00



Wash Kits

Bionix® OtoClear® Spray Wash Kit

This entry-level kit gives professionals an efficient and economic way to provide the highest level of care.

Includes:

- 20 OtoClear Tips
- 1 Spray Wash Bottle
- 1 Bionix® Ear Basin.

42546	Spray Wash Kit	\$66.00
42547	Irrigation Basins - 3pk	\$20.50



Portable WaterPik®

This convenient hand-held tool is perfect for irrigating patients' ear canals when used with OtoClear Tips. OtoClear Tips are sold separately.

Includes:

- Portable WaterPik Unit
- One-touch, high and low settings
- Automated, pulsed lavage
- 210ml reservoir
- 3 Adapter wands
- 2 year warranty

42543	Portable WaterPik	\$72.50
42541	Irrigation Tips, 40pk	\$46.50

Bionix® OtoClear® Ear Lavage System

A completely portable irrigation delivery device that can be used for procedures anywhere, in or out of the office. Perfect for high-use facilities.

Features:

- Portable WaterPik Unit
- Automated, pulsed lavage
- One-touch, high and low settings
- 210 ml self-contained water reservoir

Includes:

- 40 OtoClear Tips
- 3 Adapter wands
- 1 Bionix ear basin
- 2 year warranty

42545	Ear Lavage System	\$122.50
-------	-------------------	----------



Luer Lock Syringe

The 60cc Luer Lock Monoject Syringe is designed for the irrigation of patients' ear canals. The Luer Lock tip secures OtoClear tips safely in place.

42542	Syringe	\$3.05
-------	---------	--------



Bionix® OtoClear® Gentle Touch™ Ear Irrigation Tips

- A soft tip that fits snugly and comfortably into the ear
- Improves patient comfort and safety, virtually eliminating the risk of rupture
- Eliminates over-insertion into ear canal due to its flared tip design
- Directs water away from the tympanic membrane at a 30° angle
- Exit portals to prohibit pressure buildup and direct drainage

42541	Irrigation Tips, 40pk	\$46.50
-------	-----------------------	---------



Accessories

Mirror Go Round

Mirror Go Round allows patients to easily view their ears to see how their hearing instruments and earpieces should look when properly placed in the ear. Folds for easy transport and storage.

41157	Mirror Go Round	\$11.35
-------	-----------------	---------



Small Plastic Bags

Handy for the storage of small parts.

42325	2.5" x 4" Bags, 100pk	\$4.25
-------	-----------------------	--------



Zipper Locking Bags

Ideal for the presentation or storage of earmolds and other small parts. Bags are available in clear or with a white imprint area for writing. The Zipper Locking seal protects contents from dust.

42321	3" x 5" Bags, 25pk	\$3.15
42322	3" x 5" Bags with Imprint Area, 25pk	\$3.35
42324	4" X 6" Bags, 25pk	\$3.75



Polycon Plastic Carrying Case

Durable and compact Polycons are great for storing earpieces and hearing protection devices. Select either the standard Polycon or one modified to be compatible with AquaNot swim plugs (features a hole in the lid to aid with moisture evaporation). Polycons are sold individually, please specify color.

Note - AquaNot model is available in neon colors only.

42334	AquaNot Polycon Case w/ Chain	\$0.65
42335	Custom Imprint	\$0.85
42336-BL	Neon Blue	\$0.65
42336-GR	Neon Green	\$0.65
42336-OR	Neon Orange	\$0.65
42336-YL	Neon Yellow	\$0.65



42334



MIMUMUM ORDER OF 1,000 REQUIRED



42336-BL

42336-GR


42336-OR

42336-YL

P-42 Shipping Box

These sturdy boxes feature wire closure clips and are great for shipping impressions and hearing instruments or for the storage of small items. Boxes measure 2"W x 2"H x 3.5"L

42351	Box - 1pc	\$0.95
-------	-----------	--------



P-109 Shipping Box

P-109 boxes feature wire closure clips and are ideal for shipping impressions and hearing instruments or for the storage of small items. Boxes measure 3"W x 2"H x 4"L.

42352	Box - 1pc	\$1.05
-------	-----------	--------



Unsurpassed offering of first-class ear infection control products

With our customers' health always at the forefront of our mind, Westone offers the very best cleaners, sanitizers, disinfectants and hearing instrument care devices at a reasonable price to ensure your hearing aids always remain safe for continued use.



Equipment Cleaners

Branson Ultrasonic Cleaner

Ultrasonic energy (40 kHz) creates millions of bubbles in the cleaning solution that work their way into the small crevices of clinical equipment and earmolds to loosen dirt or debris.

Features:

- 15oz. Steel tank
- Ultra Sonic Energy: 40kHz
- Weight: 3lbs.
- Dimensions: 8.75" x 4.5" x 5"
- 5 Minute timer



40151 Ultrasonic Cleaner \$115.00

Sparkle Spa Ultrasonic Cleaner

Ultrasonic waves clean quickly and easily without the need for harsh chemicals, scrubbing, or polishing. The Sparkle Spa is great for removing dirt and grime from small clinical tools and earpieces.

Features:

- Weight: 3 lbs.
- Stainless steel tank: 1 pt. (600 ml) capacity
- Removable basket with bracelet/watch holder
- 3 programmable cleaning cycles: 3, 5 and 10 minutes
- Ultrasonic energy: 50/60 Hz
- 1 Year warranty



40153 Ultrasonic Cleaner \$49.50

GemOro Ultrasonic Cleaner

Attractive ultrasonic cleaner with an advanced digital control system, constructed in durable lightweight ABS plastic.

Features:

- Weight: 5 lbs.
- Stainless steel tank: 1.2 qt. capacity
- Removable parts basket
- Preset cycle times: 90, 180, 280, 380, or 480 seconds
- Ultrasonic energy: 50/60 Hz



40154 Ultrasonic Cleaner \$75.00

Tap and Slide® Sterilization Tray

This lift-out tray allows easy removal of items being sterilized. 31 day slide index reminds you when sterilant liquid needs to be changed. Dimensions: 10" (254 mm) L, 5" (127 mm) W, 3.5" (89 mm) D.

40171	Sterilization Tray	\$33.50
-------	--------------------	---------



Audiologist's Choice® Ultrasonic Disinfectant/Cleaner

Great in hearing healthcare dispensing offices, this product will properly disinfect probe tips, earmolds, specula, etc. An 8 oz. bottle of concentrate makes 2.67 gallons of disinfectant/cleaner.

40161	Disinfectant/Cleaner, 8oz	\$9.00
-------	---------------------------	--------



Oto-Cleans®

Designed to clean, disinfect, and deodorize a wide range of custom-fit earpieces.

54019	5ml Bottle, 3pk	\$3.95
-------	-----------------	--------



Audiologist's Choice® Earmold and ITE Disinfectant/Cleaner

Specially formulated to clean, disinfect and deodorize earpieces and ITE hearing instruments.

40141	Bottle, 2oz	\$3.25
40142	Bottle, 4oz	\$4.50



Wavicide-01

This 2% glutaraldehyde sterilant/disinfectant is premixed and requires no activator. 1 year shelf life.

40181	1 Quart	\$6.00
40182	1 Gallon	\$22.00

40181



40182

Wavicide-01 Monitor Strips

Monitor strips for testing the potency/life of Wavicide-01 sterilant/disinfectant.

40183	Monitor Strips, 50pk	\$71.00
-------	----------------------	---------



Envirocide®

This ready-to-use and versatile disinfectant/decontaminant kills TB, HIV-1, and many other pathogens. The biodegradable Envirocide cleans and disinfects without glutaraldehyde, phenols, or bleach and is safe for use on almost anything. Effective for ultrasonic use and as a pre-soak cleaner prior to sterilization.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US.

40184	24oz Spray Bottle	\$12.75
40185	1 Gallon	\$32.75

40184



40185



Sporox® II

A powerful hydrogen-peroxide based solution for the sterilization and high-level disinfection of heat-sensitive instruments (virucidal, sporicidal, tuberculocidal, bactericidal and fungicidal). Safe for instruments containing soft metals such as brass and copper

40188	1 Gallon	\$38.50
-------	----------	---------



Audiologist's Choice® Audio Wipes®

Contains a hospital-grade disinfectant that can be used to clean virtually any surface. Contain no alcohol, which can cause dryness in rubber and plastics.

40191	Pop-Up Canister, 160pk	\$7.55
40192	Mini Canister, 36pk	\$4.75
40193	Wipes Singles, 100pk	\$15.00



Covidien Webcol Alcohol Prep Pads

Large pads saturated with 70% isopropyl alcohol are ideal for use in the hearing healthcare clinic. Each pad measures 1.7" x 3.5" and is individually wrapped.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

40212	Alcohol Prep Pads, 200ct box	\$4.90
-------	------------------------------	--------



Super Sani-Cloth®

High Alcohol Formula (55%) for the disinfection of hard, non porous surfaces and equipment.

Features:

- Bactericidal, tuberculocidal, and virucidal
- Tested effective against 26 microorganisms
- Kills TB and RSV in 1 minute
- HBV, HCV, human EPA registered
- Meets CDC and OSHA guidelines
- Kills Coronavirus, Rotavirus, MRS \A, VRE and ESBL producing E. coli in 2 minutes

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

40195	Cloth, 160pk	\$8.50
-------	--------------	--------



Sani-Cloth® Plus

Low Alcohol Formula (14.85%) for the disinfection of hard, non porous surfaces and equipment.

Features:

- Bactericidal, tuberculocidal, and virucidal
- Kills RSV in 1 minute
- Kills HBV, HCV, VRE and MRSA in 3 minutes;
- Kills TB in 5 minutes
- Tested effective against 15 microorganisms
- EPA registered
- Meets CDC and OSHA guidelines

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

40196	Cloth, 160pk	\$9.50
-------	--------------	--------



Purell® Instant Hand Sanitizer

Purell Instant Hand Sanitizer is non-toxic, dye free, dermatologist tested and kills 99.9% of most common disease-causing germs.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US.

40221	Bottle, 4.25oz	\$3.50
40225	Wall Dispenser, 800ml	\$24.95
40226	Refill, 800ml	\$10.95



Purell® Instant Hand Sanitizer with Aloe

The same great product as the Instant Hand Sanitizer but enhanced with the soothing effects of Aloe.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US.

40231	Bottle, 4.25oz	\$3.50
40232	Pump, 12oz	\$8.75



Micrell® Antibacterial Lotion Soap

Contains clinically-proven chloroxylenol (PCMX), a fast-acting broad spectrum anti-microbial agent. The rich non-irritating lather cleans hands completely.

40242	Pump, 8oz	\$4.50
40245	Wall Dispenser, 800ml	\$24.95
40246	Refill, 800ml	\$8.95
40247	Refill, 1 gallon	\$15.00



GrafcO Alcohol Dispensing Pump

Refillable, alcohol dispensing bottle is ideal for use in the hearing healthcare clinic.

40251	Dispensing Bottle	\$8.95
-------	-------------------	--------



6" Tongue Depressor

6" hard wood tongue depressors for oral use, spreading ointments, mixing medications or for use as small splints.

40252	100 box	\$5.75
-------	---------	--------



AcoustiCover™

AcoustiCovers (Class II Medical Device) offer protection against cross-contamination and have passed Viral Penetration ASTM Method F 1671. The covers were found to be Acoustically Transparent by the National Institute of Standards and Technology (NIST). Designed to fit most headphones used with testing equipment.

40111	Cover, 200pk	\$33.25
40112	Cover, 500pk	\$78.00



Safe Sound Headphone Covers

An acoustically transparent disposable headphone cover for use with hearing test equipment. Available in packages of 100.

40121	100pk	\$15.00
-------	-------	---------



Procedure Masks

One-size-fits-all medical grade procedure masks aid in the control of airborne infectants.

40261	Mask, 50ct Box	\$6.98
-------	----------------	--------



TechSpray® Contact Cleaner

All purpose G3 Techspray removes oil, grease, silicone, dirt and grime from electrical contacts. Non-ozone depleting, non-flammable, and safe on most plastics.

42262	TechSpray Contact Cleaner 16oz	\$37.95
-------	--------------------------------	---------



Pro Towel

Great for use in the office or clinic, these towels feature two plies of absorbent paper on the top to soak up spills. One ply of poly on the back helps prevent spills from leaking through.

42271	White, 500ct case	\$22.50
42272	Blue, 500ct case	\$22.50



Cotton Swabs

Cotton-tipped swabs are perfect for a multitude of uses in the hearing healthcare clinic.

42281	Swabs, 100pk	\$5.15
-------	--------------	--------



Earpiece Cleaners

Convenient earpiece cleaners are the perfect size for removing cerumen and other debris from earpieces and hearing instruments.

42285	Cleaners - 50pk	\$3.25
-------	-----------------	--------



OtoCare® Kit

The OtoCare Kit provides everything required to maintain hearing instruments, all in a convenient vinyl organizer bag.

Includes:

- Standard Hearing Aid Saver
- .5 oz. bottle of Oto-Ease Lubricant
- 2 oz. bottle of Oto-Cleans
- Air blower
- Battery tester
- Telephone pad
- Wax removal tool



50111	Kit	\$20.50
50112	Organizer Bag - No Imprint	\$2.50

Nanoclean® Fine-Instrument Cleaners

Cleans moisture and debris from hearing instrument tubes, vents and sound bores - even thin tubes!

A "thread and pull" cleaning action, with the spongy nylon filament brush, wipes the surface clean and absorbs moisture. Non-metallic construction for safe use without scratching.

Total length: 10", with a 5" filament brush.

50148	20pk	\$3.75
50149	20pk / 20 boxes	\$64.95



Cleaning Gloves

Vinyl Gloves

Powder and latex-free vinyl gloves are textured for a superior grip to form an effective biological barrier.
Brand of gloves are subject to change without notice

40271	Small, 6-6.5", 100ct box	\$7.50
40272	Medium, 7-7.5", 100ct box	\$7.50
40273	Large, 8-8.5", 100ct box	\$7.50



Latex Gloves

Single-use, powder-free latex gloves are textured for a superior grip and form an effective biological barrier.
Brand of gloves are subject to change without notice

40281	Small, 100ct box	\$7.50
40282	Medium, 100ct box	\$7.50
40283	Large, 100ct box	\$7.50



Nitrile Gloves

Single-use, powder-free gloves are textured for a superior grip and form an effective biological barrier.
Brand of gloves are subject to change without notice

40286	Small, 200ct box	\$14.50
40287	Medium, 200ct box	\$14.50
40288	Large, 200ct box	\$14.50



Cerumen Removal

Cerumen Removal Tool with Brush and Magnet

Cerumen Removal Tool features a loop on one end, a brush on the other and a magnet in the handle to help remove and install hearing instrument batteries.

50121	Tool, 1pc	\$1.35
-------	-----------	--------



Cerumen Removal Tool with Brush

Cerumen Removal Tool features a wax loop on one end and a brush on the other end.

50123	Tool, 1pc	\$1.00
50124	Tool, 100pk	\$75.00



Cerumen Removal Tool

Small tool for removing built-up cerumen from sound bores and vents in earpieces and hearing instruments.

50125	Tool, 1pc	\$0.55
-------	-----------	--------



Vent Cleaning Tool

Ideal for removing hardened cerumen from tiny vent holes in earpieces and hearing instruments. The flexible tip is .025" in diameter and 1" long.

50127	Tool, 1pc	\$1.00
-------	-----------	--------



Air Blower

The Air Blower removes moisture and dirt from hearing instruments with a gentle puff of air.

50571	Air Blower, 1pc	\$3.85
50572	Display Box, 12pk	\$42.00



AcuLife® Audio-Kit™ Cleaner Tool

The Audio-Kit contains the five tools needed for fast, easy, and effective daily hearing instrument cleaning and maintenance. The tools are built into a lightweight, comfortable handle that feature a built-in battery storage compartment.

Features:

- Wax removal brush
- Wax removal pick
- Battery replacement magnet
- Tube and vent cleaner
- Battery door opener

50126	Cleaner, 1pc	\$5.95
-------	--------------	--------



Gold Wax Catcher

- Blocks cerumen before it can enter the sound bore or vents
- 302 Stainless steel wire, coated in 24kt gold
- Non-corrosive and non-allergenic
- Easy to Install

Short - .067" x .087" x .051" (1.7 mm x 2.2 mm x 1.2 mm)
 Standard - .078" x .066" x .031" (1.9 mm x 1.6 mm x .7 mm)
 Tall - .083" x .062" x .045" (2.1 mm x 1.5 mm x 1.1 mm)

50116	Tall, 1pc	\$0.35
50117	Standard, 1pc	\$0.35
50118	Short, 1pc	\$0.35



Odinell Ear Spray

Odinell Ear Spray eliminates excessive earwax from the ear canal, which helps to maintain the health of your ears and restore any hearing loss caused by excessive wax buildup. This non-aggressive and isotonic cleaning solution contains Chamomile extract, which calms and soothes the irritation of the ear canal and the ear due to daily use of hearing aids.

54059	Ear Spray, 50ml	\$5.60
-------	-----------------	--------



Dryers

PerfectDry Lux®

54054 PerfectDry Lux® \$46.00

PerfectDry Lux® UV-C

UV-C disinfection and quick drying of hearing systems

Performance, Quality, Design

PerfectDry Lux® has been designed for optimal fast drying with a real and effective disinfection of hearing systems.

The quality, the design and the size of PerfectDry Lux® make it a perfect setting for hearing systems.

- 30 minute drying cycle (standard)
- Ability to select intensive drying cycle (1.5 hours, by pressing the start button twice)
- 360° UV-C disinfection lamp destroys 99% of micro-organisms
- Forced-air fan drying system
- Quality and sturdiness
- Optimized ergonomics
- Low energy consumption
- Compact case and design
- USB or power supply (UL)
- External Dimensions: 92 mm X 48 mm X 75 mm



Forced-air fan drying system

- Temperature self-regulated at 113 °F +/-3 with a chip
- Controlled ventilation
- Optimized and homogeneous air flow
- Extraction of moisture



360° UV-C Lamp

- Shortwaves disinfection (253.7nm)
- Peripheral light 360° (2 Watts)
- Eliminates over 99% of viruses and bacteria (proven certified laboratory)
- Confirmation test on no damages on hearing aids surfaces (nano coating and varnishes)

PerfectClean®

54055 PerfectClean® \$79.00

PerfectClean® (BTE/ITE)

The First system to clean, dry and disinfect hearing aids

PerfectClean® is a CE marked and patented innovation specially developed for the safe and effective earwax removal, cleaning and drying of all types of hearing systems.

Simply place your hearing aids in the PerfectClean® and your aids will be totally cleaned, disinfected and dried. You will be able to fully enjoy the performance of your hearing system.

PerfectClean® works with a refill (about 45 cycles); an electronic safe antibacterial solution that disinfects and removes ear wax. The effect is stressed with the Vortex Power.



Benefits:

- Effectively removes earwax from hearing aid devices
- Maintains optimal performance of hearing aid devices
- Helps to prevent irritation
- Easy to Use

Hearing Aid Saver

Hearing Aid Saver is an inexpensive, moisture removal system for hearing instruments and other small moisture-sensitive items.

Includes:

- Resealable container
- Rechargeable desiccant
- Complete instructions for use

Use the standard size for BTE instruments and the mini for ITE, ITC, and CIC instruments.



50552



50551



79227

50551	Standard, 1pc	\$4.75
50552	Mini, 1pc	\$4.75
50561	Standard, Display Box - 12pk	\$50.00
50562	Mini, Display Box - 12pk	\$50.00
79227	Replacement Desiccant Pod	\$2.15\$2.15

Replacement Desiccant Pod

- 3" Diameter
- Used to protect and remove harmful moisture from hearing instruments and other small moisture-sensitive items.
- Desiccant pod can be reactivated in any microwave oven

50512	Replacement Pod, 1pc	\$4.50
-------	----------------------	--------



Green = Time to Reactivate!



AudioDry® Hearing Instrument Drying System

AudioDry Hearing Instrument Drying System is designed to safely and effectively remove moisture from hearing instruments using a gentle heating process. Also great for communication, music, and hearing protection earpieces.

- Manual on/off
- Height: 1.2" (3 cm)
- Width: 3.9" (10 cm)
- Depth: 2.8" (7 cm)
- 1 year warranty

50529	Dryer	\$32.50
-------	-------	---------



Moisture Guard™

Moisture Guard uses gentle, dry heat to remove harmful moisture and cerumen from hearing instruments.

- Manual on/off
- 1 Year warranty

50531	Dryer - Cherry	\$46.50
50533	Dryer - American Black Walnut	\$46.50



Serene® Renew Dryer & Freshener

This compact, lightweight unit is an environmentally friendly, natural, and cost effective way to dry and refresh hearing instruments, without the use of desiccants or chemicals. Renew Dryer & Freshener uses dry heat to remove moisture and condensation; accommodates up to two pairs of hearing instruments or cochlear implant processors.

- Operating cycle: 2 hr, 4 hr, or 8 hr
- UV light inhibits bacteria growth
- Automatic shut off
- Built-in battery tester
- 2.0" (5.1 cm) H x 5.8" (14.7 cm) W x 3.5" (8.9 cm) D
- Weight: 15 oz.
- 1 Year warranty

50532	Dryer	\$51.50
-------	-------	---------



Dry & Store Global II

Portable conditioning system for hearing instruments, combines the three essential elements for effective moisture removal: heat, moving air, and desiccant. The Dry & Store Global II dries cerumen, kills germs, and deodorizes hearing instruments.

Features:

- 8 hr operating cycle; automatic shut off
- Timer resets when lid is opened; short-term conditioning is possible
- UV germicidal lamp
- 3.3" (8.4 cm) H x 7.2" (18.3 cm) W x 4.9" (12.4 cm) D
- Weight: 1lb 3oz
- 1 Year warranty

50544	Dry & Store, 1pc	\$89.00
50548	Replacement UV Lamp	\$8.00



DryCaddy™

DryCaddy utilizes a molecular sieve desiccant material which has proved to great at removing tightly bound moisture from within tiny pores, openings, and tubes of hearing devices as well as other moisture-sensitive devices. Kit provides a year of protection against moisture damage.

50549 Jar and 6pk Desiccants \$18.50



Zephyr by Dry & Store®

The Zephyr uses specially formulated desiccant to remove damaging moisture, dry cerumen, and odor from hearing aid instruments and earpieces.

Features:

- 8 hr operating cycle
- Automatic shut off with on/off button
- Does not require a UV germicidal lamp
- 2.8" (7.1 cm) H x 3.6" (9.1 cm) W x 5.9" (15 cm) D
- Weight: 13.4 oz
- 1 Year warranty

50542 Zephyr, 1pc \$49.95



Dry-Brik® II

Dry-Brik II desiccant is designed for use with Dry & Store products. It works as a molecular sieve, which is a much more aggressive desiccant than silica gel beads. Minimum one year shelf life with the foil intact and stored in a climate-controlled environment. Replace approximately every two months.

50547 Dry-Brik II, 3pk \$7.00



Vacuums

Jodi® XL Vacuum Pump

The Jodi-XL is small, quiet, powerful and efficient. You're sure to appreciate the large silent desiccator! The specially designed needles range in size to efficiently clean hearing instruments, large earmolds, ITE, ITC, CIC, RIC and Thin tubes. This motor carries an industrial rating of 2000 hours of continuous operation.

- Includes:
- 2 Filtered syringe assemblies
 - 3 Needles
 - 2 Reamers
 - Reusable desiccant
 - Desiccant lid and tubing - 6"(15.2 cm) tall
 - Blue coiled tubing extends to a full length of 35"(88.90 cm)
 - Universal power supply: 12-16V, DC 2A
 - 2 pumps
- Height: 3" (7.6 cm)
 - Width: 9" (22.9 cm)
 - Depth: 6" (15.2 cm)
 - 3 Year warranty

50580	Pump, 2 pk	\$620.00
50582	XL Vacuum Pump, 3pk	\$735.00
50586	Replacement Filtered Syringe, 6pk	\$12.50



Jodi® Pro Vacuum Pump

Portable vacuum pump for hearing instruments designed for the hearing healthcare professional. The unit features a 2000-hour clinically rated motor.

- Includes:
- 2 Filtered syringe assemblies
 - 2 Reamers
 - 3 needles
 - Power supply: 12-16 V, DC .6A
 - 3 Year warranty
- Height: 2.5" (6.4 cm)
 - Width: 4" (10.2 cm)
 - Depth: 3" (7.6 cm)

50581	Pro Pump	\$225.00
50586	Replacement Filtered Syringe, 6pk	\$12.50



Jodi® Consumer

The first hearing instrument vacuum system designed for consumer use. The Jodi-Consumer is a compact, portable and affordable vacuum pump that helps keep hearing instruments clean and dramatically reduces the need for professional cleaning.

- Includes:
- 1 Filtered syringe assembly
 - 1 Standard needle
 - 1 Reamer
 - Power supply: 12 V, DC
 - 1 Year warranty for patient use
 - 90 Day warranty for professional use
- Height: 3" (7.6 cm)
 - Width: 7.25" (18.4 cm)
 - Depth: 4.5" (1.4 cm)

50583	Vacuum Pump	\$92.50
50586	Replacement Filtered Syringe, 6pk	\$12.50



Jodi® Vacuum Needle

Replacement needles for Jodi Vacuums for cleaning thin tubes (open ear).

50575	Thin Tube Needle	\$8.25
50576	CIC ITE ITC	\$2.75
50578	BTE Needle	\$2.75



Lite Vac 15 Suction System

The compact Lite Vac vacuum pump is powerful enough for use in the hearing healthcare office and inexpensive enough for home use.

- Filtered wand assembly
- 20-gauge syringe tip
- Suction unit
- Reaming needle
- AC adapter, 12 V, DC
- Height: 1.88" (4.8 cm)
- Width: 3.88" (9.9 cm)
- Depth: 5.25" (13.3 cm)
- 1 Year warranty

50591	Vacuum System	\$219.00
50592	Replacement Filter	\$10.50
50593	Replacement Needle, 5pk	\$8.25



Serene HC-200 Vacuum

Hearing instrument vacuum, designed for home use; removes cerumen, dust, moisture, and condensation. Regular use saves on costly hearing instrument repair and cleaning. Can be used on hearing instruments with sound vents. Not for clinical use.

Includes:

- 2 sizes of suction syringes with soft tips
- Washable filter
- 1 Year warranty

50594	Vacuum	\$40.00
-------	--------	---------



The most *Advanced* surplus of earmold and tubing accessories within the industry!

We have you covered with our elite tubing and tubing accessories offering that are guaranteed to enhance your listening experience, regardless of the device the tubing is plugged into.



WESTONE TUBING *Size* CHART

TUBE SIZE	ID x OD (in.)	ID x OD (mm)
○ #12 Standard	.085" x .125"	2.16 x 3.17
○ #13 Medium	.076" x .122"	1.93 x 3.10
○ #13 Thick	.076" x .130"	1.93 x 3.30
○ #13 Super Thick	.076" x .142"	1.93 x 3.60
○ #15 Standard	.059" x .116"	1.50 x 2.95
○ #16 Thin	.053" x .085"	1.35 x 2.16
○ #16 Standard	.053" x .116"	1.35 x 2.95
○ Stethoscope Tubing	.125" x .187"	2.16 x 4.75

Tubing sizes illustrated are actual size. All tubing with the exception of Dry-Tubes are made of polyvinyl chloride material.

Preformed Tubes

Quilled on one end for easy insertion into the earpiece and double bent to lie close to the head. Westone's preformed tubing is available in the six most common sizes.

25 Pack

30101	#12 Standard	\$5.50
30105	#15 Standard	\$6.50
30106	#16 Standard	\$6.50
30102	#13 Medium	\$5.50
30103	#13 Thick	\$5.50
30104	#13 Super Thick	\$5.50



Single Bend Tubes

Westone's quilled single-bend tubing is available in four sizes: Standard, Medium, Thick and Super Thick.

25 Pack

30107	#12 Standard	\$5.50
30108	#13 Medium	\$5.50
30109	#13 Thick	\$5.50
30110	#13 Super Thick	\$5.50



Free Field Tubes

Tube designed for deep canal insertion.

30171	#13 Medium, 25pk	\$7.00
-------	------------------	--------



Dry-Tubes

Specially formulated tube helps prevent moisture buildup caused by humidity or perspiration.

30224	#13 Medium - 10pk	\$9.75
30225	#13 Thick - 10pk	\$9.75



Preformed DisappEar™ Tubes

Preformed DisappEar tubes match our DisappEar earpieces and are available in two colors: A/B (pink) or C (brown). DisappEar tubing is available in four sizes



10pk

30111	#12 Standard, Color A/B - 10pk	\$9.00
30112	#12 Standard, Color C - 10pk	\$9.00
30121	#13 Medium, Color A/B - 10pk	\$9.00
30122	#13 Medium, Color C - 10pk	\$9.00
30131	#13 Thick, Color A/B - 10pk	\$9.00
30132	#13 Thick, Color C - 10pk	\$9.00
30141	#13 Super Thick, Color A/B - 10pk	\$9.00
30142	#13 Super Thick, Color C - 10pk	\$9.00

25pk

30120	#13 Medium, Color A/B - 25pk	\$12.50
30119	#13 Thick, Color A/B - 25pk	\$12.50
30163	#13 Thick, Custom Color - 25pk	\$14.00
30118	#13 Super Thick, Color A/B - 25pk	\$12.50
30164	#13 Super Thick, Custom Color - 25pk	\$14.00

Dry-Tube, 10pk

30231	#13 Medium, Dry-Tube, Color A/B - 10pk	\$17.50
30232	#13 Medium, Dry-Tube, Color C - 10pk	\$17.50
30241	#13 Thick, Dry-Tube, Color A/B - 10pk	\$17.50
30242	#13 Thick, Dry-Tube, Color C - 10pk	\$17.50

Tubing Retention System (TRS)

Designed for all soft earpiece materials. Each TRS tube uses an acrylic ring to hold it in position in the earpiece. TRS tubes are easily replaced and will not tear silicone materials.



25pk

30263	#13 Thick - 25pk	\$18.50
30264	#13 Super Thick - 25pk	\$18.50

DisappEar, 10pk

30265	#13 Thick, DisappEar, Color A/B - 10pk	\$14.75
30266	#13 Thick, DisappEar, Color C - 10pk	\$14.75
30267	#13 Super Thick, DisappEar, Color A/B - 10pk	\$14.75
30268	#13 Super Thick, DisappEar, Color C - 10pk	\$14.75

Dry-Tube, 10pk

30273	#13 Thick, Dry-Tube - 10pk	\$15.75
30275	#13 Super Thick, DisappEar, Dry-Tube, Color A/B - 10pk	\$21.50
30276	#13 Thick, DisappEar Dry-Tube, Color C - 10pk	\$21.50

Horn Tubes

Horn Tubes use the horn effect to increase response in the high-frequency range.



Normal, 6Pk

3 mm		
30302	3 mm Thick	\$18.00
30303	3 mm Medium	\$18.00
4 mm		
30304	4 mm Medium	\$18.00

Dry-Tube, 6Pk

3 mm		
30333	3 mm, Dry-Tube	\$25.50
30332	3 mm, Thick, Dry-Tube	\$19.00
4 mm		
30334	4 mm, Dry-Tube	\$25.50

DisappEar, 6Pk

3 mm		
30305	3 mm, DisappEar, Color A/B	\$25.50
30306	3 mm, DisappEar, Color C	\$15.00
30345	3 mm, DisappEar, Color A/B	\$27.00
30346	3 mm, DisappEar, Color C	\$19.00
4 mm		
30307	4 mm, DisappEar, Color A/B	\$25.50
30308	4 mm, DisappEar, Color C	\$25.50
30347	4 mm, DisappEar, Color A/B	\$19.00
30348	4 mm, DisappEar, Color C	\$19.00

Horn TRS Tubes

Horn TRS tubes use acrylic rings to hold tubes in position. Designed for soft earpiece materials.



Acrylic Ring, 6Pk

3 mm		
30313	3 mm, Acrylic Ring, Clear	\$21.50
4 mm		
30314	4 mm, Acrylic Ring, Clear	\$21.50

Dry-Tube, 6Pk

3 mm		
30353	3 mm, Medium, Dry-Tube	\$26.50
4 mm		
30354	4 mm, Medium, Dry-Tube	\$19.00

DisappEar, 6Pk

3 mm		
30325	3 mm, DisappEar, Color A/B	\$10.00
30326	3 mm, DisappEar, Color C	\$10.00
30365	3 mm, DisappEar, Color A/B	\$19.00
30366	3 mm, DisappEar, Color C	\$19.00
4 mm		
30327	4 mm, DisappEar, Color A/B	\$10.00
30328	4 mm, DisappEar, Color C	\$26.50
30367	4 mm, DisappEar, Color AB	\$29.50
30368	4 mm, DisappEar, Color C	\$19.00



Continuous Flow Adapter (CFA) Tubes

Standard size tubes fitted with a continuous flow adapter elbow allow snap in and out tubing replacement.



*Final
CLEARANCE*

DisappEar Dry-Tube

6 Pack

30215	#13 Medium, DisappEar Dry-Tube, Color A/B - 6pk	\$15.00
30216	#13 Medium, DisappEar Dry-Tube, Color C - 6pk	\$15.00
30217	#13 Thick, DisappEar Dry-Tube, Color A/B - 6pk	\$24.00
30218	#13 Thick, DisappEar Dry-Tube, Color C - 6pk	\$24.00

10 Pack

30205	#13 Medium, DisappEar, Color A/B - 10pk	\$22.00
30206	#13 Medium, DisappEar, Color C - 10pk	\$22.00
30207	#13 Thick, DisappEar, Color A/B - 10pk	\$22.00
30208	#13 Thick, DisappEar, Color C - 10pk	\$22.00
30209	#13 Super Thick, DisappEar, Color A/B - 10pk	\$22.00

Normal

10 Pack

30203	#13 Standard, 10pk	\$17.00
30201	#13 Medium, 10pk	\$17.00
30202	#13 Thick, 10pk	\$17.00
30204	Reverse Curve, 10pk	\$10.00

Dry-Tube

6 Pack

30213	#13 Medium, Dry-Tube - 6pk	\$22.00
30214	#13 Thick, Dry-Tube - 6pk	\$22.00

Continuous Flow Adapter (CFA) Seating Ring

Replacement CFA Seating Ring

30219	Seating Ring	\$0.75
-------	--------------	--------



Tube-Lock with Tube

Tubing system that uses a brass retention ring for soft earpieces.



Normal, 10 Pack

30281	#13 Thick	\$10.75
30282	#13 Super Thick	\$10.75

Dry-Tube, 10 Pack

30285	#13 Thick, Dry-Tube	\$19.25
-------	---------------------	---------

DisappEar, 10 Pack

30291	#13 Thick, DisappEar, Color A/B	\$17.00
30292	#13 Thick, DisappEar, Color C	\$17.00
30293	#13 Super Thick, DisappEar, Color A/B	\$17.00
30294	#13 Super Thick, DisappEar, Color C	\$17.00

Tube-Lock without Tube

Brass tube-lock for use with 13 thick and 13 super thick tubing.

30283	#13 Thick	\$0.65
30284	#13 Super Thick	\$0.65



30-Inch Coiled Tube

Replacement tube for the Westone Style No. 32 communication earpiece. Allows freedom of movement and concealment. Preferred by most TV anchor personnel.

Clear

30181	Left Side	\$8.75
30182	Right Side	\$8.75

DisappEar


30187	DisappEar, A/B - Left	\$10.50
30188	DisappEar, A/B - Right	\$10.50



30-Inch Coiled Tube with Continuous Flow Adapter (CFA) Elbow

Replacement tube for the Westone Style No. 32 communication earpiece ensures freedom of movement and concealment. Comes with a CFA elbow, allowing quick and easy tube replacement for right or left side earpieces.

30183	30 inch Coiled Tube - CFA Elbow	\$10.50
-------	---------------------------------	---------



Long Over the Ear (L.O.E.) Tube

L.O.E. tube allows the use of receiver-style hearing instruments when the ear is not large enough to accept a snap ring. The tube is long enough to allow the hearing healthcare professional to adjust the length.

30185	Left Side	\$3.25
30186	Right Side	\$3.25



Straight Footage Tubing

Soft and flexible, ideal for replacing old or discolored tubing. Available in the seven most common sizes. Minimum total order of 10 feet for combined tubing items. Use Quantity to indicate length in feet.

Standard

30191	#12 Standard	\$0.35/ft
30195	#15 Standard	\$0.45/ft
30196	#16 Standard	\$0.45/ft

Thin

30197	#16 Thin	\$0.45/ft
-------	----------	-----------

Medium

30192	#13 Medium	\$0.35/ft
-------	------------	-----------

Thick

30193	#13 Thick	\$0.35
-------	-----------	--------

Super Thick

30194	#13 Super Thick	\$0.35
-------	-----------------	--------

Stethoscope Tubing

30198	Stethoscope Tubing (3 foot min. order)	\$0.55/ft
-------	--	-----------



Tubing Cements

Thin Cement for Acrylic Earpieces

Designed to provide a secure bond between acrylic earpieces and acoustic tubing. By slightly softening the internal sound bore wall and external tubing surface, Thin Cement fuses the two dissimilar materials together to provide a lasting, secure bond.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

30401	8oz can with empty .5oz. bottle	\$18.50
30402	.5oz. bottle	\$3.50



Heavy Cement for Acrylic Earpieces

Heavy Cement uses the same solvent found in Thin Cement but with acrylic powder in suspension in the fluid. Heavy Cement acts as a gap filler on earpieces where the sound bore diameter has been increased due to repeated retubing.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

30403	8oz can with empty .5oz bottle	\$19.50
30404	.5oz bottle	\$4.00



Vinyl Cement

High-quality tubing cement for use on soft vinyl earpieces. By slightly softening the inner wall of the sound bore and the outer wall of the acoustic tubing, Vinyl Cement fuses the two dissimilar materials together to provide a secure bond. Not compatible with silicone materials.

THIS PRODUCT CANNOT BE SHIPPED OUTSIDE THE US

30405	8oz can with empty .5oz bottle	\$18.50
30406	.5oz bottle	\$3.50



Glass Bottle for Cement

Empty .5oz glass bottle with an applicator brush in the cap keeps a small amount of tubing cement ready for use in the office or on the road.

30455	Glass Bottle and Brush Cap	\$2.25
-------	----------------------------	--------



Un-Cure™ Debonder

Super Glue Debonder softens the adhesive bond for easier removal of tubes anchored into earpieces with cyanoacrylate (super glue). Comes in a 1oz bottle.

20262	Debonder 1oz. Bottle	\$6.50
-------	----------------------	--------



Loctite® 406™

Industrial strength cyanoacrylate (super glue). Used in our labs for anchoring tubing into silicone and SLA earpieces. Comes in a 3g tube.

20260	Loctite 712 Accelerator 1.75oz	\$13.50
20263	Loctite 406 - 3g Tube	\$6.25
20264	Loctite 770 Primer -20g Bottle	\$16.95



Tubing Accessories

Tubing Expander Plier

Stainless steel pliers with two-way lever spring allow you to quickly and easily stretch acoustic tubing before fitting it over earhooks, elbows, etc.

30421	Curved Pliers	\$24.00
30422	Straight Pliers	\$24.00



Tube-Lock Kit

Includes an easy-to-use hand tool for installing tube-lock tubes in soft earmolds plus three 13 thick and three 13 super thick tube-locks.

30295	Kit	\$8.25
30296	Tool	\$5.50



Tube Puller, Monofilament

Soft and flexible enough to easily slip through the sound bore of an acrylic earpiece, allowing quick and simple tube replacement.

30431	Monofilament	\$8.25
-------	--------------	--------



Tube Puller, Wire Loop

Rigid enough to slip easily through the sound bore of a vinyl or silicone earpiece, allowing quick and simple tube replacement.

30432	Wire Loop	\$8.25
-------	-----------	--------



Reamer Kit

Set of four reamers for cleaning out old tubing, vents and earpieces. Reamer sizes: #60 (.04"), #52 (.0635"), #44 (.086") and #35 (.1080")

30450	Reamer Kit	\$30.00
-------	------------	---------



Reamer

Hand-powered tool removes old tubing from an earpiece and prepares the sound bore for the new tube.


30451	#13 Tubing	\$7.50
30453	#16 Thin Tubing	\$7.50



Core Drill for Tubes

Quickly and easily adds a sound bore or vent to vinyl or silicone earpieces. Convenient plunger removes the core that is cut from the earpiece during the drilling process. Comes in three sizes.

30441	Small Bores .120" (3.05 mm)	\$10.50
30442	Large Bores, .148" (3.76 mm)	\$10.50
30443	Venting, .065" (1.65 mm)	\$10.50
30444	Set (Small, Large and Venting)	\$28.00



Female Adapter, Standard

Style No. 32 earpiece replacement standard female adapter.

30461	Standard - Clear	\$0.95
-------	------------------	--------



Male Adapter

Style No. 1 earpiece replacement standard male adapter.

30471	Standard - Clear	\$0.95
30472	D-90 Miniature	\$0.95



Male Adapter with 3-Inch Tube

Standard male adapter with a #13 medium tube attached

30474	#13 Medium Tube	\$1.15
30475	D-90 Miniature with #13 Medium Tube	\$1.15



Plastic Elbow

Plastic Elbow can be used as an adhesive-free alternative when tubing soft earpieces.

30476	1 piece	\$0.95
30477	#13 Medium Tube with Elbow	\$1.15



Right Angle Flanged Elbow

Polyethylene 90 degree elbow with a male insert on both ends to grasp and hold tubing. One side is .5" / 13 mm long; other side is .75" / 19 mm long

30479	1pc	\$2.10
-------	-----	--------



"Y" Adapter for Standard Tube

"Y" shaped adapter for standard #13 medium tubing allows binaural monitoring of a single acoustic signal.

30478	1pc	\$2.85
-------	-----	--------



Westone® Cell Phone Adapter Kit

Cell phone adapter kit, 6pk, size 1-6

62410	Adapter Kit - 6pk	6.75
-------	-------------------	------



Earhook, Half Moon

Standard earhooks in the 1/2 moon style.

30481	Large, Clear - 10pk	\$13.00
30482	Large, Beige - 10pk	\$13.00



Earhook, Quarter Moon

Standard earhooks in the 1/4 moon style.

30486	Small, Clear - 10pk	\$13.00
30487	Small, Beige - 10pk	\$13.00
30488	Small, Brown - 10pk	\$13.00



Receiver Snap Ring and Spring

Snap Ring and Spring for body style acrylic earpieces.

30491	1pc	\$3.50
-------	-----	--------



SwivelComm Style No. 19 Pilot Bar Insert

Replacement insert for accepting the ball and socket assembly of a Westone Style No. 19 Pilot bar.

30494	Pilot Bar Insert	\$3.15
-------	------------------	--------



SwivelComm Pilot Bar

Complete replacement bar for the Westone Style No. 19. Comes with sound transmission tube.

30495	Left side, plastic	\$10.50
30496	Right side, plastic	\$10.50
30498	Left side, black anodized aluminum	\$15.75
30499	Right side, black anodized aluminum	\$15.75
30500	Pilot Bar Universal	\$5.00



Acoustic Modification & Venting

Select-A-Vents (SAV)

Replacement Select-A-Vents are a quick and easy way to control low-frequency amplification when using custom-fit earpieces.

Contains one each:

- #1 (.125")
- #2 (.110")
- #3 (.080")
- #4 (.050")
- #5 (.025")
- #6 (solid plug)

30721	Tree	\$0.95
-------	------	--------



Mini Vent Plug (MVP)

Our smallest venting system allows you to quickly and easily change vent sizes on a custom earpiece.

Contains one each:

- #1 (.100")
- #2 (.080")
- #3 (.060")
- #4 (.040")
- #5 (.030")
- #6 (solid plug)

30722	Mini Vent Plug Tree	\$0.95
-------	---------------------	--------



Acoustic Dampers

Acoustic dampers smooth frequency responses in hearing instruments. Dampers are color coded to indicate the amount of acoustic resistance.

30731	680 ohm, White, 1pc	\$2.25
30732	1000 ohm, Brown, 1pc	\$2.25
30733	1500 ohm, Green, 1pc	\$2.25
30734	2200 ohm, Red, 1pc	\$2.25
30735	3300 ohm, Orange, 1pc	\$2.25
30736	4700 ohm, Yellow, 1pc	\$2.25



30731 30732 30733



30734 30735 30736

The most *Vast* and *Superb* selection of hearing protection options since 1959

Since our founding in 1959, Westone has always been a lead innovator within the hearing protection industry. We offer the most advanced digital, universal, and custom-fit hearing protection devices that are perfect for all ages any environment.



NEW!



ULTIMATE HIGH PERFORMANCE HEARING PROTECTION

DefendEar® DigitalX™ is the latest development in digital hearing protection and leads the trend towards obstruction free, comfortable custom noise suppression for today's shooter. DigitalX combines superb functionality with an eye catching appearance.

NRR 26 | IPIL: 37.3/44/42


IN THE BOX

- Your Custom Earmolds
- Dx Digital Modules (L&R)
- High-Impact Hard Case
- User Guide
- Cleaning Cloth
- Oto-Ease Lubricant
- Oto-Clens Cleaning Solution
- Battery Tool/Cleaning Loop
- Embroidered DefendEar Patch
- 6 - Size 13 Batteries



Water-Shield Protection

Water-Shield is a highly durable amorphous carbon polymer coating that protects your module from accidental liquid exposure.



WARNING: The purpose of hearing protection is to reduce the impact of the sound and to reduce the potential for hearing loss. Hearing protection devices, including DefendEar products, do not eliminate sound. The user is responsible for determining whether the noise reduction provided by DefendEar products will provide adequate hearing protection based on the user's unique medical history and exposure level. Your medical history, particularly any prior trauma to the ear or prior hearing loss, as well as the noise (decibel level and exposure time) the user will be exposed to should be considered prior to using DefendEar products. If you have any questions concerning the appropriate exposure levels or the level of hearing protection provided by DefendEar products, or if you have hearing loss or ear injuries, contact your local hearing professional before using.

Dx1 DefendEar DigitalX1

64001 Dx1 \$490.21 MSRP: \$651.00

Our entry level product into this new line, the DX1, features a digitally optimized program that's great for any type of general shooting environment. The DX1 offers all of the protective features you need at an affordable price, with the full color flexibility available to our higher-end models. As with the entire DEDX series, the DX1 is available in our full range of metallic colors, consisting of six different metallic color options for users to choose from. Just like the rest of our DEDX models, the new DX1 comes with an extended two year warranty and is CE EN352 certified.



Shown - Metallic Silver

Dx3 DefendEar DigitalX3

64002 Dx3 \$588.13 MSRP: \$825.00

DX3 features three digitally optimized programs or "modes" for users to choose from, which include: Clay, Game and the all new Tactical mode. The DX3 offers the same digital rocker switch that is featured with our flagship DEDX model, the DX5, and is available in our full range of metallic color choices. The DX3 also utilizes a digital "Auto-Save" feature, so your personal user settings that are in place when you turn your DX3's off will still be in place the next time you turn them back on. Just like the rest of our DEDX models, the new DX3 comes with an extended two year warranty and is CE EN352 certified.



Shown - Metallic Red

Dx5 DefendEar DigitalX5

64003 Dx5 \$728.86 MSRP: \$999.00

The new DX5 is without a doubt the most technologically advanced form of hearing protection that Westone has ever offered! Our first digital model of hearing protection to feature five different programs or modes, as well as a mute function, the DX5 features a redesign of digital controls from its previous version (the DefendEar Digital Hunter), which it now replaces. The DX5 utilizes a new "center rocker switch" that allows users to more easily toggle between the five different modes featured in each module, while also allowing users to more easily control their volume settings with the same switch. Like the rest of the DEDX line, the DX5 offers a selection of six different metallic colors to choose from, and also utilizes a digital "Auto-Save" feature so your personal settings that are in place when you turn your DX5s off will still be in place the next time you turn them back on. Finally, as is the case with all other DEDX models, the DX5 comes with an extended two year warranty and is CE EN352 certified.



Shown - Metallic Blue

DX1 Replacement Module

64111 Module, 1pc \$245.11



DX3 Replacement Module

64121 Module, 1pc \$294.07



DX5 Replacement Module

64131 Module, 1pc \$364.43



NOTE: Eartips and other products used in the ear are non-refundable due to hygiene/health concerns.

DefendEar® Digital Accessories

All doors will fit new and existing DefendEar Digital modules.

DigitalX Replacement Metallic Doors

64100	Blue Metallic Door, 1pc	\$11.95
64101	Red Metallic Door, 1pc	\$11.95
64102	Copper Metallic Door, 1pc	\$11.95
64103	Gold Metallic Door, 1pc	\$11.95
64104	Silver Metallic Door, 1pc	\$11.95
64105	Green Metallic Door, 1pc	\$11.95



**REQUIRES
SIZE 13
BATTERY**

DefendEar® Digital Battery Door

DefendEar® Digital Battery Door

78317	Black, 1pc	\$5.50
78318	Blue, 1pc	\$5.50
78319	Red, 1pc	\$5.50
78320	Yellow, 1pc	\$5.50
78321	Orange, 1pc	\$5.50
78322	Green, 1pc	\$5.50



MSRP: \$7.95 each



**REQUIRES
SIZE 312
BATTERY**

DefendEar® Digital Battery Door Kit

The DefendEar Digital Battery Door Kits provide all the necessary tools needed for a replacement or desired change in battery doors.

- Includes:
 2 DED Battery Doors
 2 DED Battery Door Hinge Pins
 1 DED Battery Door Removal Tool

DefendEar® Digital Battery Door Kit

78311	Black	\$17.45
78312	Blue	\$17.45
78313	Red	\$17.45
78314	Yellow	\$17.45
78315	Orange	\$17.45
78316	Green	\$17.45



78315

MSRP: \$24.95 each

Monitor Vault Belt Loop Carrier

Stay connected with your DefendEar Digital Modules. The Monitor Vault Belt Loop Carrier offers peace of mind, loss prevention and convenience so you can stay focused on what you love to do most.

78306 Belt Loop Carrier \$3.50



DefendEar® Patch

Replacement DefendEar® Patch

78307 DefendEar Patch \$3.00



DefendEar® Digital Battery

Replacement DefendEar Digital Module Battery Door Hinge Pin. Compatible with all DefendEar Digital Module and doors.

78308 Battery Door Hinge Pin, 1pc \$1.95



Battery Removal Tool with Wax Loop

The DefendEar Digital Battery Removal Tool w/ Wax Loop allows for the user to easily remove existing batteries with ease when being replaced. The additional wax loop will provide easy and precise earwax removal to avoid damage to internal components.

78309 Battery Removal Tool with Wax Loop, 1pc \$2.75



Replacement Door Removal Tool


Replacement DefendEar Digital Battery Door Removal Tool allows for easy door removal when replacing existing DED doors.

78310 Door Removal Tool, 1pc \$2.50



ACCESSORY CABLES FIT **ANY** DefendEar® DIGITAL CUSTOM EARMOLD!

Remove the DefendEar Digital or DigitalX module from your custom earmold and insert the cable module.



SmartPhone Cable

The DefendEar Digital SmartPhone Cable is designed to fit any DefendEar Digital Custom Earmold. The cable features an inline microphone, volume slider and a multi-function button.

TECHNICAL INFORMATION:

- Impedance: 32 ohms
- Sensitivity: 106dB

IN THE BOX:

- Smartphone Cable
- High-Impact Hard Case
- User Guide
- Charging Cable

SOLD SEPARATELY:

- DefendEar Digital Custom Earmolds

64110 Cable, 1pc \$109.00



Bluetooth Cable

Designed to fit any DefendEar Digital Custom Earmold. The cable features a microphone and a 3-button headset.

TECHNICAL INFORMATION:

- Impedance: 32 ohms
- Sensitivity: 106dB
- Bluetooth: 4.1


IN THE BOX:

- Bluetooth Cable
- High-Impact Hard Case
- User Guide
- Charging Cable

SOLD SEPARATELY:

- DefendEar Digital Custom Earmolds

64109 Cable, 1pc \$129.00



Dx5 WIRELESS NECKLOOP

Dx5 Neckloop

The neckloop is designed to interface with the Dx5 modules when using the wireless comms mode. The cable jack plugs into the users' own radio, providing full two way communication without the need to remove your custom earpiece while providing ambient awareness.

Ideal for Tactical/Practical situations. For use only with Dx5 Wireless Communications program.

Note: NeckLoop is equipped with a jack plug. Your radio may require a jack adapter.

FEATURES:

- Anti-choke quick release
- Lapel Microphone
- Small PPT switch and clip

SOLD SEPARATELY:

- DefendEar Digital Custom Earmolds
- DefendEar Digital X5 Modules

TECHNICAL INFORMATION:


- Impedance: 32 ohms
- Sensitivity: 106dB

IN THE BOX:

- Wireless Neckloop
- Lapel Microphone
- PPT switch and clip
- User Guide

WORKS ONLY WITH COMMS MODE

64113	Neckloop, 1pc	\$139.00
-------	---------------	----------



Jack fits most Motorolas radios. An adapter may be needed to accomodate other radios.

Final
CLEARANCE

CLEARANCE

Contact your Sales Associate for Pricing Details

DefendEar® Digital 1

The DefendEar Digital 1 provides excellent all-around noise suppression, wind noise reduction, and ambient amplification. The DED 1's Game Mode is ideal for intermittent gunshot environments. Battery door not included. Earpiece sold separately.

78300	Module - Left	Discounted
78301	Module - Right	Discounted
MSRP: \$507.50 each		



DefendEar® Digital 2

The DefendEar Digital 2 offers two modes for either clay or game shooting scenarios. Game Mode is optimized for suppressing intermittent shots, and clay mode is optimized for suppressing successive shots, allowing more flexible noise suppression while still amplifying ambient sounds. Battery door not included. Earpiece sold separately.

78302	Module - Left	Discounted
78303	Module - Right	Discounted
MSRP: \$607.50 each		



DefendEar® Digital Hunter

The DefendEar Digital Hunter boosts low-level sounds for detecting game moving through leaves and undergrowth while providing impulse suppression tailored for successive shots from larger caliber firearms. The DefendEar Hunter has a second mode to reduce wind noise typical in the open range and in hunting stands. Battery door not included. Earpiece sold separately.

78304	Module - Left	Discounted
78305	Module - Right	Discounted
MSRP: \$807.50 each		



DefendEar® Digital SHOOTER

Shooter ear plugs are specifically designed to fit into any ear without the need for ear impressions. The Shooter allows for clear communication, as well as enhanced ambient awareness. When a gunshot is fired, the electronic circuit instantly suppresses the loud noise, protecting your hearing!

Usable straight out of the box and featuring state-of-the-art digital technology, the new Westone Shooter offers an affordable way to experience comprehensive electronic hearing protection.

- Features:**
- NRR 30 dB
 - Output up to 86 dB
 - Three amplification levels
 - Wind noise reduction
 - Low profile
 - Universal fit
 - Comfortable
 - Hear conversations as normal
 - No gun mount interference
 - 1 Year warranty


- In the Box:**
- Right and left color-coded ear plugs
 - 5 Different sized high quality foam ear tips
 - 2 Spare wax guards
 - 1 Brush with magnet
 - 6 Size 10 batteries
 - 1 Zip pouch storage case

78346 Shooter Protection \$299.00
MSRP: \$399.00



EXAMPLE SOUNDS	
INSTANT DAMAGE	147 Gun Shot
	130 Jackhammer
THRESHOLD OF PAIN	120 Rock concert
	110 Power tools
	90 Lawn mower
SHOOTER™ PROTECTION KICKS-IN	85 Begin Danger Zone*
	80 City traffic
	50 Normal conversation
	30 Whisper
	20 Rustling leaves

30 dB BLAST PROTECTION | **86 dB** UP TO SOUND ENHANCEMENT



SAFETY NOTICE: This product is designed for people with normal hearing only! The purpose of hearing protection is to reduce the impact of the sound and to reduce the potential for hearing loss. Hearing protection devices, including DefendEar products, do not eliminate sound. The user is responsible for determining whether the noise reduction provided by DefendEar products will provide adequate hearing protection based on the user's unique medical history and exposure level. Your medical history, particularly any prior trauma to the ear or prior hearing loss, as well as the noise (decibel level and exposure time) the user will be exposed to should be considered prior to using DefendEar products. If you have any questions concerning the appropriate exposure levels or the level of hearing protection provided by DefendEar products, or if you have hearing loss or ear injuries, contact your local hearing professional before using.

Wax Guards for DefendEar® Digital Shooter

Replacement wax guards for the DefendEar Shooter with insertion/removal tool.

78347 Right - 15pk \$13.99
78348 Left - 15pk \$13.99



NOTE: Eartips and other products used in the ear are non-refundable due to hygiene/health concerns.

DefendEar® Filters

Etymotic Research® Filters for DefendEar® Concert™

Designed for use with the Westone's DefendEar Concert (Style 49) Musicians' Earplugs. Etymotic filters feature a flat frequency response that allows musicians to hear their performance accurately while providing hearing protection.

Style 49



Style 49 Filters



Single Filter - 9db

\$20.00

79311	Single - Beige
79313	Single - Clear
79315	Single - Brown
79317	Single - Red
79319	Single - Blue

MSRP: \$40.00

Filter Pair - 9db

\$40.00

79312	Pair - Beige
79314	Pair - Clear
79316	Pair - Brown
79322	Pair - Red & Blue

MSRP: \$80.00

Single Filter - 15db

\$20.00

79407	Single - Blue
79408	Single - Red
79410	Single - Clear
79413	Single - Brown
79415	Single - Beige

MSRP: \$40.00

Filter Pair - 15db

\$40.00

79409	Pair - Red & Blue
79411	Pair - Clear
79414	Pair - Brown
79416	Pair - Beige

MSRP: \$80.00

Single Filter - 25db

\$20.00

79425	Single - Beige
79427	Single - Clear
79429	Single - Brown
79424B	Single - Blue
79424R	Single - Red

MSRP: \$40.00

Filter Pair - 25db

\$40.00

79424	Pair - Red & Blue
79426	Pair - Beige
79428	Pair - Clear
79430	Pair - Brown

MSRP: \$80.00

Style 49 Solid Attenuation Adapters

This solid plug snaps into the Westone's DefendEar Concert (Style 49) Musicians' Earplugs, plugging the filter sound path and offering additional attenuation over the ER filter.

79431	Pair - Black	\$3.25
79432	Pair - Beige	\$3.25
79433	Pair - Clear	\$3.25

MSRP: \$5.00



TRU Professional Filters

TRU™ Custom Filters

Compact, low profile, and interchangeable, the TRU filters can be used with TRU Universals or TRU Customs as well as our EAS in-ear monitors. For high volume on stage environments, the TRU WM25 is our strongest attenuating filter.



Average Attenuation: 14dB | NRR: 2

WM10 - Pair		
62882	10 Green, Clear	\$24.95
62883	10 Smoke, Green Filter	\$24.95
		MSRP: \$49.95
WM10 - Single		
62862	10 Green, Clear	\$13.95
62863	10 Smoke, Green Filter	\$13.95
		MSRP: \$26.95



Average Attenuation: 18dB | NRR: 10

WM16 - Pair		
62884	16 Clear, Blue Filter	\$24.95
62885	16 Smoke, Blue Filter	\$24.95
		MSRP: \$49.95
WM16 - Single		
62864	16 Clear, Blue Filter	\$13.95
62865	16 Smoke, Blue Filter	\$13.95
		MSRP: \$26.95



Average Attenuation: 21dB | NRR: 13

WM20 - Pair		
62886	20 Clear, White Filter	\$24.95
62887	20 Smoke, White Filter	\$24.95
		MSRP: \$49.95
WM20 - Single		
62866	20 Clear, White Filter	\$13.95
62867	20 Smoke, White Filter	\$13.95
		MSRP: \$26.95



Average Attenuation: 24dB | NRR: 15

WM25 - Pair		
62888	25 Clear, Black Filter	\$24.95
62889	25 Smoke, Black Filter	\$24.95
		MSRP: \$49.95
WM25 - Single		
62868	25 Clear, Black Filter	\$13.95
62869	25 Smoke, Black Filter	\$13.95
		MSRP: \$26.95

NOTE: Eartips and other products used in the ear are non-refundable due to hygiene/health concerns.

TRU Universal Hearing Protection

TRU™ Professional

The TRU Professional offers excellent protection and comfortable long-term wear for on-the-job scenarios, as well as everyday use wherever speech or music clarity is a must.

77700	WM16 Blue	\$29.95
77701	WM25 Black	\$29.95
		MSRP: \$45.00



TRU™ Recreational

An exceptionally accurate 20 dB across the frequency range, the WR 20 is one of the most transparent hearing protection filters in today's market.

Ideal for: Concerts, Nightclubs, Sporting Events, Air Travel and Movies

77702	WR20 Pink	\$11.95
77703	WR20 Smoke	\$11.95
77704	WR20 Blue	\$11.95
		MSRP: \$19.99



TRU™ Shooter

TRU Shooter ear plugs have a switch on the end that opens and closes the filter, letting you toggle between 18dB and 5dB of hearing protection. Choose whether to also filter out additional background sound with the filters closed, or whether to allow for clear hearing with the filters open.

77724	TRU Shooter	\$15.99
		MSRP: \$19.99



TRU™ Universal POP Display

Sell more Universal Hearing protection with our TRU Universal "Point of Purchase" display. Loaded with six TRU's (2 Pink, 2 Smoke, and 2 Blue) with 20 db reduction. The TRU hearing protection display is a great way to raise hearing protection awareness and increase revenue in any practice.

79467	Display	\$37.95
62879	Display with TRU WR20	\$97.95
62876	WR20 Refill	\$60.00



TRU™ Replacement Tip

77720	Clear	\$3.25
77721	Smoke	\$3.25
77722	Blue	\$3.25
77723	Pink	\$3.25

MSRP: \$5.00



Earplugs

EarPlanes™

EarPlanes pressure-regulating earplugs reduce discomfort of travelers who fly with colds, allergies, or sinus conditions. One pair of EarPlanes is good for one round-trip flight (two flight segments).

54031	Adult, 1 pair	\$3.95
54041	Child, 1 pair	\$3.95



E-A-R® Classic™ "SOFT" Earplugs

Thermal-reactive foam earplugs start out firm for easy insertion. Then as they warm to body temperature, they soften and conform to the ear canal for more comfort and a better acoustic seal. Made from polyurethane, no latex.

54092	Earplug - 200 pair	\$32.50
54093	Earplug - 1 pair	\$0.35



E-A-R® Foam Earplugs

Bright yellow disposable earplugs expand to the size and shape of the ear canal for a great fit every time. Made from polyurethane, no latex. Each pair packaged separately in a cardboard pillow pack.

54071	Classic, Small - 1 pair	\$0.35
54072	Classic - 1 pair	\$0.35
54073	Classic Plus - 1 pair	\$0.35
54081	Classic Small, Box - 200 pair	\$32.50
54082	Classic, Box - 200 pair	\$32.50
54083	Classic Plus, Box - 200 pair	\$32.50



Migraine X Earplugs

Introducing a complete, proactive way to manage weather-related migraines, before they start. This ear device utilizes a CeramX filter specifically designed to deal with ground level changes in atmospheric pressure that reduces migraines' intensity and duration.

54042	Earplugs	\$12.99
-------	----------	---------



Filtered Foam Earplugs

Filtered Foam Earplugs have an integrated acoustic filter, are fully ventilated, provide protection from noises up to 104 dB(A) and allow users to communicate without having to remove their earplugs. This great product innovation improves safety, increases worker compliance and offers users a comfortable earplug that can be worn throughout the day. The ventilated design reduces heat build-up, eliminating discomfort and allowing users to function normally, whilst still being protected from dangerous noise levels at the same time.



70016	Black, Large	\$5.99
70018	Black, Medium	\$5.99
		MSRP: \$9.99

Ety-Plugs™ High Fidelity Universal Fit Earplugs

Great hearing protection for concerts, sporting events, airline travel and industrial use.

- Reduces sound approximately 20 dB at all frequencies
- High-fidelity hearing protection
- Lowers the overall noise level but preserve sound quality
- Each pair comes with a small carrying case and black cord
- Color cords sold separately



79444	with Case, Pair - Clear	\$10.25
-------	-------------------------	---------

Ety-Plugs™ Cords

Brightly colored cords help keep Ety-Plugs together. Each 24" (60.9 cm) long cord comes with a socket on each end that securely attaches to the Ety-Plug handle and is ideal for hanging around the neck.



79446	Pair - Blue	\$1.35
79450	Pair - Black	\$1.35

Style No. 47 Replacement Earplug Filter

One replacement filter for the Westone Style No. 47 earplugs.



42295	Replacement filter, 1pc	\$3.75
-------	-------------------------	--------

Infants & Children

EM's 4 Bubs

The world's first hearing protection earmuffs designed just for babies. The patented design is light and compact weighing only 100 grams. This sturdy design uses an interchangeable headband with a soft Velcro tab that adheres to the smaller adjustment tab.

- Headbands helps disperse pressure evenly
- The world's smallest hearing protection earmuffs
- 22dB noise reduction rating (NRR)

63000	Star's and Stripes Headband	\$7.00
63004	Blue/White	\$22.00
63005	Pink/White	\$22.00
63006	Black	\$22.00
63007	White	\$22.00



63000
Headband Only



63004



63005



63006



63007

EM's 4 Kids

The EM's 4 Kids sturdy earmuffs patented design is light and compact weighing 190 grams and are the ideal solution for children experiencing auditory sensitivities or wherever loud noise may cause distress. Perfect for all kids in a wide variety of settings wherever noise reduction is needed.

- Height 3.9" folded
- Height 6.7" unfolded
- 20dB noise reduction rating (NRR)

63001	Blue	\$19.50
63002	Mint	\$19.50
63003	Pink	\$19.50



63001



63002



63003

KidzSafe™ Earbuds

Protect children's hearing with KidzSafe™ myDesign™ Earbuds by utilizing KidzSafe™ patented volume-limited technology (<85dB). Listen safely to music and videos at moderate volume levels to help prevent noise induced hearing loss.

Includes:

- 1 Pair of Customizable Earbuds
- 47" flat tangle free cord with mic
- 16 Removable Stickers
- 1 Drawstring Bag
- 6 Silicone tips for perfect fit

79307	Pink/Purple	\$17.50
79308	Blue/Orange	\$17.50



Personalize it!



KidzSafe™ Wireless Speakers

This wireless speaker is a safe and fun way to take your music with you everywhere you go. Just like the headphones – this speaker is volume limited to ≤ 85dB to keep little ears safe.

- Bluetooth 3.0 technology with a range up to 32 feet.
- Listen for up to 5 hours on a single charge
- Easy-touch on board controls
- Built in mic, so you can answer phone calls hand-free
- Splash-resistant, boasting an IPX4 rating

Includes:

- 2 Interchangeable silicone bumpers
- 2 Carabiner Clips
- Removable suction cup to secure device
- Micro USB charging cable

63012	Gray	\$30.00
-------	------	---------

Final CLEARANCE



Tinnitus Maskers

My Butterfly Pillow

My Butterfly Pillow (Patent Pending) is the first pillow designed by side sleepers for side sleepers. Embedded beneath the ear pillow are two Night Owl speakers along with a Bluetooth receiver providing crisp, clean sound for white noise, audio books, music or even television directly into the pillow without disturbing your sleeping partner.

63014 Pillow 1.5" Riser \$109.99



www.MyButterflyPillow.com



The most *Expansive* selection of audiological consumer-based products

Westone offers an ever-growing selection of consumer-based products that are all designed to enhance or further protect customers' hearing. From ear drops and lubricants to assisted listening and hearing loss prevention devices, we have what you need to ensure your ears remain functioning properly!



Batteries



Rayovac EXTRA® Mercury-Free Zinc Air Batteries

Rayovac Extra Advanced hearing aid batteries take advantage of our new Active Core Technology; a breakthrough new cell design which maximizes internal space to hold more active ingredients, ensuring our batteries are the longest lasting on the market. This hearing aid battery is available in sizes 10, 13, 312 and 675 (for cochlear implants).

40 Count Box

52141	Size 13	\$18.50
52142	Size 312	\$18.50
52143	Size 675	\$18.50
52144	Size 10	\$18.50

60 Count Box

52161	Size 13	\$22.50
52162	Size 312	\$22.50
52163	Size 675	\$22.50
52164	Size 10	\$22.50

80 Count Box


52181	Size 13	\$25.50
52182	Size 312	\$25.50
52184	Size 10	\$25.50




52144



52161



52182



52163

Rayovac® Mercury-Free Alkaline Batteries

Provides long-lasting power and contains no mercury.

52512	AAA Batteries, 4pk	\$3.35
52521	AA Batteries, 4pk	\$3.35
52531	9 Volt Battery, 1pc	\$3.35



52512



52521



52531

Xcell™ by Rayovac® Mercury-Free Private Label

Custom imprinting provides an easy and convenient way to keep your business name in front of your customers. These Xcell Private Label batteries are mercury-free premium long-life zinc-air batteries that are manufactured in the USA.



Customize it!

40 Count Box


52196	Size 10	\$9.20
52197	Size 13	\$9.20
52198	Size 312	\$9.20
52199	Size 675	\$9.20

60 Count Box


52200	Size 10	\$13.80
52201	Size 13	\$13.80
52202	Size 312	\$13.80
52203	Size 675	\$13.80

80 Count Box


52204	Size 10	\$18.40
52205	Size 13	\$18.40
52206	Size 312	\$18.40




52204



52205



52202



52199

Xcell™ by Rayovac® Mercury-Free Unprinted

40 Count Box

52196-B	Size 10	\$7.60
52197-B	Size 13	\$7.60
52198-B	Size 312	\$7.60
52199-B	Size 675	\$7.60

60 Count Box

52200-B	Size 10	\$11.40
52201-B	Size 13	\$11.40
52202-B	Size 312	\$11.40
52203-B	Size 675	\$11.40

80 Count Box

52204-B	Size 10	\$15.20
52205-B	Size 13	\$15.20
52206-B	Size 312	\$15.20



52204-B



52205-B



52202-B



52199-B

Varta Mercury-Free Private Label Batteries

Private label batteries offer the high quality and long life of any name brand battery at a lower price. Custom imprinting provides an easy and convenient way to keep your business's name in front of your customers.

60 Count Box

52086	Size 13	\$21.50
52087	Size 312	\$21.50
52088	Size 675	\$21.50
52089	Size 10	\$21.50



Varta Mercury-Free Unprinted

Private label batteries offer the high quality and long life of any name brand battery at a lower price. Custom imprinting provides an easy and convenient way to keep your business's name in front of your customers.

60 Count Box

52086-B	Size 13	\$21.50
52087-B	Size 312	\$21.50
52088-B	Size 675	\$21.50
52089-B	Size 10	\$21.50



Duracell® Activair® Batteries

Duracell Activair batteries are different from any other batteries in the market. Color-coded tabs ensure easy size recognition by your customers.

Features:

- Highest brand recognition in the world
- Duracell Activair batteries provide a high and stable voltage level
- Comprehensive product testing and technical support
- Battery housing made of stainless steel prevents corrosion
- Superior performance in all conditions thanks to triple sealing and PTFE film
- Protective packaging

52098	Size 13, 80ct Box	\$25.75
52099	Size 312, 80ct Box	\$25.75
52100	Size 675, 60ct Box	\$20.75
52101	Size 10, 80ct Box	\$25.75



52098

PowerOne Mercury-Free Batteries

Power One Mercury Free Batteries for high-powered digital and analog hearing instruments offer more capacity than comparable hearing instrument batteries to ensure stable voltage even at high output levels. The result is maximum amplification, excellent sound quality and exceptional life span in all types of hearing instruments.

60 Count Box

52065	Size 13	\$22.50
52066	Size 312	\$22.50
52068	Size 675	\$22.50
52069	Size 10	\$22.50
52070	Implant Plus Size 675	\$39.50



52068

52066

52065

52069

PowerOne Accu plus PocketCharger

This convenient charger, the size of a credit card, fits in any pocket and allows the user to recharge their hearing instrument batteries anywhere. Batteries sold separately.

- Charges PowerOne ACCU plus Batteries: size 10, 13, and 312
- Recharges in approximately 2.5 hours
- Each charge can last 10-14 hours
- Can be used up to 18 times before recharging
- 1 Year warranty

52075	PocketCharger	\$69.95
52076	Rechargeable Batteries, Size 13 - 2pk	\$19.50
52077	Rechargeable Batteries, Size 312 - 2pk	\$19.50



Battery Testers

Universal Battery Tester

Universal Battery Tester is easy-to-use and tests all hearing instrument batteries, as well as other 1.5 V button cells. Features an easy-to-read needed meter and a convenient battery storage compartment.

52632	Tester, 1pc	\$5.25
-------	-------------	--------



Digital Battery Tester

Digital Battery Tester with bar graph. Storage compartment will hold two spare batteries. Quantity pricing discounts are available.

52629	Black, 1pc	\$4.25
52631	Maroon, 1pc	\$4.25
52630	Custom, 1pc	\$4.25



Keychain Battery Tester

Convenient key chain battery tester features a storage compartment for extra batteries and an attractive gold-plated key ring. Key chain battery testers make a great giveaway item to accompany new hearing instrument sales. Quantity pricing discounts are available.

52641	Tester, 1pc	\$3.25
-------	-------------	--------



Household Battery Tester

Household battery tester is ideal for testing the most commonly used batteries. Works with all AAA, AA, C, D, 9 V batteries, and 3 V lithium cells. Tester features a digital LCD indicator to show battery status. Not to be used with hearing aid batteries.

52635	Tester, 1pc	\$8.95
-------	-------------	--------



Battery Caddies

Battery Caddie Private Label

Custom caddies help your patients tell others about your practice. Your patients will love the "New EZ Open" feature of this caddy. The caddy will hold two hearing instrument battery of any size.

Custom imprinting is available with your business name and phone number. Minimum order of 400 or more is required.

52650	Battery Caddie, 400pk	\$275.00
-------	-----------------------	----------



Battery Holder with Chains

Battery holders provide a convenient storage place for extra batteries. Select either a two cell holder which holds all battery sizes, or the three cell which holds all sizes except the 675 battery.

52651	2 cell, White - 50pk	\$35.50
52652	3 cell, White - 50pk	\$35.50
52655	2 cell, Yellow - 50pk	\$35.50
52656	3 cell, Yellow - 50pk	\$35.50
52671	Custom Imprint (250 minimum order)	\$1.15



customization available for white holders only

Caddy LED

This lightweight 3-in-1 caddy features:

- Battery tester; tests all hearing instrument batteries
- Storage; stores 2 hearing instrument batteries
- Flashlight; features a bright LED light source

52645	Caddy, 1pc	\$4.50
-------	------------	--------



Battery Accessories

Stainless Steel Ball Chain

Replacement for lost or broken stainless steel ball chains for our battery holders

52681	Ball Chain, 1pc	\$0.50
-------	-----------------	--------



Battery Extractor Tool

Assists in removing even the smallest batteries from the tightest battery compartments

52691	Battery Extractor Tool, 1pc	\$3.85
-------	-----------------------------	--------



Magnetic Finger

Manipulate small metal objects with this powerful fingertip magnet. Magnetic Finger helps patients easily change their hearing instrument batteries and is perfect for hearing instrument adjustments in the clinic.

52695	Small/Medium	\$7.95
52696	Large/XLarge	\$7.95




Hearing Aid Pillow

Contoured pillow features a "trench" that acts as a sound channel. The user can lie comfortably on his or her side and still receive sound without feedback. Ideal for patients who have undergone recent ear surgery or those with sensitive outer ears. The pillowcase is removable for washing.

Dimensions: Length: 12"(30.5 cm), Width: 10"(25.4 cm), Depth: 3.5"(8.9 cm)

54163	Hearing Aid Pillow	\$33.50
-------	--------------------	---------



Hearing Aid Accessories

Ear Band-It® Ultra

Swimming ear protection featuring custom-designed, embossed neoprene on the inner side of the headband. Provides extra grip on the forehead, ears, and back of the head to avoid slipping during most types of water activity.

- Helps prevent swimmer's ear
- Recommended by physicians for swimmers with ear tubes
- Head band is non-reversible due to Grip-Tek™ neoprene inner lining
- Manufactured with 4 mm 100% neoprene

EAR BAND-IT *Size* CHART

SIZE	BAND LENGTH	HEAD SIZE	OPTIMUM FIT	RECOMMENDED AGE
SMALL	19.56"	16" - 20"	16.5" - 18.5"	1 - 3 YEARS
MEDIUM	20.50"	17" - 21"	18.5" - 20.5"	4 - 9 YEARS
LARGE	22.25"	19.5" - 24"	20.5" - 24"	10+ YEARS



Black

Tan

Blue

Green

Neon Orange

Pink

Neon Yellow

Small - 19.56"

\$10.50

Medium - 20.5"

\$10.50

Large - 22.25"

\$10.50

54116 Small, Blue

54126 Medium, Blue

54109 Large, Black

54117 Small, Green

54127 Medium, Green

54110 Large, Tan

54118 Small, Neon Orange

54128 Medium, Neon Orange

54136 Large, Blue

54119 Small, Pink

54129 Medium, Pink

54137 Large, Green

54120 Small, Neon Yellow

54130 Medium, Neon Yellow

54138 Large, Neon Orange

54139 Large, Pink

54140 Large, Neon Yellow

Hearing Aid Sweat Bands

Hearing Aid Sweat Bands keep moisture and dust out of BTE hearing instruments, but do not inhibit sound. Easy to install, no special tools required.



SWEAT BANDS *Size* CHART

Mini Slim - 1"	Mini - 1"	Large - 1.75"
Small Slim - 1.25"	Small - 1.25"	X Large - 2.125"
Medium Slim - 1.5"	Medium - 1.5"	5XL (FM Systems) - 3"

Mini Slim, Natural	\$15.50	Small Slim, Natural	\$15.50	Medium Slim, Natural	\$15.50
50450 Mini Slim, Natural		50460 Small Slim, Natural		50470 Medium Slim, Natural	
50451 Mini Slim, Light Brown		50461 Small Slim, Light Brown		50471 Medium Slim, Light Brown	
50452 Mini Slim, Grey		50462 Small Slim, Grey		50472 Medium Slim, Grey	

Mini, 2pk	\$15.50	Small, 2pk	\$15.50	Medium, 2pk	\$15.50
50411 Mini, White		50415 Small, White		50419 Medium, White	
50412 Mini, Natural		50416 Small, Natural		50420 Medium, Natural	
50413 Mini, Light Brown		50417 Small, Light Brown		50421 Medium, Light Brown	
50414 Mini, Medium Brown		50418 Small, Medium Brown		50422 Medium, Medium Brown	
50431 Mini, Black		50432 Small, Black		50433 Medium, Black	
50436 Mini, Grey		50437 Small, Grey		50438 Medium, Grey	

Large, 2pk	\$15.50	Extra Large, 2pk	\$15.50	5XL for FM Systems, 2pk	\$15.50
50423 Large, White		50427 Extra Large, White		50444 5XL, White	
50424 Large, Natural		50428 Extra Large, Natural		50445 5XL, Black	
50425 Large, Light Brown		50429 Extra Large, Light Brown		50446 5XL, Grey	
50426 Large, Medium Brown		50430 Extra Large, Medium Brown			
50434 Large, Black		50435 Extra Large, Black			
50439 Large, Grey		50440 Extra Large, Grey			



Protect your Hearing Instrument from Sweat, Dirt, Moisture, Loss & Wind Noise. Ear Gear is popular with both children and adults and is great for use during sports, outdoor activities, or in the workplace.

Features:

- Protects hearing instrument from dirt, moisture and loss
- Water resistant
- Double wall spandex nylon sleeve
- Acoustically transparent
- Fully washable and reusable
- Available in monaural or binaural
- Cordless - binaural only
- Standard 9 in. (22.86 cm) cord on monaural
- Standard 12 in. (30.48 cm) cord on binaural - cord optional
- 1 Year satisfaction guarantee



Ear Gear® Original

Binaural, without cord		\$18.50
55742	Camo	
55744	Black	
55747	Red/Orange	
55749	Royal Blue	
55751	Pink	
55753	Grey	
55754	Brown	
55755	Beige	

Binaural, with cord		\$20.00
55722	Camo	
55724	Black	
55727	Red/Orange	
55729	Royal Blue	
55731	Pink	
55733	Grey	
55734	Brown	
55735	Beige	

Monaural, with cord		
55702	Camo	\$15.00
55704	Black	\$15.00
55707	Red/Orange	\$15.00
55708	Pink/Purple	\$15.00
55709	Royal Blue	\$18.50
55711	Pink	\$18.50
55713	Grey	\$18.50
55714	Brown	\$18.50
55715	Beige	\$18.50

EarGear® Original for Glasses

Binaural, Beige		
57355	with cord	\$27.95
57755	without cord	\$25.70



2.5" (6.4 cm)



Ear Gear® FM Systems

Binaural, without cord

55802	Camo	\$18.50
55804	Black	\$18.50
55809	Royal Blue	\$18.50
55811	Pink	\$18.50
55813	Grey	\$18.50
55814	Brown	\$18.50
55815	Beige	\$18.50

Monaural, with cord

55762	Camo	\$15.00
55764	Black	\$18.50
55767	Red/Orange	\$15.00
55769	Royal Blue	\$15.00
55771	Pink	\$18.50
55773	Grey	\$15.00
55774	Brown	\$15.00
55775	Beige	\$15.00

Binaural, with cord

55782	Camo	\$15.00
55787	Red/Orange	\$20.00
55789	Royal Blue	\$20.00
55791	Pink	\$20.00
55793	Grey	\$20.00
55794	Brown	\$20.00

EarGear® FM System for Glasses

Binaural, Beige

57795	with cord	\$27.95
57815	without cord	\$25.70



2" (5.1 cm)



Ear Gear® Cochlear

Binaural, without cord

55982	Camo	\$15.00
55984	Black	\$21.50
55987	Red/Orange	\$21.50
55989	Royal Blue	\$21.50
55991	Pink	\$15.00
55993	Grey	\$21.50
55994	Brown	\$21.50
55995	Beige	\$21.50

Monaural, with cord

55902	Camo	\$15.00
55904	Black	\$21.50
55907	Red/Orange	\$15.00
55908	Pink/Purple	\$15.00
55909	Royal Blue	\$15.00
55911	Pink	\$15.00
55913	Grey	\$21.50
55914	Brown	\$15.00
55915	Beige	\$21.50

Binaural, with cord

55942	Camo	\$15.00
55944	Black	\$25.70
55947	Red/Orange	\$15.00
55948	Pink/Purple	\$15.00
55949	Royal Blue	\$25.70
55951	Pink	\$25.70
55953	Grey	\$25.70
55954	Brown	\$25.70
55955	Beige	\$25.70

EarGear® Cochlear for Glasses

Binaural, Beige

57955	with cord	\$30.95
57995	without cord	\$27.95



1.25" (3.2 cm)



Ear Gear® Mini Curved

Binaural, without cord		\$18.50
56082	Camo	
56084	Black	
56087	Red/Orange	
56089	Royal Blue	
56091	Pink	
56093	Grey	
56094	Brown	
56095	Beige	

Binaural, with cord		\$20.00
56062	Camo	
56064	Black	
56067	Red/Orange	
56069	Royal Blue	
56071	Pink	
56073	Grey	
56074	Brown	
56075	Beige	

Monaural, with cord

56042	Camo	\$18.50
56044	Black	\$18.50
56047	Red/Orange	\$18.50
56049	Royal Blue	\$18.50
56051	Pink	\$18.50
56053	Grey	\$18.50
56054	Brown	\$15.00
56055	Beige	\$18.50

EarGear® Mini Curved for Glasses

Binaural, Beige

57075	with cord	\$15.00
57095	without cord	\$25.70



.75" (1.9 cm)



Ear Gear® Micro

Binaural, without cord		\$18.50
56364	Black	
56371	Pink	
56372	Camo	
56373	Grey	
56374	Brown	
56375	Beige	
56377	Red/Orange	

Binaural, with cord

56344	Black	\$20.00
56349	Royal Blue	\$20.00
56351	Pink	\$20.00
56352	Camo	\$20.00
56353	Grey	\$20.00
56354	Brown	\$20.00
56355	Beige	\$20.00
56357	Red/Orange	\$15.00

Monaural, with cord

56324	Black	\$18.50
56329	Royal Blue	\$15.00
56331	Pink	\$15.00
56332	Camo	\$18.50
56333	Grey	\$18.50
56334	Brown	\$18.50
56335	Beige	\$18.50
56337	Red/Orange	\$18.50

EarGear® Micro for Glasses

Binaural, Beige

57375	without cord	\$25.70
57455	with cord	\$27.95



EarGear® ITE

Ear Gear effectively protects and secures expensive hearing instruments for both children and adults during sporting activities, outdoor events, or in the workplace. Ear Gear ITE can be used with any ITE that has an eyelet installed on the faceplate.

57798	Monaural - Beige	\$10.10
57799	Binaural - Beige	\$10.75

1.25" (3.2 cm)



EarGear® Baha

57915	Monaural with Cord - Beige	\$23.50
57935	Monaural without Cord - Beige	\$21.50
57944	Monaural without Cord - Black	\$15.00



Huggie Aids®

Designed to keep BTEs in place, Huggie Aids are easy to install and are available in sizes that fit most hearing instruments. To fit the Huggie Aid, measure the circumference of the ear (next to head) to determine appropriate size.

Mini, 2 pack	\$8.00	Large, 2 pack	\$8.00
54201A	3" (76.2 mm)	54202A	3" (76.2 mm)
54201B	3.25" (82.5 mm)	54202B	3.25" (82.5 mm)
54201C	3.5" (88.9 mm)	54202C	3.5" (88.9 mm)
54201D	3.75" (95.2 mm)	54202D	3.75" (95.2 mm)
54201E	4" (101.6 mm)	54202E	4" (101.6 mm)
54201F	4.25" (107.9 mm)	54202F	4.25" (107.9 mm)
54201G	4.5" (114.3 mm)	54202G	4.5" (114.3 mm)
54201H	4.75" (120.6 mm)	54202H	4.75" (120.6 mm)
54201I	5" (127 mm)	54202I	5" (127 mm)
54201J	5.25" (133.3 mm)	54202J	5.25" (133.3 mm)
54201K	5.5" (139.7 mm)	54202K	5.5" (139.7 mm)
54201L	5.75" (146.8 mm)	54202L	5.75" (146.8 mm)
54201M	6" (152.4 mm)	54202M	6" (152.4 mm)
54201N	6.25" (158.7 mm)	54202N	6.25" (158.7 mm)
54201O	6.5" (165.1 mm)	54202O	6.5" (165.1 mm)
54201P	6.75" (171.4 mm)	54202P	6.75" (171.4 mm)



Aid-Grip

Cushioned retention pads are designed to help hearing instruments stay in place. Folding "WINGZ"[™], on either side of the AID-GRIP, offer hearing instrument wearers additional stability and comfort, while providing excellent retention.

- Medical-grade, hypo-allergenic adhesive
- Fits most Open-Fit, RIC, and BTE instruments

42433	Aid-Grip, 30pk	\$9.95
-------	----------------	--------



3M[®] Clear Contour Strips

Clear Contour Strips help keep BTE instruments in place. Place the micro-thin, double-sided, medical-quality adhesive strips between the BTE and the head. Contour Strips are particularly useful for patients with soft ear cartilage. Measures: 3" x .875" (76 mm x 22 mm)

42431	Contour Strips - 12pk	\$2.75
42432	Contour Strips - 36pk	\$7.50



Comply Soft Wraps

Comply Soft Wraps are made of a self-adhesive sound-blocking foam that wraps around and adheres to the hearing instrument or BTE earpiece.

42441	Soft Wraps - 10pk	\$9.50
-------	-------------------	--------



Ear Suspenders

Ear Suspenders®

Ear Suspenders are used to secure personal hearing aids and help prevent their costly loss or misplacement. Headbands are great for all ages. The headbands are a great way to reduce infants pulling and removing hearing aids. They are a great accessory for children and adult sports and recreation activities. They also work great for geriatric patients.

Each Headband features two silicone sleeves, a silicone elastic band and the ability to adjust to the head size of the wearer. Each headband comes with a spare sleeve and can be altered to fit unilateral hearing aid wearers.

Adult	\$18.50	Child	\$18.50	Under 2 Years	\$18.50
58000	Adult, Light Blue	58001	Child, Light Blue	58002	Under 2, Light Blue
58003	Adult, Neon Pink	58004	Child, Neon Pink	58005	Under 2, Neon Pink
58006	Adult, Black	58007	Child, Black	58008	Under 2, Black
58009	Adult, Navy	58010	Child, Navy	58011	Under 2, Navy
58012	Adult, Light Pink	58013	Child, Light Pink	58014	Under 2, Light Pink
58015	Adult, Grey	58016	Child, Grey	58017	Under 2, Grey



Light Blue



Neon Pink



Black



Navy



Light Pink




Grey

Oto-Clips

BTE OtoClip® and BTE OtoClip II

OtoClips are designed to help prevent the loss of expensive BTE hearing instruments. The cord clips to the user's clothing and comes with an adjustable loop that goes over the BTE hearing instrument earhook. Color-coded clips identify right side (red) and left side (blue) fittings.

54211	Right - 7.5" (1 cm)	\$2.90
54212	Left - 7.5" (1 cm)	\$2.90
54231	15" (38 cm)	\$4.85



BTE OtoClip II with Bolo

This nylon cord is designed to help prevent the loss of expensive BTE hearing instruments. The cord clips to the user's clothing and comes with an adjustable loop that goes over the BTE hearing instrument earhook.

55020	Black, 15" (38 cm)	\$6.75
-------	--------------------	--------



ITE OtoClip® and ITE OtoClip II

Designed to work with ITE instruments; ITE OtoClips use a spring-loaded jeweler's clasp to secure the cord to a manufacturer-supplied eyelet in the hearing instrument faceplate. Color-coded clips identify right side (red) and left side (blue) fittings.

54251	Right - 7.5" (1 cm)	\$3.75
54252	Left - 7.5" (1 cm)	\$3.75
54271	15" (38 cm)	\$6.50



Eyelit Kit with Drill

Small eyelets for installation in the faceplate of an ITE hearing instrument. The eyelet kit contains five eyelets and a drilling tool. Perfect for use with all ITE OtoClips

42121	Kit with drill, 5pk	\$10.50
42125	Eyelit, 1pc	\$0.85



OtoClips for Glasses

Designed for use with eyeglasses. The BTE OtoClip has a much shorter cord that attaches the BTE hearing instrument to eyeglasses.

55411	BTE, 3" (76 mm)	\$3.00
55421	ITE, 3" (76 mm)	\$3.50



Color BTE and ITE OtoClip

A colorful woven nylon cord designed to help prevent the loss of expensive BTE hearing instruments. The cord clips to the user's clothing and comes with an adjustable loop that goes over the BTE hearing instrument earhook.

BTE - Monaural, 7.5" (19 cm) \$3.75

54911	Black & White
54912	Blue
54913	Green
54914	Orange
54915	Pink
54916	Purple
54917	Yellow
54918	Red
54919	Black

BTE - Binaural, 15" (38 cm) \$6.75

55011	Black & White
55012	Blue
55013	Green
55014	Orange
55015	Pink
55016	Purple
55017	Yellow
55018	Red
55019	Black

ITE - Monaural, 7.5" (19 cm) \$4.50

55111	Black & White
55112	Blue
55113	Green
55114	Orange
55115	Pink
55116	Purple
55117	Yellow
55118	Red
55119	Black

ITE - Binaural, 15" (38 cm) \$8.00

55321	Black & White
55322	Blue
55323	Green
55324	Orange
55325	Pink
55326	Purple
55327	Yellow
55328	Red
55329	Black



Accessorizing

Collar Clip

Collar Clips are an excellent cord attachment option for healthcare facilities. Ideal for users who don't need to remove their ID cards for scanning or swiping.

55630 Collar Clip \$1.00



Critter Clips®

Critter Clips bring out the child in all of us. They work just like our standard OtoClips. Each clip features a hand-painted character with a color-coordinated cord; color of cord is not an option. Available in monaural and binaural styles for BTE and ITE hearing instruments.

BTE - Monaural, 7.5" (19 cm) \$8.00 ITE - Monaural, 7.5" (19 cm) \$8.75

54311	Cool Cat, Yellow Cord	54411	Cool Cat, Yellow Cord
54312	Croc O'Dile, Green Cord	54412	Croc O'Dile, Green Cord
54313	Dogzilla, Orange Cord	54413	Dogzilla, Orange Cord
54314	Hon E. Bear, Purple Cord	54414	Hon E. Bear, Purple Cord
54315	Wabbit, Black & White Cord	54415	Wabbit, Black & White Cord

BTE - Binaural, 15" (38 cm) \$9.75 ITE - Binaural, 15" (38 cm) \$11.25

54321	Cool Cat, Yellow Cord	54421	Cool Cat, Yellow Cord
54322	Croc O'Dile, Green Cord	54422	Croc O'Dile, Green Cord
54323	Dogzilla, Orange Cord	54423	Dogzilla, Orange Cord
54324	Hon E. Bear, Purple Cord	54424	Hon E. Bear, Purple Cord
54325	Wabbit, Black & White Cord	54425	Wabbit, Black & White Cord



Cool Cat



Croc O'Dile



Dogzilla



Hon E. Bear



Wabbit

Dino Clips

Inspired by the creatures children seem to love the most, Dino Clips work just like our standard OtoClips. Each clip features a hand-painted prehistoric character with a color-coordinated cord; color of cord is not an option. Available in monaural and binaural styles for BTE and ITE hearing instruments

BTE - Monaural, 7.5" (19 cm)		\$8.00	ITE - Monaural, 7.5" (19 cm)		\$8.75
54511	Bronto, Blue Cord		54611	Bronto, Blue Cord	
54512	Rex, Orange Cord		54612	Rex, Orange Cord	
54513	Steg, Pink Cord		54613	Steg, Pink Cord	
54514	Trixie, Purple Cord		54614	Trixie, Purple Cord	

BTE - Binaural, 15" (38 cm)		\$9.75	ITE - Binaural, 15" (38 cm)		\$11.25
54521	Bronto, Blue Cord		54621	Bronto, Blue Cord	
54522	Rex, Orange Cord		54622	Rex, Orange Cord	
54523	Steg, Pink Cord		54623	Steg, Pink Cord	
54524	Trixie, Purple Cord		54624	Trixie, Purple Cord	



Bronto



Rex



Steg



Trixie

Sea Clips

Sea Clips work just like our standard OtoClips. Each clip features a hand-painted character and a color-coordinated cord; color of cord is not an option. Available in monaural and binaural styles for BTE and ITE hearing instruments.

BTE - Monaural, 7.5" (19 cm)		\$8.00	ITE - Monaural, 7.5" (19 cm)		\$8.75
55511	Orkie, Blue Cord		55611	Orkie, Blue Cord	
55512	Oscar, Blue Cord		55612	Oscar, Blue Cord	
55513	Sharkie, Blue Cord		55613	Sharkie, Blue Cord	
55514	Skipper, Blue Cord		55614	Skipper, Blue Cord	
55515	Snappie, Blue Cord		55615	Snappie, Blue Cord	

BTE - Binaural, 15" (19 cm)		\$9.75	ITE - Binaural, 15" (19 cm)		\$11.25
55521	Orkie, Blue Cord		55621	Orkie, Blue Cord	
55522	Oscar, Blue Cord		55622	Oscar, Blue Cord	
55523	Sharkie, Blue Cord		55623	Sharkie, Blue Cord	
55524	Skipper, Blue Cord		55624	Skipper, Blue Cord	
55525	Snappie, Blue Cord		55625	Snappie, Blue Cord	



Orkie



Oscar



Sharkie



Skipper



Snappie

Lubricants, Drops & Creams

Oto-Ease® Earpiece Lubricant

Inserting earpieces is made easy with the new and improved Oto-Ease. Oto-Ease is a sterile, unscented lubricant designed specifically for custom earpieces. The unique formula provides a non-greasy, water-soluble solution for patients while also helping create an effective acoustic seal.

Note: Not recommended for someone with a perforated ear drum or with patent PE tube

54011	.5oz Display Box, 12pk	\$36.50
54013	.5oz Bottle, 3pk	\$9.50
54015	.5oz Bottle, 1pc	\$3.20
54016	.5oz Bottle, 5pk	\$15.50
54017	Single Use Packets, 100pk	\$15.50
54018	Mini - 5ml Bottle, 3pk	\$8.75



54015



54017



54018

Miracell™ Pro Ear™

Miracell ProEar is a unique blend of highly concentrated natural plant extracts. Relieves itchy ears and acts as a lubricant. Also softens cerumen, so ears can clean themselves naturally.

54025	Bottle - .5oz, 1pc	\$7.50
54026	Display Box - .5oz Bottle, 12pk	\$84.00



OtoFerm Comfort Cream

Silicon-based lubricant. Improves patient comfort and provides an effective acoustic seal.

54021	Cream - 5g Tube	\$6.00
-------	-----------------	--------



ComfortEAR X™ Moisturizing Ear Solution

ComfortEAR X is physician formulated, non-drug, hypoallergenic moisturizing ear solution designed to restore essential moisture to dry, irritated ears through a rich, emollient blend of natural Oat and Aloe Extracts. ComfortEAR X is safe to use on a daily basis and safe for all ages.

54030	ComfortEAR X, 14ml	\$7.00
-------	--------------------	--------



Eargene®

Relieves the discomfort of itchy ears through a unique formula of soothing, cleansing and healing ingredients. Safe with any hearing instrument or earpiece.

54022	Bottle - 1oz, 12pk	\$52.00
54023	Display Box - .5oz, 24pk	\$78.14
54024	Bottle - 1oz, 1pc	\$6.00



Mack's® ProRinse® Earwax Removal System

Mack's ProRinse Earwax Removal System was designed for occasional use as an aid to soften, loosen and remove excessive earwax.

- Includes:
- FDA approved Carbamide Peroxide drops
 - AquaBlock earplugs
 - Steady-flow syringe and rinse tip
 - Ear wash rinse tub

54172	Earwax Removal System	\$7.95
-------	-----------------------	--------



Audiologist's Choice® Earwax Removal System

Everything a patient needs to start an effective cerumen management program.

- Includes:
- .5 oz. cerumen removal drops
 - Soft rubber ear wash syringe
 - Complete instructions

54151	Starter Kit	\$3.75
54155	Refill Drops - .5oz Bottle, 1pc	\$2.50



Audiologist's Choice® Anti-Itch Cream

Provides temporary relief from the itching associated with minor skin irritations and rashes. The 1% hydrocortisone formula is the maximum strength allowed without a prescription.

54171	Cream - 1oz Tube	\$3.95
-------	------------------	--------



Mack's® Dry-n-Clear® Ear Drying Aid

Dry-n-Clear dries excess water from the ear canal. Just a few drops of this unique moisturizing agent will help relieve the discomfort caused by water-clogged ears.

54166	1oz Bottle	\$3.95
-------	------------	--------



Mack's® Ear Dryer

Dry your ears naturally with warm, soothing air. Invented by an ENT doctor to safely and effectively dry the ear canal where bacteria and fungi can grow.

- | | |
|--|---|
| Doctor recommended for: | Includes: |
| <ul style="list-style-type: none"> • Swimming • Bathing • Water sports • Scuba • Hearing aid use and more | <ul style="list-style-type: none"> • Folding ear dryer • 4 Color- codes eartips • Cordless - Uses AA batteries (not included) • 1 Year warranty |

54164	Replacement Ear Tips, 4pk	\$3.95
54167	Dryer Only	\$44.95
54169	Value Pack	\$49.95



54169

Displays

Model Ear

An anatomically correct model of the human ear, great for display or for educating patients during consultation. Shows a cross-sectional view of the ear from the pinna to the auditory nerve.

- Includes:
- Labeled placard
 - Simulated marble base

77111	Adult - 3.75" x 5" x 6.50"	\$45.00
77112	Child - 4.75" x 5" x 6.50"	\$79.00



Clear Display Ear

An anatomically correct model of the ear that features a beige pinna and a clear concha bowl and ear canal. The model is an attractive display piece that can also be used to explain some of the procedures patients will experience.

20213	Display, Left Ear Only	\$68.50
-------	------------------------	---------



Clear Ear Model

This model is designed to aid viewing of the entire ear including the ear canal, middle ear components such as the tympanic membrane and auditory ossicles, and inner ear components such as the semicircular canals and Cochlea. All three bones that make up the auditory ossicles are represented in clear detail; malleus, incus, and stapes.

77113	2.75" x 2.25" x 2.5"	\$45.50
-------	----------------------	---------



White Silicone Display Ear

Use anatomically accurate external ear models to demonstrate to patients what hearing instrument fittings might look like in their ear. Available in white only.

74511	Left Ear	\$26.50
74512	Right Ear	\$26.50



Acrylic Ear stand with Ear

Simple yet elegant, this eye-catching clear acrylic display stand is ideal for displaying earpieces and hearing instrument fittings in waiting rooms, display windows, trade shows and information fairs. Each stand comes with a white silicone display ear.

75111	Small, 5.75" H	\$51.50
75112	Medium, 7" H	\$53.00
75113	Tall, 8" H	\$55.00



75113 75112 75111

Acrylic Head Display with Ear

This clear acrylic display stand is ideal for displaying earpieces and hearing instrument fittings in waiting rooms, display windows, tradeshow and information fairs. Each Acrylic Head Display includes a white silicone display ear. Right side only. Height 12.5" (31.7 cm)

75211	Display with Ear	\$99.00
-------	------------------	---------



Sennheiser Set 900

The Set 900 is both a mobile sound booster and an infrared stereo listening system. The Set 900 receiver offers listeners individually optimized sound in the theaters or concert halls using infrared (IR) transmission. At home, the Set 900 can be directly connected to the TV, stereo system or radio and allows listeners to enjoy cordless sound with personal hearing settings.

Features

- Easy to connect and to use
- Excellent stereo sound, individually adjustable
- 3 Selectable hearing settings on the transmitter
- Optimized speech intelligibility
- Infrared audio transmission for cordless mobility of up to 40' (12 m) within a room
- 2 Switchable microphones
- Foldable receiver with carrying case
- Amplified stereo sound when you are outside and in rooms equipped with infrared technology (e.g. theater)
- Large, ergonomic volume control and balance control
- Precisely fitting, flexible ear cushions
- Receiver with automatic on/off function
- Lithium-polymer rechargeable battery for 12 hours of use
- 2 Year warranty

61228 Set 900 \$260.00



Sennheiser Audioport A200™

A personal sound amplifier combined with the convenience of a single headset.

Features:

- 50 dB maximum acoustic amplification
- Charging compartment for a second battery
- Rechargeable battery with 12 hours operating time
- 2 Year warranty

61351 Audioport A200 \$185.00



Sennheiser EZT-3012 Neckloop

Induction neckloop (without clip) for use with the Set 830 S or Set 840 S receiver. It transmits the audio signal from the receiver to hearing aids with a telecoil ("T") switch. 3.5 mm plug.


61227 Neckloop \$45.00



Sennheiser EZT-1011 Induction Neckloop

Induction neckloop for use with the Set 250-J, Set 810-S or Set 820-S infrared receiver. It transmits the audio signal from the receiver to hearing aids with a telecoil ("T") switch. 3.5 mm plug.

61265 Induction Neckloop \$50.00



Sennheiser EZT-3011 Induction Neckloop

Induction neckloop (with clip) for use with the Set 830 S or Set 840 S receiver. It transmits the audio signal from the receiver to hearing aids with a telecoil ("T") switch. 3.5 mm plug.

61266 Induction Neckloop \$42.99



Serene® TV Direct TV95

High definition sound technology converts TV audio output into an infrared light beam for wireless transmission to receiver. This allows the user to clearly hear all dialogues and stereo sounds in a movie or TV program.

Features:

- 115 dB SPL output
- Wireless headset, 3.0 oz.
- Volume and tone controls
- Transmitter/charging stand can hold 2 headsets (headset/neckloop)
- Works with TV volume (normal, muted or off)
- Works with most TVs, plasma and LCD televisions
- Covers up to 900 sq. ft. (83.61 m²)
- 3.5 mm headset expansion jack for accessories
- 1 Year limited warranty

62176 TV Direct TV95 \$65.00



Serene® TV Direct Earbuds

Replacement earbuds for the Serene TV Direct TV95

62177 Replacement Ear Buds, 1pair \$15.50



MKE 800-TV Microphone

For use with TVs without audio out ports.
Works with the Set 810 , Set 810-S, Set 820 and Set 820-S. 2.5 mm plug

61342 Microphone \$42.99



MKE 800-TV-N Microphone

For use with TVs without audio out ports.
Works with the Set 900 , Set 840, Set 840 S, Set 830, Set 830 S, IS 410 and Set 50.

61343 Microphone \$49.99



Westone® Cell Phone Adapter Kit

Cell phone adapter kit, 6pk, size 1-6

62410 Adapter Kit \$6.75



BA-90 Battery for Set 100

BA-90 rechargeable battery is an OEM replacement battery for use with the DirectEar Set 100.

61152 Rechargeable Battery \$15.25



BA-151 Rechargeable Battery

BA-151 battery is an OEM replacement battery. Works with Sennheiser wireless systems: Set 250, Set 810 S, Set 820(S), Set 50 and Audioport A200

61251 Rechargeable Battery \$15.75



BA-300 Rechargeable Battery

BA-300 battery is an OEM replacement battery. Works with Sennheiser wireless systems: Set 830, Set 830 S, Set 840, Set 840 S, Set 900 and IS 410

61252 Rechargeable Battery \$40.00



Acousticon Personal Sound Amplification Device

The Acousticon PSAP is used to increase sound pressure levels in the ear while also lowering background noises to prevent disturbing howling or whistling due to feedback between receiver and microphone.

62465 Sound Amplification Device (PSAP) \$150.00



VitaSound TV Hearing Helper

VitaSound TV is the TV Hearing Helper. Powered by the patented Neuro-Compensator® technology, the VitaSound TV delivers crisp and natural sounds, clear voices, reduces background noise, and amplifies what's important.

Features:

- Compatible with most wired earbuds and headphones
- 33 ft. listening range
- Digital 2.4GHz RF connection
- Up to 40 dB amplification
- Connect to hearing aids with an induction neckloop
- Rechargeable lithium-ion battery
- Link up to 100 handsets
- 1 Year warranty

Includes:

- PAE-200 Handset
- PAE-200 Base unit transmitter
- Earbuds with case
- AC/DC USB adapter
- USB Cable
- Analog (RCA) AV cable
- Lanyard
- Ear tips (pair of large, regular and small) with case

62317 TV Hearing Helper - PAE-200 \$179.00



VitaSound Audio Enhance

The Personal Audio Enhancer makes life sound better. This multi-functional audio enhancement device features four sound modes (watch, talk, listen and relax) which allow users to hear better in everyday listening situations.

- Connects to TV through a wireless connection
- Powered by Neuro-Compensator® patented hearing technology
- Sleek, lightweight design
- Comes with high-end earphones

62318 Charcoal Gray \$230.00

62319 Royal Blue \$230.00



Alerting Devices

Sonic Alert® Large Display Dual Alarm Clock with Super Shaker™

Large display with dual alarm allows two different alarm times to accommodate two different sleep wake patterns.

Features:

- Equipped with dual alarm for multiple alarm setting capability
- 113 dB extra-loud alarm (with adjustable tone and volume control)
- SS12V Super Shaker bed vibrating unit (female adapter)
- Extra large display
- Hi/low dimmer switch
- 9 V battery backup (battery not included)
- Dual power capability – 110 / 220 V (requires optional transformer)
- Dual time (12/24 hour)
- 1 Year Sonic Alert blast-proof warranty

62327 Alarm Clock with Shaker - SBD375SS \$51.00



Sonic Alert® Alarm Clock with Telephone Signaler and Super Shaker™

Combines the Sonic Alert Alarm Clock/transmitter with the Sonic Alert receiver/vibrator. The alarm clock features a snooze alarm, adjustable volume, and pitch control. The powerful receiver/vibrator can be placed under a pillow or mattress. Plug a phone into the clock and the alarm will sound every time the phone rings!

Features:

- 113 dB extra-loud alarm (with adjustable tone and volume control)
- Built-in Telephone Signaler
- SS12V Super Shaker bed vibrating unit (female adapter)
- Built-in flashing alert lights
- 9V battery backup (battery not included)
- Hi/low dimmer switch
- 1 Year Sonic Alert blast-proof warranty

62330 Alarm Clock with Signaler - SBT425SS \$38.00



Sonic Alert® Super Shaker™ Bed Vibrator

The Super Shaker unit is placed under the pillow or between a mattress and box springs.

Features:

- Powerful vibrations wake up even extremely heavy sleepers
- Built-in temperature sensor protects the unit against overheating
- Cord end has female adapter
- Cord length is 6' 3" (190.5 cm)
- 1 Year Sonic Alert warranty

62332 Bed Vibrator - SS12V \$27.95
 62332W Bed Vibrator, White - SS12VW \$28.00



Sonic Alert® Alarm Clock with Super Shaker™

Combines the Sonic Alert alarm clock/transmitter with the Sonic Alert receiver/vibrator.

The alarm clock has a snooze alarm, tilting display for easy reading, large colorful buttons, adjustable volume, and pitch control. The powerful vibrator can be placed under the pillow or between the mattress and box spring.

Features:

- 113 dB extra-loud alarm (with adjustable tone and volume control)
- SS12V Super Shaker bed vibrating unit (female adapter)
- Built-in receiver (for other Sonic Alert signalers)
- Lamp connection
- Built-in flashing alert lights
- 9V battery backup (battery not included)
- 1 Year Sonic Alert blast-proof warranty



62331 Alarm Clock with Shaker - SB100SS \$64.00

Sonic Alert® Combination Alarm Clock

Both an alarm clock and a receiver for Sonic Alert transmissions.

Features:

- 113 dB extra-loud alarm (with adjustable tone and volume control)
- Adjustable tone and volume control
- Built-in receiver (for other Sonic Alert signalers)
- Lamp connection
- Built-in flashing alert lights
- 9V battery backup (battery not included)
- 1 Year Sonic Alert blast-proof warranty



62333 Alarm Clock - SB1000 \$46.00

Sonic Alert® Sonic Shaker Travel Alarm Clock

Powerful bed shaker and/or extra-loud pulsating alarm. Includes: batteries, a pillow strap with clasp, and a protective travel case.

Features:

- 90dB extra-loud alarm
- Sonic Shaker™ bed vibrating unit
- Pop-up display with convenience light
- Test button to test batteries, vibrator, and sound operation
- Easy-to-set function keys
- Easy-to-read .75" (19 mm) display
- Continuous 4 minute snooze setting
- Swivel and flush-mounted display protects settings
- 1 Year Sonic Alert blast-proof warranty



62334 Travel Alarm Clock - SBP100 \$38.50

Sonic Alert® Bedside Travel Alarm Clock with Super Shaker™

Bedside travel alarm clock with extra-loud alarm, a large display, and vibrating bed unit.

Features:

- 113 dB extra-loud alarm (with adjustable tone and volume control)
- SS12V Super Shaker™ bed vibrating unit (female adapter)
- 9V battery backup (battery not included)
- Dual power capability – 110/220 V (requires optional transformer)
- Hi/low dimmer switch
- Dual time (12/24 hour)
- 1 Year Sonic Alert blast-proof warranty

62337 Bedside Travel Clock - SB200SS \$32.00



Sonic Alert® Sonic Bomb Alarm Clock with Super Shaker™

Extra-loud 113 dB alarm is sure to wake even the soundest sleeper. Includes pulsating lights and bed vibrating unit.

Features:

- 113 dB extra-loud alarm (with adjustable tone and volume control)
- SS12V Super Shaker™ bed vibrating unit (female adapter)
- Built-in red flashing alert lights
- Hi/low dimmer switch
- 9V battery backup (battery not included)
- Dual power capability – 110/220 V (requires optional transformer)
- Dual time (12/24 hour)
- 1 Year Sonic Alert blast-proof warranty

62338 Alarm Clock with Shaker - SBB500SS \$37.50



Sonic Alert® TR75 Ring Elite

Telephone ring signaler plugs into any modular phone jack and power outlet. Features a lamp on/off button for normal use and an indicator light. Compatible with other Sonic Alert receivers.

62533 Ring Elite \$42.50



DB200 Wireless Doorbell Signal

Wireless doorbell and intercom signaler can be installed in minutes.

Features:

- Different flash codes for front and rear doors
- Optional chime
- Outlet for lamp attachment
- Sends a signal to all Sonic Alert remote receivers
- Requires 9V battery (not included)
- Telephone signaler
- 5 Year warranty

62544 Doorbell Signaler - DB200 \$84.95



Serene® CentralAlert™ Notifications System

The CA360 is an all-in-one notification system for your home. It uses a bright built-in flasher, a selectable loud yet pleasant audible sound, and a powerful bed shaker to alert you. There is no installation needed for the CA360 and all operations are wireless.

Features:

- Includes base/alarm, bed shaker, and doorbell
- Dual alarm settings and dual bed shaker connections
- Big bright icons for each alarm
- A wide selection of affordable accessories
- 2.4 GHz transmission
- 200' operating range

62346 CentralAlert CA-360 \$130.00




Serene® Vibrating Remote Receiver Serene® Vibrating Remote Receiver

The portable CA-RX vibrating remote receiver keeps you informed when you are away from the CA-360. It alerts you to all the same sounds and activities monitored by the CA360 so you will never miss anything anywhere in your home and great for hotel and other travel uses.

Features

- Lightweight and portable so you can stay informed anywhere in your home
- Adjustable loud audio sound and bright flasher to alert you
- Big and easy to see alerting icons
- "Recall" button to check and review the last activated alert
- Output jack for bed shaker to alert you at night
- Universal AC adapter powered (100-240 VAC) with battery backup (batteries not included)

62347 Vibrating Remote Receiver - CA-RX \$46.50



Serene® Audio Alarm Sensor

This Serene sensor detects the alarm sounds of virtually all your audio smoke, gas, and fire alarms, even the new alarms with T3 or T4 type alarm sounds, and sends a wireless signal to trigger the CA-360 to alert you. Install it next to your existing alarms for sound detection or connect it directly to your alarms for activation.

Features

- Mount it next to your audio alarms for alarm sound detection. It is the only sensor that will respond to all alarm sounds, including those alarms with T3 or T4 smoke alarm sounds.
- It can also be hardwired to your audio alarms that comes with connecting terminals for direct activation
- It sends a wireless signal to activate the CA360 to alert you.
- Battery powered (batteries not included) with low battery indicator

62348 Audio Alarm Sensor - CA-AX \$37.50



Serene® Baby Cry Sensor

Monitors the sound of your baby on a real time basis and it sends a wireless signal to trigger the CA360 to alert you. It stops when your baby stops making any noise. With adjustable sensitivity, there will be no more false alarms; a peace of mind at last.

Features

- Place it near your baby and when it detects a loud baby sound, it sends a wireless signal to trigger the CA360 to alert you
- It monitors the baby sound in real-time, it stops when the baby stops making noise
- Visual loudness indicators to help you gauge and set the desired sound sensitivity to minimize false alarms

62349 Baby Cry Sensor - CA-BX \$37.50



Serene® Door Sensor

install device on existing door chime or intercom door announcer for vibration detection. When someone rings your door bell, it detects vibration and sends a wireless signal to trigger the CA360 to alert. Or, you can attach the CA-DX to the inside surface of your door to detect the knocking on your door.

Features

- Detects the vibration of your existing door bell, door chimes, door intercom announcer or the knocking on your door, it then sends a signal to trigger the CA-360 to alert you
- Simple to install with adjustable sensitivity to avoid false alarms
- 2.5mm Jack for direct connection to your existing doorbell/door chime.

*With this method, electrical signal from your doorbell terminal activates CA-DX

62350 Door Sensor - CA-DX \$37.50



Serene® Storm Alert

Detects the weather warning signals from your NOAA radio and sends a signal to trigger the CA360 to warn you.

Features

- Connects directly to your NOAA radio output jack (radio and connecting cable not included)
- When warning signal is detected, it sends a wireless signal to trigger the CA-360 to alert you that there has been a weather or storm warning issued by NOAA

62351 Storm Alert - CA-NOAA \$37.50



Serene® Emergency Button

Alerts other family members in the home when you are in need. It is compact and waterproof so you can keep it with you in every part of your home.

Features

- Sends an signal to the CA-360 to notify others in your home
- Lightweight and wearable with any lanyard
- Waterproof design for installation in the shower

62352 Emergency Button - CA-SOS \$26.00



Serene® Wireless Doorbell

The CA-DB is a weather resistant and easy to install door bell for your back or side doors.

Features:

- Sends wireless signals to let you know that someone is at the door
- Weather resistant for easy outdoor mounting

62353 Wireless Doorbell - CA-DB \$26.00



Sleeping Devices

Sound Oasis® Bluetooth Sleep Sound Therapy System

Now with twice as many sounds, our most popular compact sound machine is loved for a reason – a small size that creates a robust and full sound experience. Perfect for travel or your night stand, this Bluetooth speaker and sound machine can create the perfect white noise experience at home or on the road.

63030 BST-100 Therapy System \$59.99



Sound Oasis® Smart Sleep Mask

Our modern lifestyle often prevents us from receiving natural cues that aid in a healthy sleep cycle. The illumy smart sleep mask is a comfortable and lightweight smart sleep mask that blocks out unwanted light and provides its own light to trigger your body's natural responses. Blue tones wake you gently and effectively in the morning before your chosen alarm time, and red-toned light helps you fall asleep at night. These are the same light cues used by NASA to regulate sleep on the space station.

63031 Sleep Mask \$99.00



SleepPhones®

Traditional headphones made bed friendly! Our patented "headphones in a headband" design secures ultra-thin stereo headphones inside a soft, stretchy headband. SleepPhones plug into or sync with most media devices for hours of relaxing audio without disturbing who's next to you. Color: Quiet Lavender.

Features:

- Get to sleep faster; stay asleep longer and wake up refreshed
- Reduce ambient noise, conversations, traffic sounds, even snoring
- Removable speakers and machine-washable (hypoallergenic) headband
- Designed for all night comfort and fit
- Adjustable speakers
- No ear buds or clunky ear pads



63015	Classic, Midnight Black	\$30.00
63016	Classic, Quiet Lavender	\$36.00
63017	Classic, Galaxy Blue	\$30.00

63020	Wireless, Midnight Black	\$60.00
63021	Wireless, Quiet Lavender	\$90.00
63022	Wireless, Galaxy Blue	\$90.00

CLASSIC



63017

WIRELESS



63021

PocketTalkers

Williams Sound® PockeTalker® Ultra

The PockeTalker Ultra features a lightweight, ergonomic design for portability and easy use.

Features:

- 20 40dB acoustic gain
- External tone and volume controls
- Removable microphone with extension cord
- Compatible with neckloops for telecoil
- Belt clip
- Lanyard
- Requires (2) AAA batteries
- Up to 100 hours of battery life
- 5 Year warranty

62112	with Right Earphone, EAR 008	\$105.00
62114	with Dual Earbuds, EAR 014	\$105.00
62115	with Deluxe Headset, HED 021 & 013	\$105.00
62116	with Rear-ear Headphone, HED 026	\$105.00
62117	with 18" (45 cm) Induction Neckloop, NKL 001	\$120.00
62118	Amplifier Only, PKT D1-0	\$92.00

62114



62115



62116



62117



WHILE SUPPLIES LAST!

Williams Sound® PockeTalker® PRO 1

Versatile assistive listening device that can serve as a personal listener, television, or telephone amplifier when properly equipped. Available with three choices of earpieces.

Compact enough to bring anywhere and is ideal for one-on-one conversations, small groups, and much more. The Pocketalker Pro amplifies sound clearly and easily for better understanding, and 100 hours of battery life ensures long-lasting performance.

62131	PRO 1-2 with monaural earbud, EAR 013	\$105.00
62141	PRO 1-3 with binaural over-the-head headset, HED 021	\$105.00



Williams Sound® PockeTalker® Ultra and PRO 1 Accessories

62211	Single Earbud, EAR 013	\$11.00
62215	Dual Earbud, EAR 014	\$14.00
62220	Sanitary Headphone Covers, EAR 045 - 100pk	\$40.00
62221	Replacement Earbud Pads - 10pk	\$4.00
62225	Deluxe Headset, HED 021	\$11.00
62226	Rear Wear Headphone - HED 026	\$13.50
62231	Windscreen - WND 002	\$2.95

62235	Neckloop - NKL 001	\$35.00
62241	Plug Mount Microphone, MIC 014	\$31.00
62244	Directional Mini Lapel Microphone, MIC 054	\$65.00
62245	Microphone Extension Cord, WCA 007	\$9.00
62255	Rechargeable Battery Kit - BAT KT3	\$21.00
62261	Wide Range Earphone - EAR 008	\$10.50
62265	Replacement Headset Pads, HED 021/026 - 1pair	\$1.50



62211



62215



62220



62221



62225



62226



62231



62235



62241



62245



62255



62261



62265

Williams Sound® PockeTalker® 2.0

The Williams Sound PockeTalker 2.0 personal amplifier is a slim, lightweight design that fits easily in your hand.

Features:

- Up to 63dB amplification
- Tone control
- Balance control
- Omnidirectional internal microphone
- External omnidirectional plug-in microphone
- T-coil receiver for looped environments
- 3.5 mm stereo headphone jack
- Up to 105 hours of battery life
- 5 Year warranty

Includes:

- Personal amplifier, PKT E1
- Dual stereo isolation earphones, EAR 042
- Stereo headphones, HED 024
- Plug-in microphone, MIC 014-R
- (2) AAA batteries
- Microphone extension cord, WCA 007 WC, 12'
- Belt clip
- Lanyard

62152	PockeTalker® 2.0	\$163.00
-------	------------------	----------

Accessories

62153	Rechargeable Battery Kit, BAT KT7	\$37.50
62154	Dual Stereo Isolation Earphones, EAR 042	\$16.50
62155	Stereo Headphones, HED 024	\$11.65
62241	Plug Mount Microphone, MIC 014	\$31.00
62245	Microphone Extension Cord, WCA 007	\$9.00



Comfort Duett® Personal Listener

Comfort Duett is a small, easy to use, rechargeable, assistive listening device with superior sound. Use it with earphones in place of a hearing aid, or use it with an induction neck-loop, to amplify conversations and sounds from a variety of sources like TV, radio, telephone, etc. The Duett is discreet (resembling a small MP3 player), stylish and the perfect solution for hearing impaired individuals struggling to hear properly in various situations.

Features:

- Operating time: 15-20 hours, depending on volume level
- Charge time: 1 hour charge = 1 hour use
- Operating spectrum: 275 Hz - 9000 Hz
- Max. amplification: 60dB, with max gain and bass reduction controls
- Weight including battery: 1.75 oz.
- Size: 3.375" x 1.5" x 1.125"
- Color: Metallic blue
- 2 Year limited manufacturer's warranty

62670	with Neckloop	\$160.00
62671	with Headphone and Earbud	\$125.00



Telephones & Accessories

ClearSounds® Big Button Phone

The CSC600 Amplified Big Button eco friendly phone, with an amplified speakerphone, features exceptional amplification and tone management technology.

Features:

- Up to 50 dB amplification
- Maximum 130 dB SPL amplification
- 30 dB amplification volume control
- 95 dB adjustable extra loud ringer
- Bright visual ringer
- Built-in outgoing speech amplification (8 dB)
- ClearDigital™ Full Digital Power
- Hi/Med/Lo/Off ringer volume control
- Caller ID
- Last number redial
- 5 One-touch speed-dial
- One-touch emergency dial button
- Big button dialing with backlit keypad
- Uses 50% less power than other amplified phones
- Optional bed shaker with jack
- IonArmour Natural Antimicrobial Protection
- Hearing instrument compatible
- 2 Amplified jacks (2.5 mm hands-free, 3.5 mm audio)
- 1 Year warranty



62632 Big Button Phone - WCSC600 \$140.58

ClearSounds® Bed Shaker

ClearSounds Vibrating Pad connects to any ClearSounds amplified phone, amplified ringer, or digital clock. Features a 3.5 mm plug.



62627 Bed Shaker \$18.95

ClearSounds® Amplifier

The WIL95 is designed for use with any ordinary corded phone.

Features:

- Compatible with most corded digital, analog, and VOIP phones
- 30 dB maximum amplification
- Volume and tone controls
- Boost button engages amplifier
- Powered by AC adapter (included)
- Optional 9V battery (not included)
- 1 year limited warranty



62622 Amplifier - WIL95 \$26.00

ClearSounds® Amplified Power Neckloop

One device – multiple applications: cellular phones, home/office phones, laptop computers, portable gaming devices, personal/TV listening equipment, and other audio devices. Reduces or eliminates the interference between cellular phones and hearing instruments or cochlear implant processors.

Features:

- 30 dB adjustable amplification
- Built-in microphone and quick disconnect cable
- Requires telecoil-equipped hearing instruments
- Requires (2) AAA batteries
- 1 Year warranty

62637	CLA7v2 Amplified Neckloop	\$57.50
-------	---------------------------	---------



Pocket Telephone Amplifier

Works with all telecoil hearing instruments and features 28 dB HF average full gain. The Pocket Telephone Amplifier is designed to increase receiver volume and is ideal for use both at home and while traveling. Includes (2) AAA batteries.

62464	Pocket Telephone Amplifier	\$11.50
-------	----------------------------	---------



Telephone Comfort Pads

Soft cushion helps eliminate feedback, reduce room noise, and increase the intelligibility of phone conversation for users with hearing instruments.

Black, Round

54061	Black, Round - 12pk	\$10.75
54062	Black, Round - 2pk	\$2.50

Black, Square

54063	Black, Square - 2pk	\$2.50
54065	Black, Square - 12pk	\$10.75

White, Square

54064	White, Square - 2pk	\$2.50
54066	White, Square - 12pk	\$10.75

Black, Round



Black, Square



White, Square

Cell-U-Hear Cell Phone Pads

Foam comfort pads help direct sound to the hearing instrument microphone, which helps prevent feedback. Each package contains two different size pads that fit most cell phones.

54067	Pad - 2pk	\$2.25
-------	-----------	--------



Bluetooth Devices

Clarity® HA-40 In-Line Phone Amplifier

Compact, portable in-line amplification compatible with most home and office telephones.

Features:

- 40 dB maximum amplification
- Volume control
- Tone control
- Works with digital and analog corded phones
- Boost button mutes microphone in the headset to eliminate background noise
- 9 V battery (included)
- 1 Year limited warranty



62426 Amplifier \$21.50

ClearSounds® Quattro 4.0

Connect and stream up to five Bluetooth devices to your hearing aids with the Quattro 4.0 Adaptive Bluetooth Listening System. The system enables you to stream audio from your devices (phone, tv, computer, etc) directly to your hearing aids or cochlear implants (with t-coils) via Bluetooth just by wearing the Neckloop.

Includes:

- Quattro 4.0 Amplifier
- USB charging cable with AC adapter
- Earbuds
- Lanyard and magnetic clip/charger for detachable mic



Amplifier

62462 Quattro 4.0 \$279.99

ClearSounds® QLink™ Bluetooth Transmitter

The ClearSounds QLink™ is an advanced Bluetooth adapter that allows the user to link a variety of sound devices with the Quattro Neckloop/Headset system.

Features:

- Supports TV's, iPods®, MP3 players, media players, laptops, and other audio devices
- Bluetooth profile A2DP V2.1+EDR, Class II
- Frequency range, 2.4 GHz Spectrum
- Includes stereo RCA to stereo 3.5 mm cables
- Provides up to 10 hours battery life
- 5x longer than battery life than v1.2 devices
- Charging time 2-3 hours, with full functionality
- Range of up to 33' (10 m)
- 1 Year warranty




62463 Bluetooth Transmitter \$45.00

Books & Reference Materials

Five Keys to Communication Success Book

78123	Book, Spanish	\$14.95
78124	Book, English	\$12.95

A Tool to Educate & Motivate



www.5keys.org


"At last, a book that truly tells it like it is. ...Offers constructive, realistic solutions discussed in an environment-specific format."

Robert W. Sweetow, PhD
University of California


5 Keys to Communication Success

A helpful handbook of real-life solutions for:

- New hearing aid fittings
- Tested-not-sold patients
- Family members
- Physician outreach




also available in Spanish



Hearing Loss and Hearing Aids- A Bridge to Healing Book

Written for patients with hearing loss and their families, this handbook addresses how hearing loss affects the patient and those around them. This resource is perfect for patients who want to learn more about their hearing loss.

78131	Book	\$9.99
-------	------	--------



How Hearing Loss Impacts Relationships Book

Families learn how to manage the resistance of a loved one facing hearing loss with the aid of this valuable resource. By understanding the roots of their resistance, families can help the hearing-impaired individual move from "struggling to hear" to "hearing independence".

78132	Book	\$12.95
-------	------	---------



Hearing Aids: A User's Guide Book

Potential and first-time patients can realize the full physical and psychological benefits of hearing instruments with the help of this informative text. Topics covered include: How Hearing Works, What To Expect from a Hearing Instrument, Care and Maintenance, Troubleshooting and much more.


78135 Book \$7.60



Mayo Clinic on Hearing Book

From the Mayo on Health series, this concise, complete and easy-to-understand book includes sections that explain hearing loss and the treatment options that are currently available. Many photographs and illustrations are included, as well as a glossary and list of additional resources.


78139 Book \$9.99



I'm the Boss of my Hearing Loss Book

Using colorful drawings and a positive approach, this handbook is designed for kids with hearing loss and their parents. Kids learn important concepts about hearing loss as well as how to manage challenging listening situations. This book uses lively and amusing illustrations the whole family will love.


78165 Book \$12.95



Children with Hearing Loss Book

For parents, siblings and extended family members who want a better understanding of the impact hearing loss can have on children. This book is a guide on how to address the most important educational issues through the school years, the profound emotional impact hearing loss can have on a child, as well as some of the pitfalls common among families new to a child with hearing loss. It also covers the latest technology available to children and explains how families can give hearing impaired children the proper tools that will empower them for life.


78125 Book \$13.25



Infection Control in the Audiology Clinic Book

Audiologists and students can use this resource for learning the latest infection control theories, standards and practical ready-to-use protocols. Special features in this book include infection control issues related specifically to cochlear implants and cerumen management, a complete section of frequently asked questions and a detailed index.


78183 Book \$35.00



Infection Control in the Hearing Aid Clinic Book

An outstanding reference for hearing healthcare professionals practicing in hearing instrument clinics. This book includes a review of the most common infections and infection control protocols. Infection control plans include forms useful for designing custom infection control procedures. The extensive glossary is an indispensable resource.

78182 Book \$25.00



Hear the Music by Marshall Chasin Book

Learn about hearing loss exposure for musicians and strategies to conserve their hearing. This is an excellent reference for anyone involved with music on a regular basis, rock 'n' rollers, orchestra musicians, even school band members and their teachers, as well as the hearing healthcare professionals who work with them.

78185 Book \$9.95



Cerumen Management - Book and DVD

A must-have book and DVD set for audiologists and AuD students focusing on cerumen removal. Features six chapters of information including anatomy of the ear canal, light sources, equipment, instrumentation, cerumen management and infection control. The accompanying DVD demonstrates proper cerumen removal technique.

78179 Book & DVD \$74.00



Environmental Sound Simulator

Program, verify and counsel patients on their aided-hearing performance in real-world environments without leaving your office. The simulator uses audio and full motion video to provide a realistic, virtual fitting environment. High-quality video of typical listening environments such as churches, restaurants, the outdoors and small groups, help maximize your patients hearing performance.

78201 Sound Simulator \$399.00



QuickSIN™ Speech-in-Noise Test CD

Quickly and easily measures the ability to hear noise in typical listening environments. Useful for counseling patients regarding realistic expectations. Results can be used to make specific treatment recommendations. An excellent addition to the audiometric test battery for all adult patients.

78205 Test CD \$145.00



Posters

Your Ears Hearing and Balance Chart

One of the wall charts in a series by Kate Sweeney entitled Blueprint for Health. Contains a wealth of interesting information for all ages, ranging from basic human ear anatomy to how a cricket hears. Plastic lamination on this colorful 20 x 26 wall chart provides a tear-resistant, markable surface.


77511 20" x 26" Poster \$17.50



Organs of Hearing and Balance Poster

Bright, informative, and anatomically correct, these high-quality posters can be hung in an examination area as a reference or to help explain procedures to patients. Lithographed in authentic colors and laminated with a plastic coating, each 20" x 26" poster includes eyelets for convenient hanging.

77311 20" x 26" Poster \$17.50



Middle Ear Conditions Poster

Bright, informative, and anatomically correct, these high-quality posters can be hung in an examination area as a reference or to help explain procedures to patients. Lithographed in authentic colors and laminated with a plastic coating, each 20" x 26" poster includes eyelets for convenient hanging.


77313 20" x 26" Poster \$17.50



Ear, Nose & Throat Poster

Bright, informative, and anatomically correct, these high-quality posters can be hung in an examination area as a reference or to help explain procedures to patients. Lithographed in authentic colors and laminated with a plastic coating, each 20" x 26" poster includes eyelets for convenient hanging.

77314 20" x 26" Poster \$17.50



Apparel

Westone Since 1959 T-shirt

Westone's First to Innovate 1959 shirt is a soft charcoal grey that comes in both men and women's extra-small to 2XL sizes. 60% Cotton and 40% Polyester blend. Shirts tend to run small, so consider ordering a size up!

Men's

70002	Extra Small	\$10.00
70003	Small	\$10.00
70004	Medium	\$10.00
70005	Large	\$10.00
70006	Extra Large	\$10.00
70007	Extra, Extra Large	\$10.00



Women's

70008	Small	\$10.00
70009	Medium	\$10.00
70010	Large	\$10.00
70011	Extra Large	\$10.00
70012	Extra, Extra Large	\$10.00



Westone® Wristband

Westone's silicone wristband comes in grey with Westone's logo displayed in white. Only one size available that fits most wrists.

70013 Wristband \$1.50



Monitor Vault Belt Loop Carrier

Stay connected with your DefendEar Digital Modules. The Monitor Vault Belt Loop Carrier offers peace of mind, loss prevention and convenience so you can stay focused on what you love to do most.

78306 Belt Loop Carrier \$3.50



Westone Lanyard

Lanyard with Westone logo, clips onto Plastic Polycon Cases

42337 Westone Lanyard \$3.00



Hearing Protection Device (HPD) Zipper Pouch

Zippered Pouch measures 4 x 3 inches

42300 Pouch - 1 pc \$3.50



Westone® Squeeze Pouch

Perfect for storage of custom-fit earpieces, hearing protection devices, or other small items. Custom imprinting is also available (call for details).

42339 Pouch - 1 pc \$1.75



The *Ear* and *Earpiece* Guide according to Westone

When discussing critical information about anatomical features and fittings for an earpiece or impression, it helps if the hearing healthcare professional and the lab share common terminology. You may find commonly used terms, illustrations and definitions in this guide useful when communicating with Westone about the design of your patient's earpiece.

Ear Anatomy

The anatomical landmarks described here are of primary importance to the hearing healthcare provider. Several of these hold special significance for the creation of a proper fitting earpiece.

Antihelix

A ridge of cartilage that is part of the pinna, just above the crus that the earpiece fits under

Antitragus

A ridge of cartilage opposite and slightly below the tragus

Aperture

The opening of the ear canal

Auricle or Pinna

The external ear

Concha

The bowl-shaped cavity just outside of the opening of the ear canal

Crus (ear) valley (earpiece)

The horizontal piece of cartilage located outside the ear canal that divides the upper and lower parts of the ear

Ear Canal

The 2 to 3 cm long external auditory meatus, consisting of an outer cartilaginous portion and a inner bony portion

Heel

The intertragal notch located at the bottom of the concha bowl

Helix Lock (curl)

A recessed "crevice" at the upper part of the ear, just above the valley

Helix Extension (not including lock)


An extension up to but not including the Helix lock


Tragus

A small ridge in front of the external opening of the ear canal


Canal Length

Since each earpiece is as unique as the patient's ear, discussion of canal length can be a bit confusing.





At Westone, we use the anatomical features of each individual ear to describe canal length. Short canals are defined by the center of the first bend of the ear canal, while long canals are centered at the second bend. Using this as the basis for canal length, a medium canal would be cut halfway between the 1st bend and the 2nd bend of the ear, while a medium long canal would be cut halfway between a medium and long canal.



Impressions

The ear impression is the blueprint from which any custom-fit product will be produced, so the finished product can only be as good as the impression. Any deficiencies in the impression will greatly affect the patient's satisfaction and use of the finished earpiece.

The impression materials listed below are formulated by Westone to be the finest available to our industry. When used according to proper practices and procedures, our material will help you provide the best impression possible for your patient.

NOTE: Impressions that are deficient in critical areas will require a new impression to be submitted before your order can be processed.

Impressive!

For years, we've made what our customers regard as the finest silicone ear impression material on the market. Today our material is formulated to reveal every detail and feature of the ear, and provide the best possible blueprint for any of our custom-fit products. They're still pink, and as always, the best impression material money can buy.



IMPRESSION MATERIALS

Numerous tools and techniques have been used throughout the years to make impressions of the human ear. Recent advances in materials and procedures have led to the most consistent, dimensionally stable, and easy-to-use impression materials and delivery systems to date. Even so, only strict attention to the manufacturer's instructions will ensure the best impressions possible.

There are three primary families of impression material available: two-part vinyl polysiloxane, addition-cure silicone, and methyl methacrylate polymer and monomer (powder and liquid).

VINYL POLYSILOXANE

The latest advancement in impression material is available as individual portion packs, in specially formulated cartridges, or in bulk containers. Vinyl polysiloxane is the most popular impression material in use today. Westone offers several different vinyl polysiloxane products to meet your needs:

Silicone Singles® are available as pre-measured units, making them extremely easy and convenient to use.

SiliClone®, SiliClone FIRM, and SiliClone X-act Cartridges are conveniently packaged for use with manual or electric impression guns. When propelled through the special mixing tip, the two-part material is automatically mixed to the proper consistency.

Silicast® is a two-part material packaged in bulk containers. Accurate measurement of the material is ensured by using the two different color-specific scoops provided.

Silicone Singles and Silicast can be mixed either by "Spleading" (recommended, page 12) or by hand.

To mix using the "Spleading" technique, use a metal spatula or similar tool and a SPLEAD pad. Mix the materials as quickly as possible until the material is a uniform color, then load the syringe and make the impression. It is very important not to deviate from the manufacturer's instructions.

To mix by hand, thoroughly wash your hands and then combine the material for 20 to 30 seconds or until the material is a uniform color. Work as quickly as possible to mix the material, load the syringe, and make the impression.

All vinyl polysiloxane impression materials provide accurate and dimensionally stable ear impressions.

ADDITION-CURE SILICONE

Blue Silicast® is a base and catalyst material that consists of a base putty and a tube of hardener. A scoop is provided for easy measuring. Blue Silicast can be mixed either by “Spleading” (recommended, see page 12) or by hand. As with any impression material, it is very important not to deviate from the manufacturer’s instructions.

To use the “Spleading” technique, fill the scoop with impression material. Flatten the impression material onto the SPLEAD pad and add the catalyst. Thoroughly mix the material until it’s uniform in color. Once the two components come into contact with one another, the curing process has begun. Work as quickly as possible to mix the material, load the syringe, and make the impression.

To mix by hand, thoroughly wash your hands and then fill the scoop with the recommended amount of impression material. Flatten the impression material in the palm of your hand and add catalyst. Mix the material for 20 to 30 seconds or until the material is a uniform color. While it may be tempting to speed up the cure of the impression material by using a little extra catalyst, the minute that may be saved in cure time is not worth the risk of a possible ill-fitting earpiece. Remember that once the two components have come into contact with one another, the curing process has begun. Work as quickly as possible to mix the material, load the syringe, and make the impression.

Once the impression has cured, it is dimensionally stable. Place the impression in the mailing box along with the order form and ship.

POWDER AND LIQUID

Westone’s BLEND® impression material has been available for several decades and produces an accurate impression.

When using powder and liquid impression material, remember:

- Always use pre-measured material. A slight change in the ratio of liquid to powder can produce an impression that is either too wet, which can lead to drooping, or too dry, which can cause over-expansion of the ear canal.
- Always put the liquid into the mixing bowl first, followed by the powder.
A glass mixing bowl and metal spatula are a must. Paper cups and tongue depressors absorb liquid, thereby altering the critical ratio of powder to liquid.
- The curing process starts as soon as the powder and liquid come into contact. It is imperative that you work quickly while mixing the material and loading the impression syringe.
- Powder and liquid impressions are susceptible to damage during shipping. As such, they must be carefully packaged to avoid distortion. Secure the impression to the bottom of the shipping box with an adhesive, such as Duco® Cement, and keep the order form from contacting the cured impression.
- BLEND impressions cannot be scanned.

Impression Tools

OTO-DAMS

All oto-dams work on the same basic principle. An oto-dam is manufactured to a specific diameter from a material that has an inherent resistance to reduction in diameter. As the oto-dam is placed in an ear, it is compressed by the canal walls. The resistance of the oto-dam to this compression, or what is called “radial pressure,” holds the oto-dam in place against the force of the impression material. The two types of oto-dams in common use today are vented and non-vented.

VENTED OTO-DAMS

With any conventional oto-dam (foam or cotton), once the impression material has made full contact with the ear canal wall, any air trapped behind the dam is forced to compress against the eardrum. This sensation of ear “fullness” can range from minimal to significant discomfort. A vented foam oto-dam allows for the equalization of air pressure within the ear during the impression making, curing, and removal process, which can significantly reduce discomfort.

NON-VENTED FOAM OTO-DAMS

Non-vented foam oto-dams are available in a wide variety of sizes and are soft, comfortable, and easy to use. Size selection is the most common challenge in using a foam oto-dam. Be sure to choose the appropriate size for your patient, using visual inspection to ensure that it fits snugly in the ear canal.

NON-VENTED COTTON OTO-DAMS

Non-vented cotton oto-dams are soft, comfortable, and function well. However, because of the manufacturing process, many cotton oto-dams are supplied as hard, compressed wads, so be sure to fluff the cotton oto-dam before use. The softer, more pliable the oto-dam, the more effectively and comfortably it will seal off the ear canal.

IMPRESSION SYRINGES

Historically, the syringe method has been used to take impressions of the human ear, offering clean, quick, and precise control of the impression material as it enters the patient’s ear. Syringes are available in several designs to work with specific impression materials.

IMPRESSION GUNS

Manual or electric impression guns effectively eliminate the guesswork associated with the mixing of materials by using dual-chamber cartridges and special mixing tips. The precise amounts of “A” and “B” components, propelled through the tip, ensure proper mixing prior to delivery into the ear.

Caution:

Proper bracing of the tools being used is of paramount importance to the safety of the patient. During the impression process, some patients may experience a cough or gag reflex. This is a normal response due to stimulation of the vagus nerve. To avoid injury, always brace your hand to the patient’s head by using your little finger. This will keep the tool in the same relative location to the ear in the event the patient coughs, jerks or moves suddenly.



Impression Techniques

OTOSCOPIC EXAMINATION

Before making an impression, the patient's ear must be thoroughly examined to determine if it is acceptable for the impression process. The exam not only gives you an overall view of the health of the ear, but it can also reveal concerns or contraindications to impression making such as:

- Impacted cerumen (earwax)
- Fluid discharge
- Bony growths (exostosis)
- Inflammatory conditions
- Prolapsed canal
- Perforated eardrum
- Foreign objects
- Enlarged canal
- Thick hair growth
- Cauliflower ear
- Swimmer's ear
- Surgical ear
- Foreign objects
- Enlarged canal

EARLITE USE

The earlite serves as an aid in oto-dam placement. When using earlites and otoscopes, it is always important to hold the tool properly and to use the proper bracing technique.

OTO-DAM PLACEMENT

Proper oto-dam placement is critical to ensuring the safety of your patient and ensuring a successful ear impression. The primary acoustic seal is established in the ear canal and is complicated by the fact that the ear canal actually changes dimension whenever the jaw moves.

Use the otoscopic examination to gauge the size of the oto-dam you will use for your patient. An oto-dam with too small a diameter allows the impression material to flow beyond the oto-dam and impinge on the eardrum. An oto-dam with too large a diameter prevents proper placement in the ear canal at the depth required for a good impression. It can also cause discomfort for your patient. Apply Oto-Ease® lubricant to the oto-dam to make insertion easier and more comfortable. If the size you have chosen for the oto-dam does not seal the ear canal, remove it and try a different size. Once you have selected the proper size, use the tip of the earlite to place the oto-dam into the ear canal and to make any slight adjustments. Use the otoscope to verify that the oto-dam is placed properly in the ear canal, checking to see that it is in full contact with the entire canal wall.

OPEN VERSUS CLOSED JAW IMPRESSIONS

There has been much debate regarding open versus closed jaw impression techniques.

Bite block/open jaw impressions should be taken for hearing protection products, musicians' products, and whenever acoustic feedback is the primary concern.

IMPRESSING THE EAR


Using the proper bracing technique, place the end of the mixing tip or syringe approximately 1/4 inch (6 mm) inside the ear canal opening. Gently begin to express the impression material into the ear canal, allowing the material to flow back over the syringe or mixing tip. Once the material starts to flow past the tip, start backing the tip out of the ear canal. It is critical that you fill all the landmarks in the ear. For instance, fill the concha bowl, move up into the antihelix crevice, and then follow the contour of the ear, finishing at the center of the concha.

After about 3 to 5 minutes, check to see if the impression has cured. Press an earlite tip gently into the material to see if an indentation remains. If there is no indentation and the material is not tacky, the impression is cured. It is then safe to remove.

REMOVING THE IMPRESSION

The removal string or vent tube of the oto-dam should not be used to aid in removal of the impression. Instead, gently pull the patient's pinna up and back to break the seal. Remove the impression from the antihelix area and gently rotate forward and out. Should the oto-dam remain in the ear canal, the removal string or vent tube can be used to remove it.

Examine the ear for any material remaining from the impression process. Some slight redness may be evident and is normal. In the rare case that any impression material remains and cannot be easily removed, refer the patient to the appropriate medical authority for treatment.



Note: Prior to sending the patient home, inspect the impression to make sure that the necessary canal length and all anatomical landmarks are present.


Splead Technique

“Spleading” is a technique for mixing ear impression materials that reduces the chance of contamination. We’ve found that products containing sulfur—often found in latex gloves and hand lotions—can react with vinyl polysiloxane impression materials, preventing proper setup. “Spleading” is recommended for Silicone Singles and other silicone materials. All you need is a Westone SPLEAD pad, metal spatula, silicone syringe, and Westone impression material.

Easy Steps for Spleading:

1. Place equal amounts of Part “A” and Part “B” impression material on the SPLEAD pad.
2. Mix the two materials together by blending them with the spatula. Don’t touch the materials with your hands! Continue mixing until the material is a uniform color.
3. To load the silicone syringe, remove the plunger and place the chamber of the syringe on the SPLEAD pad at a slight angle. Using the syringe chamber like a cookie cutter, press down on the material, progressively scooping up the material until all of it is loaded into the syringe.
4. Put the plunger back into the syringe to complete the loading process. Tear off and discard the used sheet of the SPLEAD pad before mixing material for the next impression.

Note: Check Westone’s website for our latest impression technique videos.



Earpiece Materials

Westone offers a range of earpiece materials, including digitally manufactured and traditional acrylic, several types of silicone, and injection molded vinyl. Our materials are available in a wide variety of both natural skin tones and vibrant colors.

TIPS & GUIDELINES for EARPIECE MATERIAL SELECTION

Many decisions must be made in order to successfully fit your patients with hearing instruments and earpieces. The following are simple guidelines based on Westone's more than 50 years of experience making custom-fit products for the ear.

Discomfort

If you suspect that discomfort might be an issue, there are a few things to consider before selecting an earpiece material. First, check your patient's skin texture. If the skin texture is firm, consider using a soft earpiece material. If your patient's skin is soft or delicate, as can often be the case with elderly patients, try using a firmer earpiece material, or a silicone earpiece with Slick-Sil AM coating.

The helix area of an earpiece can sometimes cause problems. Many earpiece styles do not require a helix to ensure a secure fit. You can reduce potential problems by ordering a 3/4 shell or half-shell style earpiece. This is a particularly good choice for patients who may experience sensitivity in the antihelix area.

Loose Fit

If loose fit is a concern, request that the earpiece be built up for a tighter fit. If you are using an acrylic material, consider changing to a softer material as well.

Allergies

Use OtoBlast clear silicone earpiece material (specify "no glaze/matte finish") if allergies are a concern. Since earpieces made from silicone can be abrasive to delicate skin, particularly in elderly patients, request our Slick-Sil AM earpiece coating to ensure ease of insertion. For minor itching with minimal skin irritation, you can order clear acrylic or clear vinyl and request that we boil the earpiece in a saline solution. Polyethylene earpiece material is available for extreme allergy situations.

Mandibular Action

Use the softest earpiece material possible for fitting a patient with severe mandibular action. Consider tapering the canal or specify a KC (hollow) style that will allow the earpiece to flex with jaw movement.

Feedback

If you suspect that feedback could be a problem, use a bite-block, deep-impression technique. Consider using a silicone material (with matte finish) if the patient's skin texture will allow.

ACRYLIC MATERIALS

Acrylic is Westone's standard hard earpiece material. It is an ideal material for most hearing healthcare needs as well as many occupation/communication applications. Acrylic is available in a wide variety of colors, ranging from traditional clear and natural skin tones to some of the most vivid colors imaginable, including our unique Swyrl color combinations.

FLEX CANAL

Flex Canal is a body temperature-reactive, semi-soft canal additive available on all standard Westone acrylic earpieces. The unique properties of this material allow the canal portion of the earpiece to remain firm at room temperature for ease of insertion and then soften at body temperature for increased comfort and seal. The flex canal option is also available in combination with the e-Compound canal option (see below).

e-COMPOUND ADDITIVE

e-Compound canal additive is designed to reduce feedback by diffusing reflected sound between the tip of the earpiece and the eardrum. e-Compound is available on all standard acrylic earpieces, but it can also be ordered in combination with the flex canal option.

ACRYLIC EARPIECE FINISHES

Westone's acrylic earpieces are available in a high gloss shiny or matte finish.

SHINY

The high gloss shiny finish looks great on more vivid colors such as our custom colors, neon colors, acrylic Glitter, and Swyrl earpieces. Not available in FIT.

MATTE

The smooth matte finish helps hide the earpiece in the ear when using clear or one of Westone's many skin tone or DisappEar color choices.

Acrylic Standard



Acrylic Custom



Acrylic Glitter



Acrylic Neon



Acrylic DisappEar



Acrylic Canal Additives



SILICONE MATERIALS

OtoBlast™

OtoBlast* is Westone's premium silicone earpiece material. It is soft, durable, and available in a myriad of bright vivid colors or simply crystal clear. OtoBlast options also include beautiful skin tones, our popular DisappEar colors, Glitter, Cat Eyes, and of course Swyrl.

*Westone OtoBlast material compatibility tested according to MIL-STD-810F by an independent NIST-NVLAP / OSHA NRTL / DSCC / NCC CB / A2LA accredited laboratory. OtoBlast has passed testing according to ISO Standard 10993: Biological Evaluation of Medical Devices Part 10, Tests for Irritation & Delayed Hypersensitivity Biological Evaluation of Medical Devices Part 5, Tests for In Vitro Cytotoxicity

W1

W-1 is our softest silicone material. It is available in water clear, six skin tones, DisappEar ("A" and "B" tints only), sixteen colors, Cat Eyes, as well as Glitter.

AquaNot™

AquaNot is the ultimate floatable swim plug material. It's 100% American made silicone and is available in thirteen bright colors. AquaNot material can also be mixed (up to three colors) to create our unique Swyrl finish.

e-COMPOUND ADDITIVE

e-Compound canal additive is designed to reduce feedback by diffusing reflected sound between the tip of the earpiece and the eardrum. e-Compound can be added to the canal portion of OtoBlast, and W-1 earpieces.

SLICK-SIL® AM (Anti-Microbial)

Slick-Sil AM is a chemically bonded coating for silicone earpieces that is proven to be 99.9% effective* against bacteria, germs, and other pathogens. As a result, Slick-Sil AM can be part of an aggressive infection control protocol. Slick-Sil AM produces a beautiful velvet finish that not only looks great but also helps with the insertion of silicone earpieces by providing a low friction surface finish. Available in a matte finish only. Not available with DisappEar skin tone additive.

*Independent testing performed October 23, 2009, in compliance with JIS Z 2801 standards.

SILICONE EARPIECE FINISHES

Westone's silicone earpieces are available in two different finishes - shiny and matte.

SHINY

Shiny is a high gloss glaze that looks great, especially with vivid colors such as our AquaNot/Swyrl, neon colors, Glitter, or Cat Eyes earpiece options.

MATTE

A smooth matte finish helps hide the earpiece in the ear when using clear or one of Westone's many skin tone color choices.

OtoBlast **Standard**



OtoBlast + W1 **Cat Eyes** *Up to 2 colors can be added to Cat Eye*



CUSTOM-FIT EARPIECE GUIDE

OtoBlast **Glitter**



OtoBlast **DisappEar™** *SlickSil and shiny finish not available with DisAppear Colors*



DisappEar™

An earpiece made with Westone's remarkable DisappEar skin tone additive virtually vanishes when inserted in the ear. Our unique manufacturing process embeds microfibers in the material that work just like camouflage, breaking up the earpiece's outline and reflection. DisappEar is available in three different colors, pink (A), beige (B), and brown (C), to blend with almost any skin tone, in your choice of hard acrylic or soft silicone materials. Not available with Slick-Sil®.

OtoBlast **Camo**



AquaNot™ **Standard**



Swyrl it all together!

Can't decide on your favorite color? You can swyrl 2-3 silicone colors together for a completely custom look!



NOTE: DisappEar, Smoke, Rose, and Transparent Brown are not available to swyrl.

W1 Standard *All W-1 earpieces are matte finish*



W1 Glitter



W1 DisappEar




SlickSil not available with DisAppear Colors

OtoBlast + W1 Canal Additives



SlickSil and shiny finish not available with e-Compound



What's growing on your patient's earpieces?

It's a fact – hearing instrument earpieces are breeding grounds for bacteria, fungi, and other germs. The solution is Slick-Sil AM (Anti-Microbial), a chemically-bonded coating for silicone earpieces, proven to be 99.9% effective against bacteria, germs, and other pathogens. Slick-Sil AM also provides a low-friction surface for easier insertion of silicone earpieces.

VINYL MATERIALS

FORMULA II (F-2)

Formula II is a soft vinyl material. F2 comes standard in light tan but is also available in light, medium, and dark brown.

FORMULA II CLEAR (F-2C)

Formula II Clear is the same material as our Formula II but with no color additive. F2C is also available in a choice of nine translucent color tints.

FORMULA II (F-2C) NEON COLOR

Formula II Neon is the same material as our Formula II. It is available in four bright colors: neon pink, neon green, neon orange, and neon yellow.

VINYL MARBLE

Vinyl Marble is the same material as our Formula II but with a white base material swirled with your choice of one of these colors: blue, pink, yellow, orange, purple, green, or red.

SUPERFLEX

Superflex is softer than Formula II but firmer than the Rx material. Superflex is available only with a semi-transparent pink tint.

Rx

Westone's Rx is the softest vinyl material available. This material is ideal for difficult fits when a silicone material might not be an option. Rx comes standard in a medium tan color but is also available in a light, medium, or dark brown.

VINYL EARPIECE FINISHES

Westone's vinyl earpieces are available in a shiny or matte finish (default is matte).


SHINY

The shiny finish looks great, especially with our neon and translucent color tints or vinyl marble color option.

MATTE

A smooth matte finish helps hide the earpiece in the ear when using clear or one of Westone's many skin tone color choices.

Formula II (F2) **Standard**



Superflex



Formula II Clear (F2) **Custom Tints**




Formula II (F2) **Neon**



Formula II (F2) **Marble**



RX



Note: Color samples depict actual earpiece colors as closely as the printing process will allow. Please be aware that neon colors are very difficult to reproduce. Neon earpieces will be significantly brighter and more vivid than can be shown here.

Earpiece Venting Guidelines

Venting typically affects frequencies below 1000 Hz. In general, the larger the vent, the more low-frequency energy is “bled off.” By using different vent diameters, varying amounts of low-frequency amplification can be reduced. Earpiece venting can range from totally open or non-occluding to very precise and interchangeable venting via a vent control insert such as a Mini Vent Plug (MVP).

Open earpiece fittings are earpieces with very large vents. These types of earpieces can be quite effective for hearing loss that does not require low-frequency amplification. They allow the patient to perceive low-pitched sounds naturally (through the open vent), while using an instrument for high-frequency amplification. Reduction or elimination of the occlusion effect is an added benefit for the patient.

VENTING OPTIONS

	VENT SIZE	INCHES	MILLIMETERS
◦	Pressure	.040"	1.0 mm
○	Standard	.063"	1.6 mm
○	Medium	.081"	2.0 mm
○	Large	.096"	2.4 mm
⊙	MVP	(all sizes)	(all sizes)
⊙	SAV	(all sizes)	(all sizes)

WESTONE TUBING *Size* CHART

	TUBE SIZE	ID x OD (in.)	ID x OD (mm)
○	#12 Standard	.085" x .125"	2.16 x 3.17
○	#13 Medium	.076" x .122"	1.93 x 3.10
○	#13 Thick	.076" x .130"	1.93 x 3.30
○	#13 Super Thick	.076" x .142"	1.93 x 3.60
○	#15 Standard	.059" x .116"	1.50 x 2.95
○	#16 Thin	.053" x .085"	1.35 x 2.16
○	#16 Standard	.053" x .116"	1.35 x 2.95
○	Stethoscope Tubing	.125" x .187"	2.16 x 4.75

Tubing sizes illustrated are actual size. All tubing with the exception of Dry-Tubes are made of polyvinyl chloride material.

VENTING GLOSSARY



Angle Vent (AV)

The vent intersects the sound bore as close as possible to the tip of the canal.



Half-External Vent

A channel from the canal tip, half way down the canal length, that transitions to an internal vent. Used when drainage or discharge from the ear is a problem.



Parallel Vent (PV)*

The sound bore and the vent travel through the earpiece side-by-side and do not intersect.




Semi-IROS Vent

The vent is cut half way from the end of the canal.



Mini Vent Plugs (MVP)

Normally drilled with a 2.8 mm diameter bore. MVPs come with each of the following vent sizes:
 #1 (.100" / 2.54 mm)
 #2 (.080" / 2.03 mm)
 #3 (.060" / 1.53 mm)
 #4 (.040" / 1.02 mm)
 #5 (.030" / .76 mm)
 #6 (solid plug).



Select-A-Vent (SAV)**

Normally drilled with a 4.5 mm diameter bore. SAVs come with each of the following vent sizes:
 #1 (.125" / 3.18 mm)
 #2 (.110" / 2.79 mm)
 #3 (.080" / 2.03 mm)
 #4 (.050" / 1.27 mm)
 #5 (.025" / .64 mm)
 #6 (solid plug).

* Vents will be parallel as a default if no vent configuration is specified. If the physical size of the earpiece prevents use of a parallel vent, an angle vent will be used.

** Defaults to MVP if SAV cannot be fitted into the earpiece.



An Easy Way to Clean Vents

For those very small vents on hearing instruments or earpieces, try cleaning them with a dental floss loop (designed for cleaning braces and bridgework). We recommend NanoClean (Westone PN 50148). The stiff end can be used for pressure vents. Larger vents can be cleaned by using a wax loop or pipe cleaner.





Hearing Instrument *Earpieces*

Westone has literally made over millions of custom-fit earpieces, but even so, we approach each earpiece as if it is the most important one we've ever made.

We make earpieces for virtually every application you will encounter. If you're not certain which version might be the best choice for your patient, call us. Our technical support staff can help you select the style and material that is the most appropriate: 1-800-525-5071

Ear Openers/Stents

Our "ear opener" earpieces keep collapsing ear canals open. Sometimes the tissue of an aging or surgical ear canal becomes very soft and flaccid, to the point of losing its definition and collapsing. When this happens, the canal can become occluded, causing a conductive hearing loss. This can be corrected by opening the ear canal, using a stent made from acrylic material. "Ear opener" earpieces are available in style numbers 2, 3, 3 Inverted, 4, 5, 6, 6 Thick, 8, and 8 Thick.





Applications: (Snap ring) Body instruments, stethoscopes (with an adapt-a-scope adapter), and communications. Works with some Plantronics® and Telex® headsets

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Adapt-a-scope, male adapter, colors, skin tones, DisappEar, Glitter, venting, initials, Slick-Sil AM, Wire Retention Adapter

Finish: Shiny, matte




Applications: (Skeleton) Amplification, telecommunication, teleprompter, law enforcement, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, venting and tube options, initials, T-Mic, Slick-Sil AM

Finish: Shiny, matte



Applications: Advanced Bionics T-Mic

Available materials: Formula 2, Formula 2 Clear, RX, Superflex, OtoBlast

Options/accessories: Colors, Skin tones, DisappEar, Glitter, initials, Slick-sil

Finish: Shiny, matte




Applications: (Semi-skeleton) Amplification, telecommunication, teleprompter, law enforcement, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, venting and tube options, initials, Slick-Sil AM

Finish: Shiny, matte




Applications: (Semi-skeleton) Amplification, telecommunication, teleprompter, law enforcement, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, venting and tube options, initials, Slick-Sil AM

Finish: Shiny, matte




Applications: (Canal) Amplification, communications, teleprompter, law enforcement, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, venting and tube options, initials, Slick-Sil AM

Finish: Shiny, matte

Style No. 4
(KC)




Applications: (Canal) For cases of pronounced mandibular action and/or severe feedback

Available materials: Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, initials, venting (no MVP/SAV), Slick-Sil AM

Finish: Shiny, matte

Style No. 5




Applications: (Canal lock) Amplification, communications, teleprompter, law enforcement, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, initials, venting, Slick-Sil AM

Finish: Shiny, matte

Style No. 6




Applications: (Shell) Amplification, communications, teleprompter, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, initials, venting and tubing options, Slick-Sil AM

Finish: Shiny, matte

Style No. 6
Thick



Applications: (Thick/Flat shell) Amplification, communications, teleprompter, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, initials, venting and tubing options, Slick-Sil AM, Wire Retention Adapter

Finish: Shiny, matte

Style No. 8




Applications: (Half-shell) Amplification, communications, teleprompter, stenotic ear opener (Order "ear opener" on the comment line of the order form, available in Acrylic only)

Available materials: Acrylic, Acrylic with Flex Canal, Acrylic with hard or soft e-Compound canal, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, initials, venting and tubing options, Slick-Sil AM

Finish: Shiny, matte

Style No. 10




Applications: (Acoustic modifier) Normally used for mild high-frequency loss; comes standard with open bore and MVP (Not guaranteed against feedback)

Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, tubing options and initials, Slick-Sil AM

Finish: Shiny, matte

Style No. 12




Applications: (Non-occluding) For high-frequency losses, Tinnitus Masker, open-fit hearing instruments, CROS and BiCROS fittings, open fittings, and cochlear implants (Not guaranteed against feedback)

Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, tubing options and initials, Venting (MVP/SAV), Slick-Sil AM

Finish: Shiny, matte

Style No. 13




Applications: (Non-occluding) For high-frequency losses, CROS and BiCROS fittings, Tinnitus Masker, open-fit hearing instruments, open fittings, and cochlear implants (Not guaranteed against feedback). Available as 13edu for Edulink and Phonak iSense.

Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, tubing options and initials, Slick-Sil AM

Finish: Shiny, matte

Style No. 15



Applications: (Free field) For high-frequency losses, Tinnitus Masker, open-fit hearing instruments, CROS and BiCROS fittings, open fittings, and cochlear implants (Not guaranteed against feedback)

Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex

Options/accessories: Colors, skin tones, DisappEar, initials

Finish: Shiny, matte


RIC, RTE & CRT *Earpieces*

The sound quality and discreet size of many RIC (Receiver In Canal), RTE (Receiver In The Ear), and CRT (Canal Receiver Technology) instruments have made them extremely popular. But despite the advancements in hearing instrument technology, many of these hearing instrument fittings require custom-fit earpieces for a variety of reasons, including:

- To provide the necessary gain before feedback
- Help retain the receiver in the ear canal
- To avoid insertion difficulties and provide improved comfort

Westone custom-fit earpieces are the solution. Select from an open recessed canal style, a canal-lock style, or even a skeleton style, based on your patient's individual needs. Westone has solutions for virtually any application, which will help maximize your fittings, improve your patient satisfaction, and reduce your hearing aid returns.

Style No. 87



Applications: (Canal) RIC earpiece for Siemens™ family of receivers. Uses the Siemens™ C-Grid™ wax protection system.

Available materials: OtoBlast, acrylic


Venting Options: Pressure vent, standard vent, medium vent, large vent, MVP, SAV, largest vent possible (space permitting)

Receiver bore options: 45 dB, 55 dB, 60 dB, 65 dB, 70 dB Siemens™ family of receivers – please specify receiver dB size and hearing aid model being used

Other Options: Skin tones, DisappEar, Glitter, Canal lock (87-5), Half-shell (87-8), Skeleton (87-2), Full-shell (87-6), initials, Slick-Sil AM

Finish: Shiny, matte

Style No. 88



Applications: (Canal) Earpiece for older power receivers. For use with Phonak®, Unitron™, Bernafon®, Oticon®, Cochlear Hybrid, and similar power receivers

Available materials: OtoBlast, acrylic


Venting Options: Pressure vent, standard vent, medium vent, MVP (space permitting)

Receiver bore options: Older power receivers listed above

Other Options: Colors, skin tones, DisappEar, Glitter, Canal lock (88-5), Half-Shell (88-8), Skeleton (88-2), Full-Shell (88-6), initials, Slick-Sil AM

Finish: Shiny, matte

Style No. 89



Applications: (Canal) For use with many manufacturers' standard size receivers and newer power size receivers, including, but not limited to, Oticon™, ReSound, Phonak™, Widex, Siemens™, Unitron™, and similar receivers. Receiver friction-fits into earpiece.

Available materials: OtoBlast, acrylic

Venting Options: Pressure vent, standard vent, medium vent, large vent, MVP, SAV, largest vent possible (space permitting)

Receiver bore options: Will drill to accommodate most manufacturers' receivers – please specify manufacturer, hearing aid model, and receiver type when ordering

Other Options: Colors, skin tones, DisappEar, Glitter, Canal lock (89-5), Half-Shell (89-8), Skeleton (89-2), Full-Shell (89-6), initials, Slick-Sil AM


Finish: Shiny, matte

Slim Tube Earpieces

The acoustic and cosmetic advantages of Slim Tube/Open-Fit hearing instruments have increased the demand for this style of fitting. However, many patients are borderline candidates because of their degree of hearing loss, or the fact that stock, generic ear tips may not fit securely in their ear canals.

Westone custom-fit earpieces are the solution. Select a style that meets your patient's needs, improve your patient's satisfaction, and reduce your returns for credit.

Slim Tubes can be purchased through your hearing aid manufacturer.



Applications: (Skeleton) Features excellent retention and fitting flexibility, largest vent possible, and semi-IROS vent


Available materials: Acrylic, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Venting Options: Largest vent possible (default), SAV, MVP

Tube Options: None (Slim tube fit)

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Canal) Cosmetically appealing, largest vent possible, and semi-IROS vent

Available materials: Acrylic, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Venting Options: Largest vent possible (default), SAV, MVP

Tube Options: None (Slim tube fit)

Options/accessories
Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Canal lock) Cosmetically appealing, features good retention and fitting flexibility, largest vent possible, and semi-IROS vent

Available materials: Acrylic, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Venting Options: Largest vent possible (default), SAV, MVP

Tube Options: None (Slim tube fit)

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Shell) Features great retention for active users, largest vent possible, and semi-IROS vent


Available materials: Acrylic, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Venting Options: Largest vent possible (default), SAV, MVP, and semi-IROS

Tube Options: None (Slim tube fit)

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Half shell) Easy to insert and remove, features largest vent possible, and semi-IROS vent


Available materials: Acrylic, Vinyl, Formula II, Formula II Clear, Superflex, Rx, W-1, OtoBlast

Venting Options: Largest vent possible (default), SAV, MVP

Tube Options: None (Slim tube fit)

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Skeleton) Non-occluding, features excellent retention, largest vent possible, and minimal canal contact

Available materials: Acrylic, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Venting Options: Largest vent possible (default), SAV, MVP

Tube Options: None (Slim tube fit)

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Skeleton) Non-occluding style earpiece features excellent retention


Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Venting Options: None

Tube Options: None (Slim tube fit)

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte




Applications: Specific to the Phonak™ Nadia S BTE hearing instrument, earpiece can be ordered in any style. Includes Phonak™ Power slim tube installed in earpiece. *Please specify old or V series.*

Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex, W-1, OtoBlast

Options/accessories: Venting, Initials, Colors, skin tones, DisappEar, Slick-Sil AM

Other Options: Canal lock (81-5), Half-Shell (81-8), Skeleton (81-2), Full-Shell (81-6)

Finish: Shiny, matte



Applications: (Canal) Specific to thin tube diameter, open-fitting hearing instruments. Use where vent must be controlled (not guaranteed against occlusion)

Available materials: Acrylic, Formula II, Formula II Clear, Rx, Superflex, W-1, OtoBlast


Options/accessories: Colors, skin tones, DisappEar, venting, initials, Slick-Sil AM

Other Options: Canal lock (85-5), Half-Shell (85-8), Skeleton (85-2), Full-Shell (85-6)

Finish: Shiny, matte

OCCUPATIONAL *Earpieces*

Westone has worked closely with pilots, doctors, television personalities, and communications professionals to develop highly effective earpieces that bring new levels of comfort and efficiency to their professions. Several of our earpieces are designed specifically to work with the most popular headsets on the market.




Applications: (Thick-shell) Comes with standard snap ring. Compatible with anesthesia stethoscopes (with an adapt-a-scope adapter and scope tubing) and audiologist listening scopes. Works with some Plantronics® StarSet® and Telex® headsets

Available materials: Acrylic, Acrylic with Flex Canal, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Options/accessories: Adapt-a-scope, male adapter, colors, skin tones, DisappEar, venting, initials, wire retention adapter, Slick-Sil AM

Finish: Shiny, matte




Applications
(Thick-shell) Works with the Plantronics MS-Series and the Telex 5x5 Pro Series

Available materials: Acrylic, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Options/accessories: Filter in vent bore, colors, skin tones, DisappEar, initials, Black Aluminum Flyer Bar

Finish: Shiny, matte

Note: Comes standard with Mini Vent Plug (MVP)



Applications
(Skeleton) Communications, IFB Teleprompter, law enforcement

Available materials: Acrylic, Acrylic with Flex Canal, Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Options/accessories: Straight tubing, collar clips, initials, venting, colors, skin tones, DisappEar, Slick-Sil AM*. Available as full-shell 32-6, half-shell 32-8, canal lock 32-5, canal 32-4, 32-12 non-occluding, 32-13 non-occluding.

Finish: Shiny, matte

* Silicone only

PERSONAL *Earpieces*

For years, Westone has been on the leading edge of what we call “recreational audiology.” We have products that allow motorcyclists to run the open road, let race car drivers communicate with their pit crews, and music lovers hear exactly what they want to hear without ear-wrenching volume.

“Sounds like fun” is a catchy little phrase that we like to use to describe “recreational audiology,” but in fact, it is very serious business.

Note: Second-bend, open-mouth, bite-block impressions are recommended for all Westone personal earpieces.



Description
(Thick-shell) Racing headset, compatible with many racing communications systems, features a 18" hardwired EPIC cable and 3.5 mm mono plug Driver: Single balanced armature, full range
Sensitivity: 114 db SPL @ 1 kHz
Impedance: 25 ohms @ 1 kHz
Frequency response: 29 Hz – 16 kHz
1-year warranty, parts and labor


Attenuation: Not rated

Applications: Designed primarily to allow race drivers to communicate with their pit crews

Available materials: Formula II, Formula II Clear

Options/accessories: Colors, carrying case

Finish: Shiny



Description: (Canal) A hollow attenuator with removal filament


Attenuation: Not rated

Applications: Designed primarily for motorcyclists to wear under a helmet to help eliminate wind noise

Available materials: Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Applications: (Thick-shell) Adapts to some cell phone headsets and button style stereo headphones to custom earpieces

Available materials: OtoBlast


Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Other Options: Drill-to-fit, available as a 7-4 canal

Finish: Shiny, matte

*Headset required with order.

Note: Comes standard with a Mini Vent Plug (MVP). Only available for classic button-style stereo earbuds




Applications: (Thick-shell) Adapts to some cell phone headsets and button style stereo headphones to custom earpieces

Available materials: Acrylic, Acrylic with Flex Canal, Formula II, Formula II Clear, Superflex, Rx

Options/accessories: Colors, skin tones, DisappEar, initials

Finish: Shiny, matte

Note: Comes standard with a Mini Vent Plug (MVP)



Description: (Canal) Custom-fit earpiece for Westone universal-fit earphones and monitors

Attenuation: Not rated

Applications: Designed for use with all Westone and most Shure® universal-fit earphones and musicians' monitors

Available materials: Formula II, Formula II Clear, Superflex, Rx, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Description: (Thick-shell) Swim plug, unsurpassed for comfort and seal; comes with left or right side identifier and nub type handle


Attenuation: NRR 27 dB (average 37 dB)

Applications: Designed for use while surface swimming or showering to prevent water from entering the ear canal. Not designed for diving or for swimming deeper than six feet.

Available materials: AquaNot

Options/accessories: Various color options, Swyrl (up to three mixed colors), cords, initials, Slick-Sil AM. *NOT AVAILABLE in glitter or clear.*

Finish: Shiny, matte



Description: Thick Shell mold designed for Etymotic Research headphones


Attenuation: Not rated

Applications: Etymotic MC series headphones, ER-6/6i headphones, ER-4 headphones, and Etymotic HF series headphones

Available materials: ww OtoBlast

Options/accessories: Colors, glitter, skin tones, DisappEar, initials, Slick-Sil AM

Finish: Shiny, matte



Description: (Thick-shell) Swim plug with a sealed membrane spanning sound bore; designed to keep cold water and cold air out of the ear canal; allows basic conversation while in the ear

Attenuation: Not rated

Applications: Designed for swimming, surfing, skiing, biking, jogging. Not designed for diving or swimming deeper than 6 feet.

Available materials: AquaNot

Options/accessories: Various color options, Swyrl (up to three mixed colors), cords, initials, Slick-Sil AM. *NOT AVAILABLE in glitter or clear.*

Finish: Shiny, matte

PASSIVE *Hearing Protection*

Westone custom-fit hearing protection earpieces provide an effective, comfortable and repeatable fit. hearing protection products are suitable for a variety of applications, whether it's for work, listening to music, shooting, motor sports or doing jobs around the house such as mowing the grass or using power tools.


Tests conducted by the United States Air Force Research Laboratory suggest that participants who wore custom-fit hearing protection experienced a greater consistency of fit and attenuation than those using generic over-the-counter products.

Note: Second-bend, open-mouth, bite-block impressions are recommended for all Westone personal earpieces.


Safe Exposure Times dB*

Instantaneous permanent damage	140+	Shotgun, rifle, jet plane takeoff
Less than one second	130	Jackhammer, heavy industry
Less than ten seconds	120	Rock concert
<hr/>		
Threshold of pain	110	Power tools, snowmobile
1.5 minutes	110	Chainsaw, motorcycle
15 minutes	100	Lawn mower
2.5 hours	90	Beginning of Danger Zone
Eight hours	85	City Traffic
<hr/>		
Prolonged exposure to noise levels 85 dB and higher can result in permanent hearing loss	80	Vacuum cleaner, hair dryer
	70	Office, sewing machine
	60	Normal conversation
	50	Refrigerator
	40	Whisper
	30	Rustling leaves
	20	Breathing
	10	Threshold of normal hearing
	0	
Common noise levels (dB), and their effect upon hearing		

*Exposure time based on NIOSH occupational noise exposure standards




- 1 Maximum (Full-shell, Style No. 40)
- 2 Motorsports/Sleeping (Canal, 40-4)
- 3 Convertible (Full-shell, Style No. 47-40, closed)
- 4 Convertible (Full-shell, Style No. 47-40, open)
- 5 Filtered (Full-shell, Style No. 47-40)
- 6 Ultra Lite (Half-shell, Style No. 42)
- 7 Hunter Passive (Full-shell, Style No. 39)



Post/Nub Type
Removal Handle



AquaNot Style
Removal Handle



Style No. 40

Hunter Passive

Description: (Thick shell) Earpiece with Y-bore containing non-linear attenuating small orifice filters, for a balance of auditory situational awareness and protection against low shot count firearm noise

Attenuation: NRR 4, IPIL Rated*

Applications: Low daily shot count hunting

Available materials: Otoblast Silicone with Slick-Sil®, an anti-microbial low friction coating


Options/accessories: Colors (except DisappEar), Swyrl (up to 3 mixed colors), Laser Engraved Initials

Finish: Matte (Shiny finish not available)

Hunter Passive IPIL Results

Peak Exposure Level (dB)	IPIL (dB)
132	15.1
150	22.8
168	26.7

Style No. 47



Convertible

Description: (Thick half-shell) Filtered attenuator allows basic verbal communication while lowering overall sound levels, features a removable cord that turns the earpiece from a filtered attenuator into a solid sound plug

Attenuation: NRR 18-28 dB, average attenuation 30-39 dB depending on material


Applications: Designed for heavy equipment operators, motor sports enthusiasts, construction workers, dental hygienists, etc.

Available materials: Formula II, Formula II Clear, Rx, Superflex, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, initials, handles, Slick-Sil AM. Available as thick shell 47-40, shell 47-6, half-shell 47-8

Finish: Shiny, matte

Style No. 47



Filtered

Description: (Thick half-shell) Filtered attenuator allows basic verbal communication while lowering overall sound levels

Attenuation: NRR 18-21 dB, average attenuation 30-32 dB depending on material


Applications: Designed for heavy equipment operators, motor sports enthusiasts, construction workers, dental hygienists, etc.

Available materials: OtoBlast, Formula II, Formula II Clear, Rx, Superflex

Options/accessories: Colors, skin tones, DisappEar, Glitter, initials, cords, handles, Slick-Sil AM. Available as thick shell 47-40, shell 47-6, half-shell 47-8 or canal 47-4

Finish: Shiny, matte

Style No. 40-4



Sleep

Description: (Canal) solid attenuator with removal filament

Attenuation: NRR 24-27 dB, average attenuation 37-39 dB depending on material


Applications: Designed to promote comfortable and uninterrupted sleep. Great for shift workers and credited with saving many a relationship!

Available materials: Formula II, Formula II Clear, Rx, Superflex, W-1, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, Glitter, initials, cords, handles, Slick-Sil AM

Finish: Shiny, matte

Style No. 42



Ultra Light

Description: (Thick-shell) Full-thickness filtered attenuator

Attenuation: NRR 10 dB, average attenuation 20 dB

Applications: Designed for low-level noise exposure

Available materials: W-1, OtoBlast, Formula II, Formula II Clear, Rx, Superflex

Options/accessories: Colors, skin tones, glitter, eyelets, cords, handles, DisappEar, initials, Slick-Sil AM. Available as a shell 42-6, half-shell 42-8 or canal 42-4

Finish: Shiny, matte

Style No. 40



Solid

Description: (Thick-shell) Solid full-thickness, The Maximum offers the highest noise reduction of any Westone DefendEar hearing protection earpiece

Attenuation: NRR 24-29 dB, average attenuation 37-40 dB depending on material

Applications: Designed for shooting enthusiasts, heavy equipment operators, construction workers, etc.

Available materials: W-1, OtoBlast, Formula II, Formula II Clear, Rx, Superflex

Options/accessories: Colors, skin tones, DisappEar, Glitter, eyelets, cords, removal handles, initials, Slick-Sil AM. Available as shell 40-6, half-shell 40-8, canal 40-4 (sleep plug)

Finish: Shiny, matte

Note: If removal handles are desired on Style No. 40, they must be requested on the order form.

Style No. 49



Musicians

Description: (Canal) Earpiece provides true high-fidelity attenuation using Etymotic Research® ER-9, ER-15 or ER-25 filters

Attenuation: 9 dB ER-9 filter, 15 dB ER-15 filter, 25 dB ER-25 filter, solid adapter (not rated)


Applications: Designed for musicians, concert attendees, music teachers, DJs, flight attendants, bartenders, waiters, airline travel, subways, dental professionals, hyperacusis/recruitment

Available materials: Formula II, Formula II Clear, OtoBlast

Options/accessories: Colors, skin tones, DisappEar, solid adapters, initials, Slick-Sil AM. Glue-in cord available. Filter colors available: clear, beige, brown, red & blue.

Finish: Shiny, matte

TRU Custom



Concert

Description: (Canal) Low profile earpiece with removal filament and a selection of four interchangeable TRU filters.

Attenuation (see chart below): WM10 filter, WM16 filter, WM20 filter, WM25 filter

Applications: Designed for musicians, sound engineers, DJs, dental professionals, concerts, nightclubs, air travel, movies, sporting events

Available material: OtoBlast only

Options/accessories: Colors, skin tones, DisappEar, initials, Slick-Sil AM, glue in cord only

Finish: Shiny, matte

Product	Avg. Attenuation	NRR
WM10 Filter	10	2
WM16 Filter	12-15	10
WM20 Filter	20	13
WM25 Filter	24	15

What are the key differences between the Style No. 49 and TRU Custom?

The TRU customs are more discreet and fit more recessed in the ear. TRU customs are recommend for patients with small and/or narrow canals. The style 49's are larger and easier to insert/remove from the ear. The filters are larger in size making it simple to interchange them in the molds.

WARNING: The purpose of hearing protection is to reduce the impact of the sound and to reduce the potential for hearing loss. Hearing protection devices, including TRU products, do not eliminate sound. The user is responsible for determining whether the noise reduction provided by TRU products will provide adequate hearing protection based on the user's unique medical history and exposure level. Your medical history, particularly any prior trauma to the ear or prior hearing loss, as well as the noise (decibel level and exposure time) the user will be exposed to should be considered prior to using TRU products. If you have any questions concerning the appropriate exposure levels or the level of hearing protection provided by TRU products, or if you have hearing loss or ear injuries, contact your local hearing professional before using.


INDUSTRIAL *Hearing Protection*

The right HPD for the job. DefendEar industrial hearing protection earpieces protect hearing in any industrial environment. This carefully designed line of products allows you to offer the perfect solution for any industrial hearing users protection need.

Situational awareness is a key component of an effective noise reduction solution. Some workers must be able to understand speech as well as detect audible warnings from equipment or facilities. Over-protection can be JUST as dangerous as being under-protected.


Whatever the requirement, Westone has an earpiece that will get the job done.

- Repeatability
- Custom-fit for maximum comfort
- Situation-appropriate attenuation
- Environmentally friendly
- Lower long term use costs compared with disposable generic fit products




DefendEar Max

Description: (Thick-shell) Solid industrial hearing protection provides maximum attenuation
Includes: Initials, Oto-Ease earpiece lubricant, Oto-Clens earpiece cleaner, zippered pouch, pop cord
Attenuation: NRR 29
Applications: Designed for high-noise environments
Available materials: OtoBlast
Options: Colors, food-processing configuration, with or without metal
Finish: Slick-Sil AM



DefendEar Convertible Comm

Description: (Thick-shell) Industrial hearing protection features a communications port that works with surveillance kit-type radios (NRR 28 of available when comm. port is closed)
Includes: Oto-Ease earpiece lubricant, Oto-Clens earpiece cleaner, zippered pouch, pop cord
Also includes an adapter to attach surveillance kit sound tube to earpiece
Attenuation: NRR 21/28
Applications: Communications/hearing-protection device, ideal for industrial environments where radio communications is required
Available materials: OtoBlast
Option: Colors, initials, food-processing configuration, with or without metal
Finish: Slick-Sil AM



DefendEar Convertible

Description: (Thick-shell) Filtered industrial hearing protection, features a removable cord that allows the user to select the amount of attenuation desired
Includes: Initials, Oto-Ease earpiece lubricant, Oto-Clens earpiece cleaner, zippered pouch, pop cord
Attenuation: NRR 21/28
Applications: Great for industrial environments with a wide-ranging amount of noise exposure
Available materials: OtoBlast
Options/accessories: Colors, initials, food-processing configuration, with or without metal
Finish: Slick-Sil AM

NOTE: If you deviate from this package, standard mold charges will apply.

DIGITAL *Hearing Protection*




DefendEar®
DIGITALX™
 CUSTOM-FIT HEARING PROTECTION

DigitalX™ is the latest development in digital hearing protection and leads the trend towards obstruction free, comfortable custom noise suppression for today's shooter. DigitalX combines superb functionality with an eye catching appearance.

DefendEar Digital Series IPIL Results*

Peak Exposure Level (dB)	IPIL (dB)
132	37.3
150	44
168	42

* See the DefendEar gunshot noise white paper on our website




DigitalX 1
Description: Features a digitally optimized program for general shooting.

Attenuation: NRR 26

Applications: General Shooting for Outdoor use only

- Features:**
- Push Button Volume Control
 - Dedicated On/Off Button
 - Digital DSP Circuitry
 - Mute
 - Low Battery Warning
 - Four stage digital volume control push button
 - Single Optimized Program
 - 6 Metallic Door Colors
 - WaterShield protection
 - Size 13 Battery (new doors)
 - Size 312 Battery (old doors)
 - NRR 26
 - IPIL: 37.3/44/42
 - 2 Year Warranty




DigitalX 3
Description: Features three optimized programs ideal for either intermittent or continuous shots in addition to ambient sound enhancement.

Attenuation: NRR 26

Applications: Outdoor use only, Recreational Shooting, Clay Target, Game Hunting and Competition Shooting

- Features:**
- Digital Volume Control Switch
 - Dedicated On/Off Switch
 - Vocal Feedback
 - Startup program indicator
 - Minimum and maximum volume
 - Program selection
 - Mute
 - Low Battery Warning
 - Volume and Program Auto-Save
 - 3 Digitally Optimized Programs
 - 6 Metallic Door Colors
 - WaterShield protection
 - Size 13 Battery (new doors)
 - Size 312 Battery (old doors)
 - NRR 26
 - IPIL: 37.3/44/42
 - 2 Year Warranty



DigitalX 5
Description: Offers five modes that detect game movement while providing impulse suppression for continuous shots in addition to ambient sound enhancement.

Attenuation: NRR 26

Applications: Outdoor use only, Indoor use only with additional protection, Tactical Situations, Recreational Shooting, Clay Target, Wingshooting, Game Hunting, Competition Shooting, Windy Conditions and Large Caliber Firearms

- Features:**
- Digital Volume Control Switch
 - Dedicated On/Off Switch
 - Gunmetal Finish
 - Vocal Feedback
 - Low Battery Warning
 - Mute
 - Volume and Program Auto-Save
 - 5 Digitally Optimized Programs
 - 6 Metallic Door Colors
 - WaterShield protection
 - Size 13 Battery (new doors)
 - Size 312 Battery (old doors)
 - NRR 26
 - IPIL: 37.3/44/42
 - 2 Year Warranty

SEE PAGE 92 FOR ACCESSORIES -



SMARTPHONE HEADSET
 PN 64110



BLUETOOTH HEADSET
 PN 64109



Dx5 NECK LOOP
 PN 64113

IN THE BOX

- Your Custom Earmolds
- Dx Digital Modules (L&R)
- High-Impact Hard Case
- User Guide
- Cleaning Cloth
- Oto-Ease Lubricant
- Oto-Clens Cleaning Solution
- Battery Tool/Cleaning Loop
- Embroidered DefendEar Patch
- 6 - Size 13 Batteries



















Note: DigitalX is designed for people with normal hearing only.

WARNING: The purpose of hearing protection is to reduce the impact of the sound and to reduce the potential for hearing loss. Hearing protection devices, including DefendEar products, do not eliminate sound. The user is responsible for determining whether the noise reduction provided by DefendEar products will provide adequate hearing protection based on the user's unique medical history and exposure level. Your medical history, particularly any prior trauma to the ear or prior hearing loss, as well as the noise (decibel level and exposure time) the user will be exposed to should be considered prior to using DefendEar products. If you have any questions concerning the appropriate exposure levels or the level of hearing protection provided by DefendEar products, or if you have hearing loss or ear injuries, contact your local hearing professional before using.

Note: Second-bend, open-mouth, bite-block impressions are recommended for all DefendEar hearing protection devices.

Hearing Protection APPLICATIONS GUIDE

PRODUCT	INDUSTRIAL & CONSTRUCTION	SHOOTING	HUNTING	AIR TRAVEL	MUSIC	MOTORCYCLE 2 WAY
Max 40						
Convertible 47-40						
Concert 49						
DefendEar Digital 1						
DefendEar Digital 2						
DefendEar Digital Hunter						
DefendEar Passive 39						
TRU Custom						

Retubing GUIDELINES

Retubing an earpiece can sometimes present more challenges than expected. By using the tools and techniques described in this section, you will be able to remove old tubing and install new tubing on your patients' earpieces.

Removing Old Tubes - Two quick and easy ways to remove old tubing from earpieces.

- 1 If you have a microwave available in your office, bring a small glass container filled with water to a rolling boil. Our powder and liquid impression mixing bowl works well for this application (Westone part #20271). While the water is still actively boiling, drop the earpiece you wish to re-tube into the container. Once the water is cool enough to take the earpiece out, remove it from the still warm water and use a small pair of needle nose pliers to remove the old tube. A sharp tug should be all that is required to pull the tubing free from the earpiece.
- 2 The second way to remove tubing from an earpiece is to use a Westone hand reamer (Westone part #30451). Cut the tubing flush with the outside surface of the earpiece, insert the hand reamer and rotate in a clockwise direction. The drilling action will remove the old tubing from inside the sound bore. After you have removed the old earpiece tubing with either method, inspect the sound bore for any remaining debris. Any debris left in the sound bore can have an effect on the high-frequency response or completely block sound transmission.



Figure 1



Figure 2

Note: Never use powered drill bits or cutting burrs when trying to remove old tubing from an earpiece. Use of powered tools can quickly alter the inside diameter of the sound bore and render the earpiece unusable.

Tube-Through

Use a small amount of the appropriate cement on a pipe cleaner to clean the sound bore along the entire length (thin cement for acrylic earpieces, vinyl cement for vinyl earpieces).




Figure 1

From the canal tip, insert a tube puller into the earpiece (monofilament for acrylic earpieces, wire for soft earpiece materials). Place the quilled end of the tube through the loop of the tube puller and pull through the earpiece.




Figure 2

When the tube is all the way through, adjust the tube by hand until you have it at the proper angle.



Figure 3

When the tube is all the way through, adjust the tube by hand until you have it at the proper angle. Pull the tubing away from one side of the sound bore and apply a small amount of cement. Capillary action will cause the cement to flow around the entire tube.



Figure 4

When the glue is dry, cut the tube flush with the end of the canal. You are now ready to refit the instrument.

Half Tubing

After cleaning the sound bore, hold the new tube against the earpiece canal and gauge approximately where you want the tube to end in the sound bore. Clip the tubing straight across. Note: It is very important that the tubing be cut at 90 degrees to the tube sidewall.



Figure 1

Pour a small amount of tubing cement (thin cement for acrylic earpieces, vinyl cement for vinyl earpieces) into a small cup or bowl.

Note: The following steps must be performed quickly. Familiarize yourself with all of the steps before you begin:

Dip the end of the tube into the cement. Tap the tube against the sides of your cement container until only a small amount of cement remains in the tube end.



Figure 2

Place the tube over the sound bore. Pinch the tube and quickly insert it into the sound bore. Pinching the tube causes the cement inside the tube to coat the wall of the sound bore, acting as a lubricant during tube insertion and then securing the tube when it dries.




Figure 3

Make sure the tube is at the proper angle and allowed to dry. The earpiece is now ready to be refitted to the hearing instrument.

Tube-Locks & TRS Tubes

Cut the old tube flush with the outside of the earpiece. Insert a tube-lock tool into the base of the earpiece and push the remaining tube and tube-lock or TRS tube out of the earpiece through the canal tip. Make sure the sound bore is free of any obstruction.




Figure 1

Insert a wire tube puller into the sound bore from the canal end.




Figure 2

Thread the quilled end of the tube-lock or TRS tube through the tube puller and then pull the tube into the earpiece.



Figure 3

Pull the quilled end of the tube and gently work the tube-lock or TRS sleeve into the earpiece.



Figure 4

Continue manipulating the earpiece and tube until the tube-lock or TRS is approximately 3 mm below the outside surface of the earpiece. Once the tube is in the proper position, trim the excess tube from the canal end and refit the instrument.





Figure 5

FIT TUBES


FIT earpieces are modeled using sophisticated three-dimensional computer software and are “printed” via an additive manufacturing process. The acrylic resin used with FIT is slightly different from conventional acrylic. The resin’s composition requires the use of cyanoacrylate (Super Glue®), rather than standard thin tubing cement, to hold the vinyl tubes in place. Retubing therefore requires a slightly different procedure.




Keeping Tubing High and Dry
In areas of high humidity, water droplets can form inside the tubing in BTE fittings. When enough of these droplets form, they may bridge the tubing walls and effectively occlude the sound channel. To your client, it might sound like the instrument is cutting out. Venting the earpiece (if possible) and using a Dry-Tube will help eliminate the problem.

Removing FIT Tubes

Tubing removal is the same for both a through tube and half-tube earpiece. To remove the old tube, we recommend holding the earpiece stationary in one hand while pulling the tube and twisting it in a circular motion with the other hand or a pair of needle nose pliers. This movement will create tension at the tip of the tube and break the glue bond.



If this action is unsuccessful, you can use a small amount of debonder (Westone part #20262) to soften the super glue at the tip of the earpiece. Once the old tube is removed, rinse any excess debonder from the surface of the earpiece.



Warning: Do not use a high-speed rotary tool or hand-held reamer when removing the old tubing as this may damage the FIT earpiece.



FIT Retubing

If the earpiece is to be tubed through, insert the tube to the desired length and position. Once in place, bend the excess tubing at the tip of the earpiece away from the opening of the sound bore. This will create a small space between the new tube and the sound bore wall.

Then, using a toothpick or similar application device, apply a small amount of super glue at the gap between the tube and the sound bore. Once the glue has been applied, let go of the tubing to close the gap, trim off the excess tubing so it is flush with the end of the canal, and allow time for the glue to dry.

FIT Half-Tubing Instructions

In cases where canal size will not permit the tube to extend to the tip of the earpiece, a half-tube will be used. To retube a half-tubed earpiece, the new tube must first be cut to the proper length and inserted into the sound bore.



Once the tube is inserted, bend it away from the opening of the sound bore. This will create a small space between the new tube and the sound bore wall. Using a toothpick or similar application device, apply a small drop of super glue to the gap between the new tube and the sound bore opening. Once the glue has been properly applied, let go of the tubing to close the gap, and allow time for the glue to dry.

SLIM TUBES




Figure 1

Remove the old tube by pulling it free of the earpiece with a firm tug. Keep this tube handy as a reference when preparing the new tube for the earpiece. (Figure 1)

Retubing with the Slim-Tube requires modifications to the tubing provided by the hearing instrument manufacturer. (Please refer to Slim-Tube Modification Chart).




Figure 2

The first step is to remove most of the "retention tail" from the tubing. All three styles of Slim-Tubes require that about 1 to 2 mm of the "retention tail" be left on the tubing. The small stub that remains from trimming the "retention tail" is used to create friction inside the sound bore of the earpiece and help hold the tube in place. You must also remove the retention threads from the end of the slim tube. Refer to the old tube or the Slim-Tube Modification Chart as a guide for these modifications. (Figure 2).




Figure 3

Next, insert the modified Slim Tube into the earpiece to verify that the tube fitting is as close as possible to its final position. Reposition if needed. (Figure 3)

Once the tube is in position, fit the earpiece and hearing instrument to your patient and visually check to see if any final adjustments need to be made.

SLIM TUBE <i>Modification</i> CHART			
Style "A"	Style "A" (modified)	Style "A" (isometric)	Phonak CROS*
Style "B" (Round)	Style "B" (modified)	Style "B" (isometric)	Oticon miniFIT Corda 2*
Style "C" (Square)	Style "C" (modified)	Style "C" (isometric)	Oticon Corda miniFIT*

* Adjust for sound bore diameter

Ordering

The Golden Rule

We treat people the way we would like to be treated. If you have a problem, we'll help you solve it. If we make a mistake, we'll correct it. Our goal is 100% customer satisfaction.

Earpiece Orders

An in-house "rush service" is available for an additional fee. Standard shipping is via First Class mail unless another shipping method is chosen. We supply First Class mailing boxes and labels for your use, free of charge.

Note: *Our postage paid labels are for mailing earpiece orders only. Please do not use these labels unless you are sending an earpiece order using a standard impression box.*

If you need your order faster than First Class, expedited shipping services are available through FedEx®, UPS®, or Priority Mail®. These services are offered at the discounted rates available to Westone. If you send us impressions via FedEx, UPS, or Priority Mail, please specify on your order form the return shipping method for the finished product. If no method is specified, we will return the order to you via the same method we received it and charge your account accordingly.

Earpiece Remake Policy

We will remake an earpiece at no additional charge if it breaks within one year due to failure of materials or workmanship. In addition, we will remake an earpiece within 90 days of the invoice date under the following circumstances:

- An earpiece we sent out is lost in the mail
- We made a manufacturing error
- An earpiece allows feedback or does not fit properly. Please note that open style earpieces are not guaranteed against feedback.

No-charge remakes are not available under the following circumstances:

- Incorrect information is supplied on the order form
- Feedback or fit problems arise after 90 days
- The earpiece is lost
- Occlusion

If a style or material change is requested on a remake, a bench fee will be charged.

Note: *Shipping charges apply on all remakes.*

Supply Orders

Orders received during normal business hours will generally be shipped the following business day. Rush orders received before 3:00 p.m. Mountain Time will be shipped the same day. Rush orders received after 3:00 p.m. Mountain Time will be shipped the following business day.

Placing Your Order

- Shop online westone.com
- Order by phone 1-800-525-5071
- Fax# 719-527-3982

ORDERING



PHONE:

Local: 719.540.9333
Toll Free: 1.800.525.5071



FAX:

1.719.527.3982



ONLINE:

westone.com/store



EMAIL:

westone@westone.com

Shipping

Note: All shipments made to a residential address will be billed \$3.20 more than the corresponding commercial rates shown in the Supplies Shipping Information chart.

Shipping and Handling

Unless a specific shipping method is requested, we will ship all supplies via FedEx Ground or UPS Ground. All supplies are shipped F.O.B. Colorado Springs, CO. The postal service is used only when specifically requested or when no other means of delivery is available.

- Estimated delivery dates are not guaranteed. Weather delays and other unforeseen circumstances may impact delivery times.
- Some items require a signature for delivery.
- For shipping related questions, please contact customer service at 1-800-525-5071.

Freight Charges

Standard freight charges are based on the dollar value of your order. Some particularly heavy items may require additional freight charges.

Back-Ordered Items

Westone reserves the right to partially ship orders, please specify if order has to ship complete. Back-ordered items, if any, will be shipped as soon as they become available. Back-ordered items will not be billed to your account until they are shipped.

Special Requests

Need something that isn't in this catalog? Call us and let us know what you are looking for. Many of the earpiece styles and new supply items in this catalog are a direct result of requests from our customers.

Earpiece Shipping Information*	One Way	Round Trip
United States Postal Service	\$0	\$0
Priority Overnight (10:30 AM)	\$19.00	\$38.00
Standard Overnight (3:00 PM)	\$16.25	\$32.50
2-Day (4:30 PM)	\$11.00	\$22.00
3-Day Saver (4:30 PM)	\$9.50	\$19.00

*Rates are for the contiguous United States only and are subject to change without notice. All shipments made to a residential address will be billed \$3.20 more than corresponding commercial rates shown in the Supplies Shipping Info chart.

Rates Effective July 15, 2019

Supplies Shipping Information*	Ground Shipping and Handling
\$0 – \$100.00	\$12.00
\$100.01 – \$200.00	\$13.00
\$200.01 – \$300.00	\$15.00
\$300.01 – \$400.00	\$16.00
\$400.01 – \$500.00	\$17.00
\$500.01 – \$700.00	\$18.00
Over \$700.00	\$20.00

*Rates are for the contiguous United States only and are subject to change without notice. All shipments made to a residential address will be billed \$3.20 more than corresponding commercial rates shown in the Supplies Shipping Info chart.

Return Policy

All returned merchandise must be accompanied by a Return Merchandise Authorization (RMA) number. Please call us to request a RMA number. RMA numbers must be obtained within 30 days from the date of the invoice.

Defective Merchandise

- May be returned within 30 days of the invoice date for credit or exchange after obtaining a RMA number.
- If you received incorrect or damaged merchandise, we will provide you with a label for return shipping.
- All paperwork, including warranties and instructions, must be returned along with the original packaging. Place the manufacturer's box into a shipping box (plastic bags are not sufficient packaging). DO NOT put any stickers, labels, tape, or marks on the original manufacturer's packaging.
- A manufacturer's warranty is included with many products. Products that become defective after the 30-day return period, and within the manufacturer's warranty period, should be sent directly to the manufacturer with a copy of your receipt enclosed. See manufacturer's warranty for further instructions.

Non-Defective Merchandise

- May be returned within 30 days of the invoice date after obtaining RMA number.*
- May be subject to a restocking fee.
- Must be returned in new and resalable condition in the original packaging. Items damaged or showing signs of use are not returnable. Items missing will be subject to additional charges.
- All paperwork, including warranties and instructions, must be included along with the original packaging. Place the manufacturer's box into a shipping box. DO NOT put any stickers, labels, tape, or marks on the original manufacturer's packaging.
- Freight charges incurred due to returned merchandise are the responsibility of the customer. We recommend that you use FedEx or UPS with insurance for returns. Westone will not be held responsible for lost or damaged packages.

Returns Procedure:

- Call us at 1-800-525-5071 and request a Return Merchandise Authorization (RMA) number or visit our website at www.westone.com. We will need the following information:
 - Purchaser/Business Name, account number and invoice number
 - Product serial number, if applicable
 - Reason for return. If defective, describe specifically what is wrong.
- Merchandise should be packaged in such a way as to protect the product and its packaging.
- Clearly mark the RMA number on the outside of the shipping box.
- Include a copy of the original invoice in the package.
- All merchandise must be shipped prepaid and fully insured to:

Westone Laboratories, Inc.
 RMA# _____
 2260 Executive Circle
 Colorado Springs, CO 80906

Credit Terms and Policies

Credit will be evaluated after submission of a Westone credit application. Please call us at 1-800-525-5071 or email customerservice@westone.com to request a credit application. Westone assigns a credit limit to each customer. All balances reflected on the monthly statement are due within 30 days of the invoice date.

Accounts that are past due or that exceed the credit limit may be placed on hold, and credit terms may be changed to a prepaid basis or suspended. Past due accounts may also be placed with a collection agency.

Note: *Westone Laboratories, Inc. does not accept liability arising out of the use, the misuse or inability to use any accessory, tool, or product described herein, or any accessory, tool, or product subsequently offered for sale by Westone Laboratories, Inc. The user assumes all risk and liability whatsoever in connection with these tools and products.*