

Smart Pill Bottle User Manual

AdhereTech's User Information for FCC

For a copy of the User Manual specific to a particular use-case, please email support@adheretech.com.

The Smart Pill Bottle can produce helpful reminder alerts to help patients stay adherent. To start using the bottle, peel & remove the entire pull-tab pouch from the Bottle. The Bottle will then light up & chime to show that it is activated & ready to use.

With normal use, the battery in the Bottle lasts up to six months from the moment when the Bottle is first used until it needs to be recharged. Battery life will vary based on the Bottle's alert settings and patient use. It takes approximately eight hours of charging for the Bottle to be fully recharged.

FCC INFORMATION FOR USER

This equipment has been tested and found to comply with the limits for a class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help

AdhereTech
11 Broadway, Suite 457
New York, NY 10004
1-800-381-9384
www.adheretech.com