

Smart Pill Bottle User Manual

AdhereTech's User Information for FCC

For a copy of the User Manual specific to a particular use-case, please email support@adheretech.com.

The Smart Pill Bottle can produce helpful reminder alerts to help patients stay adherent. To start using the bottle, peel & remove the entire pull-tab pouch from the Bottle. The Bottle will then light up & chime to show that it is activated & ready to use.

With normal use, the battery in the Bottle lasts up to six months from the moment when the Bottle is first used until it needs to be recharged. Battery life will vary based on the Bottle's alert settings and patient use. It takes approximately eight hours of charging for the Bottle to be fully recharged.

FCC COMPLIANCE STATEMENT

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions:

1) This device may not cause harmful interference, and 2) This device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

AdhereTech
11 Broadway, Suite 1115
New York, NY 10004
(646) 417-8143
www.adheretech.com