

Good reasons to take care of your health

UA-702 Manual Inflation



One of the most technologically advanced yet affordable blood pressure monitors available in the marketplace today.

The Manual Inflation monitor combines ease of use and affordability. This monitor combines the software technology of our automatic inflation monitors, but leaves cuff inflation control up to the user. The Manual Inflation monitor can assist in the management of hypertension, improve patient compliance to treatment of high blood pressure, and be used as a tool in a preventative health management program.

- Large, easy to read 16mm digital display
- Audible tone indicates proper inflation level
- Available in small, medium or large cuff sizes (cuffs can also be purchased separately)
- Large-print instruction guide with illustrations and quick start card in English, French and Spanish
- Four week blood pressure logbook included

- Professional accuracy
- Trilingual instruction guide
- Easy to use
- Easy to read digital display
- Fast measurement
- · Latex free







For questions about *LifeSource* monitors, call our toll free *LifeSource* HealthLine at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

UA-702 Manual Inflation

To obtain an accurate reading, it is important to use the correct cuff size. A cuff that is too small will result in a higher than expected reading and a cuff that is too large will result in a lower than expected reading. The UA-702 Manual Inflation unit comes in three different cuff sizes to accommodate your customer's needs.

UNITED STATES	PART NUMBER	UPC NUMBER	CUFF SIZE
MASS MARKET	UA-702V UA-702VL UA-702VS	0-93764-24702-6 0-93764-36702-1 0-93764-81702-1	Medium 7.5" - 12.2" (19-31cm) Large 11.8" - 17.7" (30-45cm) Small 5.1" - 7.9" (13-20cm)
HOME HEALTH CARE	UA-702H UA-702HL UA-702HS	0-93764-24702-6 0-93764-36702-1 0-93764-81702-1	Medium 7.5" - 12.2" (19-31cm) Large 11.8" - 17.7" (30-45cm) Small 5.1" - 7.9" (13-20cm)
CANADA	UA-702CN UA-702CNL UA-702CNS	0-93764-35702-2 0-93764-12702-1 0-93764-19702-4	Medium 7.5" - 12.2" (19-31cm) Large 11.8" - 17.7" (30-45cm) Small 5.1" - 7.9" (13-20cm)

Each unit comes with a trilingual instruction manual, quick start card, four week blood pressure logbook and warranty card. Home Health Care models also include a carrying case.

Accessories

PART #	DESCRIPTION	PART #	DESCRIPTION
11-001LF	Inflator Bulb	UA-279	Small Blood Pressure Monitor Cuff: 5.1" - 7.9" (13-20cm)
	(Latex Free)	UA-280	Medium Blood Pressure Monitor Cuff: 7.5" - 12.2" (19-31cm)
UA4-2052-B	Exhaust Valve	UA-281	Large Blood Pressure Monitor Cuff: 11.8" - 17.7" (30-45cm)
U4:2474-C	"T" Air Connector		·
U4:2452-D	"L" Air Connector		
UA-700CH	Cuff Hose		

Technical Specifications

Display Digital: 16mm character height

Pressure and pulse displayed alternately

Pressurization Manual, by inflator bulb Deflation Hand-held push button

Type Oscillometric

Accuracy Pressure: ±3mmHg or 2%, whichever is greater/Pulse: ±5%
Pressure Range 20mmHg to 280mmHg/Pulse: 40 to 200 pulses/minute
Power Source Four type AA alkaline batteries (1.5 volt). Batteries not included.

Battery Life Approx. six months with three daily measurements

 $18^{1/2} \times 11 \times 8^{1/2}$

Operating Environ. 50°F to 104°F (10°C to 40°C)
Storage Environ. 14°F to 140°F (-10°C to 60°C)

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	L x W x H (inches)	WT (lbs.)	Individual Item Cube (cubic feet)
UA-702V/CN Models	$5^{1/}_{16}$ x $3^{1/}_{2}$ x $7^{11/}_{16}$	1.05	0.08
UA-702H Models	$5^{1}/_{16} \times 3^{1}/_{2} \times 7^{11}/_{16}$	1.20	0.08
Unit Case (10 units)	LxWxH (inches)	WT (lbs.)	Case Cube (cubic feet)
UA-702V/CN Models	$18^{1/_2} \times 11 \times 8^{1/_2}$	11.70	1.00

13.45

Ask your A&D Sales Representative about point-of-purchase floor and counter displays available for this blood pressure monitor unit.



UA-702H Models

1555 McCandless Drive Milpitas, CA 95035 www.LifeSourceOnline.com







Ask your A&D Sales Representative about other *LifeSource* products.

- Blood Pressure Monitors
- Heart Rate Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about *LifeSource* monitors, call our toll-free *LifeSource* HealthLine at 1-888-726-9966 or visit our web site at www.LifeSourceOnline.com