

LAB CONTROLS

\$0.00 0 items

Home About Support Quote Request News

Shop

Home Centrifuges Eppendorf 5417C Centrifuge

SALE!



# Eppendorf 5417C Centrifuge



IN STOCK



90 90-DAY WARRANTY



LOW PRICE GUARANTEE

Item Number: 5417C Centrifuge

Manufacturer: Eppendorf

Product Category: Centrifuges

Brand Name: Eppendorf

Description: Eppendorf 5417C Centrifuge

*Have a quote from somewhere else?* [Click Here](#)

1 in stock

---

SKU: 776097b6a75c

Category: [Centrifuges](#)

---

Description

---

[Additional information](#)

---

## Description

Item Number: 5417C Centrifuge

Manufacturer: Eppendorf

Product Category: Centrifuges

Brand Name: Eppendorf

Product Overview: Eppendorf 5417C Centrifuge

---

## Overview of Eppendorf 5417C Centrifuge

---

Centrifuge 5417 C is a non-refrigerated bench-top centrifuge. It is used in research and routine laboratories in the biosciences and in the fields of medicine and industry. Centrifuge 5417 C Eppendorf Micro Test Tubes can be centrifuged at 20,800 rcf (14,000 rpm). It is possible to attain an rcf of up to 25,000 using a special 24-position rotor. For special applications, a swing-

bucket rotor is available with 8 positions for centrifuging Eppendorf Micro Test Tubes at a maximum rotational speed of 10,500 rpm (11,700 rcf). PCR tubes can be centrifuged in a rotor which can accommodate six 8-tube strips. Eppendorf 5417C Centrifuge is available from Laboratory Controls LLC. This product comes with a 90 day warranty.

## Specifications

---

Power supply	230 V / 50–60 Hz, 120 V / 50–60 Hz, 100 V / 50–60 Hz
Power requirement	400 W
Max. rotational speed	16,400 rpm
Max. centrifugal force	25,000 rcf
Max. load	30 x 2.0 ml Safe-Lock Microcentrifuge Tubes
Max. density of material to be centrifuged	1.2 g/ml
Permitted ambient temperature during operation	2 °C – 35 °C
Permitted maximum relative air humidity	75 %
Max. kinetic energy	6,100 Nm

Dimensions (H x W x D)	228 x 320 x 430 mm
Weight	17 kg
Degree of pollution	2
Overvoltage category	II

## Questions

---

Question: Is this centrifuge refrigerated? Answer: No.

Question: Does the 5417C offer continuous operation?

Answer: Yes, this can be achieved by using the arrow keys. "oo" indicated continuous operation.

## Resources

---

[Eppendorf 5417C Manual](#)

## Warranty Information

---

Laboratory Controls offers a best in the business full 90 day warranty on all refurbished items. This warranty covers your part against any material and or craftsmanship issues within one full year from purchase. Laboratory Controls will repair or replace your part at no charge. Please contact us for further information if you are in need of a repair.

## Related products



Beckman GS-6 Centrifuge

SALE!



Eppendorf 5702 Centrifuge

SALE!

Add to cart



Heraeus Labofuge 400 Centrifuge

SALE!

Contact Us

Site Info

Phone Hours

Social Media

Customer Phone  
Number  
(984) 444-9060

U.S. Customer  
Email:  
sales@laboratorycon  
trols.com  
laboratorycontrolssal  
es@gmail.com

Mail:  
6805 Dwight  
Rowland Rd - Unit  
611  
Willow Springs, NC  
27592

Order Tracking  
Policies  
Site Map  
Careers  
Privacy Policy  
Terms & Conditions

Monday - Thursday  
8am - 8:00pm (EST)  
Friday 8am - 6pm  
(EST)  
Saturday - Sunday  
8am - 5:30pm (EST)

Laboratory Controls  
Facebook  
PageLaboratory  
Controls Twitter



---

© 2017

Storefront designed by WooCommerce.