



Invacare® CareGuard™ 8000 Alternating Pressure Mattress



Yes, you can.™




Yes, you can.™

The Invacare® CareGuard™ 8000 is a full-depth (8") powered alternating pressure air mattress replacement system which reduces interface pressure for the treatment of Stage I to Stage IV pressure ulcers. The CareGuard 8000 reduces interface pressures below capillary closure through the design of 20 individual 8-inch air cell bladders that alternate every four minutes. The power unit control combines simplified operations with effective results.



Therapeutic and Portable



### User-Friendly Controls


- Adjustable pressure dial for setting patient comfort range
- Lighted ON/OFF switch for ease of use in low light
- Power unit is lightweight and compact

### Advanced Air Hose Technology

- Quick-disconnect connectors are secure yet easy to use for quick set-up or deflation
- Mattress hoses can be cross-connected to maintain inflation during transportation

### FEATURES


- Portable, lightweight mattress replacement system fits most standard medical beds
- 20 individual 8" nylon air cells reduce interface pressure below capillary closure while providing low interface pressure over the entire alternating cycle
- Mattress cells alternate every 4 minutes to reduce interface pressures
- Waterproof, anti-microbial nylon therapy cover reduces friction and shear
- Durable vinyl bottom protects the mattress cell from any abrasive surface
- CPR flag located on the mattress for quick, 20-second deflation
- 2" transport mat for patient transport or power outages
- Designed to meet EO277



### Alternating Mode

In alternating mode, the Invacare APM's 20 adjacent air cells alternately inflate and deflate at 8-minute interval, which redistributes the pressure against the skin to promote capillary circulation.

### Pressure Relief



Average pressure per point for 10 subjects

( ) Standard deviation

\* 90° side lying position

Testing completed in an independent test lab

### Specifications

Model no. CG8000	
Electrical power	115vac/60Hz
Dimensions	
Mattress	80" L x 35" W x 8" H
Power unit	10" W x 5" H x 4" D
Weight	
Mattress	16 lb.
Power unit	4 lb.
Weight capacity	250 lb.
Material	
Mattress cells	70 denier nylon with a polyurethane coating
Therapy pad	70 denier nylon top surface with a polyester backing for comfort
Safety code approval	UL544, CA106, CA117 (cover)
Limited warranty	
Power unit	1 year
Mattress	6 months

© 2002

Invacare Corporation

Form No. 00-305

Rev. 05/02



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."