

TriPlus 300 Headspace: Spare Parts Guide

This guide contains part numbers for all of the replaceable components in the TriPlus 300 Headspace. Throughout the TriPlus 300 Headspace documentation, any component with a part number can be ordered from us. Components without a part number are not available.

IMPORTANT Only components with a part number are available for purchase. The part you are looking for may be shown in an illustration, but if it does not list a part number, it is not available.

You can purchase the following column components for the TriPlus 300 Headspace. Be sure to reference the component's part number when placing an order with your local Sales/Service Representative. See [Table 1](#) for a comprehensive list of available component.

Table 1. Parts That Can Be Reordered (Sheet 1 of 3)

Component	Thermo Scientific Parts Number
TriPlus 300 Headspace — Total Flow Assembly	77010-3001
TriPlus 300 Headspace — Power Supply Assembly	77010-3002
TriPlus 300 Headspace — Oven Assembly	77010-3003
TriPlus 300 Headspace — Transfer Line 0.80 m	77010-3004
TriPlus 300 Headspace — Display	77010-3005
TriPlus 300 Headspace — Tee Union	77010-3006
TriPlus 300 Headspace — Carrier Gas Tube	77010-3007
TriPlus 300 Headspace — Pressure Sensor 60 psi “T” Tube	77010-3008
TriPlus 300 Headspace — Pressure Regulator “T” Tube	77010-3009
TriPlus 300 Headspace — Vent Tube	77010-3010
TriPlus 300 Headspace — Carrier “T” Tube	77010-3011
TriPlus 300 Headspace — Pressure Sensor 150 psi “T” Tube	77010-3012
TriPlus 300 Headspace — Zero Volume Union	77010-3031
TriPlus 300 Headspace — Vial Needle	77010-3032
TriPlus 300 Headspace — Transfer Line Needle	77010-3033
TriPlus 300 Headspace — Sample Loop 1 mL	77010-3034
TriPlus 300 Headspace — Master Board	77010-3040
TriPlus 300 Headspace — Stepper Board	77010-3041
TriPlus 300 Headspace — Support Board	77010-3044
TriPlus 300 Headspace — Carousel Board	77010-3045
TriPlus 300 Headspace — Lifter Board	77010-3046
TriPlus 300 Headspace — Valve Board	77010-3047

Table 1. Parts That Can Be Reordered (Sheet 2 of 3)

Component	Thermo Scientific Parts Number
TriPlus 300 Headspace — Temperature Control Board	77010-3048
TriPlus 300 Headspace — Pressure Control Board	77010-3049
TriPlus 300 Headspace — P2 Control Board	77010-3050
TriPlus 300 Headspace — Motor Board	77010-3051
TriPlus 300 Headspace — Power Supply Board	77010-3052
TriPlus 300 Headspace — Remote Board	77010-3053
TriPlus 300 Headspace — Main Board	77010-3054
TriPlus 300 Headspace — Solenoid Valve	77010-3055
TriPlus 300 Headspace — Pressure Sensor 150 psi	77010-3056
TriPlus 300 Headspace — Pressure Sensor 25 psi	77010-3057
TriPlus 300 Headspace — Fan	77010-3060
TriPlus 300 Headspace — Temperature Sensor	77010-3062
TriPlus 300 Headspace — Power Filter	77010-3064
TriPlus 300 Headspace — Transformer	77010-3066
TriPlus 300 Headspace — Power Supply 24 V	77010-3067
TriPlus 300 Headspace — Power Supply 12 V	77010-3068
TriPlus 300 Headspace — Hands Cable	77010-3071
TriPlus 300 Headspace — LEDs Cable	77010-3074
TriPlus 300 Headspace — Vial Movement Assembly	77010-3115
TriPlus 300 Headspace — Vial Magnet Assembly	77010-3117
TriPlus 300 Headspace — Tray 120	77010-3119
TriPlus 300 Headspace — Six Ports Valve	77010-3137
TriPlus 300 Headspace — Manifold Assembly	77010-3139
TriPlus 300 Headspace — Lifter Assembly	77010-3145
TriPlus 300 Headspace — Remote Assembly	77010-3146
TriPlus 300 Headspace — Vial Pressure Sensor Assembly	77010-3147
TriPlus 300 Headspace — Vent Valve	77010-3148
TriPlus 300 Headspace — Needle Tube	77010-3149
TriPlus 300 Headspace — Manifold Support Assembly	77010-3150
TriPlus 300 Headspace — Vent Valve Support	77010-3151
TriPlus 300 Headspace — MV Valve	77010-3160
TriPlus 300 Headspace — Heater 50 W	77010-3161
TriPlus 300 Headspace — Heater 35 W	77010-3162
TriPlus 300 Headspace — Front Panel	77010-3166
TriPlus 300 Headspace — Right Panel	77010-3167

Table 1. Parts That Can Be Reordered (Sheet 3 of 3)

Component	Thermo Scientific Parts Number
TriPlus 300 Headspace — Left Panel	77010-3168
TriPlus 300 Headspace — Standard Outfit	77010-3178
TriPlus 300 Headspace — Sample Path Kit	77010-3179
TriPlus 300 Headspace — Preventive Maintenance Kit	77010-3180
TriPlus 300 Headspace — Firmware Program	77010-3185
TriPlus 300 Headspace — Vial Flow Check	77010-3186
TriPlus 300 Headspace — M6-M5 Kit	77010-3187
TriPlus 300 Headspace — Oven Alignment Tool	77010-3188
TriPlus 300 Headspace — Grease	77010-3189
TriPlus 300 Headspace — Cooling Device	77010-3190
TriPlus 300 Headspace — Barcode Assembly	77010-3242
TriPlus 300 Headspace — Nut 5M; (Set of 20)	77010-3248
TriPlus 300 Headspace — Nut 6MB; (Set of 20)	77010-3249
TriPlus 300 Headspace — Seal 5M; 1.6 mm; Stainless Steel; (Set of 50)	77010-3250
TriPlus 300 Headspace — Seal 5M; 1.6 mm; Aluminium; (Set of 50)	77010-3251
TriPlus 300 Headspace — Seal 6MB; 1.6 mm; Stainless Steel; (Set of 20)	77010-3252
TriPlus 300 Headspace — Seal 6MB; 1.6 mm; Aluminium; (Set of 50)	77010-3253
TriPlus 300 Headspace — Seal 5M; 1 mm; Stainless Steel; (Set of 10)	77010-3254
TriPlus 300 Headspace — Seal 5M; 0.8 mm; Stainless Steel; (Set of 10)	77010-3255
TriPlus 300 Headspace — Reducer 1/8-inch; 5M	77010-3256
TriPlus 300 Headspace — BulkHead 6MB	77010-3257
TriPlus 300 Headspace — T 5M	77010-3258
TriPlus 300 Headspace — Tube 1.6 x 1mm; 1 m; Stainless Steel	77010-3259
TriPlus 300 Headspace — Tube 1 x 0.5 mm; 5 m; Stainless Steel	77010-3260