Focus on Retainers		Accutech Orthodontic Lab Inc accutechortho.com (800) 734-7855	Align Technology Inc aligntechinstitute.com (408) 470-1000	Allesee Orthodontic Appliances (AOA Labs (800) 263-5221	ClearCorrect clearcorrect. (888) 331-3323
Product Name		QuickBond BLR (Bonded Lingual Retainer)	Vivera® Retainers	Retainers	ClearCorrect Retainer
Туре	of retainer	Bonded	Clear plastic	Active/Spring, bonded, clear plastic, Hawley	Clear plastic
Mate	rial	Composite, wire	Plastic	Acrylic, composite, plastic, wire	Plastic
Turn	around time	7 business days	Varies based on type of submission	3 to 5 days	Average lab turnaround is 6 days
Initia	al cost	\$44.35	Varies based on single or subscription	Costs vary per design.	\$95
Repa rema	ir or ke cost	Zero	Varies	Varies depending on warranty.	\$95
How long is it designed to last?		Years	3 months	6 months to indefinitely	3 months when used as recommended. Can vary based on doctor-prescribed use and patient wear times.
Leng	th of warranty	90 days	N/A	90 days	3-month workmanship warranty
	Aesthetic wire			✓	
su	Clear plastic		✓	1	✓
: optio	Labeling				
Aesthetic options	Lingual wire placement	✓		✓	
Ae	Multiple colors			1	
	Personalized images			1	
Special features		Price includes bonding trays and shipping.	Vivera* retainers are manufactured using the advanced fabrication technology behind Invisalign system aligners, resulting in a precise, smooth, and comfortable fit. They can be used with Invisalign system patients or any other orthodontic patients.	All retainer options can be fabricated from digital impressions/scans provided by your office directly to AOA. This technology may allow for a more accurate fit due to the elimination of impression-taking variables. All of our acrylic designs can be customized by our acrylic artists using AOA's wide variety of colors, custom art, and decals.	ClearCorrect retainers are made with Zendura*. Dental models (the production molds) are included at no additional charge so doctors can easily make replacements in-house if the retainers are lost or wear out.

		ln ti	e e			
Focus on Retainers		DENTSPLY GAC Intl gacintl.com (800) 645-5530	DENTSPLY Raintree Essix.com (800) 883-8733	DynaFlex dynaflex (800) 489-4020	Forestadent USA forestadentusa.com (800) 721-4940	
Prod	uct Name	Hawley Retainer	MTM® Clear•Aligner	Digital Models	Track A	
Туре	of retainer	Hawley	Clear plastic	Active/spring, bonded, clear plastic, Hawley	Clear plastic	
Mate	rial	Acrylic, wire	Plastic	Acrylic, plastic, wire	Plastic	
Turn	around time	5 business days	N/A (made in-office)	3 days in laboratory per request	It takes minutes to make.	
Initia	al cost	N/A	N/A	\$48.42 plus \$8.75 per arch for digital models	Call for pricing.	
Repa rema	ir or ke cost	N/A	Less than \$10	6-month warranty, no charge	Minimal	
How desig	long is it ned to last?	Duration of treatment	18 to 24 months (varies based on selected plastic)	1 year	Depends on the care provided by the patient.	
Leng	th of warranty	90 days	N/A	6 months	N/A	
	Aesthetic wire					
tions	Clear plastic	✓	✓	✓	✓	
optic	Labeling	✓				
Aesthetic opt	Lingual wire placement					
Aes	Multiple colors	✓				
	Personalized images	✓				
Special features		Variations in design, acrylic color, and extras are also available and can be specified along with any special instructions for each customized appliance.	MTM* Clear•Aligner In-Office utilizes Essix* Brand plastics and the Hilliard Thermoplier* Pliers to fabricate clear plastic aligners that move teeth with integrated "force points." This technique can be used as a stand-alone aesthetic treatment, for mid-course correction or refinement of orthodontic cases, and in preparation for other cosmetic or therapeutic procedures.	DynaFlex is able to fabricate any traditional orthodontic appliance from digitally scanned orthodontic impressions.	Track A has an integrated separating foil, which reduces preparation time and supply inventory.	

Focus on Retainers		G&H Wire Co ghwire.com (800) 526-1026	Gold'n Braces Inc goldnbraces.com (727) 785-5491	Great Lakes Orthodontics Ltd greatlakesortho.com (716) 871-1161	Johns Dental johnsdental.com (800) 457-0504
Prod	uct Name	everStick® Ortho	ZMRW	Hawley Retainer	Retention Appliances
		everStick*ORTHO			Watasdentalcon
Туре	of retainer	Bonded	Bonded	Hawley	Active/spring, bonded, clear plastic, Hawley
Mate	rial	Composite	Wire	Acrylic, wire	Acrylic, plastic
Turn	around time	N/A	N/A	5 days	Ask for "3-day Hawleys" on your prescription
Initia	l cost	N/A	N/A	\$59.65	\$75.50 and up
Repair or remake cost		N/A	N/A	Acrylic: no charge; wire after less than 90 days: no charge	With an initial cost of \$25, remakes are 50% for up to 2 years.
How long is it designed to last?		Indefinitely	N/A	Many years if properly taken care of	Contact company
Length of warranty		N/A	N/A	Acrylic: lifetime; wire: 90 days	Contact company
	Aesthetic wire	N/A	✓	✓	✓
suc	Clear plastic	N/A		1	
c optio	Labeling	N/A		1	✓
Aesthetic options	Lingual wire placement	N/A	✓		
Ae	Multiple colors	N/A		✓	1
	Personalized images	N/A		✓	1
Special features		everStick® Ortho incorporates 1,000 glass fiber strands with a polymer/resin gel matrix that holds the individual glass fibers in a bundle, making it malleable and tacky. Completely metal free and aesthetic, it adapts directly to the teeth. No models or second visits are needed. With a 0.75-mm diameter, it flattens significantly during adaptation.	ZMRW is a 24-carat gold- plated retainer wire	There are many different patterns and colors to choose from.	The Hawley retainer can be easily adapted with the addition of simple features such as bite planes, anterior springs, QCM bows, habit options, and occlusal coverage.

Focus on Retainers		NorthStar Orthodontics Inc northstarorthodontics.coa (800) 346-0011	Ortho Organizers Inc (800) 547-2000	PAR Orthodontic Lab parortholab.com (949) 472-4788	PREAT Corp preat.com (800) 232-7732
Prod	luct Name	Hawley Bow Lasered to "C" Clasps	Hawley Spring Retainer	Hawley Retainer	eFiber
Туре	of retainer	Hawley	Active/spring	Hawley	Bonded
Mate	erial	Acrylic, wire	Plastic, wire	Acrylic, wire	N/A
Turn	around time	4 days in lab	5 days	3 days in lab	None
Initi	al cost	\$50.50	N/A	\$36.50 (base price)	\$20
Repa rema	ir or ike cost	Variable cost after warranty term	N/A	No charge for the first 90 days for manufacture defect	N/A
How desig	long is it gned to last?	3+ years	Depends on patient care	1 year+	N/A
Leng	th of warranty	60 days for repair or replacement at no charge	N/A	90 days	N/A
	Aesthetic wire	✓			N/A
tions	Clear plastic		✓		N/A
: optio	Labeling	1	✓	1	N/A
Aesthetic op	Lingual wire placement				N/A
Ae	Multiple colors	✓	✓	✓	N/A
	Personalized images	✓	✓	√	N/A
Special features		This style of retainer gives you the benefit of a wraparound in that it does not cross the occlusion but adds to the flexibility and security of retention. Clasp design is your choice, as are laser weld or solder options.	Free retainer case with every retainer.	The Hawley is primarily for retention, but a variety of springs can be added to correct individual teeth. Shown with soldered "c" and color pattern option.	eFiber makes one-visit invisible bonded ligual retainers. The translucent eFiber is impregnated with both PMMA and Bis-GMA and bonds directly to etched enamel. Easy to adjust and polish to a smooth surface. Eliminates metal, food traps, and patient discomfort—all chairside.

cus on tainers	Ribbond Inc ribbond.com (800) 624-4554	Space Maintainers Lab, part of the Appliar appliancetherapy.com (800) 423-3270	Specialty Appliances Orthodontic Laborato (800) 522-4636	Tortorici Orthodontic Lab Inc (631) 642-1900	
ıct Name	Ribbond Bondable Reinforcement Ribbon	COR Final Retainer	Fixed Lingual Retainer	Clear Fusion Hybrid Retainer	
				W The second	
of retainer	Bonded	Active/spring, clear plastic, Hawley	Bonded	Clear plastic, Hawley	
rial	Composite, plastic	Plastic, wire	Composite, wire	Acrylic, plastic	
around time	Made in office	5 days	3 days in the lab	3 working days	
l cost	\$258 for 20 3-3 retainers	\$85	\$38 to \$52	\$40 to \$65	
r or ke cost	Contact company	Contact company	No charge for 1 year	Guaranteed to fit or remake at no charge	
long is it ned to last?	The bond of the composite to the teeth dictates the mode of failure.	Forever with proper care	Indefinitely	Variable depending on patient care	
h of warranty	Contact company	90 days	1 year	3 months (negligence not included)	
Aesthetic wire	Contact company		✓		
Clear plastic	Contact company	✓		✓	
Labeling	Contact company			✓	
Lingual wire placement	Contact company		✓		
Multiple colors	Contact company			✓	
Personalized images	Contact company				
al features	Fiber reinforced composite resin retainer.	The Adjustable COR (Clear Orthodontic Retention) Retainer is easy to use, comfortable to wear, and constructed of unbreakable Biocryl. It delivers stability and permits natural bite function—especially following extensive orthodontic treatment— because the biting surfaces and edges are not covered. Designed for final retention, the COR can also be used for minor corrections.	The appliance can be constructed with scalloped wire to prevent flossing issues.	This wireless patented hybrid retainer combines the aesthetics of a clear invisible retainer and the rigidity of a Hawley appliance.	
	of retainer rial round time cost r or se cost dong is it ned to last? h of warranty Aesthetic wire Clear plastic Labeling Lingual wire placement Multiple colors Personalized images	Ribbond Bondable Reinforcement Ribbon of retainer Bonded cial Composite, plastic Made in office l cost \$258 for 20 3-3 retainers r or ce cost Contact company long is it to the teeth dictates the mode of failure. th of warranty Contact company Aesthetic wire Contact company Clear plastic Labeling Contact company Lingual wire placement Multiple colors Contact company Fiber reinforced composite	Ribbond Bondable Reinforcement Ribbon Active/spring, clear plastic, Hawley Final Composite, plastic Plastic, wire Bonded Plastic, wire Fround time Made in office Plastic, wire Contact S258 for 20 3-3 retainers Plastic From Contact company Contact company Contact company Contact company Contact company Plastic Ho dwarranty Contact company 90 days Acesthetic wire Contact company Clear plastic Contact company Lingual wire placement Contact company Multiple colors Contact company Fiber reinforced composite resin retainer. Fiber reinforced composite resin retainer. Fiber reinforced composite resin retainer. The Adjustable COR (Clear Orthodontic Retention) Retainer is easy to use, comfortable to wear, and constructed of unbreakable Biocryl. It delivers stability and permits natural bite function—especially following extensive orthodontic treatment—because the biting surfaces and edges are not covered. Designed for final retention, the COR can also be used for	Ribbond Bondable Reinforcement Ribbon Ribbond Bondable Reinforcement Ribbon Active/spring, clear plastic. Hawley Bonded Active/spring, clear plastic. Hawley Composite, wire Composite, wire Composite, wire Says 3 days in the lab Lost \$258 for 20 3-3 retainers \$85 \$38 to \$52 Contact company Contact company Contact company No charge for 1 year Indefinitely The bond of the composite to the teeth dictates the mode of failure. Active/spring, clear plastic. Bonded Active/spring, clear plastic. Composite, wire Composite, wire Composite, wire Contact company Contact company Forever with proper care Indefinitely Indefinitely Active/spring, clear plastic. Contact company Vo charge for 1 year Indefinitely Indefinitely V Contact Company V Active/spring, clear plastic, Bonded Active/spring, clear plastic, Forever with proper Contact company Vo charge for 1 year Indefinitely V Contact Company V Active/spring, clear plastic, Bonded Bonded Bonded Active/spring, clear plastic, Contact company No charge for 1 year Indefinitely V Contact Company V Contact company V Labeling Contact company Contact company V Contact company The Adjustable COR (Clear plastic) Contact company The Adjustable COR (Clear plastic) Contact company The Adjustable COR (Clear plastic) The Adjustable COR (Clear plastic	

SCAN ME!



Scan the QR code above with your smartphone to connect to our Online Buyer's Guide.

ORTHODONTIC PRODUCTS



Ask your questions about Lingual Orthodontics.

Orthodontic Products acknowledges Forestadent for its generous sponsorship of the Lingual Orthodontics expert.



Dr. Marco A. Pinto Specialist in Orthodontics and Dentofacial Orthopedics

Looking for expert advice?

OrthodonticProductsOnline.com

Focus on Retainers		TP Orthodontics In tportho.com (800) 348-8856	Tru-Tain Orthodont & Dental Supply tru-tain.com (800) 843-0904	
Product Name		FLEX-Palate®	Tru-Tain DX Retainer Plastic	
Type	of retainer	Hawley	Clear plastic	
Mate	rial	Composite	Plastic	
Turna	around time	7 business days	Immediate	
Initia	l cost	Contact company	\$1	
Repai rema	ir or ke cost	Guaranteed against breakage for life under normal use.	\$1	
How design	long is it ned to last?	Lifetime of prescribed wear time	8 to 12 months or longer	
Lengt	th of warranty	5 years on materials and craftsmanship	N/A	
	Aesthetic wire	✓		
ns	Clear plastic		✓	
coptions	Labeling	✓		
Aesthetic	Lingual wire placement			
Ae	Multiple colors			
	Personalized images			
Special features		It is 98% stronger than traditional acrylic. The modified composite material bends without breaking and is constantly self-relining for greater patient comfort.	It can be made in-office immediately after debonding orthodontic brackets.	

ntic